INTENDED USE

Calprest NG is a quantitative ELISA for detecting concentration of fecal calprotectin that is intended to aid in the diagnosis of Inflammatory Bowel Diseases (IBD), specifically Crohn's disease and ulcerative colitis, and to differentiate IBD from Irritable Bowel Syndrome (IBS) in conjunction with other clinical and laboratory findings.

SUMMARY AND EXPLANATION

Calprotectin is a 36 kDa heterodimeric complex belonging to the S100 family of proteins that derives from polymorphonuclear cells, monocytes and epithelial cells. It has direct antimicrobial effects and it plays a role within the innate immune response¹⁻². Calprotectin is found in several biological samples including plasma, serum, synovial fluid and stool in proportion to the existing level of inflammation³. The use of calprotectin as a biological marker has been studied in inflammatory gastrointestinal disorders (IBD), as well as in other pathologies like rheumatoid arthritis and colorectal cancer⁴⁻¹².

TEST PRINCIPLE

Calprest NG is an enzyme linked immuno-sorbent assay system with colorimetric detection based on the use of polyclonal and monoclonal antibodies against calprotectin. Calprotectin present in the diluted sample is bound by the antibody adsorbed to the surface of the plastic well. The enzyme-conjugate antibody binds to the captured antigen and subsequently the enzyme catalyses the conversion of the substrate to a colored product. The intensity of the color is proportional to the amount of conjugate bound, and thus to the amount of captured calprotectin. Concentration of calprotectin in the samples is calculated using the provided Calibrators.

MATERIALS PROVIDED WITH THE KIT (QUANTITY SUFFICIENT FOR 96 TESTS)

Antibody coated plate	12 x 8 wells
Enzyme-conjugate antibody (IgG, mouse)	1 x 15 ml
Substrate	1 x 15 ml
Stop solution	1 x 15 ml
Washing solution (20x)	1 x 50 ml
Diluent (10x)	1 x 20 ml
Extraction solution (2.5x)	2 x 50 ml
Calibrators	6 x 1.5 ml
Control 1 (low)	1 x 1.5 ml
Control 2 (high)	1 x 1.5 ml

COMPOSITION OF SUPPLIED REAGENTS/MATERIALS

1. Antibody coated plate: 12x8 wells coated with antibody against calprotectin. Plastic sealed bag containing a desiccant.
2. Enzyme conjugate antibody (IgG): 1 vial containing 15 ml Horse-radish peroxidase-labeled mouse anti-human calprotectin IgG antibodies in a buffer solution with Proclin 300 as a preservative. Ready to use.
3. Substrate: 1 vial containing 15 ml substrate reagent (TMB). Ready to use.
4. Stop Solution: 1 vial containing 15 ml H₂SO₄ (0.5M). Ready to use.
5. Washing solution (20x): 1 vial containing 50 ml concentrated washing solution. To be diluted with distilled water.
6. Diluent (10x): 1 vial containing 20 ml concentrated diluent solution to be diluted with distilled water.
7. Extraction solution (2.5x): 2 vials containing 50 ml concentrated extraction solution to be diluted with distilled water. This concentrated solution is irritating to eyes and skin.
8. Calibrators: 6 vials containing 1.5 ml Calprotectin solution at six known concentrations (0, 2.5, 12.5, 25, 50, 150 ng/ml). The value of each Calibrator is printed on the vial label. Ready to use.
9. Control 1: 1 vial containing 1.5 ml of control 1. Ready to use. Do not dilute. The range of values is printed on the vial label.
10. Control 2: 1 vial containing 1.5 ml of control 2. Ready to use. Do not dilute. The range of values is printed on the vial label.

MATERIALS REQUIRED BUT NOT PROVIDED

Feces sample collection

- Sample collection device (optional alternative use, i.e. **EasyCal** ref. 9062).
- Transport container.

Feces preparation

- Disposable, breakable sterile inoculation loops or wooden stick.
- Disposable polystyrene screw cap tubes, 14 ml.
- Eppendorf tubes (1-1.5 ml).
- Sensitive digital scale (40-150 mg).
- Vortex mixer.
- Shaker.
- Microcentrifuge (10,000 x g).
- Freezer (-20°C).

EQUIPMENT FOR ELISA MEASUREMENTS

1. Multi-channel pipette, 50, 200 µl.
2. Precision pipettes 5, 100 and 1,000 µl.
3. ELISA plate washer.
4. ELISA plate reader (filter 450 nm).
5. Distilled water.

PRECAUTIONS AND WARNINGS

1. For *in vitro* use only.
2. Follow universal precautions. The product does not use human origin materials.
3. All patient samples must be treated as potential biohazard, and shall be handled and disposed of according to local laboratory legislation.
4. Reagents, samples and microtiter strips should be allowed to reach room temperature (20-25°C) before starting the test.
5. **Warning:** do not interchange components from different kit batches. Satisfactory performance of the test is guaranteed only when components from the same batch of **Calprest NG** are used.
6. Unused microtiter strips should be re-sealed airtight in the foil bag with the enclosed drying pad and stored at 2-8°C.
7. Insufficient washing of the ELISA plate can lead to erroneous values of calprotectin due to incomplete removal of reagents. Routine maintenance of aspiration/wash system is strongly recommended.
8. Substrate is light sensitive, store in the dark and shake before use.
9. TMB substrate solution is irritating to the skin and mucous membranes. In case of contact, rinse eyes with plenty of water and wash skin with soap and plenty of water. Wash contaminated clothing before reuse. Do not discard the product into sewers. Store in the dark.
10. The Stop Solution contains sulfuric acid. Although diluted may cause burns must be handled with gloves, goggles and lab coat. In case of contact, rinse thoroughly with water. Use the same precautions when handling the extraction solution.
11. All reagents, except to the substrate and the concentrated washing solution, contain Proclin 300. It may cause an allergic reaction.
12. In case of external damage to packaging, make sure that all components listed in the "Materials Provided with the kit are present and undamaged". In the presence of liquid leakage from the vials, the kit should not be used.

PREPARATION OF WORKING SOLUTIONS

Extraction solution

Dilute concentrated Extraction solution by adding 1 part (50 ml) of it to 1.5 parts (75 ml) of freshly distilled water to obtain 125 ml working solution. Mix well.

Diluent

Dilute concentrated Diluent solution by adding 1 part (20 ml) of it to 9 parts (180 ml) of distilled water to obtain 200 ml working solution. Mix thoroughly.

Washing solution

Prepare the washing solution by diluting the content of the whole vial (50 ml) into 950 ml of distilled water.

STORAGE AND STABILITY OF REAGENTS AND WORKING SOLUTIONS

1. All reagents and working solutions, except washing solution, must be stored at 2-8 °C.
2. The expiration date is printed on all component labels. See table below for open-vial stability.
3. Avoid exposure to high temperature, direct sunlight or extreme humidity.
4. Unused microtiter strips should be resealed airtight in the plastic bag with the desiccant inside and stored at 2-8 °C.

Reagent stability (open-vial reagents)

The Substrate reagent is light sensitive. Store in the dark and shake before use.

Reagent	Storage conditions	Storage time
Conjugate	2-8°C	30 days
Substrate	2-8°C	90 days
Calibrators	2-8°C	30 days
Controls	2-8°C	30 days

Stability of working solutions

Reagent	Storage conditions	Storage time
Washing solution	20-25°C	30 days
Extraction solution	2-8°C	90 days
Diluent solution	2-8°C	30 days

SAMPLE COLLECTION

Random stool collection. Loose or liquid stool samples are acceptable as normalization to stool weight is part of the calculation of the result. Submission of stool samples from diapers should be avoided unless the sample submitted can be taken from a portion of the stool which is not in contact with the diaper material.

Sample requirements

1-5 g stool in a screw-top clean vial. No preservative is necessary or indicated.

Sample transport

Stool specimen should be received by the laboratory and extracted within 10 days of collection. Temperature during shipment should not exceed 30°C.

Sample Storage

Stool samples at the laboratory should be stored at 2-8 °C for up to 4 days before testing. If not immediately tested, freeze the stored samples at -20°C.

PROCEDURE**Specimen Preparation**

Thaw frozen stool samples at room temperature and ensure that all reagents reach room temperature (20-25 °C). There are two alternative stool sample preparation methods, manual and by using *EasyCal* stool extraction device. Users can choose their preferred sample preparation method.

A. EasyCal Specimen Preparation

- For the collection/extraction procedure, please refer to *EasyCal* ref. 9062 box insert.
- After the extraction procedure the fecal extract is ready to be tested either manually or automatically by placing the *EasyCal* tube (devoid of funnel and stick) directly into the ELISA instrumentation.

B. Manual Specimen Preparation

- Weigh (tare) the empty screw cap tube together with the inoculation loop or the wooden applicator stick.
- Mix well the stool and then take out approx 100 mg (between 80-120 mg) stool by means of the inoculation loop or wooden applicator stick and place into a screw-cap-tube.
- Weigh tube and loop with feces and calculate net feces weight (between 80-120 mg).
- Break off the loop handle or the wooden applicator stick, leaving the loop or wooden stick with feces and a 4-6 cm handle inside the screw cap tube.
- Add pre-diluted extraction solution (weight/volume ratio 1:50), e.g. 100 mg feces +4.9 ml diluted extraction solution as described in the table below. Close the tube.

Feces (mg)	Extraction solution (ml)	Feces (mg)	Extraction solution (ml)
120	5.9	100	4.9
115	5.6	95	4.7
110	5.4	90	4.4
105	5.2	85	4.2
		80	3.9

- Shake/mix vigorously for 30 seconds by means of a vortex.
- Homogenize 25 ± 5 minutes on a shaker or roller. The inoculation loop or the wooden stick inside the tube will act as an agitator.
- Transfer the homogenate (1 ml) to an Eppendorf tube and centrifuge for 20 minutes at 10,000 x g at Room Temperature (RT) using a bench-top centrifuge.
- Transfer 0.5 ml of the clear extract supernatant to a new Eppendorf tube. Avoid contact with the pellet as aggregates or particles can cause erroneous calprotectin values.
- The extracts may be tested immediately or stored at -20°C (up to three months) for later measurement.

ELISA Procedure

1. Ensure that all reagents reach room temperature (20-25°C).
2. Thaw frozen sample at room temperature.
3. Dilute samples 1:400 using two consecutive steps. In the first step, dilute 5 µl of extracted sample into 995 µl of diluent solution (1:200). Then take 500 µl of the first dilution and further dilute it 1:2 by adding 500 µl of diluent solution.
4. Suggested plate layout is shown below. Fit the strip holder with the required number of micro ELISA strips. Use uncoated strips to complete the strip holder if the washer requires a full plate. Calibrators and Controls must be included in each run.

	1	2	3	4	5
A	Calibrator 1	Calibrator 5	Sample	Sample	Sample
B	Calibrator 1	Calibrator 5	Sample	Sample	Sample
C	Calibrator 2	Calibrator 6	Sample	Sample	
D	Calibrator 2	Calibrator 6	Sample	Sample	
E	Calibrator 3	Control 1	Sample	Sample	
F	Calibrator 3	Control 1	Sample	Sample	
G	Calibrator 4	Control 2	Sample	Sample	
H	Calibrator 4	Control 2	Sample	Sample	

5. Add 100 µl of each Calibrator in duplicate wells (A1-B1, C1-D1, ...).
6. Add 100 µl of each control in duplicate wells (E2-F2, G2-H2, ...).
7. Mix diluted sample well before application to the plate and add 100 µl of each sample in following wells (A3, B3, C3, D3, ...).
8. Cover plate with plate cover and incubate the plate at room temperature (20-25°C) for 60 minutes.
9. At the end of the incubation time, wash the plate by adding 0.3 ml of diluted washing solution to each well. Remove as much liquid as possible. Repeat this step two more times up to a total of 3 washing steps. After the final aspiration, invert plate and tap gently on absorbent tissue to ensure complete removal of washing solution.
10. Add 100 µl conjugate to each well.
11. Cover plate with plate cover and incubate the plate at room temperature for 30 minutes.
12. Repeat washing step as above (see step 9).
13. Add 100 µl substrate solution to each well. **Note:** a multi-channel pipette is recommended in order to avoid variation in substrate development time. Avoid also the formation of bubbles during substrate pipetting.
14. Incubate the plate at room temperature for approximately 15 minutes in a dark place or wrap the plate with aluminium foil.
15. Add 100 µl of Stop Solution to all wells.
16. Read the O.D. values by means of an ELISA reader at 450 nm. The OD value of Calibrator 6 must be higher than 1.4 OD. It is possible to monitor the development of the color by reading the OD values at 620 nm. Values around 0.6 OD at 620 nm correspond to 1.8 ÷ -2.0 at 450 nm.
17. After addition of Stop Solution, plates may be stored at 4°C for 24 hours.

Assay calibration and quality control

1. A new calibration curve is used with each run.
2. The O.D. value of Calibrator S1 should be < 0.2 OD.
3. Control 1 and 2 are to be included in each run.

Reportable range

1.355 - 150 ng/ml (27.1 - 3,000 mg/kg).

Calculation of test results

Calculation of results is done by means of logit-log 4 parameter. Plot the Calibrator curve with actual calprotectin concentration of Calibrators as ng/ml and corresponding mean OD values on an x-y system. If the data processing system utilizes a log scale for X axis (concentration), the Calibrator 0 value should be entered as > 0 (e.g. 0.001). The readings of the Samples from the Calibrator curve is corrected for the dilution and converted to mg/kg by multiplying 20 (e.g. a reading of 20 ng/ml becomes 400 mg/kg). If samples are further diluted this must be compensated for during calculations. Concentrations can also be determined by use of a computer linked to the ELISA reader.

LIMITATIONS

1. False-negative results could occur in patients who have granulocytopenia due to bone marrow depression.
2. Some patients who are taking Non-Steroid-Anti-Inflammatory-Drugs (NSAIDs) will have elevations in their fecal calprotectin levels.
3. Patients with IBD fluctuate between active (inflammatory) and inactive stages of the disease. These stages must be considered when using the **Calprest NG** assay.
4. The use of proton pump inhibitors (PPIs), microscopic colitis and diverticular disease may also lead to elevated calprotectin level. Patients affected by untreated celiac disease may occasionally show elevated calprotectin value.
5. Other intestinal impairments, including many gastrointestinal infections and colorectal cancer, can result in elevated levels of calprotectin. These specimens will test positive with the **Calprest NG** assay. Therefore, a diagnosis of active IBD cannot be established solely on the basis of a positive result with the **Calprest NG** assay.
6. Fecal calprotectin is an indicator of neutrophilic presence in the stool and is not specific for IBD.

EXPECTED VALUES

An internal cut-off study was performed and the values are reported in the table below:

Calprotectin Concentration	Interpretation	Follow-Up
< 50 mg/kg	Normal	None
50 - 120 mg/kg	Borderline	Re-evaluate after 4-6 weeks
> 120 mg/kg	Abnormal	Repeat as clinically indicated

Further evaluations including asymptomatic patients, as well as patients with IBS (to differentiate from IBD) and international studies¹³⁻¹⁴ confirmed the appropriateness of such values.

CLINICAL EVALUATION

For the **Calprest NG** assay, the clinical study included 273 samples, of which 130 patients affected by Crohn's disease, Ulcerative and Intermediate Colitis and 143 negative samples from Irritable Bowel Syndrome (IBS), Recurrent Abdominal Pain (RAP) and other disease patients. The positive patient samples were diagnosed by means of clinical findings and/or confirmed with colonoscopy. Tables 1 and 2 below demonstrate the clinical performance of the **Calprest NG** assay.

Table 1 - Clinical Performance of **Calprest NG** assay (with **Calprest NG** borderline sample values considered positive with a cut-off of 50 mg/kg).

Borderline sample values considered Positive				
Clinical Diagnosis		Pos	Neg	Total
	Positive	123	14	137
Calprest NG	Negative	7	129	136
	Total	130	143	273
Sensitivity		94.6%	95% CI (89.2% - 97.8%)	
Specificity		90.2%	95% CI (84.1% - 94.5%)	
PPV*		89.8%	95% CI (83.4% - 94.3%)	
NPV**		94.9%	95% CI (89.7% - 97.9%)	

Table 2 - Clinical Performance of **Calprest NG** assay (with **Calprest NG** borderline sample values considered negative with a cut-off of 120 mg/kg).

Borderline sample values considered Negative				
Clinical Diagnosis		Pos	Neg	Total
	Positive	108	3	111
Calprest NG	Negative	22	140	162
	Total	130	143	273
Sensitivity		83.1%	95% CI (75.5% - 89.1%)	
Specificity		97.9%	95% CI (94.0% - 99.6%)	
PPV*		97.3%	95% CI (92.3% - 99.4%)	
NPV**		86.4%	95% CI (80.2% - 91.3%)	

* **PPV** - Positive Predictive Value | ** **NPV** - Negative Predictive Value | **CI** - Confidence Interval

Method Comparison

A method comparison study was performed which compared the **Calprest NG** to a comparator test using 157 clinical samples. These samples consist of clinically diagnosed IBD patients (clinical history and/or biopsy), while the negative samples were obtained from patients with IBS or other diseases. All samples were tested with the **Calprest NG** (y) and the predicate device (x) according to their corresponding package inserts. The results are summarized in Tables 3, 4 and 5.

Table 3 - Deming Regression Analysis.

	Slope (95% CI)	Y-intercept (95% CI)
$Y = 0.9232x + 0.1285$	0.9232 (0.8410 ÷ 1.005)	0.1285 (-4.643 ÷ 4.9)

Table 4 - Method comparison (with calprotectin borderline sample values considered positive with a cut-off of 50 mg/kg).

Borderline Value considered as Positive				
		Comparator test		
		Pos	Neg	Total
Calprest NG	Pos	113	0	113
	Neg	4	40	44
		117	40	157
Positive Agreement	96.6%	95% CI (91.5% - 98.7%)		
Negative Agreement	100.0%	95% CI (91.2% - 100.0%)		
Overall agreement	97.5%	95% CI (93.9% - 99.0%)		

Table 5 - Method comparison (with calprotectin borderline sample values considered negative with a cut-off of 120 mg/kg).

Borderline Value considered as Negative				
		Comparator test		
		Pos	Neg	Total
Calprest NG	Pos	73	3	76
	Neg	3	78	81
		76	81	157
Positive Agreement	96.1%	95% CI (89.0% - 98.6%)		
Negative Agreement	96.3%	95% CI (89.7% - 98.7%)		
Overall agreement	96.2%	95% CI (91.9% - 98.2%)		

PERFORMANCE CHARACTERISTICS

Matrix and aqueous linearity

3 high positive extracted stool samples, 3 low concentration extracted stool samples (matrix linearity) or a calprotectin reference material for calprotectin kit (aqueous linearity) were used in the study. The 3 high positive extracted stool samples were pooled (H) and diluted with a pool of the three low concentration extracted stool samples (L). The results of the assessment of the linearity ranges show that **Calprest NG** has both acceptable linearity and accuracy from 25.2 - 3,066 mg/kg for matrix and from 27.8 - 2,889.3 mg/kg for aqueous linearity.

Matrix Linearity					
Test Range (mg/kg)	Slope (95% CI)	Y-intercept (95% CI)	R ²	% Recovery Rate (Obtained/Theoretical)	
25.2 ÷ 3,066	1.014 (0.9868 ÷ 1.041)	9.279 (-26.54 ÷ 45.10)	0.994	86.9 ÷ 112.8	

Aqueous Linearity					
Test Range (mg/kg)	Slope (95% CI)	Y-intercept (95% CI)	R ²	% Recovery Rate (Obtained/Theoretical)	
27.8 ÷ 2,889.3	1.0322 (1.006 ÷ 1.0574)	-5.7463 (-36.40 ÷ 24.91)	0.999	93.0 ÷ 118.7	

Limits of Blank, Detection and Quantization

LoB, LoD and LoQ studies were performed according to EP17-A. The results are reported in the following table:

	Value	Concentration
	(ng/ml)	(mg/kg)
LoB	0.87	17.5
LoD	1.17	23.41
LoQ	1.35	27.1

Accuracy / Recovery

The study was carried out by testing 7 different stool samples. Each extracted stool sample was “spiked” with a constant quantity/volume of calprotectin or with an equal volume of sample diluent to compensate for volume adjustments. Each sample was assayed in triplicate. Data are shown in table below:

Calprotectin Recovery Data								
Sample		#1	#2	#3	#4	#5	#6	#7
Baseline	mg/kg	9.3	51.7	162.7	307.7	458.3	1,204.0	1,975.2
Spike Value	mg/kg	32.3	32.3	32.3	115.0	115.0	115.0	115.0
Theoretical (Base + Spike)	mg/kg	41.7	84.0	195.0	422.7	573.3	1,319.0	2,090.2
Observed (Base + Spike)	mg/kg	42.3	79.7	200.7	425.3	592.3	1,397.7	2,343.8
Recovery	%	101.6%	94.8%	102.9%	100.6%	103.3%	106.0%	112.1%

Extraction Reproducibility

In order to determine the extraction-to-extraction reproducibility, 4 samples (2 positive, 1 around cut-off and 1 negative) were extracted 10 times each, and each extract tested in duplicate. The results are reported below:

Extracted stool sample #				
	1	2	3	4
Mean (mg/kg)	29.6	52.7	227.0	1,973.3
SD	4.3	5.8	15.5	117.3
CV (%)	14.6	11.1	6.8	5.9

Single-Site Precision Evaluation study

The study has been carried out by extracting 7 different stool samples (6 positive and 1 around cut-off) and assaying each extract in 2 replicates during 2 separate runs per day for 10 days by 3 different operators. The CV% values should be equal and lower than 20%. Data are presented in the table below:

Single-Site Precision Evaluation study: Results												
ID #	N	Mean (mg/kg)	Within Run		Between Run		Between Day		Between Operator		Total	
			SD	CV (%)	SD	CV (%)	SD	CV (%)	SD	CV (%)	SD	CV (%)
1	120	170.88	3.40	2.0%	6.77	4.0%	5.51	3.2%	0.80	0.5%	9.40	5.5%
2	120	741.60	17.38	2.3%	43.12	5.8%	55.37	7.5%	0.00	0.0%	72.30	9.7%
3	120	248.74	6.44	2.6%	11.98	4.8%	12.21	4.9%	0.00	0.0%	18.28	7.3%
4	120	42.60	1.28	3.0%	4.30	10.1%	0.00	0.0%	0.55	1.3%	4.52	10.6%
5	120	1,193.10	52.91	4.4%	73.33	6.1%	118.78	10.0%	60.32	5.1%	161.01	13.5%
6	120	1,006.17	41.69	4.1%	39.67	3.9%	55.80	5.5%	28.24	2.8%	84.98	8.4%
7	120	2,267.29	223.08	9.8%	83.22	3.7%	155.52	6.9%	58.33	2.6%	290.31	12.8%

Multisite Precision Evaluation Study

The study has been carried out in 3 centers by testing 8 different stool sample extracts (6 positive, 1 around cut-off and 1 negative). Each center has tested each extract in 5 replicates over 5 days. The CV% values should be equal and lower than 20%.

ID Sample	N	Mean (mg/kg)	Repeatability		Precision		Reproducibility	
			SD	CV (%)	SD	CV (%)	SD	CV (%)
Sample A	75	433.80	13.43	3.1%	63.99	14.8%	64.55	14.9%
Sample B	75	36.54	3.69	10.1%	5.11	14.0%	5.11	14.0%
Sample C	75	1,180.64	71.08	6.0%	156.16	13.2%	156.16	13.2%
Sample D	75	682.76	35.88	5.3%	49.57	7.3%	78.95	11.6%
Sample E	75	1,629.44	96.29	5.9%	214.76	13.2%	226.84	13.9%
Sample F	75	62.08	4.19	6.8%	5.64	9.1%	8.68	14.0%
Sample G	75	356.84	12.00	3.4%	30.07	8.4%	51.93	14.6%
Sample H	75	2,152.20	173.58	8.1%	308.13	14.3%	312.19	14.5%

ID #	Repeatability 95% CIs				Precision 95% CIs				Reproducibility 95% CIs			
	SD		CV (%)		SD		CV (%)		SD		CV (%)	
A	11.39	to 16.35	2.6%	to 3.8%	46.39	to 103.09	10.7%	to 23.8%	47.26	to 101.80	10.9%	to 23.5%
B	3.13	to 4.49	8.6%	to 12.3%	4.17	to 6.61	11.4%	to 18.1%	4.17	to 6.61	11.4%	to 18.1%
C	60.33	to 86.54	5.1%	to 7.3%	118.0	to 230.93	10.0%	to 19.6%	119.47	to 225.50	10.1%	to 19.1%
D	30.45	to 43.68	4.5%	to 6.4%	40.42	to 64.14	5.9%	to 9.4%	49.28	to 193.64	7.2%	to 28.4%
E	81.72	to 117.23	5.0%	to 7.2%	162.28	to 317.59	10.0%	to 19.5%	168.95	to 345.24	10.4%	to 21.2%
F	3.56	to 5.10	5.7%	to 8.2%	4.63	to 7.22	7.5%	to 11.6%	5.42	to 21.29	8.7%	to 34.3%
G	10.19	to 14.61	2.9%	to 4.1%	22.56	to 45.07	6.3%	to 12.6%	31.11	to 149.23	8.7%	to 41.8%
H	147.32	to 211.32	6.8%	to 9.8%	240.6	to 428.66	11.2%	to 19.9%	242.64	to 437.93	11.3%	to 20.3%

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

REFERENCES / BIBLIOGRAFIA / BIBLIOGRAPHIE / BIBLIOGRAFÍA / LITERATUR

1. Fagerhol M. K. et al: Calprotectin (The L1 leukocyte protein) in: Smith V. L. and Dedman J. R. (eds.): Stimulus response coupling: The role of intracellular calcium-binding proteins. CRC Press, Boca Raton 1990, p. 187-210.
2. Johne, B. et al: Functional and clinical aspects of the myelomonocyte protein calprotectin. Mol Pathol, 1997. 50(3): p. 113-23.
3. Roseth, A. G. et al: Assessment of the neutrophil dominating protein calprotectin in feces. A methodologic study. Scand J Gastroenterol, 1992. 27(9): p. 793-8.
4. Roseth, A. G. et al: Assessment of disease activity in ulcerative colitis by fecal calprotectin, a novel granulocyte marker protein. Digestion, 1997. 58(2): p. 176-80.
5. Roseth A. G. et al: Correlation between fecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.
6. Ton H. et al: B. Improved assay for fecal calprotectin. Clinica Chimica Acta 2000; 292:41-54.
7. Limburg P. J. et al: Fecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. Am J Gastroenterol 2000; 95:2831-2837.
8. Tibble J. et al: A simple method for assessing intestinal inflammation in Crohn's disease. Gut 2000; 47:506-513.
9. Burri E, Beglinger C. Faecal calprotectin -- a useful tool in the management of inflammatory bowel disease. Swiss Med Wkly. 2012;142:w13557. Published 2012 Apr 5
10. Tibble J. A. et al: High prevalence of NSAID enteropathy as shown by a simple faecal test. Gut 1999; 45:362-366.
11. Montalto M. et al: Faecal calprotectin concentrations in untreated coeliac patients. Scand J. Gastr., 2007; 42: 957-961.
12. Calcaterra V. et al: Ann Nutr Metab 2018; 73:177–183 “Serum Calprotectin Level in Children: Marker of Obesity and its Metabolic Complications”.
13. Carroccio A., et al. Diagnostic accuracy of fecal calprotectin assay in distinguishing organic causes of chronic diarrhea from Irritable Bowel Syndrome: A prospective study in adults and children. Clin Chem 2003; 49:861-867.
14. Whitehead Ws et al.: Between-assay variability of fecal calprotectin enzyme-linked immunosorbent assay kits. Ann Clin Biochem 2012; 1-9.

Legend / Legenda / Légende / Leyenda / Legende	
 Consult instructions for use / Consultare le istruzioni per l'uso / Consulter la notice d'utilisation / Consultense las instrucciones de uso / Gebrauchsanweisung beachten	 Biological risk / Rischio biologico / Risque biologique / Riesgo biologico / Biogefährdung
 Use by / Data di scadenza / Date de péremption / Fecha de caducidad / Verwendbar bis	 Manufacturer / Fabbicante / Fabricant / Fabricante / Hersteller
 European Conformity / Conformità agli standard europei / Conformité aux normes européennes / Conformidad europea / Europäische Konformität	 In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Producto sanitario para diagnóstico in vitro / In-vitro-Diagnostikum
 Temperature limitation / Limiti di temperatura / Limites de température / Limitación de temperatura / Temperaturbegrenzung	 Sufficient for / Sufficiente per / Suffisant pour / Válido para / Ausreichend für
 Catalogue number / Numero di Catalogo / Référence du catalogue / Número de catálogo / Bestellnummer	 Solid phase / Fase solida / Phase solide / Fase sólida / Mikrotiterplatte
 Batch number / Numero di lotto / Numéro de lot / Número de lote / Chargennummer	 Calibrator / Calibratore / Calibrateur / Calibrador / Kalibrator
 Diluent / Diluente / Diluant / Diluyente / Probenverdünnungspuffer	 Substrate / Sustrato / Substrat / Substrat / Substrat
 Control / Controllo / Contrôle / Control / Kontrolle	 Washing solution / Soluzione di lavaggio / Solution de lavage / Solución de lavado / Waschlösung
 Stop solution / Soluzione di arresto / Solution d'arrêt / Solución de parada / Stopplösung	 Extraction solution / Soluzione di estrazione / Solution d'extraction / Solución de extracción / Extraktionpuffer
 Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	

Calprest NG

Code 9069, 96 test.

Codice 9069, 96 test.

Code 9069, 96 test.

Código 9069, 96 test.

Artikel 9069, 96 test.

Eurospital SpA

Via Flavia 122, 34147 Trieste, Italia
Tel. +39 040 8997.1 - Fax +39 040 280944
www.eurospital.com - info@eurospital.it