Helori CTX IgG

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Bibliografia

1. Vaira D, Holton J, Miglioli M, et al. Peptic ulcer disease and *Helicobacter pylori* infection. *Current opinion in Gastroenterology* 1994; 10: 98-104.
2. Vaira D, Miglioli M, Mulè P, et al. Prevalence of peptic ulcer in *Helicobacter pylori* positive blood donors. *Gut* 1994; 35: 309-12.
3. Figura N, Bugnoli M, Cusi MG, et al. Pathogenic mechanisms of *Helicobacter pylori*: production of cytotoxin. In: Malfertheiner P, Ditschuneil H (ed): *Helicobacter pylori, gastritis and peptic ulcer*: Springer, Berlin, 1990, p.86-95
4. Weel JF, Van der Hulst RW, Gerrits Y, et al. The interrelationship between cytotoxin-associated gene A, vacuolating cytotoxin, and *Helicobacter pylori*-related diseases. *J Infect Dis*, 173:5, 1996 May, 1171-5

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Helori CTX IgG

Enzyme immunoassay for the determination of IgG antibodies against CagA (Helicobacter pylori) in human serum.
For in vitro diagnostic use only

Application field

Helori-CTX is used for the determination of IgG antibodies directed against the cytotoxin associated protein (CagA) of Helicobacter pylori in the human serum. These antibodies have a relevant importance in the gastric pathologies. Helori CTX IgG represents a support to the clinical data for the diagnosis of this disease since the link between H.p. infection and ulcer development has been proved. Recent studies, moreover, show that H.p. strains expressing cytotoxic activity are more pathogen than those without the protein.

Materials provided with the kit

Solid phase	12x8 wells
Conjugate (IgG, goat antibody) (2000x)	1x300 µl
Substrate	1x14 ml
Positive control (human)	1x1.5 ml
Negative control (human)	1x1.5 ml
Calibrator	1x1.5 ml
Sample/conjugate diluent (10x)	1x30 ml
Washing solution (20x)	1x50 ml
Stop solution	1x15 ml

Composition of supplied reagents/materials

1. Solid phase

12 strips of 8 wells coated with purified antigen, in a plastic bag containing a desiccant.

2. Conjugate (IgG)

One vial containing 300 µl concentrate horse-radish peroxidase-labelled goat anti-human IgG antibodies (2000x) in a buffer solution with a preservative.

3. Substrate

One vial containing 14 ml TMB solution in a buffer solution with a preservative.

4. Positive control

One vial containing 1.5 ml positive control serum (human) with a preservative. Ready to use.

5. Negative control

One vial containing 1.5 ml negative control serum (human) with a preservative. Ready to use.

6. Calibrator

One vial containing 1.5 ml positive control serum (human) in a buffer solution with a preservative. The value of calibrator Units is printed on the vial label. Ready to use.

7. Sample/conjugate diluent

One vial containing 30 ml concentrated buffer (10x) solution with a preservative.

8. Washing solution

One vial containing 50 ml concentrate solution (20x) with detergent and preservative.

9. Stop solution

One vial containing 15 ml H₂SO₄ 1N.

Materials required but not provided

1. Microplate EIA reader (450 nm)
2. Precision pipettes (10, 20, 100, 1000 µl)
3. Incubator at 37°C.
4. Distilled water

Principle of the procedure

Helori CTX IgG is a sandwich type enzyme-immunoassay. The purified antigen is attached to microtitration wells. Antigen-specific antibodies in the patient's serum bind to the antigens. Non-specific antibodies are removed by washing.

Horse-radish peroxidase labelled goat-anti-human IgG are added and bind to human antibodies in the well. The excess conjugate is washed away; then a chromogenic substrate is added. After an appropriate incubation period, the Optical Density of the substrate is measured by means of a EIA reader. The higher the amount of specific antibodies present in the sample, the greater the O.D. values.

Performance criteria and limits of the procedure

1. Intra-assay

Four reference sera were tested 4 times.

The results are reported below:

Serum	1	2	3	4
Mean OD Value	0.144	0.253	0.488	1.005
C.V.%	10.1	3.9	4.3	2.6

2. Inter-assay variation

Three reference sera were tested in duplicate in three different runs. The results are reported below:

Serum	1	2	3
Mean Value (units)	7.9	14.8	31.8
C.V.%	7.2	9.3	6.5

Sensitivity and specificity

Sensitivity: 94.1%

Specificity: 97.9%

Preparation of reagents

Washing solution

Prepare the washing solution diluting the content of the whole vial with distilled water to a final volume of 1000 ml.

Sample/conjugate diluent

To be diluted 1:10 with distilled water.

Conjugate

Prepare the conjugate solution diluting 1:2000 with conjugate diluent. First dilute 1:50 (20 µl concentrated conjugate in 980 µl sample/conjugate diluent). Then transfer 300 µl (1:50 dilution) in 12 ml of sample/conjugate diluent. The diluted conjugate is stable for 24 hours.

Storage

All reagents must be stored at 2-8°C and used before the expiration date printed on the vial label.

Sample collection

Helori CTX IgG requires serum. It is important to preserve the integrity of the blood sample from its collection till the end of the test. Perform phlebotomy using a non-traumatic venipuncture technique. Allow blood to clot for at least 20-30 minutes, until the clot starts to retract. Before centrifugation, remove the clot from the collection tube walls. It is not necessary to add preservatives to the serum to keep integrity of the serum sample. Samples should be stored at 2-8°C and assayed

within 24 hours after collection. If the test cannot be run within 24 hours, serum sample should be frozen. Avoid repeated freezing and thawing of serum. Sample must be diluted 1:100 (10 µl serum in 0.99 ml sample diluent) before use. Diluted samples can be stored for 24 hours at 4°C or for a week at -20°C.

Procedure

Define the appropriate amount of patients sera needed to complete the test. Prepare the samples by diluting each serum 1:100 (10 µl serum in 0.99 ml sample diluent) and mix well. Calibrator, Positive and Negative controls provided with the kit are ready to use and are not to be diluted.

1. Pipette 100 µl sample/conjugate diluent into wells A1-B1.
2. Pipette 100 µl of calibrator into wells C1-D1.
3. Pipette 100 µl positive control, negative control and diluted samples into the remaining wells.
4. Cover the wells in order to prevent evaporation and incubate for 60 minutes at 37°C.
5. After incubation remove the liquid completely.
6. Wash 3 times using 200 µl/well of washing solution.
7. Remove the liquid as before.
8. Pipette 100 µl anti-IgG diluted conjugate in each well.
9. Incubate covered wells for 60 minutes at 37°C.
10. Repeat washing steps as above.
11. Pipette 100 µl chromogenic substrate in each well.
12. Incubate covered wells for 30 minutes at room temperature in the dark.
13. Stop the reaction by adding 100 µl stop solution. The colour of the reagent solution will turn from blue to yellow.
14. Read the results at 450 nm.

Calculation of results

Determine the average O.D. values of Blank, Controls and patients samples. Values must be compared with mean values of calibrator. Subtract the Blank O.D. (A1-B1) from all the obtained results. The concentration of each sample is expressed

as Units obtained from the following formula:

$$\text{Unit\`a} = \frac{\text{Sample O.D.} \times \text{Calibrator value (Units)}}{\text{Calibrator O.D.}}$$

O.D. value of Calibrator must be > 0.600. The value of calibrator Units is printed on the vial label. Do not use calibrators from different batches.

Quality control

O.D. of negative control must be always <0.100.

O.D. of positive control must be always >0.800.

Reference limits

Good laboratory practice indicates that each laboratory should establish its own range of normal values using a statistically significant number of samples without any clinical evidence of present or previous pathology connected with *Helicobacter pylori* infection. The negative range is usually indicated as Mean plus 2 Standard Deviations, while the lower limit of the positive is indicated by Mean plus 3 Standard Deviations. The range between 2 and 3 Standard Deviations refers to the borderline samples. The borderline values need further investigation and should re-tested. Therefore, each laboratory should establish its own cut-off value. Preliminary analysis of Italian samples provided with the following values:

Negative	< 5 Units
Borderline	5 - 7.5 Units
Positive	> 7.5 Units

The values reported above are just an indication and must not replace the routine method used in each laboratory.

Precautions

1. Concentration of Thymerosal used as preservative in the reagent solutions is lower than 0,05% (w/w). Avoid any unnecessary contact.

The stop solution may cause irritation of the skin, of the eyes and of the respiratory tract in sensitive individuals. Handle wearing suitable gloves.

2. Allow all the reagents to reach room temperature.

3. All sera used in this kit were found to be not reactive for Hepatitis B surface antigen (HBsAg), for anti-HCV and for anti-HIV 1-2. However, this material should be handled as potentially infectious.

4. Do not mix components from different boxes.

Bibliography

1. Vaira D, Holton J, Miglioli M, et al. Peptic ulcer disease and *Helicobacter pylori* infection. Current opinion in Gastroenterology 1994; 10: 98-104.
2. Vaira D, Miglioli M, Mulè P, et al. Prevalence of peptic ulcer in *Helicobacter pylori* positive blood donors. Gut 1994; 35: 309-12.
3. Figura N, Bugnoli M, Cusi MG, et al. Pathogenic mechanisms of *Helicobacter pylori*: production of cytotoxin. In: Malfertheiner P, Ditschuneil H (ed): *Helicobacter pylori*, gastritis and peptic ulcer: Springer, Berlin, 1990, p.86-95
4. Weel JF, Van der Hulst RW, Gerrits Y, et al. The interrelationship between cytotoxin-associated gene A, vacuolating cytotoxin, and *Helicobacter pylori*-related diseases. J Infect Dis, 173:5, 1996 May, 1171-5.

Helori CTX IgG

code 8723, 96 tests

Date of preparation

31.10.2004

Manufactured by:
ravo Diagnostika GmbH
Oltmannstr. 5
79100 Freiburg
Germany

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Eurospital SpA