Rapid immunochromatographic assay for quantitative detection of calprotectin from human fecal specimens to aid in the diagnosis of IBD and to differentiate IBD from IBS

REF 9191NT

INTENDED PURPOSE

CalFast[®] NeXT is an immunochromatographic rapid non-automated assay for the quantitative detection of human calprotectin in feces. Calprotectin is extracted with the dedicated pre-analytical processing device **EasyCal[®]** and the test cassette is read with the proprietary **CalFast[®] Reader**. **CalFast NeXT** can be used as an in vitro diagnostic to aid in the diagnosis in generic population of Inflammatory Bowel Diseases (IBD), specifically Crohn's disease and ulcerative colitis, and to differentiate IBD from Irritable Bowel Syndrome (IBS) in conjunction with other clinical and laboratory findings. For professional use only.

SUMMARY AND EXPLANATION

Calprotectin is a 36 kDa heterodimeric complex belonging to the S100 family of proteins that derives from polymorphonuclear cells, monocytes and epithelial cells. It has direct antimicrobial effects and it plays a role within the innate immune response^{1,2}. Calprotectin is found in several biological samples including plasma, serum, synovial fluid and stool in proportion to the existing level of inflammation^{1,2}. The use of calprotectin as a biological marker has been studied in inflammatory gastrointestinal disorders (IBD), as well as in other pathologies like rheumatoid arthritis and colorectal cancer^{1,3,10}.

TEST PRINCIPLE

CalFast NeXT is an immunochromatographic assay which combines a mixture of anti-calprotectin polyclonal and monoclonal antibodies. Calprotectin present in the diluted sample is bound by the gold particle-conjugated antibody. The so formed complex moves along the reaction strip and binds to the anti-calprotectin antibodies present in the reaction area (T) while unbound antibodies will bind in the control area (C). The concentration of calprotectin in the sample is calculated by means of the **CalFast Reader**.

MATERIALS PROVIDED (QUANTITY SUFFICIENT FOR 20 TESTS)

- Cassettes: 20 cassettes single-packed in a sealed aluminum foil containing a desiccant. Cassettes are disposable.
- Diluent: 1 vial containing 50 ml (yellow cap). Ready to use.
- **EasyCal** (ref. 79333): 20 devices containing a ready-to-use extraction solution. Devices are disposable.
- Positive Control: 1 vial containing 250 µl of Positive Control, ready to use (red cap). Do not dilute. The range of values is printed on the vial label.
- Negative Control: 1 vial containing 250 µl of Negative Control, ready to use (transparent cap). Do not dilute. The range of values is printed on the vial label.
- Caps: 20 white caps for **EasyCal**.
- Calibration strip: 1 paper strip with a printed out QR code containing reference calibration curve. Lot-dependent.

COMPOSITION OF THE REAGENTS

- Cassette: adsorbed mixture of monoclonal and polyclonal antibodies and conjugated gold particles (< 10 µg), detergents (< 100 µl), BSA (< 500 µg), sugars and salts (< 500 µg), sodium azide (< 0.02% v/v, final concentration after resuspension).
- Diluent: salts (< 200 mM), detergents (< 1% v/v) and preservatives (< 0.05% v/v) in dd water.
- Extraction solution in **EasyCal**: salts (< 500 mM), BSA (< 5% w/w), urea (< 3 M) and preservatives (< 0.05% v/v) in dd water.
- Positive Control: solution of a polyclonal antibody (< 1 mg/ml) in sample diluent.
- Negative Control: solution of non-reactive salts (< 200 mM), detergents (< 1% v/v) and preservatives (< 0.05% v/v) in dd water.

MATERIALS REQUIRED BUT NOT PROVIDED

- Plastic tubes for fecal extract dilution
- Vortex (or similar)
- Rolling shaker
- Precision pipettes with disposable tips (10-200 µl, 1 ml)
- Centrifuge for microtubes
- Eppendorf tubes (or similar) for extract dilution
- Timer
- **CalFast Reader** (ref. 9190)
- **CalFast NeXT** quality control set (optional, ref. 9191NQ, lot dependent, 3 x Positive Control and 3 x Negative Control)

STORAGE AND EXPIRY DATE OF REAGENTS

All reagents must be stored at 2-8 °C and used before expiration date printed on the label.

STABILITY AFTER FIRST OPENING

The Sample Diluent and quality controls are stable for 60 days after first opening, if stored at 2-8 °C in their original container.

TRANSPORTATION STABILITY

An accelerate stability test proved that all components of the kit are stable after storage for 72 hours at 37 °C.

WARNING AND PRECAUTIONS

- All patient samples must be treated as potential biohazard, and shall be handled and disposed according to local laboratory legislation. It is advised to handle stool as potentially infectious material.
 - **EasyCal** and Cassette are for single use. Unused reagents shall be stored at 2-8 °C and used before expiration date printed on the label.
 - Wait until stool samples reach room temperature (18-25 °C) before starting the sampling with **EasyCal**.
 - Wait until all reagents and extracted stool samples reach room temperature before starting the test.
 - Do not expose the kit to heat source.
 - In order to obtain reliable quantitative results, it is mandatory to centrifuge extracts before proceeding to their dilution.
 - When using supplemental quality controls (ref. 9191NQ), ensure the lot is the same as the **CalFast NeXT** kit in use.
 - **EasyCal**: EUH208 Contain Mixture of: 5-Chloro-2-methyl-2H-isotiazol-3-one [EC no. 247-500-7]; 2-Methyl-2H-isotiazol-3-one [EC no. 220-239-6] (3:1). May produce an allergic reaction. EUH210 Material safety data sheet available on request.
 - Positive and Negative controls: EUH208 Contain Mixture of: 5-Chloro-2-methyl-2H-isotiazol-3-one [EC no. 247-500-7]; 2-Methyl-2H-isotiazol-3-one [EC no. 220-239-6] (3:1). May produce an allergic reaction. EUH210 Material safety data sheet available on request.
 - Diluent: EUH208 Contain Mixture of: 5-Chloro-2-methyl-2H-isotiazol-3-one [EC no. 247-500-7]; 2-Methyl-2H-isotiazol-3-one [EC no. 220-239-6] (3:1). May produce an allergic reaction. EUH210 Material safety data sheet available on request.
 - In case of external damage to packaging, make sure that all components listed in the “Materials Provided” with the kit are present and undamaged. In the presence of liquid leakage from the vials, the kit should not be used.
 - Do not interchange components from different product lots.
 - The waste materials generated by the use of the kit must be disposed as dangerous special waste and treated in accordance with local laws.
 - Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
 - The product does not use human origin materials.
 - Any serious incident shall be reported to manufacturer and competent authority of the State in which user and/or patient is established.
-

PROCEDURE FOR STOOL COLLECTION

The stool specimen to be analyzed with **CalFast NeXT** must be collected by using **EasyCal** included in the kit package.

A. Sampling procedure for feces of normal consistency

1. Collect a stool sample in a sterile container.
2. Pull out the shaped stick of the pre-analytical processing by turning the screw-cap counter-clockwise.
3. Dip the shaped stick in the previously collected stool sample.
4. Repeatedly rotate the stick in the stool sample until all grooves are filled with fecal material.
5. Before placing the stick into the device, remove the stool excess by rotating the stick on the internal wall of the stool container.
6. Fully insert the stick, with fecal material into the tube containing the extraction solution through the filter (blue). While turning the screw cap clockwise until completely closed, keep the filter (blue) by the fingers in order to avoid its opening.

B. Sampling procedure for liquid feces or anyhow having a consistence which does not allow the use of the EasyCal device

1. It is necessary to have a laboratory pipette.
2. By using the pipette, aspirate 56 µl of liquid stool.
3. Pull out the shaped stick and the filter (blue) by turning them clockwise.
4. Pipette the liquid stool sample into the tube containing the extraction solution.
5. Place the shaped stick and the filter back into the tube containing the extraction solution only. Turn the filter (blue) counter-clockwise until completely closed.

C. Sampling procedure in case of absence of the EasyCal device

1. Should the laboratory temporarily not have the **EasyCal**, collect a stool sample in a sterile container and keep it at 2-8 °C for no more than 2 days. Should the storage period exceed 2 days at 2-8 °C, the sample must be frozen at -20 °C for maximum 1 year.
 2. Once the device is available again, proceed to the sampling procedure described above. If the sample was frozen, allow it to reach room temperature (18-25 °C) before to collect the stool sample by means of the device.
 3. Proceed from point 2 of the “Sampling procedure for feces of normal consistency”.
-

STORAGE OF SPECIMENS

- Stool specimens should be received by the laboratory and extracted within 3 days of collection. Temperature during shipment should not exceed 30 °C, use of ice pack is highly recommended.
- Stool specimens at the laboratory can be stored into their original container at 2-8 °C up to 2 days, if not immediately tested. For longer storage, it is necessary to store them frozen at -20 °C. Frozen stool specimens are stable up to 1 year at -20 °C and can be subjected to no more than 3 freeze/thaw cycles.
- The **EasyCal** containing the collected stool specimen should be processed immediately. If there is a delay in continuing with the laboratory procedure, the **EasyCal** containing the stool specimen can be stored at room temperature (18-25 °C) up to 72 hours, otherwise it must be stored at 2-8 °C (up to 14 days). If the **EasyCal** device containing the specimen has to be delivered to the laboratory it is advisable to store it in a thermal container with cooling packs for a proper preservation.

EXTRACTION PROCEDURE

1. Bring the specimens to room temperature (18-25 °C).
2. Shake each **EasyCal** device by vortex mixer for 30-60 seconds in order to properly homogenate the content. If the grooves of the devices are visually empty of any stool material, proceed to step 3, otherwise repeat the vortexing for other 30-60 seconds. After a maximum vortexing time of 120 seconds, proceed to step 3 regardless if there is some residual stool sample in the grooves.
 - 2.1 Alternatively, after a brief vortexing, capped devices can be left on a rolling shaker for 20 minutes.
3. After the completion of the previous steps, remove the capped stick (white) together with the funnel (blue) from devices by rotating the funnel clockwise. Close the tube with the provided additional screw cap.
4. Centrifuge the sample for 10 minutes at 1,000 - 3,000 x g in order to remove the possible residuals of fecal material. After centrifugation, remove the device from the centrifuge.
5. The fecal extract is now ready to be tested.

STORAGE OF THE EXTRACTED SAMPLES

After the completion of the aforementioned extraction procedure, the fecal extract is stable under the following conditions:

- Up to 72 hours at room temperature (18-25 °C).
- Up to 14 days at 2-8 °C.
- Up to 90 days at -20 °C. Frozen extracted samples can be subjected to up to 4 freeze/thaw cycles.

SETTING UP CALFAST READER BEFORE RUNNING THE TEST

1. Check for the **CalFast Reader** being properly connected to a PC via the USB cable. Turn on the computer, double-click on the icon Client Reader and enter User Name and Password.
2. Verify that the **CalFast NeXT** lot information present in the software correspond to the ones present on the product label (lot number and expiry date). If not, please do as follows:
 - 2.1 Place the two-dimensional barcode included in the package into the slot of the instrument.
 - 2.2 Check that the barcode is properly placed inside the **CalFast Reader** correctly.
 - 2.3 Click on this icon  in the **CalFast Reader** Client software for reading the barcode. All information regarding the **CalFast NeXT** lot will be automatically uploaded into the **CalFast Reader**.
 - 2.4 Save the info by clicking on .
 - 2.5 Exit by clicking on this icon .
3. This procedure must be repeated every time a new lot of **CalFast NeXT** is used. The uploaded information will be stored into the software until a new lot information will be read.

For more detailed information, please refer to **CalFast Reader** User Manual. Additional information:

- Batch code: indicates the lot currently in use.
- Kit name: indicates the product in use (**CalFast**, **CalFast XT** or **CalFast NeXT**).
- Check Strip: to be used before starting a test (optional). When activated, the image of the cassette appears and a check is made on the presence of the Control line. If the Control line is present and correctly positioned, the Assay Result window in the Client Reader software will display "READABLE". On the contrary, if control line is missing or positioning of the cassette is incorrect, it will show "UNREADABLE".

SAMPLES DATA INPUT

1. To start the **CalFast NeXT** test reading procedure, at least one of the following data must be entered: Family Name, Name, Date of Birth or Patient code.
2. Once completed, the green LED in the Client Reader software lights up, indicating that the cassette can be placed inside the slot of the **CalFast Reader**.
3. Place the cassette into the slot of the **CalFast Reader**: Check if the position is correct by clicking on "Check the strip" (optional).

TEST PROCEDURE

1. Ensure that all reagents reach room temperature (18-25 °C).
2. If there are frozen samples to be analyzed, thaw them at room temperature.
3. Dilute stool extracts 1:300 (e.g., add 5 µl of extracted sample into 1,495 µl of diluent solution). Mix by vortex (10 seconds).
4. Open the sealed aluminum foil containing the cassette. Take it and place onto a horizontal surface.
5. Pipette either 100 µl of diluted fecal extract, 100 µl of Positive Control (ready-to-use) or 100 µl of Negative Control (ready-to-use), directly into the round well of the test cassette.
6. The use of Positive and/or Negative Control is not mandatory. Additional positive and negative lot dependent controls can be purchased separately (ref. 9191NQ).
7. Incubate at room temperature for 20 minutes.

TEST READING

1. If the cassette is readable, start the test by pressing the icon  in the Client Reader software. While reading the cassette, the yellow LED turns on and the icon  is replaced by .
2. The Check Strip box shows two different images: one is the “real” image of the cassette, while the captured image is used to perform the calculations for the determination of calprotectin concentration.
3. Captured image: it is the image used for calculations. The red and blue lines define the reading range for the calculation.
4. **CalFast Reader** performs the reading and reports the result in mg/kg or µg/g.
5. The concentration value is automatically saved and associated to the specific sample information.

RESULT VALIDATION

- To consider a test result valid, the Control line (C) must be always present. This line is just a functional control of the test and cannot be used for the interpretation of the Test line (T).
- The amount of calprotectin in the stool sample (T) is calculated by the **CalFast Reader**.
- If only the line T appears, the test is not valid and must be repeated by using a new cassette.
- If neither the line (C) and the line (T) show up, the result is not valid and the test must be repeated by using a new cassette. Note: The usability of the test cassette can be assessed via the “Check the Strip” function of the **CalFast Reader** Client software.
- The results obtained with the Positive and Negative Controls must fall into the range of concentration described on the Controls labels. Additional controls lot dependent can be purchased separately (ref. 9191NQ).

METROLOGICAL TRACEABILITY

In order to guarantee the metrological traceability of **CalFast NeXT**, concentration values recorded for the calibration points and the controls are compared to the reference range values internally defined.

RESULTS OBTAINED BY CALFAST READER

The **CalFast Reader** uses a lot specific calibration curve to calculate the calprotectin concentration. Calprotectin values obtained by **CalFast Reader** will be displayed as reported in table below:

Calprotectin concentration detected	Calprotectin concentration displayed by the CalFast Reader
0 - 50 µg/g	< 50 µg/g
51 - 2,100 µg/g	51 - 2,100 µg/g
> 2,101 µg/g	> 2,100 µg/g

REFERENCE RANGE AND INTERPRETATION OF RESULTS

Eurospital performed an internal cutoff study and the values are reported in the table below:

Calprotectin concentration	Interpretation	Follow-Up
< 70 µg/g	Normal	None
70 - 140 µg/g	Borderline	Re-evaluate after 4-6 weeks
> 140 µg/g	Abnormal	Repeat as clinically indicated

Samples having a calprotectin concentration higher than 140 µg/g must be considered positive to the **CalFast NeXT** test hence being submitted to further clinical investigation to assess the pathological cause for the high calprotectin concentration. In healthy adult patients, the average value of calprotectin is around 25 µg/g. This value, if obtained by means of **CalFast Reader**, will be reported as < 50 µg/g. While in adult and pediatric samples calprotectin values are similar, they may be significantly higher in newborns^{2, 4, 7, 11, 12, 13}.

LIMITATIONS

- Patients who are taking Non-Steroid-Anti-Inflammatory-Drugs (NSAIDs) medicaments can result having higher than normal fecal calprotectin levels^{9,18,24}.
- Patients with IBD fluctuate between active (inflammatory) and inactive stages of the disease²⁵. These stages must be considered when using the **CalFast NeXT**.
- The use of Proton Pump Inhibitors (PPIs)²⁴, as well as presence of microscopic colitis^{21,22} and diverticular disease²³, may also lead to elevated calprotectin level. Patients affected by untreated celiac disease may occasionally show elevated calprotectin values^{19,20}.
- Other intestinal impairments, including many gastrointestinal infections^{15,16,17} and colorectal cancer^{13,14,22}, can result in elevated levels of calprotectin. These specimens will test positive with the **CalFast NeXT** assay. Therefore, a diagnosis of active IBD cannot be established solely on the basis of a positive result with the **CalFast NeXT** assay.
- Fecal calprotectin is an indicator of neutrophilic presence in the stool and is not specific for IBD²⁶.

ANALYTICAL PERFORMANCES

Linearity

Linearity study was performed according to CLSI EP06 2nd Edition. Linearity ranges assessment of **CalFast NeXT** showed acceptable results for both linearity and accuracy for calprotectin concentrations in the range from 40 to 2,100 µg/g.

Limits of Blank, Detection and Quantitation

LoB, LoD and LoQ studies were performed according to CLSI EP17-A. The results are reported in the following table:

LoB	17.7 µg/g
LoD	29.4 µg/g
LoQ	35.0 µg/g

Single-Site Precision Evaluation study

The study has been carried out by extracting 6 different stool samples (4 positive, 1 in the cutoff area and 1 negative) and assaying each extract in 3 replicates during 2 separate runs per day for 20 days. The obtained CV% values were equal or lower than 20%. Results are presented in the following table:

ID #	N	Mean (µg/g)	Within Run		Between Run		Between Day		Total	
			SD (µg/g)	CV	SD (µg/g)	CV	SD (µg/g)	CV	SD (µg/g)	CV
1	120	56.0	6.9	12.4%	3.1	5.5%	2.9	5.1%	7.6	13.6%
2	120	75.9	10.5	13.8%	2.9	3.8%	3.0	3.9%	10.8	14.3%
3	120	238.0	29.5	12.4%	7.2	3.0%	9.4	4.0%	30.4	12.8%
4	120	446.6	58.9	13.2%	12.8	2.9%	19.3	4.3%	60.4	13.5%
5	120	783.5	87.8	11.2%	14.8	1.9%	20.4	2.6%	89.2	11.4%
6	120	1559.2	300.0	19.2%	0.0	0.0%	28.9	1.9%	290.9	18.7%

Lot-to-Lot Precision Evaluation study

The study has been carried out by extracting 6 different stool samples (4 positive, 1 in the cutoff area and 1 negative) and assaying each extract in 5 replicates during 1 run per day for 5 days using 3 different lots of **CalFast NeXT**. The obtained CV% values were equal or lower than 20%. Results are presented in the following table:

ID #	N	Mean (µg/g)	Repeatability		Within Lot		Between Lot		Reproducibility	
			SD (µg/g)	CV	SD (µg/g)	CV	SD (µg/g)	CV	SD (µg/g)	CV
A	75	54.4	6.2	11.4%	7.1	13.1%	5.9	10.9%	9.2	17.0%
B	75	66.4	8.3	12.6%	9.0	13.6%	7.6	11.4%	11.7	17.7%
C	75	277.6	34.3	12.4%	39.0	14.1%	0.0	0.0%	39.0	14.1%
D	75	409.8	44.8	10.9%	52.4	12.8%	16.1	3.9%	54.8	13.4%
E	75	764.7	109.6	14.3%	127.5	16.7%	0.0	0.0%	127.5	16.7%
F	75	1326.3	182.6	13.8%	184.7	13.9%	41.1	3.1%	189.2	14.3%

Interference and Cross-Reactive study

An interference and cross-reactive study was conducted to determine whether potentially interfering or cross-reacting substances might affect the performance of **CalFast NeXT**. Tested molecules are listed in the table below and, at the tested concentrations, do not affect calprotectin concentrations in human stool samples when measured with **CalFast NeXT**.

Interfering/cross-reacting substance	Tested quantity / 50 mg stool
Hemoglobin	5.83 mg
Vancomycin	0.67 mg
Ciprofloxacin	0.50 mg
Prevacid	0.02 mg
Azathioprine	0.07 mg
Prednisone	0.01 mg
5-aminosalicylic acid	1.33 mg
Vitamin A	8 IU
Vitamin D	1.1 IU

Vitamin C	0.05 mg
Vitamin E	0.10 mg
Barium sulfate	0.25 mg
Loperamide HCl	0.34 µg
Bismuth subsalicylate	0.04 mg
Metronidazole	0.63 mg
Mucin	0.17 mg
Aluminum hydroxide + Magnesium hydroxide	0.21 mg
Palmitic acid	0.07 mg
Stearic acid	0.13 mg
Polyethylene glycol 3350	3.95 mg
Omeprazole	0.03 mg
S100A12 protein	1.05 µg
Simethicone	0.03 mg
Cimetidine	0.03 mg
Calcium carbonate	0.03 mg

CLINICAL PERFORMANCES

For the **CalFast NeXT** assay, the clinical study included 81 samples, of which 42 patients affected by Crohn's disease, Ulcerative, Intermediate and Microscopic Colitis and 39 negative samples from Irritable Bowel Syndrome (IBS), Recurrent Abdominal Pain (RAP) and other disease patients. The positive patient samples were diagnosed by means of clinical findings and/or confirmed with colonoscopy. In order to support the clinical information, the samples used for this study have been previously tested by **Calprest NG** IVD. Tables below demonstrate the clinical performance of the **CalFast NeXT** assay.

Borderline sample values considered as Positive				
		CalFast NeXT (cutoff = 70 µg/g)		
		Positive	Negative	Tot
Clinical diagnosis	IBD	42	0	42
	NON IBD	12	27	39
	Total	54	27	81
Sensitivity		100.0%	95% CI* (91.6% - 100.0%)	
Specificity		69.2%	95% CI* (53.6% - 81.4%)	
PPV**		76.5%	95% CI* (67.0% - 83.9%)	
NPV***		100.0%	95% CI* (100.0% - 100.0%)	

*CI = Confidence Interval (Mercaldo-Wald)

**PPV = Positive Predictive Value

***NPV = Negative Predictive Value

Borderline sample values considered as Negative				
		CalFast NeXT (cutoff = 140 µg/g)		
		Positive	Negative	Tot
Clinical diagnosis	IBD	40	2	42
	NON IBD	4	35	39
	Total	44	37	81
Sensitivity		95.2%	95% CI* (84.2% - 98.7%)	
Specificity		89.7%	95% CI* (76.4% - 95.9%)	
PPV**		90.9%	95% CI* (79.8% - 96.2%)	
NPV***		94.6%	95% CI* (81.8% - 98.5%)	

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Limit of Blank (LoB), Limit of Detection (LoD) und Limit of Quantification (LoQ)

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
***NPV = Negativer Vorhersagewert

REFERENCES / BIBLIOGRAFIA / BIBLIOGRAPHIE / BIBLIOGRAFÍA / LITERATUR

1. Fagerhol M. K. et al: Calprotectin (The L1 leukocyte protein) in: Smith V. L. and Dedman J. R. (eds.): Stimulus response coupling: The role of intracellular calcium-binding proteins. CRC Press, Boca Raton 1990, p. 187-210.
2. Johne, B. et al: Functional and clinical aspects of the myelomonocyte protein calprotectin. *Mol Pathol*, 1997. 50(3): p. 113-23.
3. Roseth, A. G. et al: Assessment of disease activity in ulcerative colitis by fecal calprotectin, a novel granulocyte marker protein. *Digestion*, 1997. 58(2): p. 176-80.
4. Roseth A. G. et al: Correlation between fecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. *Scand J Gastroenterol* 1999; 34:50-54.
5. Ton H. et al: B. Improved assay for fecal calprotectin. *Clinica Chimica Acta* 2000; 292:41-54.
6. Limburg P. J. et al: Fecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. *Am J Gastroenterol* 2000; 95:2831-2837.
7. Tibble J. et al: A simple method for assessing intestinal inflammation in Crohn's disease. *Gut* 2000; 47:506-513.
8. Burri E, Beglinger C. Faecal calprotectin -- a useful tool in the management of inflammatory bowel disease. *Swiss Med Wkly*. 2012;142:w13557. Published 2012 Apr 5.
9. Tibble J. A. et al: High prevalence of NSAID enteropathy as shown by a simple faecal test. *Gut* 1999; 45:362-366.
10. Calcaterra V. et al: *Ann Nutr Metab* 2018; 73:177–183 “Serum Calprotectin Level in Children: Marker of Obesity and its Metabolic Complications”.
11. Campeotto F. et al.: Fecal calprotectin: cutoff values for identifying intestinal distress in preterm infants. *J Pediatr Gastroenterol Nutr*. 2009 Apr; 48(4):507-10.
12. Campeotto et al.: High faecal calprotectin concentrations in newborn infants. *Arch Dis Child Fetal Neonatal Ed* 2004; 89:F353–F355.
13. Johne B. et al.: A new fecal calprotectin test for colorectal neoplasia. Clinical results and comparison with previous method. *Scand J Gastroenterol* 2001; 36:291-296.
14. Luley K, Noack F, Lehnert H, Homann N. Local calprotectin production in colorectal cancer and polyps—active neutrophil recruitment in carcinogenesis. *Int J Colorectal Dis* 2011;26:603-7.
15. Shastri YM, Bergis D, Povse N, Schäfer V, Shastri S, Weindel M, et al. Prospective multicenter study evaluating fecal calprotectin in adult acute bacterial diarrhea. *Am J Med* 2008;121:1099-106.
16. Nielsen HL, Engberg J, Ejlersten T, Nielsen H. Evaluation of fecal calprotectin in *Campylobacter concisus* and *Campylobacter jejuni/coli* gastroenteritis. *Scand J Gastroenterol* 2013;48:633-5.
17. Chen CC, Huang JL, Chang CJ, Kong MS. Fecal calprotectin as a correlative marker in clinical severity of infectious diarrhea and usefulness in evaluating bacterial or viral pathogens in children. *J Pediatr Gastroenterol Nutr* 2012;55:541-7.
18. Maiden L, Thjodleifsson B, Theodors A, Gonzalez J, Bjarnason I. A quantitative analysis of NSAID-induced small bowel pathology by capsule enteroscopy. *Gastroenterology* 2005;128:1172-8.
19. Ertekin V, Selimoğlu MA, Turgut A, Bakan N. Fecal calprotectin concentration in celiac disease. *J Clin Gastroenterol* 2010;44:544-6.
20. Balamtekin N, Baysoy G, Uslu N, Orhan D, Akçören Z, Özen H, et al. Fecal calprotectin concentration is increased in children with celiac disease: relation with histopathological findings. *Turk J Gastroenterol* 2012;23:503-8.
21. Lissette Batista et al.: Usefulness of fecal calprotectin as a biomarker of microscopic colitis in a cohort of patients with chronic watery diarrhoea of functional characteristics. *Dig Liver Dis*. 2019 Dec;51(12):1646-1651.
22. Gian P Caviglia, et al. Fecal calprotectin: beyond intestinal organic diseases. *Panminerva Med* 2018 Mar;60(1):29-34.
23. Davor Štimac, et al. What's New in Diagnosing Diverticular Disease . *J Gastrointestin Liver Dis*, 2019 Vol. 28 Suppl 4: 17-21.
24. David Lundgren, et al. Proton pump inhibitor use is associated with elevated faecal calprotectin levels. A cross-sectional study on subjects referred for colonoscopy. *Scand J Gastroenterol*. 2019 Feb;54(2):152-157.
25. Jost Langhorst, et al. Noninvasive markers in the assessment of intestinal inflammation in inflammatory bowel diseases: performance of fecal lactoferrin, calprotectin, and PMN-elastase, CRP, and clinical indices. *Am J Gastroenterol*. 2008 Jan;103(1):162-9.
26. Ganga W Pathirana, et al. Faecal Calprotectin. *Clin Biochem Rev* 39 (3) 2018.

Legend / Legenda / Légende / Leyenda / Legende

	Consult electronic instructions for use / Consultare le istruzioni per l'uso elettroniche / Consulter la notice d'utilisation électronique / Consultarse las instrucciones electrónicas de uso / Elektronischen Gebrauchsanweisung beachten		Manufacturer / Fabbricante / Fabricant / Fabricante / Hersteller
	Use by / Data di scadenza / Date de péremption / Fecha de caducidad / Verwendbar bis		<i>In vitro</i> diagnostic medical device / Dispositivo medico-diagnostico <i>in vitro</i> / Dispositif médical de diagnostic <i>in vitro</i> / Producto sanitario para diagnóstico <i>in vitro</i> / <i>In-vitro</i> -Diagnostikum
	European Conformity / Conformità agli standard europei / Conformité aux normes européennes / Conformidad europea / Europäische Konformität		Sufficient for / Sufficiente per / Suffisant pour / Válido para / Ausreichend für
	Temperature limitation / Limiti di temperatura / Limites de température / Limitación de temperatura / Temperaturbegrenzung		Disposable / Monouso / Usage unique / Desechable / Einmalige Verwendung
	Catalogue number / Numero di Catalogo / Référence du catalogue / Número de catálogo / Bestell Nummer		Stool collection and extraction device / Dispositivo raccolta ed estrazione feci / Dispositifs de collecte et extraction des selles / Dispositivo de recogida y extracción de heces / Stuhlröhrchen zur Probenextraktion
	Cassette / Cassetta / Cassette / Cartucho / Kassetten		Diluent / Diluente / Diluant / Diluyente / Probenverdünnungspuffer
	Batch number / Numero di lotto / Numéro de lot / Número de lote / Chargennummer		Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negativkontrolle
	Caps / Tappi / Casquettes / Tapas / Kappen		Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
	Two-dimensional barcode / Codice a barre bidimensionale / Code à barres bidimensionnel / Código de barras bidimensional / Zweidimensionaler Barcode		

English
Italiano
Français
Español
Deutsch

Calfast NeXT

Ref. 9191NT, 20 tests.

Ref. 9191NT, 20 test.

Ref. 9191NT, 20 tests.

Ref. 9191NT, 20 ensayos.

Ref. 9191NT, 20 tests.

Eurospital SpA

Via Flavia 122, 34147 Trieste, Italia
Tel. +39 040 8997.1 - Fax +39 040 280944
www.eurospital.com - info@eurospital.it