FOB-Test Turbo

Immunoturbidimetric assay for quantitative detection of occult blood in feces

REF 9297

English

Italiano

Français

Español

Deutsch

INTENDED PURPOSE

FOB-Test Turbo is a quantitative latex turbidimetric assay for detecting concentration of occult blood (or hemoglobin, Hb) in feces. For professional use only.

TEST PRINCIPLE

Hemoglobin is a globular protein found in red blood cells that ensures the transport of oxygen in the blood. Since hemoglobin is sensitive to degradation in the gastrointestinal tract (stomach and small intestine), its presence in the feces is due to the bleeding of specific colorectal lesions¹. Intact hemoglobin alone can therefore be detected quickly and non-invasively in fecal samples, and can be used as a marker for colorectal cancer^{2,3}, thus directing the patient to further clinical investigations⁴.

FOB-Test Turbo is a particle-enhanced turbidimetric immunoassay (PETIA) for the quantitative detection of hemoglobin in human stool samples. The test was designed to be used on several automated analyzers. It is easy to use, widely applicable and the results are obtained in a very short period of time.

FOB-Test Turbo is based on the in vitro agglutination of antibodies anti-hemoglobin coated polystyrene latex particles with the hemoglobin antigen present in the analyzed solution. If there is no hemoglobin present in the sample, the mixture will keep its smooth appearance. With the presence of hemoglobin, the turbidity of the medium will increase after the aggregation of the antigen with the antibodies on the latex particles. Such turbidity is measured as an increase in absorbance at the determinate wavelength and is proportional to the quantity of hemoglobin present in the sample. Compared to normal values, a high concentration of hemoglobin in the feces indicates a conspicuous bleeding of the lower gastrointestinal tract, attributable to colorectal cancer.

MATERIALS PROVIDED (QUANTITY SUFFICIENT FOR AT LEAST 100 TESTS)

Component	Quantity	Cap color	Composition
FOB-Test Turbo Reagent 1	1 bottle x 22 ml		Enhancer reagent, ready to use
FOB-Test Turbo Reagent 2	1 bottle x 8 ml		Solution with latex particles coated with antibodies against human hemoglobin, ready to use
FOB-Test Turbo Calibrator 0	1 vial x 1 ml	Yellow	Hemoglobin solution (0 ng/ml), ready to use
FOB-Test Turbo Calibrator 1	1 vial x 1 ml	Pink	Hemoglobin solution (50 ng/ml), ready to use
FOB-Test Turbo Calibrator 2	1 vial x 1 ml	Red	Hemoglobin solution (100 ng/ml), ready to use
FOB-Test Turbo Calibrator 3	1 vial x 1 ml	Violet	Hemoglobin solution (250 ng/ml), ready to use
FOB-Test Turbo Calibrator 4	1 vial x 1 ml	Green	Hemoglobin solution (500 ng/ml), ready to use
FOB-Test Turbo Calibrator 5	1 vial x 1 ml	Blue	Hemoglobin solution (1,000 ng/ml), ready to use
FOB-Test Turbo Control 1	2 vials x 1 ml	White	Human hemoglobin in buffer solution ready to use, do not dilute The range of values is printed on the label
FOB-Test Turbo Control 2	2 vials x 1 ml	Brown	Human hemoglobin in buffer solution ready to use, do not dilute The range of values is printed on the label

MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen container
- Transport container
- Disposable gloves and laboratory equipment
- Stool sample collection device (i.e. **EasyCal FOB** ref. 9293)
- Vortex mixer
- Centrifuge for 1,000 - 3,000 x g

Equipment for turbidimetric measurements

Any turbidimetric automated analyzer; for the instrument-specific protocol, see additional Application note.

STORAGE AND STABILITY OF REAGENTS AND WORKING SOLUTIONS

- All reagents and working solutions must be stored at 2-8°C.
- The expiration date is printed on all component labels.
- Avoid exposure to high temperature, direct sunlight or extreme humidity.

PRECAUTIONS AND WARNINGS

- Patient specimens and samples should be handled in accordance with Good Laboratory Practice (GLP) using appropriate precautions and anyway as if capable of transmitting infections.
- For *in vitro* diagnostic use.
- A trained person in turbidimetric technique and autoanalyzer use is required.
- Do not exchange components from different lots. Good performance can only be guaranteed by using components from the same lot of **FOB-Test Turbo**.
- Do not mix screw caps from different reagent vials.
- The analyzer should be ready before performing any assay.
- Components must not be used after the expiry date printed on the labels.
- If the outer packaging is damaged, ensure that all components listed in the “Materials Provided” section are present and intact. In case of leakage of liquid from the bottles, the kit must not be used.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The solutions should be discarded in a proper container after testing following local regulations.

SPECIMEN COLLECTION

Collect sufficient and representative quantity of human stool samples. These specimens should be collected in clean and dry single-use containers (no preservatives or transport media). Loose or liquid stool specimens are acceptable as normalization to stool weight is part of the calculation of the result. Collection of stool specimen from diapers should be avoided unless the specimen can be taken from a portion which is not in contact with the diaper material.

Specimen requirements

1-5 g stool in a screw-cap clean vial. No preservative is necessary or indicated.

PROCEDURE

Sample preparation

Thaw frozen stool specimen and ensure that all reagents for extraction procedure reach room temperature (20-25°C).

EasyCal FOB sample preparation

1. For the collection/extraction procedure, please refer to **EasyCal FOB** ref. 9293 package insert.
2. After the extraction procedure the fecal extract is ready to be tested. Insert the **EasyCal FOB** tube (devoid of funnel and stick) directly into the turbidimetric analyzer or transfer the supernatant to an adaptor cup (not provided) for the analytical procedure.

STORAGE OF THE FECAL EXTRACT

The fecal extract is stable at 4°C for up to 7 days if protected from direct light.

ANALYTICAL PROCEDURE

1. Ensure that all reagents reach room temperature (20-25°C).
2. Thaw frozen sample at room temperature.
3. **FOB-Test Turbo** Reagent 1 and **FOB-Test Turbo** Reagent 2 are ready to use solutions filled in bottles that can be loaded directly into the Eu-Turbo Analyzer (ref. 9194). Shall the test be used on a different turbidimetric analyzer, the content of the bottles has to be transferred into bottles proprietary of the used system.
4. **FOB-Test Turbo** Calibrator vials (0 to 5) and **FOB-Test Turbo** Control vials (1 and 2) are ready to use.
5. Setting up the calibration curve:
 - 5.1 For the calibration of the curve use only **FOB-Test Turbo** Calibrator vials. The concentration is indicated on the label of each vial.
 - 5.2 Once set up, calibration curve is stable up to 30 days on board, although it is highly recommended to calibrate the system every week.
 - 5.3 Recalibrate the system when reagent lot is changed or when the controls are out of the assigned range given on the control label or certificate of analysis provided with the kit.
6. Quality control:
 - 6.1 For quality control use **FOB-Test Turbo** Control vials. The concentration is indicated on the label of each vial.
 - 6.2 **FOB-Test Turbo** Controls should be assayed every day before running patient fecal samples extract to validate the calibration curve.
 - 6.3 If the obtained results are out of the control lot-specific indicated range, the equipment, the reagents or the technique must be reviewed.
7. Running the test:
 - 7.1 The methodology can be applied to the most common turbidimetric analyzer. For the list of validated turbidimetric analyzers and the relative application notes, please contact Eurospital at the following email address:

marketing.diagnostics@eurospital.it

INTERPRETATION OF TEST RESULTS

- Negative results, hemoglobin values below 50 ng/ml indicate a normal hemoglobin level in the fecal sample analyzed.
- Positive results, hemoglobin values above 200 ng/ml indicate an abnormal presence of hemoglobin in the fecal sample analyzed. Adult individuals with these values are identified as patients with gastrointestinal bleeding and therefore require further diagnostic investigations.
- Haemoglobin values between 50 and 200 ng/ml of stool are considered as an abnormal presence of hemoglobin, being indicative of mild bleeding, therefore monitoring and follow-up of the patient is recommended.

The following chart summarizes the interpretation of results:

Hemoglobin concentration	Interpretation	Follow-Up
< 50 ng/ml	Normal	None
50 - 200 ng/ml	Mild bleeding	Monitoring and follow-up of the patient
> 200 ng/ml	Abnormal	Repeat as clinically indicated

PERFORMANCE CHARACTERISTICS

FOB-Test Turbo cutoff value:

- Hemoglobin concentration values lower than 50 ng/ml determined from stool samples are to be considered normal values, therefore not indicative of a colorectal bleeding.
- Hemoglobin concentration values higher than 200 ng/ml determined from stool samples are indicative of a conspicuous colorectal bleeding.

LIMITATIONS

- Negative results do not rule out bleeding caused by polyps or colorectal cancer, as they can manifest as intermittent bleeding. Also, blood may not be evenly distributed in stool samples.
- This test may be less sensitive to detect upper gastrointestinal bleeding because blood degrades as it passes through the gastrointestinal tract.
- Samples collected during the menstrual period or samples from patients with bleeding haemorrhoids, blood in the urine or constipation are not valid.
- **FOB-Test Turbo** should be used for the determination of fecal hemoglobin only from human stool samples. The use of the test on samples different than human stool has not been validated yet. The quality of the test depends also on the quality of the analyzed sample; proper fecal specimens must be obtained.
- Positive results indicate an abnormal presence of hemoglobin in fecal samples. A positive result should be followed up by additional clinical investigations, such as colonoscopies, in order to properly identify the nature of bleeding and establish its extent.

INTERFERENCES AND CROSS REACTIVITY

In an evaluation to determine cross reactivity, no cross reactivity was found against other fecal markers occasionally present in feces, such as bovine serum albumin, ascorbic acid and hemin.

CLINICAL PERFORMANCE

A method comparison study was performed between Eurospital **FOB-Test Turbo** and a predicate. In the study 100 stool samples were tested with both methods according to the manufacturer's instruction for use. Results are summarized in the following tables:

		Predicate		
		Positive	Negative	Total
FOB-Test Turbo ref. 9297	Positive	20	2	22
	Negative	0	78	78
	Total	20	80	100

Positive percent agreement	100%
Negative percent agreement	97.5%
Overall percent agreement	98%

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REFERENCES / BIBLIOGRAFIA / BIBLIOGRAPHIE / BIBLIOGRAFÍA / LITERATUR

1. Garcia et al: Fecal hemoglobin concentration as a measure of risk to tailor colorectal cancer screening: are we there yet? *Eur J Cancer Prev.* 2015;24(4):321–7.
2. Yeasmin F. et al: A comparative study of chemical and immunological method of fecal occult blood test in the diagnosis of occult lower gastrointestinal bleeding. *Bangladesh Med Res.* 2013.
3. Towler B.P. et al: Screening for colorectal cancer using the faecal occult blood test, Hemoccult. *Cochrane Database Syst Rev.* 2000.
4. Ransohoff DF and Lang CA: Screening for colorectal cancer with the Fecal Occult Blood Test: a background paper. *Ann Intern Med.* 1997; 126: 811-822.

English**Italiano****Français****Español****Deutsch**

Legend / Legenda / Légende / Leyenda / Legende

	Consult electronic instructions for use / Consultare le istruzioni per l'uso elettroniche / Consulter la notice d'utilisation électronique / Consultense las instrucciones electrónicas de uso / Elektronischen Gebrauchsanweisung beachten		Sufficient for / Sufficiente per / Suffisant pour / Válido para / Ausreichend für
	Use by / Data di scadenza / Date de péremption / Fecha de caducidad / Verwendbar bis		Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller
	European Conformity / Conformità agli standard europei / Conformité aux normes européennes / Conformidad europea / Europäische Konformität		In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Producto sanitario para diagnóstico in vitro / In-vitro-Diagnostikum
	Temperature limitation / Limiti di temperatura / Limites de température / Limitación de temperatura / Temperaturbegrenzung		Reagent / Reagente / Réactif / Reactivo / Reagens
	Catalogue number / Numero di Catalogo / Référence du catalogue / Número de catálogo / Bestell Nummer		Control / Controllo / Contrôle / Control / Kontrolle
	Batch number / Numero di lotto / Numéro de lot / Número de lote / Chargennummer		Calibrator / Calibratore / Calibrateur / Calibrador / Kalibrator

English
Italiano
Français
Español
Deutsch

**ANNEX INFORMATION OF THE VALIDATED COMPATIBLE TURBIDIMETRIC INSTRUMENTS
INFORMAZIONI DEGLI STRUMENTI TURBIDIMETRICI COMPATIBILI VALIDATI
ANNEXE INFORMATIONS SUR LES INSTRUMENTS TURBIDIMÉTRIQUES COMPATIBLES VALIDÉS
ANEXO INFORMACIÓN DE LOS INSTRUMENTOS TURBIDIMÉTRICOS COMPATIBLES VALIDADOS
ANHANG INFORMATIONEN DER VALIDIERTEN KOMPATIBELN TURBIDIMETRISCHEN INSTRUMENTE**

Manufacturer	Model
Eurospital	Eu-Turbo Analyzer
Abbott	Alinity c-series; Architect c4000 / c8000 / c16000
Awareness	Chemwell-T
Beckman Coulter	AU480/AU680
Biosystems	A15
DyaSis	Respons 910
Mindray	BS-200/200E
Ortho Clinical Diagnostics	Vitros 5600
Roche	Cobas c301/c501/c111
Tecom	TC220
Tokyo Boeki	Biolis i24/i50
Erba	XL-180/XL-200
Binding Site	Optilite
Siemens	Atellica
Furuno	CA800
Biobase	BK-Mini

FOB-Test Turbo

Ref. 9297, 100 tests.

Ref. 9297, 100 test.

Ref. 9297, 100 tests.

Ref. 9297, 100 ensayos.

Ref. 9297, 100 tests.

Eurospital SpA

Via Flavia 122, 34147 Trieste, Italia
Tel. +39 040 8997.1 - Fax +39 040 280944
www.eurospital.com - info@eurospital.it