Elaprest Turbo

Immunoturbidimetric assay for quantitative detection of pancreatic elastase from fecal specimens

REF 9298

INTENDED PURPOSE

Elaprest Turbo is a quantitative latex turbidimetric assay for specifically detecting the concentration of fecal human pancreatic elastase E1 that is intended to aid in the evaluation of pancreatic exocrine functionality. For professional use only.

TEST PRINCIPLE

Pancreatic elastase E1 is an enzyme produced and secreted by the pancreas which, thanks to its proteolytic activity, promotes the digestion of food ingested as it passes through the gastrointestinal tract. Unlike other pancreatic enzymes, pancreatic elastase is not significantly degraded during intestinal transit, therefore determining the concentration of this marker in the stool allows to evaluate the exocrine functionality of the pancreas¹. The use of pancreatic elastase as a biological marker in fecal samples makes it possible to distinguish patients with pancreatic exocrine insufficiency (PEI) from those with similar symptoms (abdominal pain, vomiting and diarrhoea), attributable, however, to other pathologies such as peptic ulcer, Crohn's disease and ulcerative colitis²⁻⁴.

Elaprest Turbo is a particle-enhanced turbidimetric immunoassay (PETIA) for the quantitative detection of pancreatic elastase in human stool samples. The test was designed to be used on several automated analyzers. It is easy to use, widely applicable and the results are obtained in a very short period of time. **Elaprest Turbo** is based on the in vitro agglutination of antibodies anti-pancreatic elastase coated polystyrene latex particles with the pancreatic elastase antigen present in the analyzed solution. If there is no pancreatic elastase present in the sample, the mixture will keep its smooth appearance. With the presence of pancreatic elastase, the turbidity of the medium will increase after the aggregation of the antigen with the antibodies on the latex particles. Such turbidity is measured as an increase in absorbance at the determinate wavelength and is proportional to the quantity of pancreatic elastase present in the sample. Compared to normal values, a low concentration of pancreatic elastase in the feces indicates a reduced function of the pancreas, a direct consequence of PEI.

MATERIALS PROVIDED (QUANTITY SUFFICIENT FOR AT LEAST 100 TESTS)

Component	Quantity	Cap color	Composition
Elaprest Turbo Reagent 1	1 bottle x 24 ml		Enhancer reagent, ready to use
Elaprest Turbo Reagent 2	1 bottle x 4 ml		Solution with latex particles coated with antibodies against human pancreatic elastase, ready to use
Elaprest Turbo Calibrator 0	1 vial x 1 ml	Yellow	Pancreatic elastase solution (0 µg/g), ready to use
Elaprest Turbo Calibrator 1	1 vial x 1 ml	Pink	Pancreatic elastase solution (25 µg/g), ready to use
Elaprest Turbo Calibrator 2	1 vial x 1 ml	Red	Pancreatic elastase solution (50 µg/g), ready to use
Elaprest Turbo Calibrator 3	1 vial x 1 ml	Violet	Pancreatic elastase solution (100 µg/g), ready to use
Elaprest Turbo Calibrator 4	1 vial x 1 ml	Green	Pancreatic elastase solution (200 µg/g), ready to use
Elaprest Turbo Calibrator 5	1 vial x 1 ml	Blue	Pancreatic elastase solution (400 µg/g), ready to use
Elaprest Turbo Control 1	2 vials x 1 ml	White	Human recombinant pancreatic elastase in buffer solution ready to use, do not dilute The range of values is printed on the label
Elaprest Turbo Control 2	2 vials x 1 ml	Brown	Human recombinant pancreatic elastase in buffer solution ready to use, do not dilute The range of values is printed on the label

MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen container.
- Transport container.
- Disposable gloves and laboratory equipment.
- Stool sample collection device (*i.e.* **EasyCal Turbo** ref. 9302).
- Disposable, breakable sterile inoculation loops or wooden stick (if using fecal manual extraction method described on page 2).
- Disposable screw cap tubes, 14 ml (if using fecal manual extraction method described on page 2).
- Extraction solution 1x (**Calprest Turbo** Extraction kit ref. 9299, if using fecal manual extraction method described on page 2).
- Microcentrifuge tubes (1-1.5 ml).
- Sensitive digital scale (40-150 mg) (if using fecal manual extraction method described on page 2).
- Vortex mixer.
- Shaker or roller (if using fecal manual extraction method described on page 2).
- Microcentrifuge (10,000 x g) (if using fecal manual extraction method described on page 2).
- Precision pipettes.

Equipment for turbidimetric measurements

Any turbidimetric automated analyzer; for the instrument-specific protocol, see additional Application note.

STORAGE AND STABILITY OF REAGENTS AND WORKING SOLUTIONS

- All reagents and working solutions must be stored at 2-8°C.
- The expiration date is printed on all component labels.
- Avoid exposure to high temperature, direct sunlight or extreme humidity.

PRECAUTIONS AND WARNINGS

- Patient specimens and samples should be handled in accordance with Good Laboratory Practice (GLP) using appropriate precautions and anyway as if capable of transmitting infections.
- For *in vitro* diagnostic use.
- The kit does not contain materials of human origin.
- A trained person in turbidimetric technique and autoanalyzer use is required.
- Do not exchange components from different lots. Good performance can only be guaranteed by using components from the same lot of **Elaprest Turbo**.
- Do not mix screw caps from different reagent vials.
- The analyzer should be ready before performing any assay.
- Components must not be used after the expiry date printed on the labels.
- If the outer packaging is damaged, ensure that all components listed in the “Materials Provided” section are present and intact. In case of leakage of liquid from the bottles, the kit must not be used.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The solutions should be discarded in a proper container after testing following local regulations.
- EUH210 - Safety data sheet available upon request.

SPECIMEN COLLECTION

Collect sufficient and representative quantity of human stool samples. These specimens should be collected in clean and dry single-use containers (no preservatives or transport media). Collection of stool specimen from diapers should be avoided unless the specimen can be taken from a portion which is not in contact with the diaper material.

Specimen requirements

1-5 g stool in a screw-cap clean vial. No preservative is necessary or indicated.

Specimen transport

Stool specimens should be received by the laboratory and extracted within 3 days of collection. Temperature during shipment should not exceed 30°C, use of ice packs is highly recommended.

Specimen storage

Stool samples at the laboratory shall be stored at 2-8°C for up to 2 days before testing, if not tested, freeze the stool samples at -20°C for up to 1 year. Avoid more than 3 freeze / thaw cycles.

PROCEDURE

Sample preparation

Thaw frozen stool specimen and ensure that all reagents for extraction procedure reach room temperature (20-25°C). There are two alternative stool sample preparation methods, manual or by using **EasyCal Turbo** stool extraction device. Users can choose their preferred sample preparation method.

A. **EasyCal Turbo** sample preparation

1. For the collection/extraction procedure, please refer to **EasyCal Turbo** ref. 9302 instruction for use.
2. After the extraction procedure the fecal extract is ready to be tested. Insert the **EasyCal Turbo** tube (devoid of funnel and stick) directly into the turbidimetric analyzer or transfer the supernatant to an adaptor cup (not provided) for the analytical procedure.
3. If the extract is not analyzed immediately, it can be stored according to the instructions for use of **EasyCal Turbo** (Ref. 9302).

B. **Manual** sample preparation

Use a screw cap centrifugation tube for each sample to be tested. Label centrifugation tube with patient name or number.

1. Weigh (tare) the empty screw cap tube together with the inoculation loop or the wooden applicator stick.
2. Mix well the stool by means of the inoculation loop or wooden applicator stick in order to make the sample as homogenous as possible. Using the inoculation loop (or wooden applicator stick) take out approximately 50 mg (between 40-60 mg) of stool and place into the screw-cap-tube (weighted in point 1).
3. Weigh the tube with loop and feces inside and calculate net feces weight (between 40-60 mg).
4. Break off the loop handle or the wooden applicator stick, leaving the loop or wooden stick with feces and a 4-6 cm handle inside the screw cap tube.
5. Add ready to use extraction solution (**Calprest Turbo** Extraction kit ref. 9299, with final weight/volume ratio 1:100), *i.e.* 50 mg feces + 4.95 ml extraction solution as described in the table below. Fecal sample concentration in the extraction solution will be 10 mg/ml. Close the tube.

Feces (mg)	Extraction solution (ml)	Feces (mg)	Extraction solution (ml)
60	5.94	45	4.45
55	5.45	40	3.96
50	4.95		

6. Shake/mix vigorously for 30 seconds by means of a vortex.
7. Homogenize 25 ± 5 minutes on a shaker or roller. The inoculation loop or the wooden stick inside the tube will act as an agitator.
8. Transfer the homogenate (1 ml) to a microcentrifuge tube and centrifuge for 20 minutes at 10,000 x g at room temperature.
9. Transfer 0.5 ml of the clear extract supernatant to a new microcentrifuge tube. Avoid contact with the pellet as aggregates or particles may lead to erroneous pancreatic elastase values.
10. The extracts can be tested immediately or stored at the following conditions:
 - Up to 72 hours at room temperature
 - Up to 14 days at 2-8°C
 - Up to 90 days at -20°C. Frozen samples can be subjected to up to 4 freeze/thaw cycles

ANALYTICAL PROCEDURE

1. Ensure that all reagents reach room temperature (20-25°C).
2. Thaw frozen sample at room temperature.
3. **Elaprest Turbo** Reagent 1 and **Elaprest Turbo** Reagent 2 are ready to use solutions filled in bottles that can be loaded directly into the Eu-Turbo Analyzer (ref. 9194). Shall the test be used on a different turbidimetric analyzer, the content of the bottles has to be transferred into bottles proprietary of the used system.
4. **Elaprest Turbo** Calibrator vials (0 to 5) and **Elaprest Turbo** Control vials (1 and 2) are ready to use.
5. **Setting up the calibration curve:**
 - 5.1 For the calibration of the curve use only **Elaprest Turbo** Calibrator vials. The concentration is indicated on the label of each vial.
 - 5.2 Once set up, calibration curve is stable up to 30 days on board, although it is highly recommended to calibrate the system every week.
 - 5.3 Recalibrate the system when reagent lot is changed or when the controls are out of the assigned range given on the control label or certificate of analysis provided with the kit.
6. **Quality control:**
 - 6.1 For quality control use **Elaprest Turbo** Control vials. The concentration is indicated on the label of each vial.
 - 6.2 **Elaprest Turbo** Controls should be assayed every day before running patient fecal samples extract to validate the calibration curve.
 - 6.3 If the obtained results are out of the control lot-specific indicated range, the equipment, the reagents or the technique must be reviewed.
7. **Running the test:**
 - 7.1 The methodology can be applied to the most common turbidimetric analyzer. For the list of validated turbidimetric analyzers and the relative application notes, please contact Eurospital at the following email address:

marketing.diagnostics@eurospital.it

INTERPRETATION OF TEST RESULTS

- **Negative results**, pancreatic elastase values higher than 200 µg/g indicate a normal level of human pancreatic elastase in the analyzed fecal sample.
- **Positive results**, pancreatic elastase values lower than 100 µg/g indicate an abnormal concentration of human pancreatic elastase in the analyzed fecal sample. Adult individuals with these values are identified as patients with reduced pancreatic function and therefore require further diagnostic investigations.
- **Borderline results**, pancreatic elastase values in the gray zone (100-200 µg/g), must be re-analyzed after 2-6 weeks.

The following chart summarizes the interpretation of results:

Pancreatic elastase concentration	Interpretation	Follow-Up
> 200 µg/g	Normal	None
100 - 200 µg/g	Borderline	Re-evaluate after 2-6 weeks
< 100 µg/g	Abnormal	Repeat as clinically indicated

PERFORMANCE CHARACTERISTICS

Elaprest Turbo cutoff value:

- Pancreatic elastase concentration values higher than 200 µg/g determined from stool samples are to be considered normal values, therefore not indicative of reduced pancreatic function.
- Pancreatic elastase concentration values lower than 100 µg/g determined from stool samples are indicative of reduced pancreatic function.

LIMITATIONS

Elaprest Turbo should be used for the determination of pancreatic elastase only from human stool samples. The use of the test on samples different than human stool has not been validated yet. The kit is not intended to be used for different matrixes such as blood, serum, plasma, urine, cerebrospinal fluid, oral fluid, synovial fluid or empyema fluid. The quality of the test depends also on the quality of the analyzed sample; proper fecal specimens must be obtained.

INTERFERENCES AND CROSS REACTIVITY

An evaluation study was performed to determine the interferences and the cross reactivity of **Elaprest Turbo**. No interferences were founded against BSA, Hemin, Ascorbic acid and no cross reactivity was founded against porcine pancreatic elastase.

CLINICAL PERFORMANCE

A method comparison study was performed between Eurospital **Elaprest Turbo** and a predicate (ScheBo® Pancreatic Elastase 1™ Stool Test) with 93% diagnostic specificity and 93% diagnostic sensibility. In the study 40 samples from donors with gastrointestinal disorders were tested with both methods according to the manufacturer's instruction for use. Results are summarized in the following tables:

		Predicate		
		Positive	Negative	Total
Elaprest Turbo ref. 9298	Positive	10	2	12
	Negative	0	28	28
	Total	10	30	40

Positive percent agreement	100%
Negative percent agreement	93%
Overall percent agreement	95%

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REFERENCES / BIBLIOGRAFIA / BIBLIOGRAPHIE / BIBLIOGRAFÍA / LITERATUR

1. A. Szegoleit & D. Linder: Studies on the Sterol-Binding Capacity of Human Pancreatic Elastase 1; Gastroenterology; 1991
2. Andriulli et al.: Exocrine Pancreatic Insufficiency, as Assessed by Fecal Elastase-1 Levels, in Diabetic Patients: An Estimate of Prevalence in Prospective Studies; Journal of Diabetes and Metabolism; 2014
3. Dominguez-Muñoz, J.E. et al.: Potential for Screening for Pancreatic Exocrine Insufficiency Using the Fecal Elastase-1 Test; Digestive Diseases and Sciences; 2017.
4. R. Vanga et al.: Diagnostic Performance of Measurement of Fecal Elastase-1 in Detection of Exocrine Pancreatic Insufficiency: Systematic Review and Meta-analysis; Clinical Gastroenterology and Hepatology; 2018

Legend / Legenda / Légende / Leyenda / Legende	
 Consult electronic instructions for use / Consultare le istruzioni per l'uso elettroniche / Consulter la notice d'utilisation électronique / Consultense las instrucciones electrónicas de uso / Elektronischen Gebrauchsanweisung beachten	 Sufficient for / Sufficiente per / Suffisant pour / Válido para / Ausreichend für
 Use by / Data di scadenza / Date de péremption / Fecha de caducidad / Verwendbar bis	 Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller
 European Conformity / Conformità agli standard europei / Conformité aux normes européennes / Conformidad europea / Europäische Konformität	 In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Producto sanitario para diagnóstico in vitro / In-vitro-Diagnostikum
 Temperature limitation / Limiti di temperatura / Limites de température / Limitación de temperatura / Temperaturbegrenzung	 Reagent / Reagente / Réactif / Reactivo / Reagens
 Catalogue number / Numero di Catalogo / Référence du catalogue / Número de catálogo / Bestel Nummer	 Control / Controllo / Contrôle / Control / Kontrolle
 Keep dry / Conservare in luogo asciutto / Garder au sec / Almacenar en lugar seco / Bleib trocken	 Calibrator / Calibratore / Calibrateur / Calibrador / Kalibrator
 Keep out of the sunlight / Tenere lontano dalla luce del sole / Gardez hors de la lumière du soleil / Manténgase fuera de la luz del sol / Vor dem Sonnenlicht schützen	 Batch number / Numero di lotto / Numéro de lot / Número de lote / Chargennummer

**ANNEX INFORMATION OF THE VALIDATED COMPATIBLE TURBIDIMETRIC INSTRUMENTS
 INFORMAZIONI DEGLI STRUMENTI TURBIDIMETRICI COMPATIBILI VALIDATI
 ANNEXE INFORMATIONS SUR LES INSTRUMENTS TURBIDIMÉTRIQUES COMPATIBLES VALIDÉS
 ANEXO INFORMACIÓN DE LOS INSTRUMENTOS TURBIDIMÉTRICOS COMPATIBLES VALIDADOS
 ANHANG INFORMATIONEN DER VALIDIERTEN KOMPATIBLEN TURBIDIMETRISCHEN INSTRUMENTE**

Manufacturer	Model
Eurospital	Eu-Turbo Analyzer
Abbott	Architect c4000 / c8000 / c16000
Abbott	Alinity c-series
Awareness	Chemwell-T
Binding Site	Optlite
Mindray	BS-200
Roche	cobas c111
Tecom	TC220
Tokyo Boeki	Biolis i24/i50

Elaprest Turbo

- Ref. 9298**, 100 tests.
Ref. 9298, 100 test.
Ref. 9298, 100 tests.
Ref. 9298, 100 ensayos.
Ref. 9298, 100 tests.

Eurospital SpA
 Via Flavia 122, 34147 Trieste, Italia
 Tel. +39 040 8997.1 - Fax +39 040 280944
 www.eurospital.com - info@eurospital.it