Immunoturbidimetric assay for quantitative detection of calprotectin from fecal specimens

REF 9300

INTENDED PURPOSE

Calprest® Turbo is a quantitative latex turbidimetric assay for detecting concentration of fecal calprotectin that is intended to aid in the diagnosis of Inflammatory Bowel Diseases (IBD), (e.g. Crohn's disease and ulcerative colitis), and to differentiate IBD from Irritable Bowel Syndrome (IBS) in conjunction with other clinical and laboratory findings.

TEST PRINCIPLE

Calprotectin is a 36 kDa heterodimeric complex belonging to the S100 family of proteins that derives from polymorphonuclear cells, monocytes and epithelial cells. It has direct antimicrobial effects and it plays a role within the innate immune response¹⁻². Calprotectin is found in several biological samples including plasma, serum, synovial fluid and stool in proportion to the existing level of inflammation³. The use of calprotectin as a biological marker has been studied in inflammatory gastrointestinal disorders (IBD), as well as in other pathologies like rheumatoid arthritis and colorectal cancer⁴⁻¹².

Calprest Turbo is a particle-enhanced turbidimetric immunoassay (PETIA) for the quantitative detection of calprotectin in human stool samples. The test was designed to be used on several automated analyzers. It is easy to use, widely applicable and the results are obtained in a very short period of time.

Calprest Turbo is based on the in vitro agglutination of antibodies anti-calprotectin coated polystyrene latex particles with the calprotectin antigen present in the analyzed solution. If there is no calprotectin present in the sample, the mixture will keep its smooth appearance. With the presence of calprotectin, the turbidity of the medium will increase after the aggregation of the antigen with the antibodies on the latex particles. Such turbidity is measured as an increase in absorbance at the determinate wavelength and is proportional to the quantity of calprotectin present in the sample.

MATERIALS PROVIDED (QUANTITY SUFFICIENT FOR AT LEAST 200 TESTS)

Component	Quantity	Cap color	Composition
Calprest Turbo Reagent 1	2 bottles x 27 ml		Enhancer reagent, ready to use
Calprest Turbo Reagent 2	1 bottle x 8 ml		Solution with latex particles coated with antibodies against human calprotectin, ready to use
Calprest Turbo Calibrator 0	1 vial x 1 ml	Yellow	Calprotectin solution (0 µg/g), ready to use
Calprest Turbo Calibrator 1	1 vial x 1 ml	Pink	Calprotectin solution (50 µg/g), ready to use
Calprest Turbo Calibrator 2	1 vial x 1 ml	Red	Calprotectin solution (100 µg/g), ready to use
Calprest Turbo Calibrator 3	1 vial x 1 ml	Violet	Calprotectin solution (250 µg/g), ready to use
Calprest Turbo Calibrator 4	1 vial x 1 ml	Green	Calprotectin solution (750 µg/g), ready to use
Calprest Turbo Calibrator 5	1 vial x 1 ml	Blue	Calprotectin solution (1,500 µg/g), ready to use
Calprest Turbo Control 1	2 vials x 1 ml	White	Human recombinant calprotectin in buffer solution ready to use, do not diluite The range of values is printed on the label
Calprest Turbo Control 2	2 vials x 1 ml	Brown	Human recombinant calprotectin in buffer solution ready to use, do not diluite The range of values is printed on the label

MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen container.
- Transport container.
- Disposable gloves and laboratory equipment.
- Stool sample collection device (i.e. **EasyCal Turbo** ref. 9302).
- Disposable, breakable sterile inoculation loops or wooden stick (if using fecal manual extraction method described on page 2).
- Disposable screw cap tubes, 14 ml (if using fecal manual extraction method described on page 2).
- Extraction solution 1x (**Calprest Turbo** Extraction kit ref. 9299, if using fecal manual extraction method described on page 2).
- Microcentrifuge tubes (1-1.5 ml).
- Sensitive digital scale (40-150 mg) (if using fecal manual extraction method described on page 2).
- Vortex mixer.
- Shaker or roller (if using fecal manual extraction method described on page 2).
- Microcentrifuge (10,000 x g) (if using fecal manual extraction method described on page 2).
- Precision pipettes.

Equipment for turbidimetric measurements

Any turbidimetric automated analyzer; for the instrument-specific protocol, see additional Application note.

STORAGE AND STABILITY OF REAGENTS AND WORKING SOLUTIONS

- All reagents and working solutions must be stored at 2-8°C.
- The expiration date is printed on all component labels.
- Avoid exposure to high temperature, direct sunlight or extreme humidity.

PRECAUTIONS AND WARNINGS

- Patient specimens and samples should be handled in accordance with Good Laboratory Practice (GLP) using appropriate precautions and anyway as if capable of transmitting infections.
- For *in vitro* diagnostic use.
- The kit does not contain materials of human origin.
- A trained person in turbidimetric technique and autoanalyzer use is required.
- Do not exchange components from different lots. Good performance can only be guaranteed by using components from the same lot of **Calprest Turbo**.
- Do not mix screw caps from different reagent vials.
- The analyzer should be ready before performing any assay.
- If the result exceeds the measurement range, use the ready to use Extraction solution 1x (ref. 9299) to dilute the sample and repeat the assay again.
- Components must not be used after the expiry date printed on the labels.
- If the outer packaging is damaged, ensure that all components listed in the “Materials Provided” section are present and intact. In case of leakage of liquid from the bottles, the kit must not be used.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The solutions should be discarded in a proper container after testing following local regulations.

SPECIMEN COLLECTION

Random stool collection. Collect sufficient quantity of human stool samples. These specimens should be collected in clean and dry single-use containers (no preservatives or transport media). Loose or liquid stool specimens are acceptable as normalization to stool weight is part of the calculation of the result. Collection of stool specimen from diapers should be avoided unless the specimen can be taken from a portion which is not in contact with the diaper material.

Specimen requirements

1-5 g stool in a screw-cap clean vial. No preservative is necessary or indicated.

Specimen transport

Stool specimen should be received by the laboratory and extracted within 10 days of collection. Temperature during shipment should not exceed 30°C.

Specimen Storage

Stool samples at the laboratory should be stored at 2-8°C for up to 7 days before testing, if not tested, freeze the stool samples at -20°C.

PROCEDURE

Sample Preparation

Thaw frozen stool specimen and ensure that all reagents for extraction procedure reach room temperature (20-25°C). There are two alternative stool sample preparation methods, manual or by using **EasyCal Turbo** stool extraction device. Users can choose their preferred sample preparation method.

A. **EasyCal Turbo** Sample Preparation

1. For the collection/extraction procedure, please refer to **EasyCal Turbo** ref. 9302 instruction for use.
2. After the extraction procedure the fecal extract is ready to be tested. Insert the **EasyCal Turbo** tube (devoid of funnel and stick) directly into the turbidimetric analyzer or transfer the supernatant to an adaptor cup (not provided) for the analytical procedure.
3. If the extract is not analyzed immediately, it can be stored according to the instructions for use of **EasyCal Turbo** (Ref. 9302).

B. **Manual Sample Preparation**

Use a screw cap centrifugation tube for each sample to be tested. Label centrifugation tube with patient name or number.

1. Weigh (tare) the empty screw cap tube together with the inoculation loop or the wooden applicator stick.
2. Mix well the stool by means of the inoculation loop or wooden applicator stick in order to make the sample as homogenous as possible. Using the inoculation loop (or wooden applicator stick) take out approximately 50 mg (between 40-60 mg) of stool and place into the screw-cap-tube (weighted in point 1).
3. Weigh the tube with loop and feces inside and calculate net feces weight (between 40-60 mg).
4. Break off the loop handle or the wooden applicator stick, leaving the loop or wooden stick with feces and a 4-6 cm handle inside the screw cap tube.
5. Add ready to use extraction solution (**Calprest Turbo** Extraction kit ref. 9299, with final weight/volume ratio 1:100), *i.e.* 50 mg feces + 4.95 ml extraction solution as described in the table below. Fecal sample concentration in the extraction solution will be 10 mg/ml. Close the tube.

Feces (mg)	Extraction solution (ml)	Feces (mg)	Extraction solution (ml)
60	5.94	45	4.45
55	5.45	40	3.96
50	4.95		

- Shake/mix vigorously for 30 seconds by means of a vortex.
- Homogenize 25 ± 5 minutes on a shaker or roller. The inoculation loop or the wooden stick inside the tube will act as an agitator.
- Transfer the homogenate (1 ml) to a microcentrifuge tube and centrifuge for 20 minutes at $10,000 \times g$ at room temperature.
- Transfer 0.5 ml of the clear extract supernatant to a new microcentrifuge tube. Avoid contact with the pellet as aggregates or particles may lead to erroneous calprotectin values.
- The extracts can be tested immediately or stored at the following conditions:
 - Up to 72 hours at room temperature
 - Up to 14 days at $2-8^{\circ}\text{C}$
 - Up to 90 days at -20°C . Frozen samples can be subjected to up to 4 freeze/thaw cycles

ANALYTICAL PROCEDURE

- Ensure that all reagents reach room temperature ($20-25^{\circ}\text{C}$).
- Thaw frozen sample at room temperature.
- Calprest Turbo** Reagent 1 and **Calprest Turbo** Reagent 2 are ready to use solutions filled in bottles that can be loaded directly into the Euro-Turbo Analyzer (ref. 9194). Shall the test be used on a different turbidimetric analyzer, the content of the bottles has to be transferred into bottles proprietary of the used system.
- Calprest Turbo** Calibrator vials (Calibrator 0 to 5) and **Calprest Turbo** Control vials (Control 1 and Control 2) are ready to use.
- Setting up the calibration curve:
 - For the calibration of the curve use only **Calprest Turbo** Calibrator vials. The concentration is indicated on the label of each vial.
 - Once set up, calibration curve is stable up to 30 days on board, although it is highly recommended to calibrate the system every week.
 - Recalibrate the system when reagent lot is changed or when the controls are out of the assigned range given on the control label or certificate of analysis provided with the kit.
- Quality control:
 - For quality control use **Calprest Turbo** Control 1 and **Calprest Turbo** Control 2 vials. The concentration is indicated on the label of each vial.
 - Calprest Turbo** Controls should be assayed every day before running patient fecal samples extract to validate the calibration curve.
 - If the obtained results are out of the control lot-specific indicated range, the equipment, the reagents or the technique must be reviewed.
- Running the test:
 - The methodology can be applied to the most common turbidimetric analyzer. For the list of validated turbidimetric analyzers and the relative application notes, please contact Eurospital at the following email address:

marketing.diagnostics@eurospital.it

INTERPRETATION OF TEST RESULTS

- Negative results, calprotectin values lower than $50 \mu\text{g/g}$, indicate a normal presence of calprotectin in the analyzed stool samples.
- Positive results, calprotectin values higher than $120 \mu\text{g/g}$, indicate an abnormal presence of calprotectin in the analyzed stool samples. Adult patients with such results are identified with gastrointestinal inflammatory pathology that will require further investigation to diagnose the exact pathology.
- Borderline results, calprotectin values falling in the grey zone ($50-120 \mu\text{g/g}$), should be re-evaluated after 4-6 weeks.

The following chart summarizes the interpretation of results:

Calprotectin Concentration	Interpretation	Follow-Up
$< 50 \mu\text{g/g}$	Normal	None
$50 - 120 \mu\text{g/g}$	Borderline	Re-evaluate after 4-6 weeks
$> 120 \mu\text{g/g}$	Abnormal	Repeat as clinically indicated

Further evaluations including asymptomatic patients, as well as patients with IBS (to differentiate from IBD) and international studies confirmed the appropriateness of such values.

PERFORMANCE CHARACTERISTICS

Calprest Turbo cutoff value

- Calprotectin concentration values lower than $50 \mu\text{g/g}$ determined from stool samples are to be considered normal values, therefore not indicative of an inflammation of the gastrointestinal tract.
- Calprotectin concentration values higher than $120 \mu\text{g/g}$ determined from stool samples are indicative of an inflammation of the gastrointestinal tract.

PROZONE

Studies have been made up to a fecal calprotectin concentration of 8,000 µg/g and no false negative results have been observed. Studies using higher concentrations have not been carried out yet.

LIMITATIONS

- False negative results could occur in patients who have granulocytopenia due to bone marrow depression.
- Some patients who are taking Non-Steroid-Anti-Inflammatory-Drugs (NSAIDs) will have elevations in their fecal calprotectin levels.
- Patients with IBD fluctuate between active (inflammatory) and inactive stages of the disease. These stages must be considered when using the **Calprest Turbo** assay.
- The use of Proton Pump Inhibitors (PPIs), microscopic colitis and diverticular disease may also lead to elevated calprotectin level. Patients affected by untreated celiac disease may occasionally show elevated calprotectin value.
- Other intestinal impairments, including many gastrointestinal infections and colorectal cancer, can result in elevated levels of calprotectin. These specimens will test positive with the **Calprest Turbo** assay. Therefore, a diagnosis of active IBD cannot be established solely on the basis of a positive result with the **Calprest Turbo** assay.
- **Calprest Turbo** should be used for the determination of fecal calprotectin only from human stool samples. The kit is not intended to be used for different matrixes such as blood, serum, plasma, urine, cerebrospinal fluid, oral fluid, synovial fluid or empyema fluid. The use of the test on samples different than human stool has not been validated yet. The quality of the test depends also on the quality of the analyzed sample; proper fecal specimens must be obtained.
- Positive results indicate an abnormal presence of calprotectin in fecal samples. A positive result should be followed up by additional clinical investigations, such as colonoscopies and biopsies, in order to properly identify the nature of the inflammation and establish its extent.
- Fecal calprotectin is an indicator of neutrophilic presence in the stool and is not specific for IBD.
- If symptoms or situation still persist, pathology evaluation should be carried out with other alternative methods. Negative results do not exclude IBD with inflammation, some diseases such as coeliac sprue and/or microscopic colitis polyps that mainly involve mononuclear inflammation.
- Neonatal fecal calprotectin levels have been reported higher than those in normal children with a mean of 167 µg/g (range 22-860 µg/g).

INTERFERENCES AND CROSS REACTIVITY

In an evaluation to determine cross reactivity, no cross reactivity was found against other fecal markers occasionally present in feces, such as bovine serum albumin, ascorbic acid and hemin.

PERFORMANCE CHARACTERISTICS

A method comparison study was performed between Eurospital **Calprest Turbo** and the predicate Eurospital **Calprest NG** with 97.9% diagnostic specificity and 83.1% diagnostic sensibility (when borderlines are considered as positive).

Method Comparison

One hundred and three (103) human stool samples were extracted in parallel with both **EasyCal Turbo** and **EasyCal**. The extracted stool samples were tested in duplicate for calprotectin concentration determination using the **EasyCal Turbo** with **Calprest Turbo** and **EasyCal** with **Calprest® NG** assays, according to their respective package inserts. Obtained results are summarized in Table 1 and Table 2, as well as represented in figure 1.

Method comparison EasyCal Turbo with Calprest Turbo (y) versus EasyCal with Calprest NG (x)	
Passing-Bablok regression	EasyCal Turbo with Calprest Turbo = -7.407 + 1.138 EasyCal with Calprest NG
Slope (95% CI)	1.138 (1.021 to 1.246)
Y-intercept (95% CI)	-7.407 (-14.37 to -1.100)
Correlation - r	0.924
Bias at 120 µg/g (95% CI)	7.7% (-1.2% to 14.2%)

Table 1. Quantitative comparison between the results obtained using **EasyCal Turbo** with **Calprest Turbo** and **EasyCal** with **Calprest NG** assays. Passing-Bablok regression analysis was used.


Figure 1. Regression analysis using Passing-Bablok fit comparing the results obtained using *EasyCal Turbo* with *Calprest Turbo* and *EasyCal* with *Calprest NG* assays.

Qualitative agreement of <i>EasyCal Turbo</i> with <i>Calprest Turbo</i> versus <i>EasyCal</i> with <i>Calprest NG</i>					
		<i>EasyCal Turbo</i> with <i>Calprest Turbo</i>			
		BORD	NEG	POS	Total
<i>EasyCal</i> with <i>Calprest NG</i>	BORD	15	1	4	20
	NEG	1	12	0	13
	POS	2	0	68	70
	Total	18	13	72	103
Borderline as positive					
Negative agreement (95% CI)		92.3% (66.7 - 98.6%)			
Positive agreement (95% CI)		98.9% (94.0 - 99.8%)			
Total agreement (95% CI)		98.1% (93.2 - 99.5%)			
Borderline as negative					
Negative agreement (95% CI)		87.9% (72.7 - 95.2%)			
Positive agreement (95% CI)		97.1% (90.2 - 99.2%)			
Total agreement (95% CI)		94.2% (87.9 - 97.3%)			

Table 2. Qualitative comparison between the results obtained using *EasyCal Turbo* with *Calprest Turbo* and *EasyCal* with *Calprest NG* assays analyzing the data using the borderline range as positive and negative separately.

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
REFERENCES / BIBLIOGRAFIA / BIBLIOGRAPHIE / BIBLIOGRAFÍA / LITERATUR

1. Fagerhol M. K. et al: Calprotectin (The L1 leukocyte protein) in: Smith V. L. and Dedman J. R. (eds.): Stimulus response coupling: The role of intracellular calcium-binding proteins. CRC Press, Boca Raton 1990, p. 187-210.
2. Johne, B. et al: Functional and clinical aspects of the myelomonocyte protein calprotectin. *Mol Pathol*, 1997. 50(3): p. 113-23.
3. Roseth, A. G. et al: Assessment of the neutrophil dominating protein calprotectin in feces. A methodologic study. *Scand J Gastroenterol*, 1992. 27(9): p. 793-8.
4. Roseth, A. G. et al: Assessment of disease activity in ulcerative colitis by fecal calprotectin, a novel granulocyte marker protein. *Digestion*, 1997. 58(2): p. 176-80.
5. Roseth A. G. et al: Correlation between fecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. *Scand J Gastroenterol* 1999; 34:50-54.
6. Ton H. et al: B. Improved assay for fecal calprotectin. *Clinica Chimica Acta* 2000; 292:41-54.
7. Limburg P. J. et al: Fecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. *Am J Gastroenterol* 2000; 95:2831-2837.
8. Tibble J. et al: A simple method for assessing intestinal inflammation in Crohn's disease. *Gut* 2000; 47:506-513.
9. Burri E, Beglinger C. Faecal calprotectin -- a useful tool in the management of inflammatory bowel disease. *Swiss Med Wkly*. 2012;142:w13557. Published 2012 Apr 5.
10. Tibble J. A. et al: High prevalence of NSAID enteropathy as shown by a simple faecal test. *Gut* 1999; 45:362-366.
11. Montalto M. et al: Faecal calprotectin concentrations in untreated coeliac patients. *Scand J. Gastr.*, 2007; 42: 957-961.
12. Calcaterra V. et al: *Ann Nutr Metab* 2018; 73:177–183 “Serum Calprotectin Level in Children: Marker of Obesity and its Metabolic Complications”.

Legend / Legenda / Légende / Leyenda / Legende

	Consult electronic instructions for use / Consultare le istruzioni per l'uso elettroniche / Consulter la notice d'utilisation électronique / Consultense las instrucciones electrónicas de uso / Elektronischen Gebrauchsanweisung beachten		Sufficient for / Sufficiente per / Suffisant pour / Válido para / Ausreichend für
	Use by / Data di scadenza / Date de péremption / Fecha de caducidad / Verwendbar bis		Manufacturer / Fabbricante / Fabricant / Fabricante / Hersteller
	European Conformity / Conformità agli standard europei / Conformité aux normes européennes / Conformidad europea / Europäische Konformität		In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Producto sanitario para diagnóstico in vitro / In-vitro-Diagnostikum
	Temperature limitation / Limiti di temperatura / Limites de température / Limitación de temperatura / Temperaturbegrenzung		Reagent / Reagente / Réactif / Reactivo / Reagens
	Catalogue number / Numero di Catalogo / Référence du catalogue / Número de catálogo / Bestell Inummer		Control / Controllo / Contrôle / Control / Kontrolle
	Keep dry / Conservare in luogo asciutto / Garder au sec / Almacenar en lugar seco / Bleib trocken		Calibrator / Calibratore / Calibrateur / Calibrador / Kalibrator
	Keep out of the sunlight / Tenere lontano dalla luce del sole / Gardez hors de la lumière du soleil / Manténgase fuera de la luz del sol / Vor dem Sonnenlicht schützen		Batch number / Numero di lotto / Numéro de lot / Número de lote / Chargennummer

English
Italiano
Français
Español
Deutsch

**ANNEX INFORMATION OF THE VALIDATED COMPATIBLE TURBIDIMETRIC INSTRUMENTS
INFORMAZIONI DEGLI STRUMENTI TURBIDIMETRICI COMPATIBILI VALIDATI
ANNEXE INFORMATIONS SUR LES INSTRUMENTS TURBIDIMÉTRIQUES COMPATIBLES VALIDÉS
ANEXO INFORMACIÓN DE LOS INSTRUMENTOS TURBIDIMÉTRICOS COMPATIBLES VALIDADOS
ANHANG INFORMATIONEN DER VALIDIERTEN KOMPATIBELN TURBIDIMETRISCHEN INSTRUMENTE**

Manufacturer	Model
Eurospital	Eu-Turbo Analyzer
Abbott	Alinity; Architect c1000 / c4000 / c8000
Awareness	Chemwell-T
Beckman Coulter	AU480 / AU680
Biosystems	A15
DyaSis	Respons 910
Mindray	BS-200 / BS-200E
Ortho	Vitros 5600
Roche	cobas c111, c501, c702
Siemens	ADVIA 1800/2400; Atellica
Tokyo Boeki	Biolis i24/i50
Binding site	Optilite
Erba	XL-180 / XL-200
Tecom	TC220

Calprest Turbo

Ref. 9300, 200 tests.

Ref. 9300, 200 test.

Ref. 9300, 200 tests.

Ref. 9300, 200 ensayos.

Ref. 9300, 200 tests.

Eurospital SpA

Via Flavia 122, 34147 Trieste, Italia
Tel. +39 040 8997.1 - Fax +39 040 280944
www.eurospital.com - info@eurospital.it