DiabeGen - 1° step

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- genetics consortium families*. *Diabetes*, 2008. **57**(4): p. 1084-92.
9. Noble, J.A. and A.M. Valdes, *Genetics of the HLA region in the prediction of type 1 diabetes*. *Curr Diab Rep*, 2011. **11**(6): p. 533-42.
10. Van der Auwera, B., et al., *DRB1*0403 protects against IDDM in Caucasians with the high-risk heterozygous DQA1*0301-DQB1*0302/DQA1*0501-DQB1*0201 genotype*. *Belgian Diabetes Registry*. *Diabetes*, 1995. **44**(5): p. 527-30.
11. Awata, T., et al., *Genetic analysis of HLA class II alleles and susceptibility to type 1 (insulin-dependent) diabetes mellitus in Japanese subjects*. *Diabetologia*, 1992. **35**(5): p. 419-24.
12. Rewers, A., et al., *Ethnic differences in the associations between the HLA-DRB1*04 subtypes and type 1 diabetes*. *Ann NY Acad Sci*, 2003. **1005**: p. 301-9.
13. Thomson, G., et al., *Relative predispositional effects of HLA class II DRB1-DQB1 haplotypes and genotypes on type 1 diabetes: a meta-analysis*. *Tissue Antigens*, 2007. **70**(2): p. 110-27.
14. Nepom, G.T., *A unified hypothesis for the complex genetics of HLA associations with IDDM*. *Diabetes*, 1990. **39**(10): p. 1153-7.
15. Pugliese, A., et al., *HLA-DQB1*0602 is associated with dominant protection from diabetes even among islet cell antibody-positive first-degree relatives of patients with IDDM*. *Diabetes*, 1995. **44**(6): p. 608-13.
16. Baisch, J.M., et al., *Analysis of HLA-DQ genotypes and susceptibility in insulin-dependent diabetes mellitus*. *N Engl J Med*, 1990. **322**(26): p. 1836-41.
17. Sacks, D.B., et al., *Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus*. *Diabetes Care*, 2011. **34**(6): p. e61-99.
18. Hitman, G.A., *The message for MODY*. *Diabet Med*, 2011. **28**(9): p. 1009.
19. Atkinson, M.A. and E.H. Leiter, *The NOD mouse model of type 1 diabetes: as good as it gets?* *Nat Med*, 1999. **5**(6): p. 601-4.
20. Nanto-Salonen, K., et al., *Nasal insulin to prevent type 1 diabetes in children with HLA genotypes and autoantibodies conferring increased risk of disease: a double-blind, randomised controlled trial*. *Lancet*, 2008. **372**(9651): p. 1746-55.
21. Todd, J.A., M. Knip, and C. Mathieu, *Strategies for the prevention of autoimmune type 1 diabetes*. *Diabet Med*, 2011. **28**(10): p. 1141-3.

-
22. Barker, J.M., *Clinical review: Type 1 diabetes-associated autoimmunity: natural history, genetic associations, and screening*. J Clin Endocrinol Metab, 2006. **91**(4): p. 1210-7.

DiabeGen - 1° step

Cod 9192, 12 pazienti

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Rev. 0

Eurospital 

Eurospital SpA

Via Flavia 122 - 34147 Trieste - Italia

Tel. +39 040 89971 - Fax +39 040 280944

www.eurospital.com - info@eurospital.com

DiabeGen – 1^o step

For in vitro diagnostic use only.

Intended use

Molecular biology kit for the detection of Class II HLA alleles predisposing to Type I Diabetes Mellitus (or Insulin Dependent Diabetes Mellitus T1D) in human DNA from peripheral blood by means of RT-PCR.

DiabeGen 1^o step MUST NOT be used for tissue typing

Principle of the method

Type I Diabetes is a complex and multifactorial disease, is the most common form of diabetes among children and young adults in the Caucasian population with prevalence of 0.4% [1]. This incidence varies geographically and by age: prevalence in China and Venezuela is 0.1 / 100,000 cases, while it reaches 36/100.000 in Sardinia and Finland [1]. The disease is characterized by the subtotal (> 80%) or total destruction of the β -cells in the pancreas islets caused by a lymphocytic insulinitis (autoantibodies) [2].

Although more than 40 gene loci [3] have been associated with Type I diabetes (T1D), the HLA region, due to the number of genes and their polymorphisms, contributes from 30 up to 50% to the genetic susceptibility to this disease [4, 5]. The highest probability to develop T1D is associated with the presence of the DRB1*04-DQA1*03-DQB1*03:02 and DRB1*03-DQA1*05-DQB1*02 haplotypes in heterozygosis (serologically known as DR4-DQ8 and DR3-DQ2 [2,4-7]), followed by homozygote DR4-DQ8 and homozygote DR3-DQ2. The estimated risk of developing T1D in children in the general population having the HLA-DR3-DQ2/DR4-DQ8 genotype is approx. 1/20 cases. This condition becomes 1/5 cases if other cases of T1D are present in the family [5]. Although the DR4-DQ8/DR3-DQ2 genotype is associated with the highest risk in the predisposition to T1D, it is possible to differentiate the risk associated with HLA DR/DQ genotypes, which presence leads to either high, medium or absence of risk (protection).

It is well known that in the Caucasian population there is a hierarchy in the risk stratification for T1D which is primarily influenced by the allelic variants of DRB1*04 when it is associated with DQB1*03:02 (DQ8) [8, 9]. A high risk in the onset of the disease is associated with the presence of DRB1*04:05 and DRB1*04:01, followed by an medium risk associated with DRB1*04:02 and DRB1*04:04 [4-6].

Conversely, the presence of DRB1*04:03 or DRB1*04:06 confers protection against the disease [8, 10-13]. For the detection of DRB1*04 allele variants (high resolution), please refer to the DiabeGen II^o step kit package insert.

In addition to DRB1*04 allelic variants, the risk stratification in the onset of T1D is also influenced by HLA-DQB1 profile. Several studies [4, 9, 14] have demonstrated that the presence of DQB1*06:02 may negatively influence the predisposition to Type I diabetes, even if autoantibodies associated with diabetes Type I and / or risk alleles are present [15]. Therefore, the presence of DQB1*06:02 allele seems to confer protection against T1D [14, 16], even when the subject carries a at-risk haplotype (DR4-DQ8 or DR3-DQ2) on the other chromosome.

DiabeGen 1^o step allows to identify the alleles predisposing T1D, namely DR4B1*04, DRB1*03, DQB1*03:02, DQB1*02 and the presence of the protective allele DQB1*06:02.

Mix 6 (status DR4) contains primers and probes for the identification of all DRB1* alleles (DRB1*01, 03, 07, 08, 09, 10, 11, 12, 13, 14, 15 and 16) with the only exception of DRB1*04.

Mix 6 in association with Mix 1 is used to define the status of DRB1*04, and in association with may give three combinations:

1. Mix 1 and 6 are positive: the subject is heterozygote for DRB1*04 (1 copy on one chromosome)
2. Mix 1 is positive, Mix 6 is negative: the subject carries homozygote DRB1*04 (2 copies, one for both chromosomes).
3. Mix 1 is negative, Mix 6 is positive: this will occur in all subjects not carrying DRB1*04 (e.g. Case 5 reported in the table "DiabeGen - 1^o step Combination Profile").

The information above is useful for the users of the DiabeGen Step II kit, since the Eurospital DiabeGen 1^o step and II^o step system is able to provide a unambiguous information about the status of the DRB1*04 allele only if the sample has one copy of such allele (heterozygote DRB1*04). If DRB1*04 is present on both chromosomes (homozygosis, Mix 1 positive and Mix 6 negative), high resolution identification of this allele should be done by another technique, e.g. sequencing.

Preventive Strategies

Differently than other genetic HLA-linked diseases, e.g. celiac

disease, the knowledge of the genetic pattern of a subject potentially at risk of developing T1D does not influence the physician's approach to the disorder. As matter of facts, while in case of Celiac Disease the physician may decide whether the patient should or should not be undergo a gluten free diet, a similar approach is not possible in case of T1D because trigger or triggers of the disease are still unknown.

The scientific community has different opinions about the use of genetic markers in the diagnosis of T1D. As an example, the authors of "Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes (Diabetes Care, 2011. 34(6): 61-99) [17] state that HLA-DR/DQ typing in T1D might be useful only in patients already positive to T1D auto-antibodies to define the presence or absence of specific haplotypes, because some haplotypes predispose to the disease, others provide significant delay or even protection.

The authors, therefore, recommend that genetic testing should not be recommended for routine clinical diagnosis of Type I diabetes. On the other hand, some authors just affirm the opposite, namely that knowing the HLA pattern of a subject potentially at risk of T1D allows, firstly, to differentiate it from other types of diabetes such as Maturity-Onset Diabetes of the Young (MODY) and Type II diabetes [18] (American Diabetes Association, <http://www.diabetes.org>), secondly, to inform patients about the very high probability to develop the disease (1/5 cases) if other people in their family had been affected by diabetes.

Studies on animal models [19] have also shown how to prevent T1D autoimmunity is easier before the onset of the disease, rather than trying to prevent it when the autoimmune pattern is already known. If this were true for human beings too, the proposed DiabeGen Step II genetic test could represent an effective clinical tool to identify the patients at risk, mainly children and adolescents, before the onset of the disease this allowing the physicians to plan their clinical follow-up properly.

In agreement with this approach, a recent publication [20] states that genetic screening of T1D can avoid unnecessary suffering to children/adolescents, as only 4% of the subjects enrolled in this clinical trial, with known HLA pattern, were affected by typical diabetic ketoacidosis against 15-20% of children diagnosed in the general population [20, 21]. This aspect, namely how an early preventive diagnosis can prevent from hospitalization due to ketoacidosis and at least partially preserve patient β cells [22], is shared and also reported by the authors reported above [17], for whom, the genetic test for T1D has a moderate value.

Description of the DiabeGen I^o step kit

The DiabeGen I^o step kit detection system requires extracted and purified DNA samples. Each sample is then amplified by Real Time PCR reactions, according to the scheme shown in Table 1.

Supplied Materials

(quantity sufficient for 12 tests)

| | |
|------------------------------|-------------|
| Reagent A | 2 x 0,4 ml |
| Reagent B | 2 x 1,5 ml |
| 8 tube strips with Mix 1 - 7 | 2 x 1 strip |
| TE 1X Buffer | 2 x 1.5 ml |

Composition of supplied materials/reagents

1. Buffer TE 1X

2 vials of 2 ml, each containing 1,5 ml buffer (Tris/EDTA) for resuspension of primers and lyophilized probes. Ready-to-use

2. Reagent A

2 vials (blue cap) each containing 0,4 ml RT PCR Taq Polymerase. Ready-to-use

3. Reagent B

2 vials (green cap) each containing 1,5 ml reaction buffer for Master Mix. Ready-to-use


4. Strip with mix of primers and probes (lyophilised) for mixes 1-7

- **Tube 1** (Mix 1), marked with a GREEN dot.
- **Tube 7**, EMPTY and marked with a White dot (see Table 1).
- **Mix 1 to 5** contain primer and probes with different specificities allowing simultaneous amplification of aspecific HLA allele to the predisposition to T1D and human β -Globin, which acts as an internal control and guarantees the reliability of the obtained results.
- **Mix 6** (status DR4), contains primers and probes for identification of all DRB1* alleles (DRB1*01, 03, 7, 08, 09, 10, 11, 12, 13 14, 15 and 16) with the only exception of DRB1*04.

Notice: Mix 7 contains primers and probes for the amplification only of human β -Globin, which is used as No Template Control (NTC) to detect possible contaminations by human genomic DNA in the Master Mix. Due to the sensitivity of the RT system and according to the scheme shown in Table 1, Mix 7 is located in the last tube of the 8 well strip: this means that the seventh well marked with a white dot, is EMPTY.

The mixes appear as a lyophilised pink pellet (primers and probes) laying on the bottom of the tubes. Once dissolved, the

solution has a light pink colour. Each strip allows testing 6 samples.

| Table 1 Strip with GREEN DOT | | |
|---------------------------------|-----------------------------------|------------------------------------|
| Mix | HLA class II Allele (FAM channel) | Internal Control (VIC/HEX channel) |
| 1 green dot | DRB1*04 | β globin |
| 2 | DRB1*03 | β globin |
| 3 | DQB1*03:02 | β globin |
| 4 | DQB1*02 | β globin |
| 5 | DQB1*06:02 | β globin |
| 6 | STATUS DR4 | β globin |
| White dot | EMPTY | |
| 7 | NO TARGET | β globin |

Required but not supplied instruments and materials.

1. Laboratory and equipments to run Real Time PCR.
2. Eppendorf[®] like sterile 1,5 ml tubes.
3. Adjustment-volume pipettes (10, 20, 200 and 1000 μ l).
4. Multi-channel (8 tips) adjustment-volume (5 to 50 μ l) pipettes.
5. Disposable sterile tips with filter for 10,20,200,1000 μ l pipettes.
6. 96 well microplates for Real Time PCR with optical grade caps and adhesive film or 8 tube strips with optical grade caps.
7. Sealing device
8. Vortex
9. Centrifuge for Eppendorf[®]-like-tubes
10. Centrifuge for 96 well microplate.
11. Centrifuge for 8 tube strips.
12. Thermal Cycler Real Time.

Performance criteria

The method allows defining the presence/absence of alleles listed in Table 1, which, as reported in the literature, are involved in the genetic predisposition or protection to T1D.

Sensitivity and specificity

25 DNA samples coming from HLA Reference Centres have been tested. Results obtained by using DiabeGen I[®] step kit have shown 100% sensitivity and specificity for the detected alleles.

Storage

All reagents must be stored at -2-8° C.

Stability after first opening

The mixes with primers and probes, once dissolved, are stable for 60 days if stored at -2-8° C. Other reagents, if used according to the advices reported in "General Warnings", are stable until the expiry date reported on the label.

Stability during transportation

An accelerate stability study showed that all reagents were stable after storage at 37°C for 72 hours.

Sampling

The DiabeGen I[®] step requires EDTA whole blood collected in adequate Vacutainer[®] tubes (eg.K3 and Molecular Biology). For sample storage, refer to the indication provided by the collection tube manufacturer.

Sample preparation

The sample is genomic DNA extracted from EDTA whole blood. DNA can be extracted by means of the most common DNA extraction techniques normally used by molecular biology laboratories. For optimal results, DNA concentration must be between 10 and 100 μ g/ml, which corresponds to Ct values of 22-25 for human β -Globin. The acceptability limit of the test allows the use of samples with Ct values in the range of 21-28 for human β -Globin falling. For Ct > 28 the test result is not guaranteed. If a sample dilution is required, use sterile water only.

DiabeGen I[®]Step has been optimised by using reagents included in the Eurospital Eu-Gen Estrazione[®] kit (code 9132).The use of other extraction kits requires that each laboratory should check its own procedures in terms of DNA concentration, stability and storage features.

Sample storage

Purified DNA, if not immediately submitted to amplification, can be stored at 2-8° C for 3 (three) months, or at -20°C for longer times. In this case, DNA should be stored in aliquots in order to avoid repeated freezing and thawing procedures which might damage the sample. For stability and storage conditions of the sample, refer to instructions included in the relevant extraction kits.

Re-use of DNA stored at 2-8° C for three months has been tested only by using reagents included in the Eurospital Eu-Gen Estrazione[®] Kit and appropriate tubes. Data on DiabeGen I[®] step performances when using DNA samples either stored

for longer periods and/or extracted by different methods than those reported above are not available.

DNA Amplification Procedure

Lyophilised Mix Resuspension

1. Place the 8 tube strip marked with GREEN dot in a support (not supplied) having the hinge at the far left (Figure 1); doing so the GREEN dot will be located on the first tube at the top (Mix 1, DRB1*04 allele). Each tube strip contains material for 6 amplifications, therefore for 6 samples.

2. Open all tube caps. Transfer 40 μ l from the TE 1X vial (resuspension buffer) in all tubes containing the lyophilised mix making sure to use a new pipette tip every time (Figure 1). Mix a couple of time in order to resuspend the lyophilised primers and probes contained in each tube.


NOTE: the tube marked with a WHITE dot MUST not be open as it is empty.

3. Make sure that resuspended primer solution lays on the bottom of each tube. On the contrary, close the tube caps and shortly centrifuge the strips.

Dispensing Mixes for Amplification

For each samples to be tested, take 5 μ l of resuspended primers and probes solution from the strip and put them column wise in a 96 wells plate (or into a 8 well strip for Real Time PCR, if such a support is used) as shown in **Figure 1**.

Figure 1

Preparation of Master Mix

1. Due to the high viscosity of the solutions involved in the preparation of Master Mix, before using Reagents A and B vials it is advisable to shake and then shortly centrifuge them to eliminate possible drops inside the caps. Carefully transfer the required volumes to prepare the mix.

2. For each sample, transfer 60 μ l of Reagent A (blue cap vial) and 20 μ l of Reagent B (green cap vial) in an Eppendorf-like sterile tube. Shake and shortly spin the tubes containing the Master Mix. (**Fig.2**)

3. The method allows to assay up to a maximum of 12 samples if a 96 well RT thermal cycler is used.

4. Take 10 μ l of each Master Mix and transfer it into the well containing Mix 7 (corresponding to No Template Control NTC) as it contains primers and probes for human β -Globin only (**Figure 2**, Mix 7, yellow highlighted).

5. **Note:** It is advisable to close wells containing Mix 7, until the plate is fully dispensed in order to minimise possible contaminations (aerosol and environmental contaminations, etc.).

6. For each sample to be tested, transfer 20 μ l of DNA in each Master Mix tube. Shake and shortly centrifuge all the tubes. The Master Mix with DNA is ready.


7. For each sample to be tested, transfer 10 μ l of Master Mix with DNA in each well containing the primers and probes (**Figure 3**, Mix 1 to 6, yellow highlighted).

To properly mix the solutions (primer and master with DNA), proceed as indicated below:

- If Vortex and microplate centrifuges are available: close the plate by means of the adhesive film and seal it by using the sealing device (not provided), shortly Vortex (5-10 sec), low speed if possible, then shortly centrifuge for removing possible bubbles and drops.
- Other cases: it is advisable to repeatedly pipette the content of each well of the plate by means of an 8 tip multichannel pipette, then close the plate by means of the adhesive film and seal each well by using the sealing device (not provided).

The same procedure must be used also

Figure 2


| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 |
|---|-------|---|---|---|---|---|---|---|
| A | Mix 1 | | | | | | | |
| B | Mix 2 | | | | | | | |
| C | Mix 3 | | | | | | | |
| D | Mix 4 | | | | | | | |
| E | Mix 5 | | | | | | | |
| F | Mix 6 | | | | | | | |
| G | | | | | | | | |
| H | Mix 7 | | | | | | | |

Transfer 10 µl
in Mix 7
ONLY

Sample 1

Mix 7 = Negative Control
No Template Control (NTC)

Figure 3

| | 1 | 2 | 3 | 4 | 5 |
|---|-------|---|---|---|---|
| A | Mix 1 | | | | |
| B | Mix 2 | | | | |
| C | Mix 3 | | | | |
| D | Mix 4 | | | | |
| E | Mix 5 | | | | |
| F | Mix 6 | | | | |
| G | | | | | |
| H | Mix 7 | | | | |

Sample 1

when optical strip and/or caps are used.

Real Time PCR Amplification

A 96 wells RT thermal cycler allows to run 12 samples each time, one per each column of the plate. If RT PCR 8 tube strips with optical grade caps are used, 12 strips, one for each sample, will be used.

Procedure

1. Start the software and follow the program specifications of the relevant thermal cycler. DiabeGen I^ostep kit has been validated with the following thermal cyclers: AB 7300, 7500 and 7500 FAST, Bio-Rad CFX96 Touch and Connect. As an example, when using Applied Biosystems thermal cyclers, set

the *fluorophors* as described below:

| Detector name | Reporter | Quencher | Colour |
|---------------------------|----------|----------|--------|
| DiabeGen Allele | FAM | NONE | RED |
| DiabeGen β -globina | VIC * | NONE | BLACK |

(*) NOTE for MIX 7: select VIC channel on Applied Biosystems thermal cyclers and HEX channel on Bio-Rad thermal cyclers

- If requested, set *ROX* as *Passive Reference*.
- NOTE: wells containing Mix 7 are used as NTC: please mark them on the plate accordingly.
- Set the following thermal protocol

| Temperature | Time | Cycles |
|-------------|--------|--------|
| 95°C | 10 min | 1 |
| 95°C | 15 sec | 37 |
| 60°C | 90 sec | |

- Set optical reading of the instrument on Step at 60°C (box highlighted in yellow in the table above).
- IMPORTANT NOTE for Bio-RadCFX96 Touch and Connect users only**

The ramp of the thermal protocol described in table above must be set on 1.5°C/sec for each step.

- Set the reaction volume at 15 μ l.
- Start up the reaction.

Validation of the expected results.

Once the reaction protocol has been completed, set the following parameters for reading results:

- Select plate/ all samples.
- Set the following values:

| Applied Biosystem 7300, 7500 and 7500 Fast | |
|--|------|
| Baseline: Start cycle-End cycle | 7-17 |
| Threshold DiabeGen Allele (FAM) | 0,2 |
| Threshold DiabeGen β -Globin (VIC) | |

| Bio-Rad CFX96 Touch and Connect | |
|--|------|
| Baseline: Start cycle-End cycle | 7-17 |
| Threshold DiabeGen Allele (FAM) | 50 |
| Threshold DiabeGen β -Globin (HEX) | 80 |

- Analyse the run.

Test validation

A test is valid, if the following conditions are satisfied:

- For each sample, check that mix 7 (negative control or NTC) is negative to DiabeGen β -Globin reading filter (VIC or HEX channel). Due to the sensitivity of the detection method (Real Time PCR), positive signals showing Ct values higher than 34 (Ct \geq 34) might be tolerated, as they might result from environmental contaminations or out-of-control events. These signals do not affect the final result. Conversely, the presence of a fluorescent signal showing a Ct value lower than 34 (Ct < 34) indicates a reagent contamination which might affect the result obtained with that Master Mix. In this case, the test must be repeated.

- For each sample, check that all Mix from 1 to 6 show a positive values for VIC or HEX channel with Ct values in the 21 ÷ 28 range. The presence of the reaction internal control (human β -Globin) is very important in case of negative sample amplification in a given mix. If the Ct values obtained for the human β -Globin are higher than the upper reference value (Ct=28), some diabetes specific alleles could have not been amplified, or amplified beyond the deemed acceptable limit of 35 cycles. In this case, it is recommended to repeat the test by increasing the DNA concentration. If the Ct value obtained for human β -Globin is lower than 21, it is recommended to repeat the test by diluting the DNA sample. Please note that, theoretically, a 1:10 dilution of the DNA sample should provide an increased Ct value (approx. 3,3 units).

- Read Ct values for each well (FAM channel, DiabeGen Allele). Values lower than 17 (Ct < 17) for FAM and VIC/HEX channels are not possible (check components and the amplification profile) and should not be considered.

- Read Ct values for each well (FAM channel, DiabeGen Allele). Ct values between 18-35 indicate a positive result for a given allele in that specific mix.

Results interpretation

- If the acceptability limits reported in the previous section "Test validation" of the test are met, report the positive results for the FAM channel in the provided form "DiabeGen I^o step-Results".

- The combination profile obtained for each sample must be crosschecked with those present in the table "DiabeGen I^o step- Combination profile". Report the relevant combination in the last column of the form "DiabeGen I^o step- Results".

- Finally, by using the form "DiabeGen I^o step - analysis of results" each reported combination is associated with the presence/absence of predisposing alleles and/or DQB1*06:02 allele. In the same form, the column "Possible additional steps

and comments" indicate possible further steps to better characterize the sample (e.g. high resolution typing of DRB1*04 by means of the DiabeGen II^o Step kit)

4. As an example, if a sample belongs to Combination 1, the sample will be positive for DRB1*04, DRB1*03, DQB1*03:02 and DQB1*02 alleles. These allele define the DR4-DQ8/DR3-DQ2 genotype showing the highest predisposing risk to T1D. The result indicates that predisposing alleles to T1D are present, while the protective allele DQB1*06:02 is not present. As an additional step, it is advised to test the sample by means of the DiabeGen II^o Step kit to check the status of the DRB1*04 allele at high resolution. If the subject should carry either DRB1*04:03 or 04:06 allele, the patient would be anyhow protected against T1D [8, 10, 11]. If, on the contrary, if the sample were DRB1*04:05 positive, the patient would have an increased predisposing risk to the disease [8].

Precautions

- 1.** Wait for the reagent to reach room temperature before starting the test.
- 2.** Strictly follow the procedure. In case of accident or disease, seek immediately for medical help and show the physician the container, the label and, if available, the Material Safety Data Sheet.
- 3.** Warning: due to the viscosity of solutions, in order to obtain the most appropriate test performance it is mandatory to mix 2-3 times the reaction mixes by the 8 tip multi-channel pipette. New pipette tips must be used for each new well line.
- 4.** Use Vortex and plate/strip centrifuges if available.
- 5.** Warning: it is very important not to use reagents belonging to different batches.
- 6.** Always wear double gloves when running the test. If any doubt of possible contamination with reagents and/or sample arises, immediately change gloves.
- 7.** If not working under a hood, always wear a mask to avoid test tube contamination.
- 8.** It is recommended to perform each step of the test (sample and reagent preparation, test procedure and reading) in separate and defined rooms in order to minimise test tube contamination.
- 9.** Before running the test for the first time, it is also recommended to HLA type any user who might be working with the system to be able to identify any possible contamination during the test.
- 10.** Once the amplification is completed, please dispose of any tube and/or plate containing the amplified materials without ever opening them.
- 11.** The Mixes contain DNA specific sequences, which

possible hazard has never been fully tested.

12. The waste materials generated by the use of the kit should be disposed of according to National Directives and Laws regulating the treatment of wastes coming from chemicals biological. In particular, liquid wastes generated by sample extraction procedures and other waste materials (e.g. pipette tips) should also be treated as potentially infectious material and inactivated before their disposal. Avoid contact between waste sample extraction materials and bleach.

Limitations of the test.

- 1.** DiabeGen I^o step kit requires the use of EDTA whole blood collected in adequate Vacutainer[®] tubes (eg. K3 and Molecular Biology). For sample storage refer to the indication provided by the tubes manufacturer.
- 2.** DiabeGen I^o step has been optimised by using the extraction reagents included in the Eurospital[®] Eu-Gen Estrazione[®] kit (code 9132). The usage of other extraction kits requires that the laboratory should check out its own procedures in terms of DNA concentration, stability and storage features.
- 3.** DNA should be tested as soon as possible after the extraction.
- 4.** Purified DNA, if not use immediately tested, can be stored at 2-8°C for three (3) months, or frozen at -20°C for longer times. In this case, it is recommended to aliquot the DNA sample in order to avoid possible problems due to continuous thawing and freezing steps which might lead to degradation.
- 5.** DNA concentration of each sample must be in the range between 10 and 100 µg/ml.
- 6.** DiabeGen I^o step I kit has been validated on the following RT thermal cyclers: AB 7300, 7500 and 7500 FAST, Bio-Rad CFX96 Touch and Connect.

General advices

- 1.** Always use sterile tips with filter.
- 2.** Take any precaution to avoid reagents contamination.
- 3.** Expose lyophilised and resuspended mixes to light just for the time needed to perform the test.
- 4.** Place unused strips straight back in the provided aluminium minigrip sachet and make sure a desiccant inside is inside.
- 5.** Do not touch pipette tips by hands.
- 6.** Always wear protection gloves to avoid reagent contamination.
- 7.** Use optical grade tubes, caps and/or adhesive films only.
- 8.** Avoid mixing reagents from different batches.
- 9.** DiabeGen I^o step should be used ONLY according to the

purposes reported in the paragraph "Intended Use" and stated in the paragraph "Performance Criteria".

Troubleshooting

No signal

1. Wrong setting of the reading filters.
2. Dispensing mistakes or lack of reagents.
3. Sample Inhibitory effects: insufficient DNA extraction and purification processes.
4. Wrong storage of the kit.
5. Check thermal cycler performance.

False positive results.

DNA samples, mostly extracted by means of magnetic beads, might rarely interfere with the kit reagents (Reagent A master mix). This kind of interference creates positive signals (FAM channel) in wells that should be negative, giving a liner graph (rather than the typical exponential profile expected for Real Time PCR positive results). The graph, when crossing the threshold, gives Ct values. If this problem occurs, check the raw fluorescence signal recorded by the instrument (in the Applied Biosystems thermal cyclers it is named "component" in) for the well/sample which FAM and possible *passive reference* (ROX) are altered. If this scenario occurs, the result is to be considered as false positive.

Signal strength too low

1. Deterioration of the detection system due to incorrect handling and/or storage of the reagents.
2. Very low DNA concentration of DNA in the tested sample.
3. Incorrect mix/mixes resuspension.
4. Air bubbles trapped inside the amplification tubes.

Literature

1. Karvonen, M., et al., *Incidence of childhood type 1 diabetes worldwide*. *Diabetes Mondiale (DiaMond) Project Group*. *Diabetes Care*, 2000. **23**(10): p. 1516-26.
2. R. Arcari, N.G., A. Lezo, D. Boscolo, P.Cavallo Perin, *Eziopatogenesi del Diabete Mellito di Tipo 1*. *Caleidoscopio Italiano*, 1998.
3. Barrett, J.C., et al., *Genome-wide association study and meta-analysis find that over 40 loci affect risk of type 1 diabetes*. *Nat Genet*, 2009. **41**(6): p. 703-7.
4. Noble, J.A. and H.A. Erlich, *Genetics of type 1 diabetes*. *Cold Spring Harb Perspect Med*, 2012. **2**(1): p. a007732.
5. Steck, A.K. and M.J. Rewers, *Genetics of type 1 diabetes*. *Clin Chem*, 2011. **57**(2): p. 176-85.
6. Pociot, F. and M.F. McDermott, *Genetics of type 1 diabetes mellitus*. *Genes Immun*, 2002. **3**(5): p. 235-49.
7. Rani, R., A. Sood, and R. Goswami, *Molecular basis of*

predisposition to develop type 1 diabetes mellitus in North Indians. *Tissue Antigens*, 2004. **64**(2): p. 145-55.

8. Erlich, H., et al., *HLA DR-DQ haplotypes and genotypes and type 1 diabetes risk: analysis of the type 1 diabetes genetics consortium families*. *Diabetes*, 2008. **57**(4): p. 1084-92.
9. Noble, J.A. and A.M. Valdes, *Genetics of the HLA region in the prediction of type 1 diabetes*. *Curr Diab Rep*, 2011. **11**(6): p. 533-42.
10. Van der Auwera, B., et al., *DRB1*0403 protects against IDDM in Caucasians with the high-risk heterozygous DQA1*0301-DQB1*0302/DQA1*0501-DQB1*0201 genotype*. *Belgian Diabetes Registry*. *Diabetes*, 1995. **44**(5): p. 527-30.
11. Awata, T., et al., *Genetic analysis of HLA class II alleles and susceptibility to type 1 (insulin-dependent) diabetes mellitus in Japanese subjects*. *Diabetologia*, 1992. **35**(5): p. 419-24.
12. Rewers, A., et al., *Ethnic differences in the associations between the HLA-DRB1*04 subtypes and type 1 diabetes*. *Ann NY Acad Sci*, 2003. **1005**: p. 301-9.
13. Thomson, G., et al., *Relative predispositional effects of HLA class II DRB1-DQB1 haplotypes and genotypes on type 1 diabetes: a meta-analysis*. *Tissue Antigens*, 2007. **70**(2): p. 110-27.
14. Nepom, G.T., *A unified hypothesis for the complex genetics of HLA associations with IDDM*. *Diabetes*, 1990. **39**(10): p. 1153-7.
15. Pugliese, A., et al., *HLA-DQB1*0602 is associated with dominant protection from diabetes even among islet cell antibody-positive first-degree relatives of patients with IDDM*. *Diabetes*, 1995. **44**(6): p. 608-13.
16. Baisch, J.M., et al., *Analysis of HLA-DQ genotypes and susceptibility in insulin-dependent diabetes mellitus*. *N Engl J Med*, 1990. **322**(26): p. 1836-41.
17. Sacks, D.B., et al., *Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus*. *Diabetes Care*, 2011. **34**(6): p. e61-99.
18. Hitman, G.A., *The message for MODY*. *Diabet Med*, 2011. **28**(9): p. 1009.
19. Atkinson, M.A. and E.H. Leiter, *The NOD mouse model of type 1 diabetes: as good as it gets?* *Nat Med*, 1999. **5**(6): p. 601-4.
20. Nanto-Salonen, K., et al., *Nasal insulin to prevent type 1 diabetes in children with HLA genotypes and autoantibodies conferring increased risk of disease: a double-blind, randomised controlled trial*. *Lancet*, 2008. **372**(9651): p. 1746-55.

-
21. Todd, J.A., M. Knip, and C. Mathieu, *Strategies for the prevention of autoimmune type 1 diabetes*. Diabet Med, 2011. **28**(10): p. 1141-3.
 22. Barker, J.M., *Clinical review: Type 1 diabetes-associated autoimmunity: natural history, genetic associations, and screening*. J Clin Endocrinol Metab, 2006. **91**(4): p. 1210-7.

DiabeGen - 1° step

Code 9192, 12 patients

Date of preparation: 2014.09.15

Rev. 0

Eurospital 

Eurospital SpA

Via Flavia 122 - 34147 Trieste - Italy

Tel. +39 040 89971 - Fax +39 040 280944

www.eurospital.com - info@eurospital.com