Istruzioni d'uso

ver. 5 - 04/04/2022

Biologia Molecolare

CF DEL

REF

AC033

CE

25 TEST

CND

W0106010101

IVD

NUCLEAR LASER MEDICINE S.r.l.

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IVD

NUCLEAR LASER MEDICINE S.r.l.

HEADQUARTER: Via Cascina Conighetto – BUSINESS OFFICES: Viale delle Industrie, 3 - 20049 SETTALA (MI), ITALY

Phone: (+39) 02/95.24.51 - Fax (+39) 02/95.24.52.37 - (+39) 02/95.24.52.38 WEB: www.nlm.it - E-MAIL: segreteria@nlm.it

ISO 9001 Quality Management System Organization certified and ISO 13485 Medical Sector Quality Management System certified (IMQ certification body - CSQ certification).

INTENDED USE

The **AC033** device provides reagents for the identification of genetic predisposition to Cystic Fibrosis, by detection of allelic variants in the gene sequence coding for CFTR protein, through the amplification of target sequences, reverse-hybridization and color development.

Amplification enzyme assigned: **orange** cap NLM code BA092, followed by size indication, which could change depending on enzyme quantity supplied (i.e.: BA092/50, BA092/60, etc).

This device has been validated with:

- a) manual column-based extraction (NLM code AA1001) starting from fresh blood with the most common automatic systems for nucleic acids purification.
- b) PTC 100 (MJ Research); 2720 (Applied Biosystems), Mastercycler (Eppendorf); T100 and Mycycler (BioRad) thermal cycler instruments.

⇒ c) automatic detection instruments: Tecan ProfiBlot 48 and Dynex Dynablot Heat

INTRODUCTION

Cystic Fibrosis (CF) is the most common congenital, chronic disease that is transmitted with an autosomal-recessive mechanism more frequent in the Caucasian population: disease occurs in 1/2500-2700 among newborns alive.

Cystic Fibrosis is due to an anomaly in the CFTR protein (Cystic Fibrosis Transmembrane Conductance Regulator) placed in the apical membrane of epithelial cells. Its function is related to the regulation of electrolytic exchange.

In 1989 the gene coding for CFTR protein was mapped on 7q chromosome (long-arm).

Alterations in the protein structure lead to problems in the salt-transport pathways with the productions of dehydrated secretions.

Sweat is very rich in Na and Cl, mucus is dense and sticky and tends to obstruct ducts.

Cystic Fibrosis affects lots of organs and apparatus: respiratory apparatus, pancreas, liver and gut. In males reproductive apparatus is frequently affected.

The disease could manifest very early, in neonatal age or in the first weeks/months of life with different levels of severity and, in some cases, related to particular genetic mutations.

The disease may seldom manifest in the adolescent or adult life with mild clinical profile.

REFERENCES

- *M.P. Audrezet et al Hum. Mutat. 2004; 23:343-357*
- *T. Dork et al Hum Genet 2000; 106:259-268*
- *Lerer et al Hum Mutat 1999; Mutation in Brief n° 231*
- *F. Chevalier-Porst et al. 2005; Mutation in Brief n° 806*
- *Cremonesi L. et al. Hum. Mutat. 1992; 1 (4):314-9*
- *Bonizzato A. et al. Hum. Genet. 1995; 95(4): 397-402*
- *Bombieri C. et al Hum. Genet. 1998; 103(6): 718-22*
- *D'Apice M.R. et al. BMC Med. Genet. 2004; 5 (1): 8*
- *Chillon M. et al. N. Engl. J. Med. 1995; 332: 1475-80*
- *Cuppens H. et al. J. Clin. Invest. 1998; 101: 487-96*
- *Bobadilla JL et al. Hum. Mutat. 2002; 19:575-606*

PRINCIPLE OF THE PROCEDURE

The test is based on the reverse-hybridization principle, where specific oligonucleotide probes immobilized as parallel lines on membrane-based strips hybridize with biotinylated PCR products. The exact match between probes and amplified product generates a signal exploiting the bond between biotin and streptavidin conjugated with alkaline phosphatase and a subsequent color developer.

The analysis consists of three steps:

1. DNA extraction from non-coagulated blood (see the procedure in the following *Extraction* paragraph).
2. PCR of target DNA sequences using biotinylated primers.
3. Hybridization of the amplification products with the oligonucleotide probes and colorimetric detection.

The device allows simultaneous detection of the following macrodeletions: CFTR dele2ins182, CFTR dele1, CFTR dele22,23,24, CFTR dele22,23, CFTR dele2,3, CFTR dele17a,17b,18, CFTR dele14b-17b.

PRODUCT COMPOSITION

⇒ BOX A (store at -25/-15°C)

Reagents	NLM Code	Quantity	N° Vial
PCR Mix Solution containing oligonucleotides	KA902-n	1600 µl	1

⇒ BOX B (store at +2/+8°C)

Reagents		NLM Code	Quantity	N° Vial	
Strips Nitrocellulose membrane coated with oligonucleotides		-	KC007-n	25	1
DNAT Denaturation solution containing NaOH	 DANGER	H314, P280, P301+P330+P331, P303+P361+P353, P305+P351+P338	KA107-n	1,3 ml	1
Hybridization (*) Saline Solution containing preservatives	 WARNING	H319, EUH208, P280, P305+P351+P338, P337+P313	KA102-n	60 ml	1
Stringent Wash Solution Saline solution containing detergent and preservatives	-	-	KA928-n	175 ml	1
Conjugate Solution containing streptavidin labeled with alkaline phosphatase, stabilizers and preservatives	-	-	KA574-n	60 ml	1
Wash Solution B Saline Solution containing preservatives	-	EUH208	KA105-n	175 ml	1
Color developer Solution containing 5-Bromo-4-chloro-3-indolyl-phosphate and 4-Nitroblue tetrazolium	-	-	AA564-n	60 ml	1

Transparent interpretation card	-	-		1	
Collector sheet For strip storage, not suitable for results interpretation and/or as medical report	-	-		1	

(*) The exemption from the requirements for labeling and packaging is determined by quantity and severity of danger reagent. The criteria are: **amounts <125 ml and categories not serious (ref. Regulation (EC) N ° 1272/2008 - list 1.5.2.1.1).**

⇒ The components volume show above refer to the standard size of the kit. Reduced packages of the kit are available on request for evaluation and / or demonstration of the product.

STABILITY AND STORAGE

- All opened or unopened reagents are stable up to the expiry date indicated on the label when stored at the correct temperature (**BOX A -25/-15°C; BOX B +2/+8°C**).
- At the end of each assay store reagents at the appropriate temperature.
- The kit should be kept isolated from any source of contaminating DNA, especially amplified products.
- Developed dry strips should be stored in the dark at room temperature (+15/+30°C).
- The vial containing the DNAT should be closed immediately after use.

PRECAUTION

- Only professional and opportunely trained personnel should use this kit. Handle this product according to established good laboratory practices and universal precautions; wear personal protective apparel.
- When handling DNA samples or amplification products, use aerosol-resistant pipette tips to avoid pipettes contamination.
- Discard used materials as bio hazardous waste.
- Plastic trays can be reused maximum twice, if washed with water at the end of detection.
- Do not eat, drink, smoke or apply cosmetics in areas where reagents or specimens are handled.
- If skin or mucous membrane exposure occurs, immediately wash the area with copious amount of water. Seek medical advice.
- Do not use components beyond the expiration date.
- Test strips stored at +2/+8°C in the dark are stable until the expiry date.
- Do not mix reagents from different lots.
- Reagents have to be preserved separated from possible contaminants (as DNA samples and amplification products).
- It is recommended to work in three separated areas:
 - Area 1: pre-PCR (samples handling and extraction).
 - Area 2: Master Mix preparation.
 - Area 3: post-PCR (amplification and detection).
- Don't use the device if the box is damaged; contact the supplier.
- **It is advisable to have constant and uniform laboratory temperature, avoid to place the instruments near heating/cooling sources that may compromise the correct working.**
- **Use the transparent interpretation card of the specific lot.**
- **Hazard statements**
 - **H314:** Causes severe skin burns and eye damage.
 - **H319:** Causes serious eye irritation.
- **Precautionary statements**
 - **P280:** Wear protective gloves/protective clothing/eye protection/face protection.
 - **P301+P330+P331:** IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
 - **P303+P361+P353:** IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
 - **P305+P351+P338:** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
 - **P337+P313:** If eye irritation persists: Get medical advice/attention.
- **Supplemental hazard informations**
 - **EUH208:** It contains CMIT/MIT. May cause allergic reaction (skin).

MATERIALS REQUIRED BUT NOT SUPPLIED

AREA 1

- DNA Extraction Kit
- Vertical downflow air box
- Dedicated adjustable volume pipettes set
- Aerosol-resistant pipettes tips and 1,5 ml tubes

AREA 2

- Vertical downflow air box
- Dedicated adjustable volume pipettes set
- Aerosol-resistant pipettes tips
- 0,2 ml PCR tubes
- DNA polymerase with **orange cap**

⇒ AREA 3

- Validated thermal cycler or equivalent instrument*
- Agarose gel electrophoresis equipment (optional)
- Dedicated adjustable volume pipettes set and tips
- Clean tweezers to handle strips
- Plastic supports (trays) for manual or automatic detection

Manual strip processing:

- Water bath with inclined lid, shaking platform (50-100 rpm) and adjustable temperature ($45 \pm 0,5$ °C)
- Thermometer
- Vacuum aspiration apparatus
- Deionized or distilled water
- Timer

Automatic strip processing:

- Tecan ProfiBlot 48, Dynex Dynablot Heat or equivalent instrument.

* **WARNING:** if thermal cyclers different from those mentioned above are used, verify that the ramp rate (heating/cooling) is $\leq 3^\circ\text{C}/\text{sec}$. For T100 (BioRad) set default ramp rate ($4^\circ\text{C}/\text{sec}$).

PROCEDURE

DNA EXTRACTION

Manual purification: NLM code AA1001.

For extraction with automatic systems follow the specific manufacturer's IFUs.

Any other type of DNA extraction with similar results of concentration and purification can be used.

Use EDTA or citrate only as anticoagulant, do not use heparin.

- It is possible to use fresh blood, blood stored at +2/+8°C for no more than 2 days or stored at -25/-15°C for longer storage periods. The use of frozen blood is recommended since freezing and thawing enhance blood cells lysis.

AMPLIFICATION

Keep all the reagents for the amplification on ice during the execution of the whole procedure.

Avoid repeated freeze-thaw cycles of the PCR mix

Master Mix preparation

Reagents	Volume/sample
PCR Mix	46,75 µl
DNA polymerase (5U/µl), orange cap	0,25 µl

- Prepare the Master Mix for the number of samples extracted + 1 more volume (if "n" ≤ 10) or + 2 more volumes (if "n" >10).
- Mix gently and dispense **47 µl** of Master Mix in the previously marked 0,2 ml test tubes.
- Add **3 µl** of extracted DNA to each tube and mix pipetting up and down (**DNA concentration to be used should be 30-40 ng/µl**).

Amplification program

- Preheat the thermal cycler lid.
- Put reaction tubes in the thermal cycler and run the following amplification program:

Temperature	Time	Cycles
94 °C	30 sec	1
94 °C	30 sec	35
54 °C	30 sec	
72 °C	40 sec	
72°C	10 min	1
10°C	∞	-

Keep PCR products at + 2/+8°C if detection is done in the same day, or store them at -25/-15°C.

Optional: analyze amplification products by 2% agarose gel containing a DNA intercalating agent.

Macrodeletion traditional nomenclature	Macrodeletion New nomenclature*	Macrodeletion Amplicon lenght (bp)	Wild type Amplicon lenght (bp)
CFTR dele2ins182	c.54-5811_164+2186del8108ins182 ⁽¹⁾	274	364
CFTR dele1	c.4_53+69delins299 ⁽²⁾	490	310
CFTR dele22,23,24	c.3964-3890_4443+3143del9454ins5 ⁽¹⁾	175	202
CFTR dele22,23	c.3964-78_4242+577del ⁽³⁾	226	259
CFTR dele2,3	c.54-5940_273+10250del21kb ⁽³⁾	245	306
CFTR dele17a,17b,18	c.(?_2989)_(3468_?)del ⁽²⁾	180	222
CFTR dele14b-17b	c.(?_2620)_(3367_?)del ⁽²⁾	440	516

*Human Genome Variation Society; reference sequence: NM_000492.3

⁽¹⁾ Taulan et al. (2007), BMC Medical Genetics 2007, 8:22 doi:10.1186/1471-2350-8-22

⁽²⁾ Cystic Fibrosis Mutation Database (CFTR1); available at <http://www.genet.sickkids.on.ca/app>

⁽³⁾ CFTR2 - The Clinical and Functional TRanslation of CFTR; available at <https://www.cftr2.org/>

AUTOMATIC DETECTION

⇒ Detection can be performed by using automatic instrument (ProfiBlot T48, Dynablot or similar instrument): select the **FIBROSI CISTICA** program. In this case **25 µl PCR product + 25 µl DNAT** and 1,5 ml of each reagent should be used for each sample.

MANUAL DETECTION

- Adjust the water level of the shaking bath to approx. 2/3 of the height of the Plastic Tray.
 - Heat the water bath to **41°C** and check the temperature ($41 \pm 0,5^\circ\text{C}$) using a calibrated thermometer.
 - Pre-warm the Hybridization and Stringent Wash Solution to **41°C**; all crystals should be completely dissolved.
 - Allow DNAT, Conjugate, Wash Solution B and Color developer to reach room temperature.
 - Remove one strip for each sample using clean tweezers (touch strips with gloves only). Label strips with a pencil (do not use ballpoint pen, markers, etc).
- Warning:** do not allow the strips to dry during the whole procedure.

Hybridization (41°C)

- Pipette **15 µl** of **DNAT** into each lane of the Plastic Tray
Warning: do not use DNAT if is not blue.
- For each sample add **15 µl** of **amplification** product to the corresponding drop of DNAT. Mix thoroughly with a pipette (solution will remain blue).
- Incubate **5 min** at room temperature.
- Add **1 ml** of **Hybridization** (pre-warmed to 41°C) to each sample and mix gently: blue color disappears.
- Place strips into the respective lanes with marker lines facing up. Submerge completely.
- Incubate for **45 min** at **41°C** in the shaking water bath (about 50 rpm). Keep the water bath closed to avoid temperature variations.

Stringent Wash solution (41°C)

- After hybridization step remove the tray from the water bath; angle the tray upwards and aspirate the solution from each lane using a pipette or a vacuum aspiration apparatus, without touching the strip surface to avoid its damage.
- Add **1 ml** of **Stringent Wash solution** (pre-warmed at 41°C) to each lane and rinse briefly.
- Remove the solution and add **1 ml** of **Stringent Wash solution** (pre-warmed at 41°C).
- Incubate **10 min** at **41°C** in the shaking water bath (about 50 rpm).
Warning: put Stringent Wash solution back in the water bath during incubation
- Remove the solution and add **1 ml** of **Stringent Wash solution** (pre-warmed at 41°C).
- Incubate **10 min** at **41°C** in the shaking water bath (about 50 rpm).
- Remove the solution from each lane.

Color development (room temperature)

All the next steps have to be performed at room temperature with shaking. Since water takes too much time to reach room temperature, place a tray (for example of polystyrene) on the lid of the water bath to isolate the Plastic Trays from the hot water and use the agitation function of the water bath only.

- Add **1 ml** of **Conjugate**.
- Incubate **20 minutes** while shaking (about 50 rpm) at room temperature.
- Remove the solution, add **1 ml** of **Wash solution B** and rinse briefly.
- Remove the solution and add **1 ml** of **Wash solution B**.
- Incubate for **5 min** while shaking (about 50 rpm) at room temperature.
- Remove the solution and add **1 ml** of **Wash solution B**.
- Incubate for **5 min** while shaking (about 50 rpm) at room temperature.
- Remove the solution and add **1 ml** of **Color developer**.
- Incubate for **20 min** in the dark, while shaking (about 50 rpm) at room temperature.
- Remove color developer and wash strips with distilled water.
- Let strips dry in the dark on absorbent paper before proceeding with the interpretation of results.


RESULTS INTERPRETATION

The strip scheme below shows the position of the different probes on the strip. A line is considered positive when a purple band appears at the end of the detection. To identify the genotype of a sample analyze each strip with the help of the transparent interpretation card. In order to correctly align strip and transparent interpretation card, colored marker lines on the top and on the bottom have to be used.

For each mutation one of the following patterns should be obtained:

- Only wild type line Normal genotype (**homozygous wild type**)
- Wild Type and Mutant lines **Heterozygous** genotype
- Only mutant line line **Homozygous mutant**

Strip scheme

Warning:

- the color intensity could be different between different bands on the same strip but this doesn't affect the interpretation of the results.
- The detection control line (the upper one) indicates the correct reaction of both conjugate and color developer: it **always** has to stain positive.

TROUBLESHOOTING

FALSE NEGATIVE OR LOW SIGNALS

- A wrong amount of amplification product has been used. PCR products concentration could be too low due to an inefficient amplification.
- Low signals (except Control Line) and false negative results could be caused by too high temperature during hybridization and/or Stringent Wash Solution steps.

Check carefully water bath temperature

- Hybridization and Stringent Wash Solution have not been well pre-warmed or mixed. Discard remaining reagents and repeat tests with new solutions.

FALSE POSITIVE SIGNALS

- Non-specific signals could be caused by too low temperature during hybridization and/or Stringent Wash Solution step

Check carefully water bath temperature

- Hybridization solution and Stringent Wash Solution have not been well pre-warmed or mixed. Discard remaining reagents and repeat tests with new solutions.
- Water level of the water bath is too low: be sure that level reaches approximately 2/3 of the height of the Plastic Tray
- Check carefully the water bath lid during the Hybridization or Stringent Wash Solution, keep it closed.
- Keep Stringent Wash Solution in the water bath during the incubation steps at 41°C

SINGLE LINE UNHOMOGENEOUS COLORATION

Shaking rate during detection steps is too low. Check that strips are totally submerged.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

The analytical sensitivity of the kit "CF DEL", expressed as the minimum quantity of the target that can be detected (probability > 99%), is 5 ng/μl of DNA.

Diagnostic sensitivity

The diagnostic sensitivity of the kit "CF DEL" was assessed by analyzing 54 DNA samples. Since the condition of homozygosity for these macrodeletions is not vital, the specificity has been studied through the use of specific synthetic samples.

All samples were correctly genotyped, resulting in a diagnostic sensitivity of 100%.

Diagnostic specificity

The diagnostic specificity of the kit "CF DEL" was determined by analyzing 54 DNA samples. No false positive results were obtained. In agreement with these data, the kit specificity is 100%.

Reproducibility

INTRA-ASSAY

The intra-assay reproducibility was evaluated by using 1 sample in 20 replicates in the first run and the same sample in 14 replicates in a second run.

Two different operators tested the replicates using two batch lots.

INTER-ASSAY

The inter-assay reproducibility was evaluated by using 4 samples in 3 replicates. Two different operators tested the replicates in different runs and days, using three lots of reagents. For strips development both manual and automatic system were used.

Tests were performed as described in the instructions for use and the reproducibility was 100%.

WARNING

As every system based on amplification and detection of nucleic acids, it is possible to have unexpected results in the presence of unknown genetic mutation in the primer and/or probe specific target regions.