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
ver. 7 - 06/04/2022

Molecular Biology

SCREENING TEST FOR THROMBOPHILIC DISEASE 7 MUTATIONS

	AC034	
	25 TEST	
GMDN	59910	

NUCLEAR LASER MEDICINE S.r.l.

HEADQUARTER: Via Cascina Conighetto – BUSINESS OFFICES: Viale delle Industrie, 3 - 20049 SETTALA (MI), ITALY

Phone: (+39) 02/95.24.51 - Fax (+39) 02/95.24.52.37 - (+39) 02/95.24.52.38 WEB: www.nlm.it - E-MAIL: segreteria@nlm.it

ISO 9001 Quality Management System Organization certified and ISO 13485 Medical Sector Quality Management System certified (IMQ certification body - CSQ certification).

INTENDED USE

The **AC034** device provides reagents for the identification of the main gene mutations involved in thrombophilic diseases (Factor V R506Q, H1299R, Y1702C; Factor II G20210A; MTHFR C677T, A1298C; PAI-1 4G/5G), through the amplification of target sequences, reverse-hybridization and colour development.

Amplification enzyme assigned: **orange** cap NLM code BA092, followed by size indication, which could change depending on enzyme quantity supplied (i.e.: BA092/50, BA092/100, etc).

This device has been validated with:

- a) manual column-based extraction (NLM code AA1001) starting from fresh blood or dried blood-spot and with the most common automatic systems for nucleic acids purification.
- b) PTC 100 (MJ Research); 2720 (Applied Biosystems), Mastercycler (Eppendorf); T100, C1000 and Mycycler (BioRad) thermal cycler instruments.

INTRODUCTION

Venous thrombosis is among the most common cardiovascular diseases, affecting about 1 in every 1000 people each year. Thrombotic risk is determined both by environmental factors (age, surgery, pregnancy, oral contraception) and genetic predisposition. While most of known genetic defects within the blood coagulation cascade (protein S, protein C, antithrombin III) is rare, a G>A mutation at nucleotide position 1691 in the gene for coagulation factor V (FV Leiden – R506Q) is found frequently (20-60%) in patients affected by thrombosis. Other mutations in Factor V gene, as HR2 aploptype (H1299R: A>G substitution at 4070 nucleotide), can cause thrombosis particularly in addition to R506Q mutation. Other mutations, as Y1702C (A>G substitution at 5279 nucleotide), can lead to level deficiency of Factor V in the blood.

Two further point mutations in methylenetetrahydrofolate reductase gene (MTHFR C677T and A1298C) and 3'-untranslated region of the prothrombin gene (Factor II G20210A) have been reported to be associated respectively to a decreased enzymatic activity of the MTHFR protein and a higher plasma prothrombin level. Each of these variants, in association with Factor V Leiden R506Q, can lead to increased risk of venous thrombosis.

Another candidate marker for inherited susceptibility to cardiovascular disease is PAI-1 gene (Plasminogen Activator Inhibitor - type 1), a physiologic regulator of fibrinolysis. A common deletion polymorphism that results in a sequence of 4G instead of 5G in the promoter region is associated with an increased risk of venous thromboembolism.

REFERENCES

- Alhenc-Gelas M. et al. Thromb Haemost 1999 Feb; 81(2): 193-7
- Bernardi F. et al. Blood 1997 Aug 15; 90(4): 1552-7
- Bertina et al. Nature 1994(369): 64-67
- Castaman G. et al. Br J Haematol 1997 Nov; 99(2): 257-61
- Castoldi E. et al. Blood 2000 Aug 15; 96(4): 1443-8
- Castoldi E. et al. Haematologica 2001 Jun; 86(6): 629-33
- Cattaneo et al. Art, Thromb and Vasc Biology 1997 (17): 1662-1666
- De Stefano et al. Blood 1998 (91): 3562-3565
- Doggen et al. Circulation 1998 (97): 1037-1041
- Faioni E.M. et al. Blood 1999 Nov 1; 94(9): 3062-6
- Finan et al. Am. J. Hematol. 2002 (71-4): 300-305
- Gemmati et al. Haematologica 1999 (84): 824-828
- Lunghi B. et al. Thromb Haemost 1996 Jan; 75(1): 45-8
- Poort et al. Blood 1996 (88): 3698-3703
- Wald et al. BMJ 2002 (325): 1202-1208

PRINCIPLE OF THE PROCEDURE

The test is based on the reverse-hybridization principle, where specific oligonucleotide probes immobilized as parallel lines on membrane-based strips hybridize with biotinylated PCR products. The exact match between probes and amplified product generates a signal exploiting the bond between biotin and streptavidin conjugated with alkaline phosphatase and a subsequent color developer.

The analysis consists of three steps:

1. DNA extraction from non-coagulated blood see the procedure in the following *Extraction* paragraph).
2. PCR of target DNA sequences using biotinylated primers
3. Hybridization of the amplification products with the oligonucleotide probes and colorimetric detection

The device allows simultaneous detection of the 7 mutations in the following order:

FACTOR V

1. R506Q (Factor V Leiden - G1691A)
2. H1299R (HR2 aptotype)
3. Y1702C

FACTOR II (PROTHROMBIN)

4. G20210A

MTHFR

5. C677T
6. A1298C

PAI-1

7. 4G/5G

PRODUCT COMPOSITION

⇒ **BOX A (store at -25/-15°C)**

Reagents	NLM Code	Quantity	N° Vial
PCR Mix Solution containing oligonucleotides	KA899-n	1600 µl	1

⇒ **BOX B (store at +2/+8°C)**

Reagents			NLM Code	Quantity	N° Vial
Strips Nitrocellulose membrane coated with oligonucleotides		-	KC035-n	25	1
DNAT Denaturation solution containing NaOH	 DANGER	H314, P280, P301+P330+P331, P303+P361+P353, P305+P351+P338	KA107-n	1 ml	1
Hybridization (*) Saline Solution containing preservatives	 WARNING	H319, EUH208, P280, P305+P351+P338, P337+P313	KA102-n	60 ml	1
Wash Solution A Saline solution containing detergent ad preservatives	-	-	KA103-n	170 ml	1
Conjugate Solution containing streptavidin labeled with alkaline phosphatase, stabilizers and preservatives	-	-	KA574-n	60 ml	1
Wash Solution B Saline Solution containing preservatives	-	EUH208	KA105-n	170 ml	1
Color developer Solution containing 5-Bromo-4-chloro-3-indolyl-phosphate and 4-Nitroblue tetrazolium	-	-	AA564-n	60 ml	1
Transparent interpretation card	-	-		1	
Collector sheet For strip storage, not suitable for results interpretation and/or as medical report	-	-		1	

(*) The exemption from the requirements for labeling and packaging is determined by quantity and severity of danger reagent. The criteria are: **amounts <125 ml and categories not serious (ref. Regulation (EC) N ° 1272/2008 - list 1.5.2.1.1).**

STABILITY AND STORAGE

- All opened or unopened reagents are stable up to the expiry date indicated on the label when stored at the correct temperature (**BOX A -25/-15°C; BOX B +2/+8°C**)
- At the end of each assay store reagents at the appropriate temperature
- The kit should be kept isolated from any source of contaminating DNA, especially amplified products
- Developed dry strips should be stored in the dark at room temperature (+15/+30°C)
- The vial containing the DNAT should be closed immediately after use

PRECAUTION

- Only professional and opportunely trained personnel should use this kit. Handle this product according to established good laboratory practices and universal precautions; wear personal protective apparel
- When handling DNA samples or amplification products, use aerosol-resistant pipette tips to avoid pipettes contamination
- Discard used materials as bio hazardous waste
- Plastic trays can be reused maximum twice, if washed with water at the end of detection
- Do not eat, drink, smoke or apply cosmetics in areas where reagents or specimens are handled
- If skin or mucous membrane exposure occurs, immediately wash the area with copious amount of water. Seek medical advice.
- Do not use components beyond the expiration date
- Test strips stored at +2/+8°C in the dark are stable until the expiry date
- Do not mix reagents from different lots
- Reagents have to be preserved separated from possible contaminants (as DNA samples and amplification products)
- It is recommended to work in three separated areas:
 - Area 1: pre-PCR (samples handling and DNA extraction)
 - Area 2: Master Mix preparation
 - Area 3: post-PCR (amplification and detection)
- Don't use the device if the box is damaged; contact the supplier
- **It is advisable to have constant and uniform laboratory temperature, avoid to place the instruments near heating/cooling sources that may compromise the correct working**
- **Use the transparent interpretation card of the specific lot**
- **Hazard statements**
 - **H314:** Causes severe skin burns and eye damage.
 - **H319:** Causes serious eye irritation.
- **Precautionary statements**
 - **P280:** Wear protective gloves/protective clothing/eye protection/face protection.
 - **P301+P330+P331:** IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
 - **P303+P361+P353:** IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
 - **P305+P351+P338:** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
 - **P337+P313:** If eye irritation persists: Get medical advice/attention.
- **Supplemental hazard informations**
 - **EUH208:** It contains CMIT/MIT. May cause allergic reaction (skin).

MATERIALS REQUIRED BUT NOT SUPPLIED

AREA 1

- DNA Extraction Kit
- Vertical downflow air box
- Dedicated adjustable volume pipettes set
- Aerosol-resistant pipettes tips and disposable 1,5 ml tubes

AREA 2

- Vertical downflow air box
- Dedicated adjustable volume pipettes set
- Aerosol-resistant pipettes tips
- 0,2 ml PCR tubes
- DNA polymerase with orange cap

AREA 3

- Thermal Cycler or equivalent instrument*
- Agarose gel electrophoresis equipment (optional)
- Water bath with inclined lid, shaking platform (50-100 rpm) and adjustable temperature ($45^{\circ}\text{C} \pm 0,5^{\circ}\text{C}$)
- Plastic supports (trays) for manual or automatic detection
- Thermometer
- Vacuum aspiration apparatus
- Dedicated adjustable volume pipettes set and tips
- Deionized or distilled water
- Timer
- Clean tweezers to handle strips

* **WARNING:** if thermal cyclers different from those mentioned above are used, verify that the ramp rate (heating/cooling) is $\leq 3^{\circ}\text{C}/\text{sec}$. For T100 (BioRad) the ramp rate is $2,5^{\circ}\text{C}/\text{sec}$.

PROCEDURE

DNA EXTRACTION

Manual purification: NLM code AA1001.

For extraction with automatic systems follow the specific manufacturer's IFUs.

Any other type of DNA extraction with similar results of concentration and purification can be used.

Use EDTA or citrate only as anticoagulant, do not use heparin. It is possible to use fresh blood, blood stored at +2/+8°C for no more than 2 days or stored at -25/-15°C for longer storage periods. The use of frozen blood is recommended since freezing and thawing enhance blood cells lysis.

AMPLIFICATION

Keep all the reagents for the amplification on ice during the execution of the whole procedure.

Avoid repeated freeze-thaw cycles of the amplification mix.

Master Mix preparation

Reagents	Volume/sample
PCR Mix	47,25 µl
DNA polymerase (5U/µl), orange cap	0,25 µl

- Prepare the Master Mix for the number of samples extracted + 1 more volume (if "n" ≤ 10) or + 2 more volumes (if "n" >10).
- Mix gently and dispense **47,5 µl** of Master Mix in the previously marked 0,2 ml test tubes. Discard the remaining Master Mix.
- Add **2,5 µl** of extracted DNA to each tube and mix pipetting up and down (**DNA concentration to be used should be 30-40 ng/µl**).

Amplification program

- Preheat the thermal cycler lid.
- Put reaction tubes in the thermal cycler and run the following amplification program:

Temperature	Time	Cycles
95°C	30 sec	1
95°C	15 sec	35
60°C	30 sec	
72°C	45 sec	
72°C	2 min	1
10°C	∞	-

Keep PCR products at + 2/+8°C if detection is done in the same day, or store them at -25/-15°C.

Optional: analyze amplification products by 2% agarose gel containing an intercalating DNA agent.

Mutation	Amplicon lenght (bp)
Factor V R506Q	146
Factor V H1299R	174
Factor V Y1702C	337
Factor II G20210A	207
MTHFR C677T	182
MTHFR A1298C	175
PAI-1 4G/5G	177

AUTOMATIC DETECTION

Detection can be performed by using automatic instrument (ProfiBlot T30/T48 or similar instrument): select “**COAGUL 5**” program. In this case **30 µl PCR product + 30 µl DNAT** and 1,5 ml of each reagent should be used for each sample.

MANUAL DETECTION

- Adjust the water level of the shaking bath to approx. 2/3 of the height of the Plastic Tray
- Heat the water bath to **45°C** and check the temperature ($45^{\circ}\text{C} \pm 0,5^{\circ}\text{C}$) using a calibrated thermometer
- Pre-warm the Hybridization and Wash Solution A to **45°C**; all crystals should be completely dissolved
- Allow DNAT, Conjugate, Wash Solution B and Colour developer to reach room temperature
- Remove one strip for each sample using clean tweezers (touch strips with gloves only). Label strips with a pencil (do not use ballpoint pen, markers, etc)
Warning: do not allow the strips to dry during the whole procedure

Hybridization (45°C)

- Pipette **20 µl** of **DNAT** into each lane of the Plastic Tray
Warning: do not use DNAT if is not blue
- For each sample add **20 µl** of **amplification** product to the corresponding drop of DNAT. Mix thoroughly with a pipette (solution will remain blue).
- Incubate **5 min** at room temperature
- Add **1 ml** of **Hybridization** (pre-warmed to 45°C) to each sample and mix gently: blue colour disappears
- Place strips into the respective lanes with marker lines facing up. Submerge completely.
- Incubate for **30 min** at **45°C** in the shaking water bath (about 50 rpm). Keep the water bath closed to avoid temperature variations

Wash Solution A (45°C)

- After hybridization step remove the tray from the water bath; angle the tray upwards and aspirate the solution from each lane using a pipette or a vacuum aspiration apparatus, without touching the strip surface to avoid its damage
- Add **1 ml** of **Wash Solution A** (pre-warmed at 45°C) to each lane and rinse briefly
- Remove the solution and add **1 ml** of **Wash Solution A** (pre-warmed at 45°C)
- Incubate **10 min** at **45°C** in the shaking water bath (about 50 rpm)
Warning: put Wash solution A back in the water bath during incubation
- Remove the solution and add **1 ml** of **Wash solution A** (pre-warmed at 45°C)
- Incubate **10 min** at **45°C** in the shaking water bath (about 50 rpm)
- Remove the solution from each lane

⇒ Color development (to be performed inside the instrument with the thermostat off)

All the next steps have to be performed with thermostat switched off and shaking. Switched off the water thermostat and use the agitation function of the water bath only. During incubation submerge completely the tray as in previous incubations.

TROUBLESHOOTING

FALSE NEGATIVE OR LOW SIGNALS

- A wrong amount of amplification product has been used. PCR products concentration could be too low due to an inefficient amplification
- Low signals (except Control Line) and false negative results could be caused by too high temperature during hybridization and/or Wash Solution A steps

Check carefully water bath temperature

- Hybridization and Wash Solution A have not been well pre-warmed or mixed. Discard remaining reagents and repeat tests with new solutions

FALSE POSITIVE SIGNALS

- Non-specific signals could be caused by too low temperature during hybridization and/or Wash Solution A step

Check carefully water bath temperature

- Hybridization solution and Wash Solution A have not been well pre-warmed or mixed. Discard remaining reagents and repeat tests with new solutions
- Water level of the water bath is too low: be sure that level reaches approximately 2/3 of the height of the Plastic Tray
- Check carefully the water bath lid during the Hybridization or Wash Solution A, keep it closed
- Keep Wash Solution A in the water bath during the incubation steps at 45°C

SINGLE LINE UNHOMOGENEOUS COLORATION

Shaking rate during detection steps is too low. Check that strips are totally submerged.

PERFORMANCE CHARACTERISTICS

Analytical Sensitivity

The analytical sensitivity of the kit " SCREENING TEST FOR THROMBOPHILIC DISEASE, 7 MUTATIONS ", expressed as the minimum quantity of the target that can be detected, is 5 ng/μl of DNA.

Diagnostic sensitivity

The diagnostic sensitivity of the "SCREENING TEST FOR THROMBOPHILIC DISEASE, 7 MUTATIONS" was assessed by analyzing 76 DNA sample positive for the mutations investigated by the kit. For rare mutations, in case of lack of human samples, a specific synthetic DNA was used.

All samples were correctly genotyped, resulting in a diagnostic sensitivity of 100%.

Sample used

Mutation	wild type	heterozygous	mutant	Nr samples
Fattore V R506Q	53	9	14	76
Fattore V H1299R	34	26	16	76
Fattore V Y1702C	68	4	4	76
Fattore II G20210A	58	9	9	76
MTHFR C667T	29	22	25	76
MTHFR A1298C	18	43	15	76
PAI-1 4G/5G	30(5G)	22(4G/5G)	24(4G)	76

Diagnostic specificity

The specificity was determined by analyzing 76 samples. No false positive results were obtained. In agreement with these data, the kit specificity is 100%.

Reproducibility

INTRA-ASSAY

The intra-assay reproducibility was evaluated using 1 sample in 31 replicates in the same run and 5 samples with different genotypes in 2 replicates in a second run. Two different operators tested the replicates using two batch lots.


INTER-ASSAY

The inter-assay reproducibility was tested using at least 3 samples in 4 replicates. Two different operators tested several replicates in different runs and days, using two lots of reagents. The strips development was made in automatic and manual.

Tests were performed as described in the instruction for use and the reproducibility was 100%.

WARNING

As every system based on amplification and detection of nucleic acids, it is possible to have unexpected results in the presence of unknown genetic mutation in the primer and/or probe specific target regions.

 Lot number	 Catalogue number	 Temperature limitation	 Use by	 Consult instruction for use	 Modified text compared with the previous version
 European Conformity	 Manufacturer	 <i>In vitro</i> diagnostic device	 Sufficient for	 Keep away from sunlight	