Istruzioni d'uso

ver. 10 - 14/04/2022

Biologia Molecolare

ALFA GLOBINA TEST

REF

AC099

CE

10 TEST

CND

W0106010111

IVD

NUCLEAR LASER MEDICINE S.r.l.

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IVD

NUCLEAR LASER MEDICINE S.r.l.

HEADQUARTER: Via Cascina Conighetto – BUSINESS OFFICES: Viale delle Industrie, 3 - 20049 SETTALA (MI), ITALY

Phone: (+39) 02/95.24.51 - Fax (+39) 02/95.24.52.37 - (+39) 02/95.24.52.38 WEB: www.nlm.it - E-MAIL: segreteria@nlm.it

ISO 9001 Quality Management System Organization certified and ISO 13485 Medical Sector Quality Management System certified (IMQ certification body - CSQ certification).

INTENDED USE

The Alpha Globin Test device provides reagents for the identification of 22 alterations of alpha-globin gene, involved in alpha type thalassemia, through the amplification of the target sequences, reverse-hybridization and colour development.

In the device are included the **Enzyme mix-AC099** containing:

- the amplification enzyme
- Uracil-DNA Glycosylase (UNG) to prevent the contamination by amplification products.

This device has been validated with:

- a) manual column-based extraction (NLM code AA1001) and with the most common automatic systems for nucleic acids purification.
- b) PTC100 (MJ Research); GeneAmp® 9700, AB2720, Veriti e SimpliAmp (Applied Biosystems); Mastercycler ep gradient S (Eppendorf); T100 (BioRad) thermal cycler instruments.
- ⇒ c) automatic detection instruments: Tecan ProfiBlot 48 and Dynex Dynablot Heat

INTRODUCTION

Alpha-thalassemia is one of the most common haemoglobinopathies in the world and it is caused by the absence or reduced synthesis of alpha-globin chains. Originally, α -thalassemias were only endemic in malarial regions but recent global migration of people from these areas resulted in a world-wide occurrence of these diseases.

The presence on the short arm of chromosome 16 of two copies of alpha-globin genes per haploid genome makes the pathology of alpha-thalassemia far more complicated than that of beta-thalassemia.

There are four major types of alpha-thalassemia:

- α^+ thalassemia ($-\alpha/\alpha$) in which one of the two linked α -globin genes on a chromosome is inactivated
- α^0 thalassemia ($--/\alpha$) in which two α -globin genes are inactivated
- haemoglobin H disease ($--/\alpha$) in which three α -globin genes are absent with consequent chronic haemolytic anaemia of variable severity, splenomegaly and a shortened red cells survival
- Hb Bart's hydrops fetalis ($--/--$) is caused by the absence of alpha chain production and fetuses die either *in utero* or shortly after birth.

Although individuals with three functional α genes ($-\alpha/\alpha$) are clinically and haematologically silent and carriers with α -thalassemia trait ($-\alpha/\alpha$ or $--/\alpha$) only result in very mild hypochromic microcytic anaemia, couples with these genotypes are at risk of having a hydrop baby or with HbH disease.

For suitable clinical patients management, the achieved result using this kit should always be associated with patient clinical evaluations and other laboratory parameters.

REFERENCES

- Fortina, Dianzani, Serra, Gottardi, Saglio, Farinasso, Piga, Gabutti, Camaschella, *British Journal of Haematology* 1991; 78: 529-534
- Oron-Karni, Filon, Oppenheim and Rund, *American Journal of Hematology* 1998; 58: 306-310
- Liu, Old, Miles, Fisher, Weatherall and Clegg, *British Journal of Haematology* 2000; 108: 295-299
- Samuel S. Chong, Corinne D. Boehm, Douglas R. Higgs, and Garry R. Cutting, *Blood*, 1 Jan 2000 Vol 95
- Arnold S.-C. Tan, Thuan C. Quah, Poh S. Low and Samuel S. Chong, *Blood* 2001; 98: 250-251
- Neishabury, Abbasi-Moheb, Poorfathollah, *Arch Irrn Med* 2001; 4 (4):160-164
- Foglietta, Bianco, Maggio, and Giambona, *American Journal of Hematology* 2003; 74:191–195
- David H. K. Chui, Suthat Fucharoen and Vivian Chan, *Blood* 2003; 101: 791-800

PRINCIPLE OF THE PROCEDURE

The test is based on the reverse-hybridization principle, where specific oligonucleotide probes immobilized as parallel lines on membrane-based strips hybridize with biotinylated PCR products. The exact match between probes and amplified product generates a signal exploiting the bond between biotin and streptavidin conjugated with alkaline phosphatase and a subsequent colour developer.

The analysis consists of three steps:

1. DNA extraction from non-coagulated blood (see the procedure in the following *Extraction* paragraph);
2. PCR of target DNA sequences using biotinylated primers;
3. Hybridization of the amplification products with the oligonucleotide probes and colorimetric detection;

The device allows simultaneous detection of the following variations:

Strip A

	HGVS NOMENCLATURE*
-3.7 (type I)	NG_000006.1:g.34164_37967del3804
-4.2	-
-20.5	NG_000006.1:g.15164_37864del22701
--MED	NG_000006.1:g.24664_41064del16401
--SEA	NG_000006.1:g.26264_45564del19301
--FIL	NG_000006.1:g.11684_43534del31851
--CAL	NG_000006.1:g.8464_40664del32201
α1 cd 14 G>A	HBA1:c.44G>A
α1 cd 59 Hb Adana	HBA1:c.179G>A
α1 cd 47 Hasharon	HBA2:c.142G>C

Strip B

	HGVS NOMENCLATURE*
anti 3.7	-
α2 init cd	HBA2:c.2T>C
α2 cd 19 -G	HBA2:c.56delG
α2 IVS1 -5nt	HBA2:c.95+2_95+6delTGAGG
α2 cd 59 Hb Adana	HBA2:c.179G>A
α2 cd 125 Hb Quong Sze	HBA2:c.377T>C
α2 cd 142 Hb Constant Spring	HBA2:c.427T>C
α2 cd 142 Hb Icaria	HBA2:c.427T>A
α2 cd 142 Hb Pakse	HBA2:c.429A>T
α2 cd 142 Hb Koya Dora	HBA2:c.428A>C
α2 poly-A1	HBA2:c.*+94A>G
α2 poly-A2	HBA2:c.*+92A>G

*Human Genome Variation Society; reference sequence: NG_000006.1

PRODUCT COMPOSITION

⇒ BOX A (store at -25/-15°C)

Reagents	NLM Code	Quantity	N° Vial
PCR Mix A1 (transparent cap) Solution containing oligonucleotides	KA862-n	440 µl	1
PCR Mix A2 (yellow cap) Solution containing oligonucleotides	KA883-n	440 µl	1
PCR Mix A3 (violet cap) Solution containing oligonucleotides	KA939-n	440 µl	1
PCR Mix B (blue cap) Solution containing oligonucleotides	KA891-n	440 µl	1
Enzyme Mix-AC099 (white cap) Solution containing enzymes	KA1026-n	50 µl	1

⇒ BOX B (store at +2/+8°C)

Reagents		NLM Code	Quantity	N° Vial	
Strips A Strips B Nitrocellulose membranes coated with oligonucleotides		-	KC036-n KC027-n	10 10	1 1
DNAT Denaturation solution containing NaOH	 DANGER	H314, P280, P301+P330+P331, P303+P361+P353, P305+P351+P338	KA107-n	0,8 ml	1
Hybridization (*) Saline Solution containing preservatives	 WARNING	H319, EUH208, P280, P305+P351+P338, P337+P313	KA102-n	48 ml	1
Stringent Wash Solution Saline solution containing detergent ad preservatives	-	-	KA928-n	136 ml	1
Conjugate Solution containing streptavidin labeled with alkaline phosphatase, stabilizers and preservatives	-	-	KA574-n	48 ml	1
Wash Solution B Saline Solution containing preservatives	-	EUH208	KA105-n	136 ml	1
Color developer Solution containing 5-Bromo-4-chloro-3-indolyl-phosphate and 4-Nitroblue tetrazolium	-	-	AA564-n	48 ml	1
Transparent interpretation card	-	-		2	
Collector sheet For strip storage, not suitable for results interpretation and/or as medical report	-	-		1	

(*) The exemption from the requirements for labeling and packaging is determined by quantity and severity of danger reagent. The criteria are: **amounts <125 ml and categories not serious (ref. Regulation (EC) N° 1272/2008 - list 1.5.2.1.1).**

⇒ The components volume show above refer to the standard size of the kit. Reduced packages of the kit are available on request for evaluation and / or demonstration of the product.

STABILITY AND STORAGE

- All opened or unopened reagents are stable up to the expiry date indicated on the label when stored at the correct temperature (**BOX A -25/-15°C; BOX B +2/+8°C**).
- At the end of each assay store reagents at the appropriate temperature.
- The kit should be kept isolated from any source of contaminating DNA, especially amplified products.
- Developed dry strips should be stored in the dark at room temperature (+15/+30°C).
- The vial containing the DNAT should be closed immediately after use.

PRECAUTIONS

- Only professional and opportunely trained personnel should use this kit. Handle this product according to established good laboratory practices and universal precautions; wear personal protective apparel.
- When handling DNA samples or amplification products, use aerosol-resistant pipette tips to avoid pipettes contamination.
- Discard used materials as bio hazardous waste.
- Plastic trays can be reused maximum twice, if washed with water at the end of detection.
- Do not eat, drink, smoke or apply cosmetics in areas where reagents or specimens are handled.
- If skin or mucous membrane exposure occurs, immediately wash the area with copious amount of water. Seek medical advice.
- Do not use components beyond the expiration date.
- Test strips stored at +2/+8°C in the dark are stable until the expiry date.
- Do not mix reagents from different lots.
- Reagents have to be preserved separated from possible contaminants (as DNA samples and amplification products).
- It is recommended to work in three separated areas:
 - Area 1: samples handling and extraction; PCR set-up with automatic systems
 - Area 2: Master Mix preparation (manual).
 - Area 3: post-PCR (amplification and detection).
- Don't use the device if the box is damaged; contact the supplier.
- **It is advisable to have constant and uniform laboratory temperature, avoid to place the instruments near heating/cooling sources that may compromise the correct working.**
- **Use the transparent interpretation card of the specific lot.**
- **Hazard statements**
 - **H314:** Causes severe skin burns and eye damage.
 - **H319:** Causes serious eye irritation.
- **Precautionary statements**
 - **P280:** Wear protective gloves/protective clothing/eye protection/face protection.
 - **P301+P330+P331:** IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
 - **P303+P361+P353:** IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
 - **P305+P351+P338:** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
 - **P337+P313:** If eye irritation persists: Get medical advice/attention.

⇒ • **Supplemental hazard informations**

- **EUH208:** It contains CMIT/MIT. May cause allergic reaction (skin).

MATERIALS REQUIRED BUT NOT SUPPLIED

AREA 1

- DNA Extraction Kit
- Automatic system for DNA purification and PCR set-up and disposables
- Vertical downflow air box
- Dedicated adjustable volume pipettes set
- Aerosol-resistant pipettes tips and disposable 1,5 ml tubes

AREA 2

- Vertical downflow air box
- Dedicated adjustable volume pipettes set
- Aerosol-resistant pipettes tips
- 0,2 ml PCR tubes

⇒ **AREA 3**

- Validated thermal cycler or equivalent instrument*
- Agarose gel electrophoresis equipment (optional)
- Dedicated adjustable volume pipettes set and tips
- Clean tweezers to handle strips
- Plastic supports (trays) for manual or automatic detection

Manual strip processing:

- Water bath with inclined lid, shaking platform (50-100 rpm) and adjustable temperature ($45 \pm 0,5$ °C)
- Thermometer
- Vacuum aspiration apparatus
- Deionized or distilled water
- Timer

Automatic strip processing:

- Tecan ProfiBlot 48, Dynex Dynablot Heat or equivalent instrument.

*** WARNING**

- ✓ If thermal cyclers different from those mentioned above are used, verify that the ramp rate (heating/cooling) is $\leq 3^{\circ}\text{C}/\text{sec}$. Specifically:
 - for T100 (BioRad) and SimplyAmp (Applied Biosystems) the ramp rate is $2,5^{\circ}\text{C}/\text{sec}$;
 - for Mastercycler ep gradient set as 40% (heating ramp rate) and 50% (cooling ramp rate)
- ✓ Contact NLM for informations about compatibility between thermal cycler and DNA extraction

PROCEDURE

DNA EXTRACTION

Manual purification: NLM code AA1001.

For extraction with automatic systems follow the specific manufacturer's IFUs.

Any other type of DNA extraction with similar results of concentration and purification can be used.

Use EDTA or citrate only as anticoagulant, do not use heparin. It is possible to use fresh blood, blood stored at $+2/+8^{\circ}\text{C}$ for no more than 2 days or stored at $-25/-15^{\circ}\text{C}$ for longer storage periods.

The use of frozen blood is recommended since freezing and thawing enhance blood cells lysis.

AMPLIFICATION

Keep all the reagents for the amplification on ice during the execution of the whole procedure. Avoid repeated freeze-thaw cycles of the amplification and Enzyme mix.

Master Mix preparation (manual)

For each of the four mixes (Mix A1, Mix A2, Mix A3, Mix B) proceed as follows:

Reagents	Volume/sample
PCR Mix A1, A2, A3 and B	26,3 µl
Enzyme mix-AC099	0,7 µl

- Prepare the Master Mix for the number of samples extracted + 1 more volume.
- Mix gently and dispense **27 µl** of Master Mix in the previously marked 0,2 ml test tubes. Discard the remaining Master Mix.
- Add **3 µl** of DNA to each tube and mix pipetting up and down (**DNA concentration to be used should be 10-40 ng/µl**).

Amplification program

- Preheat the thermal cycler lid.
- Put reaction tubes in the thermal cycler and run the following amplification program:

Temperature	Time	Cycles
37°C	10 min	1
97°C	3 min	1
97°C	40 sec	15
66°C*	50 sec	
72°C	1min 30sec	
97°C	40 sec	25
56°C	50 sec	
72°C	1 min 30sec	
72°C	3 min	1
10°C	∞	


* for PTC 100 MJ Research, GeneAmp® 9700 (Applied Biosystems) with sample tray only, Mastercycler ep gradient S set up the annealing temperature to 65°C.

Keep PCR products at + 2/+8°C if detection is done in the same day, or store them at -25/-15°C.

Optional: analyze amplification products by 2% agarose gel containing a DNA intercalating agent.

Mix	Amplicon	Amplicon length (bp)
A1	-3.7	1600
	Alfa 1 Gene	925
	A1 PCR Control	273
A2	-4.2	1000
	-20.5	750
	--MED	650
	--SEA	250
	--FIL	700
	A2 PCR Control	192
A3	--CAL	1800
	A3 PCR Control	141
B	anti 3.7	1900
	Alfa 2 Gene	926
	B PCR Control	176

Examples

AUTOMATIC DETECTION

Detection can be performed by using automatic instrument (ProfiBlot T48, Dynablot or similar instrument): select the **“COAGUL 5”** program. In this case **15 µl A1 PCR product + 15 µl A2 PCR product + 15 µl A3 PCR product + 45 µl DNAT** and **15 µl B PCR product + 15 µl DNAT** and 1,5 ml of each reagent should be used for each sample.

MANUAL DETECTION

- Adjust the water level of the shaking bath to approx. 2/3 of the height of the Typing Tray.
 - Heat the water bath to **45°C** and check the temperature ($45 \pm 0,5^\circ\text{C}$) using a calibrated thermometer.
 - Pre-warm the Hybridization and Stringent Wash Solution to **45°C**; all crystals should be completely dissolved.
 - Allow DNAT, Conjugate, Wash Solution B and Colour developer to reach room temperature.
 - Remove one strip for each sample using clean tweezers (touch strips with gloves only). Label strips with a pencil (do not use ballpoint pen, markers, etc).
- Warning:** do not allow the strips to dry during the whole procedure.

Hybridization (45°C)

- Pipette **45 µl of DNAT** for strip **A** and **15 µl of DNAT** for strip **B** into each lane of the Typing Tray
- Warning:** do not use DNAT if is not blue.
- For each sample add **15 µl of amplification product A1 + 15 µl of amplification product A2 + 15 µl of amplification product A3** for **strip A** and **15 µl of amplification product B** for **strip B** to the corresponding drop of DNAT. Mix thoroughly with a pipette (solution will remain blue).
 - Incubate **5 min** at room temperature.
 - Add **1 ml of Hybridization** (pre-warmed to 45°C) to each sample and mix gently: blue colour disappears.
 - Place strips into the respective lanes with marker lines facing up. Submerge completely.
 - Incubate for **30 min** at **45°C** in the shaking water bath (about 50 rpm). Keep the water bath closed to avoid temperature variations.

Stringent Wash solution (45°C)

- After hybridization step remove the tray from the water bath; angle the tray upwards and aspirate the solution from each lane using a pipette or a vacuum aspiration apparatus, without touching the strip surface to avoid its damage.
- Add **1 ml** of **Stringent Wash solution** (pre-warmed to 45°C) to each lane and rinse briefly.
- Remove the solution and add **1 ml** of **Stringent Wash solution** (pre-warmed to 45°C).
- Incubate **10 min** at **45°C** in the shaking water bath (about 50 rpm).
Warning: put Stringent Wash solution back in the water bath during incubation
- Remove the solution and add **1 ml** of **Stringent Wash solution** (pre-warmed to 45°C).
- Incubate **10 min** at **45 °C** in the shaking water bath (about 50 rpm).
- Remove the solution from each lane.

Colour development (room temperature)


All the next steps have to be performed at room temperature with shaking. Since water takes too much time to reach room temperature, place a tray (for example of polystyrene) on the lid of the water bath to isolate the Typing Trays from the hot water and use the agitation function of the water bath only.

Add **1 ml** of **Conjugate**.

- Incubate **20 minutes** while shaking (about 50 rpm) at room temperature.
- Remove the solution, add **1 ml** of **Wash solution B** and rinse briefly.
- Remove the solution and add **1 ml** of **Wash solution B**.
- Incubate for **5 min** while shaking (about 50 rpm) at room temperature.
- Remove the solution and add **1 ml** of **Wash solution B**.
- Incubate for **5 min** while shaking (about 50 rpm) at room temperature.
- Remove the solution and add **1 ml** of **Colour developer**.
- Incubate for **20 min** in the dark, while shaking (about 50 rpm) at room temperature.
- Remove colour developer and wash strips with distilled water.
- Let strips dry in the dark on absorbent paper before proceeding with the interpretation of results.

RESULTS INTERPRETATION

The strip scheme below shows the position of the different probes on the strip:

A line is considered positive when a purple band appears at the end of the detection. To identify the genotype of a sample analyze each strip with the help of the transparent interpretation card. In order to correctly align strip and transparent interpretation card, coloured marker lines on the top and on the bottom have to be used.

The detection control line (the upper one) indicates the correct reaction of both conjugate and colour developer: it **always** has to stain positive.

→ **The PCR control line always has to stain positive.**

- For each point mutation one of the following patterns should be obtained:
- Only wild type line Normal genotype (**homozygous wild type**)
- Wild type and mutant lines **Heterozygous** genotype
- Only mutant line **Homozygous mutant** genotype


Warning:

- the colour intensity could be different between different bands on the same strip but this doesn't affect the interpretation of the results.
- Since some of the point mutations are located within a few nucleotides on the α -globin gene, they have a common wild type probe. Samples that are compound heterozygous for two of these mutations (e.g. Hb Constant Spring + Hb Pakse) don't show a positive signal on the common wild type probe.
- In composed -3.7/Hasharon case the 16 probe may be weak or absent. However, reference the sample as a double heterozygous.
- This kit does not allow to distinguish between the heterozygous and the homozygous mutant state of the anti 3.7 gene triplication.
- In macro deletions presence, the wild type probes of the two strips behave differently depending on the macro deletions zigosity status:

Deletion	Homozygous mutant*	Heterozygous
-3.7	WT signals on strip B absent	All WT signals present
-4.2		
-20.5	WT signals 14 and 16 on strip A and all WT on strip B absent	
--MED	All WT signals absent	
--SEA		
--FIL		
--CAL		

* Its possible achieve same result in presence of following macro deletion associated with alpha1 or alpha1/alpha2 deletions not detected by this kit. It is recommended to always compare the strip result with the other laboratory parameters for the same sample.

Examples


Normal

-3.7 heterozygous

-3.7 homozygous

-4.2/IVS1-5nt

Hb Constant Spring/Pakse

Hb Constant Spring homozygous

-3.7/--MED

Negative Ctr/PCR failure

TROUBLESHOOTING

FALSE NEGATIVE OR LOW SIGNALS

- Insufficient amplification product due to inefficient PCR, related to:
 - wrong reagent volumes in PCR setup.
 - low concentration and/or purity DNA; verify that the A260/A230 ratio is about 2 ($1.8 \leq \text{value} \leq 2.2$).
 - thermal cycler failure.
- Absence and/or weak signal (except Control Line) could be caused by temperature higher than indicated range for hybridization and Stringent Wash steps.

Check carefully water bath temperature

- Hybridization and Stringent Wash Solution have not been well pre-warmed or mixed. Discard remaining reagents and repeat tests with new solutions.

FALSE POSITIVE SIGNALS

- Non-specific signals could be caused by too low temperature during hybridization and/or Stringent Wash Solution step

Check carefully water bath temperature

- Hybridization solution and Stringent Wash Solution have not been well pre-warmed or mixed. Discard remaining reagents and repeat tests with new solutions.
- Water level of the water bath is too low: be sure that level reaches approximately 2/3 of the height of the Typing Tray
- Check carefully the water bath lid during the Hybridization or Stringent Wash Solution, keep it closed.
- Keep Stringent Wash Solution in the water bath during the incubation steps at **45°C**

SINGLE LINE UNHOMOGENEOUS COLORATION

Shaking rate during detection steps is too low. Check that strips are totally submerged.

WARNING

--CAL samples amplified with efficient thermal cyclers could show a diffuse background on the strip. This doesn't affect the interpretation of the results.

As every system based on amplification and detection of nucleic acids, it is possible to have unexpected results in the presence of unknown genetic mutation in the primer and/or probe specific target regions.

PERFORMANCE CHARACTERISTICS

Analytical Sensitivity

The analytical sensitivity of the kit "ALPHA GLOBIN TEST", expressed as the minimum quantity of the target that can be detected, is equal to 10 ng/μl DNA.

Diagnostic specificity

The specificity was determined by analyzing 65 samples. No false positive results were obtained.

In agreement with these data the kit specificity is 100%.

Reproducibility

INTRA-ASSAY:

The intra-assay reproducibility was assessed by analyzing at least 3 samples in 3 replicates.

INTER-ASSAY:

The reproducibility was tested by analyzing 3 samples in 3 independent seats between them. The revelation was made in automatic and manual.

Tests were performed as described in the method and the reproducibility was 100%.

Diagnostic sensitivity

The diagnostic sensitivity of the kit "ALPHA GLOBIN TEST" was assessed by analyzing 65 DNA sample positive for the mutations investigated by the kit. We used 58 samples of DNA obtained from whole blood and 7 synthetic sample.

All samples were correctly genotyped, resulting in a diagnostic sensitivity of 100%.

Samples used

Nr samples	Mutation	HGVS nomenclature
5	-3,7	NG_000006.1:g.34164_37967del3804
4	-4,2	-
4	-20,5	NG_000006.1:g.15164_37864del22701
4	--MED	NG_000006.1:g.24664_41064del16401
5	--SEA	NG_000006.1:g.26264_45564del19301
3	--FIL	NG_000006.1:g.11684_43534del31851
4	--CAL	NG_000006.1:g.8464_40664del32201
2	α1 cd 14 G>A	HBA1:c.44G>A
1	α1 cd 59 Hb Adana	HBA1:c.179G>A
2	α1 cd 47 Hb Hasharon	HBA2:c.142G>C
5	anti 3.7	-
4	α2 init cd	HBA2:c.2T>C
1	α2 cd 19-G	HBA2:c.56delG
4	α2 IVS1-5nt	HBA2:c.95+2_95+6delTGAGG
1	α2 cd 59 Hb Adana	HBA2:c.179G>A
1	α2 cd 125 Hb Quong Sze	HBA2:c.377T>C
4	α2 cd 142 Hb Constant Spring	HBA2:c.427T>C
3	α2 cd 142 Hb Icaria	HBA2:c.427T>A
1	α2 cd 142 Hb Pakse	HBA2:c.429A>T
1	α2 cd 142 Hb Koya Dora	HBA2:c.428A>C
4	α2 poly A-1	HBA2:c.*+94A>G
2	α2 poly A-2	HBA2:c.*+92A>G