Istruzioni d'uso

ver. 6 - 21/04/2022

Biologia Molecolare

BETA GLOBINA PLUS TEST

REF

AC104

CE

20 TEST

CND

W0106010111

IVD

NUCLEAR LASER MEDICINE S.r.l.

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BIBLIOGRAFIA

- Cao A. & Moi P. (2002). Regulation of the globin genes. *Pediatric Research*, 51, 4, pp. 415-421
- Esposito G., Grosso M., Gottardi E., Izzo P., Camaschella C. & Salvatore F. (1994). A unique origin for the Sicilian ($\delta\beta$)⁰-thalassemia in 33 unrelated families and its rapid diagnostic characterization by PCR analysis. *Human Genetics*, 93, pp. 691-693.
- Galanello R. & Origa R. (2010). β -thalassemia. *Orphanet Journal of Rare Diseases*.
- Thein S.L., Menzel S., Lathrop M. & Garner C. (2009). Control of fetal hemoglobin: new insights emerging from genomics and clinical implications. *Human Molecular Genetics*, 18, pp. 216-223
- Weatherall D.J. (2001). Phenotype-genotype relationships in monogenic disease: lessons from the thalassaemias. *Nature Reviews*, 2, pp. 245-255
- Weatherall D.J. & Clegg J.B. (2001). The Inherited haemoglobin disorders: an increasing global health problem. *Bulletin of the World Health Organization*, 2001,79(8)
- Wood, W.G. (2001) Hereditary persistence and fetal hemoglobin and $\delta\beta$ thalassemia. In *Disorders of Hemoglobin* (Steinberg, M.H., Forget, B.G., Higgs, D.R., and Nagel, R.L., eds.). pp. 356-388. Cambridge University Press, Cambridge.

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
 **NUCLEAR LASER MEDICINE S.r.l.**

HEADQUARTER: Via Cascina Conighetto – BUSINESS OFFICES: Viale delle Industrie, 3 - 20049 SETTALA (MI), ITALY

Phone: (+39) 02/95.24.51 - Fax (+39) 02/95.24.52.37 - (+39) 02/95.24.52.38 WEB: www.nlm.it - E-MAIL: segreteria@nlm.it

ISO 9001 Quality Management System Organization certified and ISO 13485 Medical Sector Quality Management System certified (IMQ certification body - CSQ certification).

INTENDED USE

The **AC104** device provides reagents for the identification of the α -Globin Gene Triplications (anti 3.7) and 14 mutations/deletions of beta-globin gene (**HBB**), through the amplification of the target sequences, reverse-hybridization and color development.

Amplification enzyme assigned: **orange** cap NLM code BA092, followed by size indication, which could change depending on enzyme quantity supplied (i.e.: BA092/50, BA092/60, etc).

This device has been validated with:

- a) manual column-based extraction (NLM code AA1001) starting from fresh blood, NLM code AA340 (QIAamp DNA Mini Kit, Qiagen) for DNA extraction from chorionic villus samples and with the most common automatic systems for nucleic acids purification.
- ⇒ b) PTC 100 (MJ Research); Veriti®, 2720 and GeneAmp® 9700 (Applied Biosystems); Mastercycler ep gradient S (Eppendorf); T100 (BioRad) thermal cycler instruments
- ⇒ c) automatic detection instruments: Tecan ProfiBlot 48 and Dynex Dynablot Heat

INTRODUCTION

The hemoglobinopathies are a heterogeneous group of monogenic diseases spread worldwide.

They are commonly divided into three subgroups partially overlapped: the structural variants (characterized by the production of abnormal globin chains), thalassemias (due to an impaired synthesis of one or more globin chains) and the conditions of hereditary persistence of fetal haemoglobin (HPFH) in adult life (Weatherall & Clegg, 2001; Galanello & Origa, 2010).

The subgroup of structural variants include, for example, the HbE (widespread in India and in Asia) and Hb Lepore (due to unequal crossing-over events which involves the δ - and β -globin gene and determines the formation of a hybrid globin chain).

The β -thalassemias are due to alterations in β -globin gene (HBB) located on chromosome 11 (single base changes, small deletions or insertions) that determine a reduced (beta +) or absent (beta 0) synthesis of the beta chains of the hemoglobin tetramer.

Concerning the clinical and haematological features, β -thalassemias are classified according to severity: thalassemia minor (carrier of beta thalassemia), thalassemia intermedia and thalassemia major.

The presence of α -globin gene triplication is an important modulator of severity in individuals carrying a single mutation in the β -globin gene: the excessive production of α -globin chains results in a greater clinical severity due to the increased imbalance globin chain production (Cao & Moi, 2002; Weatherall, 2001)

HPFH syndroms are characterized by a high activity of γ -genes and decreased activity of the β -globin gene, so that the normal levels of HbA2 globin synthesis are balanced (Wood, 2001).

Forms of HPFH may be caused by large deletions within the β -globin cluster or point mutations in the proximal or distant promoter regions of the fetal globin genes (Cao & Moi, 2000). An example of HPFH deletions in the β -globin cluster is the Sicilian $\delta\beta^0$ deletion (Esposito G. et al., 1994). A similar but less marked effect on HbF increase is found in the Hb Lepore (Weatherall & Clegg, 2001)

Traditionally, the diagnosis of the rearrangements of the gene was performed by Southern blot and the use of certain restriction enzymes, with a considerable use of time and costs. The assays that are based on the PCR, however are more rapid, less expensive and more sensitive.

REFERENCES

- Cao A. & Moi P. (2002). Regulation of the globin genes. *Pediatric Research*, 51, 4, pp. 415-421
- Esposito G., Grosso M., Gottardi E., Izzo P., Camaschella C. & Salvatore F. (1994). A unique origin for the Sicilian ($\delta\beta$)⁰-thalassemia in 33 unrelated families and its rapid diagnostic characterization by PCR analysis. *Human Genetics*, 93, pp. 691-693.
- Galanello R. & Origa R. (2010). β -thalassemia. *Orphanet Journal of Rare Diseases*.
- Thein S.L., Menzel S., Lathrop M. & Garner C. (2009). Control of fetal hemoglobin: new insights emerging from genomics and clinical implications. *Human Molecular Genetics*, 18, pp. 216-223
- Weatherall D.J. (2001). Phenotype-genotype relationships in monogenic disease: lessons from the thalassaemias. *Nature Reviews*, 2, pp. 245-255
- Weatherall D.J. & Clegg J.B. (2001). The Inherited haemoglobin disorders: an increasing global health problem. *Bulletin of the World Health Organization*, 2001,79(8)
- Wood, W.G. (2001) Hereditary persistence and fetal hemoglobin and $\delta\beta$ thalassemia. In *Disorders of Hemoglobin* (Steinberg, M.H., Forget, B.G., Higgs, D.R., and Nagel, R.L., eds.) pp. 356-388. Cambridge University Press, Cambridge.

PRINCIPLE OF THE PROCEDURE

The test is based on the reverse-hybridization principle, where specific oligonucleotide probes immobilized as parallel lines on membrane-based strips hybridize with biotinylated PCR products. The exact match between probes and amplified product generates a signal exploiting the bond between biotin and streptavidin conjugated with alkaline phosphatase and a subsequent color developer.

The analysis consists of three steps:

1. DNA extraction from non-coagulated blood or chorionic villus samples using extraction kit (see the procedure in the following *DNA Extraction* paragraph);
2. Multiplex PCR of target DNA sequences using biotinylated primers;
3. Hybridization of the amplification products with the oligonucleotide probes and colorimetric detection.

The device allows simultaneous detection of the following variations:

	Mutation	HGVS nomenclature	Category
1.	anti 3.7	-	
2.	-86 C>A	HBB:c.-136C>A	β^+
3.	-29 A>G	HBB:c.-79A>G	β^+
4.	cap+1 A>C	HBB:c.-50A>C	β^+
5.	cap+33 C>G	HBB:c.-18C>G	β^+
6.	cod 7 A>G HbG-San José	HBB:c.23A>G	Hb variants
7.	cod 19 C>A HbD-Ouled Rabah	HBB:c.60C>A	Hb variants
8.	cod 26 G>A HbE	HBB:c.79G>A	Hb variants
9.	cod 47 G>A HbG-Copenaghen	HBB:c.142G>A	Hb variants
10.	cod 104 G>C Hb Camperdown	HBB:c.315G>C	Hb variants
11.	cod 121 G>C HbD-Punjab	HBB:c.364G>C	Hb variants
12.	cod 121 G>A Hb O-Arab	HBB:c.364G>A	Hb variants
13.	cod 126 T>G Hb Neapolis	HBB:c.380T>G	Hb variants
14.	Hb Lepore-BW	NG_000007.3:g.63632_71046del	β^+
15.	Sicilian ($\delta\beta$) ⁰ deletion	NG_000007.3:g.64336_77738del13403	β^0

PRODUCT COMPOSITION

⇒ BOX A (store at -25/-15° C)

Reagents	NLM Code	Quantity	N° Vial
PCR Mix Solution containing oligonucleotides	KA888-n	800 µl	1

⇒ BOX B (store at +2/+8° C)

Reagents	NLM Code	Quantity	N° Vial
Strips Nitrocellulose membrane coated with oligonucleotides	-	20	1
DNAT Denaturation solution containing NaOH	 DANGER	H314, P280, P301+P330+P331, P303+P361+P353, P305+P351+P338	KA107-n 0,3 ml 1
Hybridization (*) Saline Solution containing preservatives	 WARNING	H319, EUH208, P280, P305+P351+P338, P337+P313	KA102-n 48 ml 1
Stringent Wash Solution Saline solution containing detergent and preservatives	-	-	KA928-n 136 ml 1
Conjugate Solution containing streptavidin labeled with alkaline phosphatase, stabilizers and preservatives	-	-	KA574-n 48 ml 1
Wash Solution B Saline Solution containing preservatives	-	EUH208	KA105-n 136 ml 1
Color developer Solution containing 5-Bromo-4-chloro-3-indolyl-phosphate and 4-Nitroblue tetrazolium	-	-	AA564-n 48 ml 1
Transparent interpretation card	-	-	1
Collector sheet For strip storage, not suitable for results interpretation and/or as medical report	-	-	1

(*) The exemption from the requirements for labeling and packaging is determined by the quantity and severity of danger reagent. The criteria are: **amounts <125 ml and categories not serious (ref. Regulation (EC) N° 1272/2008 - list 1.5.2.1.1).**

The components volume show above refer to the standard size of the kit. Reduced packages of the kit are available on request for evaluation and / or demonstration of the product.

STABILITY AND STORAGE

- All opened or unopened reagents are stable up to the expiry date indicated on the label when stored at the correct temperature (**BOX A -25/-15° C; BOX B +2/+8° C**).
- At the end of each assay store reagents at the appropriate temperature.
- The kit should be kept isolated from any source of contaminating DNA, especially amplified products.
- Developed dry strips should be stored in the dark at room temperature (+15/+30° C).
- The vial containing the DNAT should be closed immediately after use.

PRECAUTIONS

- Only professional and opportunely trained personnel should use this kit. Handle this product according to established good laboratory practices and universal precautions; wear personal protective apparel.
- When handling DNA samples or amplification products, use aerosol-resistant pipette tips to avoid pipettes contamination.
- Discard used materials as bio hazardous waste.
- Plastic trays can be reused maximum twice, if washed with water at the end of detection.
- Do not eat, drink, smoke or apply cosmetics in areas where reagents or specimens are handled.
- If skin or mucous membrane exposure occurs, immediately wash the area with copious amount of water. Seek medical advice.
- Do not use components beyond the expiration date.
- Test strips stored at +2/+8° C in the dark are stable until the expiry date.
- Do not mix reagents from different lots.
- Reagents have to be preserved separated from possible contaminants (as DNA samples and amplification products).
- It is recommended to work in three separated areas:
 - Area 1: pre-PCR (samples handling and extraction).
 - Area 2: Master Mix preparation.
 - Area 3: post-PCR (amplification and detection).
- Don't use the device if the box is damaged; contact the supplier.
- **It is advisable to have constant and uniform laboratory temperature, avoid to place the instruments near heating/cooling sources that may compromise the correct working.**
- **Use the transparent interpretation card of the specific lot.**
- **Hazard statements**
 - **H314:** Causes severe skin burns and eye damage.
 - **H319:** Causes serious eye irritation.
- **Precautionary statements**
 - **P280:** Wear protective gloves/protective clothing/eye protection/face protection.
 - **P301+P330+P331:** IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
 - **P303+P361+P353:** IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

- **P305+P351+P338:** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
- **P337+P313:** If eye irritation persists: Get medical advice/attention.

- **Supplemental hazard informations**

- **EUH208:** It contains CMIT/MIT. May cause allergic reaction (skin).

MATERIALS REQUIRED BUT NOT SUPPLIED

AREA 1

- DNA Extraction Kit
- Vertical downflow air box
- Dedicated adjustable volume pipettes set
- Disposable aerosol-resistant pipettes tips and 1,5 ml tubes

AREA 2

- Vertical downflow air box
- Dedicated adjustable volume pipettes set
- Disposable aerosol-resistant pipettes tips
- 0,2 ml PCR tubes
- DNA polymerase with **orange cap**

AREA 3

- Validated thermal cycler or equivalent instrument*
- Agarose gel electrophoresis equipment (optional)
- Dedicated adjustable volume pipettes set and tips
- Clean tweezers to handle strips
- Plastic supports (trays) for manual or automatic detection

Manual strip processing:

- Water bath with inclined lid, shaking platform (50-100 rpm) and adjustable temperature (45 ± 0,5 °C)
- Thermometer
- Vacuum aspiration apparatus
- Deionized or distilled water
- Timer

Automatic strip processing:

- Tecan ProfiBlot 48, Dynex Dynablot Heat or equivalent instrument.

* **WARNING:** if thermal cyclers different from those mentioned above are used, verify that the ramp rate (heating/cooling) is ≤ 3° C/sec. For T100 (BioRad) set 4° C/sec as ramp rate.

PROCEDURE

DNA EXTRACTION

Manual purification: NLM code AA1001* for DNA extraction from blood samples and AA340** (QIAamp DNA Mini Kit, Qiagen) for DNA extraction from chorionic villus samples.

* AA1001: Use EDTA or citrate only as anticoagulant, do not use heparin. It is possible to use fresh blood, blood stored at +2/+8°C for no more than 2 days or stored at -25/-15°C for longer storage periods. The use of frozen blood is recommended since freezing and thawing enhance blood cells lysis.

** AA340: Add about 80µl of sterile PBS to 10-25 mg of chorionic villus biopsy, contained in a 1.5 ml eppendorf. Mix well with the tip in order to homogenize the sample, and, if possible, to partially degrade the biopsy. Then add 100 µl of ATL lysis buffer and proceed to the next step. Depending on the amount of starting biopsy resuspend in a final volume of 50-200 µl of AE elution buffer. Keep biopsies -25/-15°C.

For extraction with automatic systems follow the specific manufacturer's IFUs.

Any other type of DNA extraction with similar results of concentration and purification can be used.

AMPLIFICATION

Keep all the reagents for the amplification on ice during the execution of the whole procedure. Avoid repeated freeze-thaw cycles of the PCR Mix.

Master Mix preparation

Reagents	Volume/sample
PCR Mix	26,6 µl
DNA polymerase (5U/µl), orange cap	0,4 µl

- Prepare the Master Mix for the number of samples extracted + 1 more volume (if "n" ≤ 10) or + 2 more volumes (if "n" >10).
- Mix gently and dispense **27 µl** of Master Mix in the previously marked 0,2 ml test tubes. Discard the remaining Master Mix.
- Add **3 µl** of extracted DNA to each tube and mix pipetting up and down (**DNA concentration to be used should be 10-40 ng/µl**).

Note: Always dilute 1:3 sample obtained with Maxwell (Promega) and Magcore (RBCBioscience) before use.

Amplification program

- Preheat the thermal cycler lid.
- Put reaction tubes in the thermal cycler and run the following amplification program:

Temperatura	Tempo	Cicli
97° C	3 min	1
97° C	40 sec	15
66° C*	50 sec	
72° C	1min 30 sec	
97° C	40 sec	25
56° C	50 sec	
72° C	1min 30 sec	
72° C	3 min	1
10° C	∞	

* for PTC 100 MJ Research, AB9700 (with sample tray only) and EPGradient S set up the annealing temperature to 65° C.

It's also possible to amplify the samples using the amplification program of ALPHA GLOBIN TEST (NLM code AC099).

Keep PCR products at + 2/+8° C if detection is done in the same day, or store them at -25/-15° C.


Optional: analyze amplification products by 2% agarose gel containing a DNA intercalating agent.

Mutation	Amplicon Length (bp)
PCR control	182
-86 C>A, -29 A>G, cap+1 A>C, cap+33 C>G, HbG-San José, HbD-Ouled Rabah, HbE	343
HbG-Copenaghen, Hb Camperdown, Hb Lepore-BW	413
HbD-Punjab, Hb O-Arab, Hb Neapolis	506
Amplicons 343pb +413pb	770
Anti 3.7	1700
Sicilian ($\delta\beta$) ⁰	499

Example:

Agarose gel 2%

- 1) Wt
- 2) Anti 3.7
- 3) Hb Lepore-BW eterozygous
- 4) Siciliana($\delta\beta$)⁰ eterozygous
- 5) Hb Lepore-BW homozygous
- 6) Siciliana($\delta\beta$)⁰ homozygous

Warning:

You can also observe the presence of amplified approximately 770 and 1700 bp. The presence or not of these bands does not affect the results.

The amplicon of 1700 bp does not necessarily correspond to anti 3.7.

The amplicon of 182 pb could be not visible with gel

AUTOMATIC DETECTION

⇒ Detection can be performed by using automatic instrument (ProfiBlot T48, Dynablot or similar instrument): select the **COAGUL 5** program. In this case **8 µl PCR product + 8 µl DNAT** and 1,5 ml of each reagent should be used for each sample.

MANUAL DETECTION

- Adjust the water level of the shaking bath to approx. 2/3 of the height of the Plastic Tray.
- Heat the water bath to **45° C** and check the temperature ($45 \pm 0,5^\circ \text{C}$) using a calibrated thermometer.
- Pre-warm the Hybridization and Stringent Wash Solution to **45° C**; all crystals should be completely dissolved.
- Allow DNAT, Conjugate, Wash Solution B and Color developer to reach room temperature.
- Remove one strip for each sample using clean tweezers (touch strips with gloves only). Label strips with a pencil (do not use ballpoint pen, markers, etc).
Warning: do not allow the strips to dry during the whole procedure.

Hybridization (45° C)

- Pipette **8 µl** of **DNAT** into each lane of the Plastic Tray
Warning: do not use DNAT if is not blue.
- For each sample add **8 µl** of **amplification** product to the corresponding drop of DNAT. Mix thoroughly with a pipette (solution will remain blue).
- Incubate **5 min** at room temperature.
- Add **1 ml** of **Hybridization** (pre-warmed to 45° C) to each sample and mix gently: blue color disappears.
- Place strips into the respective lanes with marker lines facing up. Submerge completely.
- Incubate for **30 min** at **45° C** in the shaking water bath (about 50 rpm). Keep the water bath closed to avoid temperature variations.

Stringent Wash solution (45° C)

- After hybridization step remove the tray from the water bath; angle the tray upwards and aspirate the solution from each lane using a pipette or a vacuum aspiration apparatus, without touching the strip surface to avoid its damage.
- Add **1 ml** of **Stringent Wash solution** (pre-warmed at 45° C) to each lane and rinse briefly.
- Remove the solution and add **1 ml** of **Stringent Wash solution** (pre-warmed at 45° C).
- Incubate **10 min** at **45° C** in the shaking water bath (about 50 rpm).
Warning: put Stringent Wash solution back in the water bath during incubation
- Remove the solution and add **1 ml** of **Stringent Wash solution** (pre-warmed at 45° C).
- Incubate **10 min** at **45° C** in the shaking water bath (about 50 rpm).
- Remove the solution from each lane.

Color development (room temperature)

All the next steps have to be performed at room temperature with shaking. Since water takes too much time to reach room temperature, place a tray (for example of polystyrene) on the lid of the water bath to isolate the Plastic Trays from the hot water and use the agitation function of the water bath only.

- Add **1 ml** of **Conjugate**.
- Incubate **20 minutes** while shaking (about 50 rpm) at room temperature.
- Remove the solution, add **1 ml** of **Wash solution B** and rinse briefly.
- Remove the solution and add **1 ml** of **Wash solution B**.
- Incubate for **5 min** while shaking (about 50 rpm) at room temperature.
- Remove the solution and add **1 ml** of **Wash solution B**.
- Incubate for **5 min** while shaking (about 50 rpm) at room temperature.
- Remove the solution and add **1 ml** of **Color developer**.
- Incubate for **20 min** in the dark, while shaking (about 50 rpm) at room temperature.
- Remove color developer and wash strips with distilled water.
- Let strips dry in the dark on absorbent paper before proceeding with the interpretation of results.


RESULTS INTERPRETATION

The following strip scheme shows the position of the different probes on the strip.
A line is considered positive when a purple band appears at the end of the detection.
To identify the genotype of a sample analyze each strip with the help of the transparent interpretation card. In order to correctly align strip and transparent interpretation card, colored marker lines on the top and on the bottom have to be used.

Warning:

- the color intensity could be different between different bands on the same strip but this doesn't affect the interpretation of the results.
- The Conjugate control line (the upper one) indicates the correct reaction of both conjugate and color developer.
- A positive reaction of the PCR control line (the second one) indicates the PCR reaction was successful
- Both Conjugate control and PCR control lines always have to stain positive**

Strip scheme

For each mutation one of the following patterns should be obtained:

- Only wild type line **Normal genotype (wild type)**
- Wild Type and Mutant lines **Heterozygous genotype**
- Only mutant line **Homozygous mutant genotype**

Warning:

- HbD-Punjab and HbO-Arab, since they affect the same codon, share the wt probe (25-W. Cod.121).

In case of double- heterozygous sample for these two mutations, there is no signal on the corresponding wt.

- For single and double gene deletions, several wild type probes discriminate between the heterozygous and the homozygous mutant state:


Deletion	Heterozigous	Homozygous mutant
Hb Lepore-BW	All WT signals present	WT signals from 16 to 23 absent
Sicilian ($\delta\beta$) ⁰	All WT signals present	All Wt signals absent

- This kit does not allow to distinguish between the heterozygous and the homozygous mutant state of the anti 3.7 gene triplication.
- Some mutations, because are very close each other on beta-globin gene, share the same wild type probe with mutation investigate with the kit “Beta Globin Test” (AC091).


Beta globin test (AC091)	Beta globin plus test (AC104)
-92 C>T, -87C>G	-86 C>A
-30 T>A	-29 A>G
codon 5–CT, codon 6 HbC, codon 6 HbS, codon 6–A, codon 8 –AA, codon 8/9 +G	San José
IVS2.1 G>A	Hb Camperdown

In case of homozygous results for one of these mutations (positive signal on the probe changed and no signal on the corresponding wt), we suggest to test the sample with the kit “Globina Beta test”, as the sample may be a double heterozygote.


Examples


Hb Lepore heterozygous


Hb Lepore homozygous


Hb Lepore/cap +33C>G


Hb Lepore/Hb O-Arab


Sicilian ($\delta\beta$)⁰ heterozygous

Sicilian ($\delta\beta$)⁰ homozygous

Sicilian ($\delta\beta$)⁰ /cap +33C>G

Sicilian ($\delta\beta$)⁰ / Hb O-Arab

TROUBLESHOOTING

FALSE NEGATIVE OR LOW SIGNALS

- A wrong amount of amplification product has been used. PCR products concentration could be too low due to an inefficient amplification.
- Weak signals only on the last two wild type probes may result from the use of not pure extracts (OD 260/230 <2.0). In this case, repeat amplification using the samples more diluted.
- Low signals (except Control Line) and false negative results could be caused by too high temperature during hybridization and/or Stringent Wash Solution steps.

Check carefully water bath temperature

- Hybridization and Stringent Wash Solution have not been well pre-warmed or mixed. Discard remaining reagents and repeat tests with new solutions.
- The missing signal on a wt probe can be determined by the presence (in homozygous or double heterozygous) of mutations not investigated by the kit but located on β -globin gene within a few nucleotide of the mutations analyzed. Changes in the nucleotide sequence prevent the amplicon to hybridize to the wt probe.
- The absence of signal on the PCR control band may indicate that the sample was not amplified (absence of signal on all the line of HBB gene) or that an incorrect combination mix / strip occurred: it has been revealed a sample amplified with BETA GLOBIN TEST (NLM code AC091) on the strips of the BETA GLOBIN TEST PLUS (NLM code AC104). In this case, a signal is present on the bands of the HBB gene. Repeat the detection step.

FALSE POSITIVE SIGNALS

- Non-specific signals could be caused by too low temperature during hybridization and/or Stringent Wash Solution step

Check carefully water bath temperature

- Hybridization and Stringent Wash Solution have not been well pre-warmed or mixed. Discard remaining reagents and repeat tests with new solutions.
- Water level of the water bath is too low: be sure that level reaches approximately 2/3 of the height of the Plastic Tray.
- Check carefully the water bath lid during the Hybridization or Stringent Wash Solution, keep it closed.
- Keep Stringent Wash Solution in the water bath during the incubation steps at **45° C**
- In case of a sample homozygous for Sicilian ($\delta\beta$)⁰ deletion, you may see a weak nonspecific signal of wt probes.

SINGLE LINE UNHOMOGENEOUS COLORATION

- Shaking rate during detection steps is too low. Check that strips are totally submerged.

PERFORMANCE CHARACTERISTICS

Diagnostic sensitivity

The diagnostic sensitivity of the kit "BETA GLOBIN TEST PLUS" was assessed by analyzing 66 DNA samples positive for the mutations investigated by the kit. We used 59 samples of DNA obtained from whole blood or chorionic villus and 7 synthetic samples.

All the samples were correctly genotyped, resulting in a diagnostic sensitivity of 100%.

Samples used:

Nr samples	Mutation	HGVS nomenclature
8	anti 3.7	-
2	-86 C>A	HBB:c.-136C>A
4	-29 A>G	HBB:c.-79A>G
3	cap+1 A>C	HBB:c.-50A>C
3	cap+33 C>G	HBB:c.-18C>G
3	cod 7 A>G HbG-San José	HBB:c.23A>G
4	cod 19 C>A HbD-Ouled Rabah	HBB:c.60C>A
5	cod 26 G>A HbE	HBB:c.79G>A
7	cod 47 G>A HbG-Copenaghen	HBB:c.142G>A
3	cod 104 G>C Hb Camperdown	HBB:c.315G>C
6	cod 121 G>C HbD-Punjab	HBB:c.364G>C
3	cod 121 G>A Hb O-Arab	HBB:c.364G>A
4	cod 126 T>G Hb Neapolis	HBB:c.380T>G
6	Hb Lepore-BW	NG_000007.3:g.63632_71046del
5	Sicilian ($\delta\beta$) ⁰ deletion	NG_000007.3:g.64336_77738del13403

Analytical Sensitivity

The analytical sensitivity of the kit "BETA GLOBIN TEST PLUS", expressed as the minimum quantity of the target that can be detected, is equal to 10 ng/μl DNA.

Diagnostic specificity

The specificity was determined by analyzing 68 samples. No false positive results were obtained. In agreement with these data the specificity of the kit is 100%.

Reproducibility:

INTRA-ASSAY

The intra-assay reproducibility was assessed by analyzing at least 3 samples in 2 replicates.

INTER-ASSAY

The reproducibility was assessed by analyzing 4 samples by three different operators. For the tests 4 different lot of the kit were used. Amplification was performed using different thermalcycler.

The tests were performed as described in the method and the reproducibility was 100 %

WARNING

As every system based on amplification and detection of nucleic acids, it is possible to have unexpected results in the presence of unknown genetic mutation in the primer and/or probe specific target regions.