



ColoScape™ Colorectal Cancer Mutation Detection Kit

Instructions For Use

Catalog # DC-30-0024E, DC-30-0006E

Instruction Version: Rev. 1

Date of Revision: February 2017



DiaCarta Inc.
2600 Hilltop Drive
Richmond, CA 94806
TEL: (510) 878-6662
FAX: (510) 735-8636
E-MAIL: information@diacarta.com



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany

Table of Contents

1. Intended Use.....	3
2. Gene Mutations and Colorectal Cancer (CRC).....	4
3. QClamp® Technology for Mutation Detection	4
4. Reagents and Instruments.....	6
• Package Contents	6
• Materials Required but Not Provided in the Kit	7
• Instruments	7
• Handling and Storage	7
• General Considerations.....	8
• Warnings and Precautions	8
5. Instructions for Use.....	9
• DNA Isolation.....	9
• Preparation of Reagents.....	9
• Preparation of Assay Mixes	9
• Suggested Run Layout (96-well plate, 384-well plate, tube strips, or tubes)	10
• Instrument Set-Up	11
• Roche Light cycler 96 and Roche Light cycler 480, BioRad CFX 384.....	11
• Assessment of Real-Time PCR Results	11
• Data Analysis for Light Cycler 480	11
• Judging Validity of Sample Data Based on Internal Control Results	12
• Scoring Mutational Status	12
• Differentiating KRAS c12/KRAS c13 Mutational Status.....	13
6. Assay Performance Characteristics	13
• Analytical Performance	13
• Analytical Accuracy and Comparison to Reference Method	14
• Precision and Reproducibility.....	14
• Analytic Sensitivity (LOD)	17
• Analytic Specificity.....	18
• Cut-Offs	19
• Limit of Blank	19
• Interfering Substances	19
• Multiple Freeze-Thaw Cycles	19
• Shelf-Life.....	19
• Clinical performance of the assay	19
7. Symbols Used in Packaging	20
8. References	21

Intended Use

ColoScape™ is a real-time PCR based *in vitro* diagnostic assay for qualitative detection of colorectal cancer associated biomarkers including *APC* (codons 1309, 1367, 1450), *KRAS* (codons 12 and 13), *BRAF* (codon 600) and *CTNNB1* (codons 41 and 45). The detection kit identifies the presence or absence of mutations in the targeted regions but does not specify the exact nature of the mutation. The detection kits are designed to detect any mutation at or near the stated codon site without specifying the exact nucleotide change. ColoScape™ Colorectal Cancer Mutation Detection kit is for use in diagnostic procedures.

Table 1. List of Mutations and Cosmic Identities Found in Coloscape targeted genes

KRAS			
Exon	Amino Acid Change	Nucleotide change	Cosmic No.
2	G12>A	c.35G>C	522
	G12>R	c.34G>C	518
	G12>D	c.35G>A	521
	G12>C	c.34G>T	516
	G12>S	c.34G>A	517
	G12>V	c.35G>T	520
	G13>D	c.38G>A	532
	G13>C	c.37G>T	527
	G13>R	c.37G>C	529
APC			
Exon	Amino Acid Change	Nucleotide change	Cosmic No.
15	E1309fs*	c.3921_3925delAAAAG	COSM18764
	Q1367*	c.4099C>T	COSM13121
	R1450*	C.4348C>T	COSM13127
CTNNB1			
Exon	Amino Acid Change	Nucleotide change	Cosmic No.
3	p.T41A	c.121A>G	COSM5664
	p.T41I	c.122C>T	
	p.S45P	c.133T>C	COSM5663
	P.S45F	c.134C>T	
	P.S45del	C133-135delTCT	
BRAF			
Exon	Amino Acid Change	Nucleotide change	Cosmic No.
15	p.V600E	c.1799T>A	COSM476
	p.V600K	c.1798_1799GT>AA	COSM473
	p.V600R	c.1798_1799GT>AG	COSM474
	p.V600D	c.1799_1800TG>AT	COSM477

Table 1 shows a list of mutations commonly found in the targeted gene that can be detected by the kit. The kit is to be used by trained laboratory professionals within a laboratory environment.

Gene Mutations and Colorectal Cancer (CRC)

Complex signal pathways are involved in the colorectal cancer pathogenesis such as the WNT and RAS /RAF/MAPK pathways. Genetic and epigenetic changes in the pathway components have been studied extensively in relation to their roles in the initiation and development of CRC. KRAS mutations are found in several cancers including colorectal, lung, thyroid, and pancreatic cancers and cholangiocarcinoma. More than 90% KRAS mutations are located within codons 12 and 13 of exon 2, which may lead to abnormal growth signaling by the p21– ras protein. These alterations in cell growth and division may trigger cancer development as signaling is excessive. KRAS mutations have also been detected in many colorectal cancer patients.

The B-type Raf Kinase (BRAF) protein is a serine/threonine kinase that has important roles in regulating the MAP kinase/ERK signaling pathways, affecting cellular proliferation, differentiation, and programmed cell death. A BRAF mutation is commonly found in many human cancers including melanoma, colorectal cancer, lung cancer, and papillary thyroid carcinoma. The most common mutations in BRAF occur in codon 600, where an amino acid substitution in the activation segment of the kinase domain creates a constitutively active form of the protein. The V600E and V600K mutations are found in high frequencies in human cancer V600E 70–90% and V600K 10–15%. BRAF mutations are generally found in tumors that are wild-type for KRAS.


The adenomatous polyposis coli (APC) gene is a key tumor suppressor gene and APC mutation has been found in most colon cancers. The gene encodes a multi-domain protein that binds to various proteins, including β -catenin, axin, CtBP, Asefs, IQGAP1, EB1 and microtubules. Most (~60%) cancer-linked APC mutations occur in a region referred to as the mutation cluster region (MCR) and result in C-terminal truncation of the protein. Mutations in the tumor suppressor gene APC result in the accumulation of catenin which activates the Wnt signaling pathway, leading to tumorigenesis. APC also plays roles in other fundamental cellular processes including cell adhesion and migration, organization of the actin and microtubule networks, spindle formation and chromosome segregation. Mutations in APC cause deregulation of these cellular processes, leading to the initiation and expansion of colon cancer. APC has been used as a biomarker for early colon cancer detection.

The β -catenin gene (*CTNNT1*) is also an important component of the Wnt pathway. Mutations in the serine or threonine phosphorylation sites in the regulatory domain (exon 3, codon 29–48) of the gene leads to accumulation of the gene product (β -catenin) which activates the Wnt pathway.

QClamp® Technology for Mutation Detection

The QClamp® ColoScape™ mutation detection assay is based on xenonucleic acid (XNA) mediated PCR clamping technology. XNA is a synthetic DNA analog in which the phosphodiester backbone has been replaced by a repeat formed by units of (2-aminoethyl)-glycine. XNAs hybridize tightly to complementary DNA target sequences only if the sequence is a complete match. Binding of XNA to its target sequence blocks strand elongation by DNA polymerase. When there is a mutation in the target site, and therefore a mismatch, the XNA:DNA duplex is unstable, allowing strand elongation by DNA polymerase. Addition of an XNA, whose sequence with a complete match to wild-type DNA, to a PCR reaction, blocks amplification of wild-type DNA allowing selective amplification of mutant DNA. XNA oligomers are not recognized by DNA polymerases and cannot be utilized as primers in subsequent real-time PCR reactions.

Figure 1. Principle of the QClamp® Coloscape™ Mutation Test in targeted genes



Reagents and Instruments

Package Contents

Table 2. Package Contents.

Vial No.	Name of Component	Description	Volume, 24-test kit	Volume, 6-test kit
1	APC c1309 and c1367 Primer/Probe Mix	APC c1309 and c1367 Primers and probe	1 X 62 uL	1 X 15 uL
2	APC c1309 and c1367 XNA	APC c1309 and c1367 XNA	1 X 28 uL	1 X 7 uL
3	BCT c41 Primer/probe Mix	BCT c41 Primers and probe	1 X 62 uL	1 X 15 uL
4	BCT c41 XNA	BCT c41 XNA	1 X 28 uL	1 X 7 uL
5	APC c1450 Primer/Probe Mix	APC c1450 Primer and probe	1 X 62 uL	1 X 15 uL
6	APC C1450 XNA	APC C1450 XNA	1 X 28 uL	1 X 7 uL
7	BCT c45 Primer/probe Mix	BCT c45 Primers and probe	1 X 62 uL	1 X 15 uL
8	BCT c45 XNA	BCT c45 XNA	1 X 28 uL	1 X 7 uL
9	KRAS c12 Primer/Probe mix	KRAS c12 Primer/probe	1 X 62 uL	1 X 15 uL
10	KRAS c12 XNA	XNA for KRAS c 12	1 X 28 uL	1 X 7 uL
11	KRAS c13 Primer/Probe mix	KRAS c13 Primer/probe	1 X 62 uL	1 X 15 uL
12	KRAS c13 XNA	XNA for KRAS c 13	1 X 28 uL	1 X 7 uL
13	BRAF c600 Primer/Probe Mix	BRAF V600 Primers and probe	1 X 62 uL	1 X 15 uL
14	BRAF c600 XNA	BRAF V600 XNA	1 X 28 uL	1 X 7 uL
15	2X Coloscape™ qPCR Master Mix	PCR Reaction Premix	1013 uL	266 uL
16	ColoScape™ Negative Control	Wild-type DNA	1 X 56 uL	1 X 28 uL
17	ColoScape™ Positive Control	APC c1309, c1367, c1450, BCT c41, BCT c45, KRAS c12, KRAS c13 BRAF c600 mutant templates	1 X 56 uL	1 X 28 uL
18	Non template control	Nuclease-Free Water	1 X 76 uL	1 X 56 uL

Materials Required but Not Provided in the Kit

Reagents for DNA Isolation

QIAamp DSP DNA FFPE Tissue Kit (QIAGEN, cat. no. 60404) or equivalent

Consumables

- 0.2 ml DNase-free PCR tubes or plates
- Nuclease-free, low-binding micro centrifuge tubes
- Nuclease-free pipet tips with aerosol barriers

Equipment

- Permanent marker
- Real time PCR instrument
- Dedicated pipettes* (adjustable) for sample preparation
- Dedicated pipettes* (adjustable) for PCR master mix preparation
- Dedicated pipettes* (adjustable) for dispensing of template DNA
- Micro centrifuge
- Bench top centrifuge* with rotor for 1.5 ml tubes
- Vortexer
- PCR rack
- Reagent reservoir
- Distilled water

* Prior to use ensure that instruments have been maintained and calibrated according to the manufacturer's recommendations.

Instruments

The assays have been developed and validated on the instruments shown in the table below. Instrument platforms not listed in the table should be validated by the individual labs. Guidance for validation can be obtained from DiaCarta upon request.

Table 3. List of Instruments Validated with This Kit.

Company	Model
Roche	Light cycler 96
Roche	Light cycler 480
Bio-rad	CFX384

Handling and Storage

This kit is shipped on dry ice. If any component of the kit is not frozen on arrival, the outer packaging has been opened during transit, or the shipment does not contain a packaging note or the reagents, please contact DiaCarta or the local distributors as soon as possible.

The kit should be stored at –20 °C immediately upon receipt, in a constant-temperature freezer and protected from light. When stored under the specified storage conditions, the kit is stable until the stated expiration date. It is recommended to store the PCR reagents (Box 1 and 2) in a pre-amplification area and the controls (Box 3) in a postamplification (DNA template-handling) area.

The kit can undergo up to 6 freeze-thaw cycles without affecting performance.

All reagents must be thawed at ambient temperature for a minimum of 30 minutes before use. Do not exceed 2 hours at ambient temperature. The primer and probe mixes contain fluorophore labeled probes and should be protected from light.

Attention should be paid to expiration dates and storage conditions printed in the box and labels of all components. Do not use expired or incorrectly stored components.

General Considerations

Effective use of real-time PCR tests requires good laboratory practices, including maintenance of equipment that is dedicated to molecular biology. Use nuclease-free lab ware (pipettes, pipette tips, reaction vials) and wear gloves when performing the assay. Use aerosol-resistant pipette tips for all pipetting steps to avoid cross contamination of the samples and reagents.

Prepare the assay mixes in designated pre-amplification areas using only equipment dedicated to this application. Add template DNA in a separate area (preferably a separate room). Use extreme caution to prevent DNase contamination that could result in degradation of the template DNA, or PCR carryover contamination, which could result in a false positive signal.

Reagents supplied are formulated specifically for use with this kit. Make no substitutions in order to ensure optimal performance of the kit. Further dilution of the reagents or alteration of incubation times and temperatures may result in erroneous or discordant data.

Warnings and Precautions

- Use extreme caution to prevent contamination of PCR reactions with the positive and wild type DNA controls provided.
- Minimize exposure of the 2X PCR Master mix to room temperature for optimal amplification.
- Avoid overexposing the primer-probe mixes to light for optimal fluorescent signal.
- Use of non-recommended reagent volumes may result in a loss of performance and may also decrease the reliability of the test results.
- Use of non-recommended volumes and concentrations of the target DNA sample may result in a loss of performance and may also decrease the reliability of the test results.
- Use of non-recommended consumables with instruments may adversely affect test results.
- Do not re-use any remaining reagents after PCR amplification is completed.
- Additional validation testing by user may be necessary when using non-recommended instruments.
- Perform all experiments under proper sterile conditions using aseptic techniques.
- Perform all procedures using universal precautions.
- Wear personal protective apparel, including disposable gloves, throughout the assay procedure.
- Do not eat, drink, smoke, or apply cosmetics in areas where reagents or specimens are handled.
- Dispose of hazardous or biologically-contaminated materials according to the practices of your institution.
- Discard all materials in a safe and acceptable manner, in compliance with all legal requirements.
- Dissolve reagents completely, then mix thoroughly by vortexing.
- If exposure to skin or mucous membranes occurs, immediately wash the area with large amounts of water. Seek medical advice immediately.
- Do not use components beyond the expiration date printed on the kit boxes.
- Do not mix reagents from different lots.
- Return all components to the appropriate storage condition after preparing the working reagents.
- Do not interchange vial or bottle caps, as cross-contamination may occur.
- Keep all the materials on ice when in use.
- Do not leave components out at room temperature for more than 2 hours.
- Reagents supplied are formulated specifically for use with this kit. Make no substitutions in order to ensure optimal performance of the kit. Further dilution of the reagents or alteration of incubation times and temperatures may result in erroneous or discordant data.

Instructions for Use

DNA Isolation

Human genomic DNA must be extracted from fixed paraffin-embedded tissue, frozen tissue or plasma prior to use. Several methods exist for DNA isolation. For consistency, we recommend using a commercial kit, such as Qiagen DNA extraction kit (QIAamp DNA FFPE Tissue Kit, cat No. 56404, for paraffin embedded specimens; DNeasy Blood & Tissue kit, cat. No. 69504 or 69506, for tissue and blood specimens, QIAamp Circulating Nucleic Acid Kit, cat. No. 55114 for plasma). Follow the genomic DNA isolation procedure according to manufacturer's protocol. Sufficient amounts of DNA can be isolated from FFPE blocks or fresh frozen sections as well as plasma (approx. 2–10 µg).

This ColoScape™ assay requires a total of 22.5 – 35 ng of DNA per sample (2.5–5 ng/reaction). After DNA isolation, measure the concentration using fluorometric analysis (i.e. Qubit) and dilute to 1.25–2.5 ng/µl. If using spectrophotometric analysis, make sure the A260/A230 value is greater than 2.0 and A260/A280 value between 1.8 and 2.0.

Preparation of Reagents

A 24-test kit contains enough control material for 3 runs. Thaw all primer and probe mixes, XNAs, Positive Control, WT Negative Control, Nuclease-Free Water and 2X PCR mastermix provided. Thaw all reaction mixes at room temperature for a minimum of 30 minutes. Vortex all components except the PCR Master Mix and Primer and probe Mix for 5 seconds and perform a quick spin. The PCR Master Mix and Primer/probe mix should be mixed gently by inverting the tube a few times. Prior to use, ensure that any precipitate in the PCR Master Mix is re-suspended by pipetting up and down multiple times. Do not leave kit components at room temperature for more than 2 hours. The PCR reactions are set up in a total volume of 10 µl/reaction.

Table 4 shows the component volumes for each 10 ul reaction.

Table 4. QClamp Taqman Assay Components and Reaction Volume.

Components	Volume/Reaction
2X PCR Master mix	5 µl
Primer and probe Mix	2 µl
XNA	1 µl
DNA sample or Controls	2 µl
Total volume	10 µl

For accuracy, 2x PCR Mastermix, primers and XNA should be pre-mixed into assay mixes as described in *Table 5* below

Preparation of Assay Mixes

Assay mixes should be prepared just prior to use. Label a micro centrifuge tube (not provided) for each reaction mix, as shown in *Table 5*. For each control and mutation detection reaction, prepare sufficient working assay mixes for the DNA samples, one Positive Control, one Nuclease-Free Water for No-Template Control (NTC), and one WT Negative Control, according to the volumes in *Table 5*. Include reagents for 1 extra sample to allow sufficient overage for the PCR set up. The assay mixes contain all of the components needed for PCR except the sample.

Each sample requires one reaction for each mutation site detected by the kit and an Internal Control assay. The External Control assay uses ACTB as reference gene to determine if an appropriate level of amplifiable DNA is present in the sample.

Table 5. Preparation of Assay Mixes.

	Volume of 2X PCR Master Mix	Volume of Primer and probe Mix	Volume of XNA*
APC 1309 and 1367 Mix	5 µl x (n+1)	2 µl x (n+1)	1 µl x (n+1)
APC 1450 Mix	5 µl x (n+1)	2 µl x (n+1)	1 µl x (n+1)
BCT 41 Mix	5 µl x (n+1)	2 µl x (n+1)	1 µl x (n+1)
BCT 45 Mix	5 µl x (n+1)	2 µl x (n+1)	1 µl x (n+1)
KRAS11 Mix	5 µl x (n+1)	2 µl x (n+1)	1 µl x (n+1)
KRAS 13 Mix	5 µl x (n+1)	2 µl x (n+1)	1 µl x (n+1)
BRAF Mix	5 µl x (n+1)	2 µl x (n+1)	1 µl x (n+1)

n = number of reactions (DNA samples plus 3 controls). Prepare enough for 1 extra sample (n + 1) to allow for sufficient coverage for the PCR set.

* For the Internal Control Mix, use 1 µl of water provided in the kit in place of XNA. For accuracy, do not pipette less than 10 µl.

A reaction mix containing all reagents except for the DNA samples or controls should be prepared for the total number of samples and controls to be tested in one run. The ColoScape™ Positive Control (PC), ColoScape™ Negative Control (NC), and a No Template Control (NTC) should be included in each run. Negative Controls use wild-type DNA as the template. Wild-type DNA should have no mutations, therefore the XNA probes will bind strongly, blocking the polymerase from making amplicons. However, the External Control Mix with the Negative Control should make amplicons efficiently, providing another way to monitor performance of the primers, polymerase, and sample.

A set of positive controls must also be run with each reaction mix, every time the assay is run. The Positive Control contains one mutant template for each reaction mix. Positive controls contain mutations; therefore XNA probes will not bind, allowing amplification of the mutant template. Positive controls must show the appropriate values for the reaction to be valid. A set of no template control (tube NTC) is run with each reaction mix every time the assay is run. Nuclease-Free Water is used in the place of template. The NTC serves as a negative control and assesses potential contamination during assay set-up.

Suggested Run Layout (96-well plate, 384-well plate, tube strips, or tubes)

Add 8 ul of the appropriate assay mix to the plate or tubes. Add 2 ul of template.

Table 6. Suggested Plate Layout.

	1	2	3	4	5	6
A	NTC APC1309, 1367 Mix	PC APC 1309, 1367 Mix	CC APC 1309, 1367 Mix	S1 APC 1309, 1367 Mix	S2 APC 1309, 1367 Mix	S3 APC 1309, 1367 Mix
B	NTC APC 1450 Mix	PC APC 1450 Mix	CC APC 1450 Mix	S1 APC 1450 Mix	S1 APC 1450 Mix	S3 APC 1450 Mix
C	NTC BCT41 Mix	PC BCT41 Mix	CC BCT41 Mix	S1 BCT41 Mix	S2 BCT41 Mix	S3 BCT41 Mix
D	NTC BCT45 Mix	PC BCT45 Mix	CC BCT45 Mix	S1 BCT45 Mix	S2 BCT45 Mix	S3 BCT45 Mix
E	NTC KRAS 12 Mix	PC KRAS 12 Mix	CC KRAS 12 Mix	S1 KRAS 12 Mix	S2 KRAS 12 Mix	S3 KRAS 12 Mix
F	NTC KRAS 13 Mix	PC KRAS 13 Mix	CC KRAS 13 Mix	S1 KRAS 13 Mix	S2 KRAS 13 Mix	S3 KRAS 13 Mix
G	NTC BRAF 600 Mix	PC BRAF 600 Mix	CC BRAF 600 Mix	S1 BRAF 600 Mix	S2 BRAF 600 Mix	S3 BRAF 600 Mix

PC: Positive Control, NTC: No-Template Control (water), CC: Negative Control (Wild-type DNA), S1-3: Samples 1-3.

Table 6 is a suggested plate set-up for a single experiment analyzing 3 unknown samples. Please disregard any assay mixes listed below that are not part of your kit.

When all reagents have been added to the plate, tightly seal the plate to prevent evaporation. Spin at 1000rpm for 1 minute to collect all the reagents. Place in the real-time PCR instrument immediately.

Instrument Set-Up

Roche Light cycler 96 and Roche Light cycler 480, BioRad CFX 384

1. Use FAM/HEX as the Detector
2. Setup the cycling parameters as shown in the table below
3. Start the run.

Table 7. Roche Light Cycler Platforms Cycling Parameters.

Step	Temperature (°C)	Time (Seconds)	Cycles	Data Collection
Preincubation	95	300	1	OFF
Denaturation	95	20	X50	OFF
XNA Annealing	70	40		OFF
Primer Annealing	64	30		OFF
Extension	72	30		FAM and HEX

Assessment of Real-Time PCR Results

The real-time PCR instrument generates a cycle threshold (Ct) value for each sample. Ct is the cycle number at which a signal is detected above the set threshold for fluorescence. The lower the cycle number at which signal rises above background, the stronger the PCR reaction it represents (**please see MIQE Guidelines under References for more information).

Negative Controls use wild-type DNA as the template. Wild-type DNA contains none of the targeted mutations, therefore the XNA probes will bind strongly, inhibiting the target amplification. In contrast, the Internal Control Assay assessed with HEX channel should make amplicons efficiently for the Negative Control, providing another way to monitor performance of the primers, polymerase, and sample DNA quality/quantity.

The Positive Control contains one mutant template for each assay mix. Positive controls contain mutations; therefore XNA probes will not bind, allowing amplification of the mutant template. Positive controls must show the appropriate values in both HEX and FAM channels for the run to be valid.

A set of no-template control (NTC) is run with each assay mix every time the assay is run. Nuclease-Free Water is used in the place of template. The NTC serves as a negative control and assesses potential contamination during assay set-up.

Data Analysis for Light Cycler 480

For the Light Cycler 480, open the LightCycler480 SW 1.5.1.61 and select Abs Quant/2nd Derivative Max algorithm to analyze the run file data.

No-Template Controls

Verify that no amplification is observed in the no-template controls (NTC) for each of the reaction mixes. Cq should be Undetermined.

For each control or sample, calculate the difference in Cq value between the mutation assay and the External Control Assay as follows:

$$\text{Cq difference } (\Delta\text{Cq}) = \text{Mutation Assay Cq} - \text{External Control Assay Cq}$$

Negative and Positive Controls

For the assay to be valid, the Negative Control and Positive Control must meet the criteria in *Table 8*.

Table 8. Acceptable Values for Positive Controls and Negative Controls.

Assay	Positive Control	Negative Control
Internal Control	$25 < \text{Cq} < 31$	$25 < \text{Cq} < 31$
APC 1309/1367	$\Delta\text{Cq} \leq 3$	$\Delta\text{Cq} > 17$
APC 1450	$\Delta\text{Cq} \leq 4.5$	$\Delta\text{Cq} > 7$
BCT 41	$\Delta\text{Cq} \leq -0.5$	$\Delta\text{Cq} > 7.3$
BCT 45	$\Delta\text{Cq} \leq 1.5$	$\Delta\text{Cq} > 9.2$
KRAS12	$\Delta\text{Cq} \leq 5$	$\Delta\text{Cq} > 10.4$
KRAS 13	$\Delta\text{Cq} \leq 4.5$	$\Delta\text{Cq} > 9$
BRAF V600	$\Delta\text{Cq} \leq 4$	$\Delta\text{Cq} > 8.3$

Judging Validity of Sample Data Based on Internal Control Results

The Cq value of the Internal Control Mix serve as an indication of the purity and concentration of DNA in each well. Thus, the validity of the test can be decided by the Cq value of the Internal Control mix. Cq values of any sample with Internal Control Mix should be in the range of $25 < \text{Cq} < 31$. If the Cq values fall outside this range, the test results should be considered invalid. The experiment should be repeated following the recommendations in Table 9.

Table 9. Acceptable Internal Control Cq Ranges for Samples.

Validity	Cq Value of Internal Control Mix	Descriptions and Recommendations
Optimal	$25 < \text{Cq} < 31$	The amplification and amount of DNA sample were optimal.
Invalid	$\text{Cq} < 25$	Possibility of a false positive is high. Repeat the PCR reaction with less DNA.
Invalid	$\text{Cq} \geq 31$	Not enough DNA or DNA not pure. The amplification is not optimal. Check DNA amount and purity. Repeat the experiment with more DNA or a new DNA prep may be required

Scoring Mutational Status

If a Cq value is Undetermined, assign a Cq of 50 and proceed to analysis.

The table below should be used to determine mutational status.

Table 10. Scoring Mutational Status for Roche LC480 system

Mutation	APC c1309/1367	APC c1450	BCT c41	BCT c45	KRAS c12	KRAS c13	BRAF c600
Positive:	<16.4	< 7.6	< 8.40	< 9.5	< 10.8	< 10.9	< 8.5
Negative	≥16.4	≥ 7.6	≥ 8.40	≥ 9.5	≥ 10.8	≥ 10.9	≥ 8.5

Table 11. Scoring Mutational Status for BioRad CFX384 system

Mutation	APC c1309/1367	APC c1450	BCT c41	BCT c45	KRAS c12	KRAS c13	BRAF c600
Positive:	<10	< 7	< 8.1	< 10.5	< 10.8	< 9.4	< 8.2
Negative	≥10	≥ 7	≥ 8.1	≥ 10.5	≥ 10.8	≥ 9.4	≥ 8.2

Differentiating KRAS c12/KRAS c13 Mutational Status

The KRAS c12 reaction mix detects both KRAS c12 and KRAS c13 mutations, whereas the KRAS c13 reaction mix detects only KRAS c13 mutations. Therefore, in order to differentiate between KRAS c12 and KRAS c13 Mutations a combination of results from the 2 mixes should be used, as described in *Table 12* below.

Table 12. Interpretation of G12/G13 Mutational Status

Reaction Mix	Result Based on Tables 10, 11	Mutational Status
KRAS c12 Reaction Mix KRAS c13 Reaction Mix	Positive Negative	G12 Mutation
KRAS c12 Reaction Mix KRAS c13 Reaction Mix	Positive Positive	G13 Mutation
KRAS c12 Reaction Mix KRAS c13 Reaction Mix	Negative Positive	G13 Mutation

Assay Performance Characteristics

The performance characteristics of this product were established on the, Roche LightCycler 96, Roche LightCycler 480 and BioRad CFX 384 real-time PCR instruments.

Analytical Performance

The specific performance characteristics of the ColoScape™ kit were determined by studies involving APC, Beta Catenin (CTNNB1), KRAS and BRAF-defined genomic DNA reference samples and FFPE samples of the cell lines with defined mutations obtained from Horizon Diagnostics (Cambridge, England) and ATCC. These samples have been characterized genetically as containing heterozygous or homozygous mutations in the coding sequence of the respective target regions. These single nucleotide polymorphisms in the target regions of the four genes in the assay have been confirmed by genomic DNA sequencing. Additional samples consisted of formalin-fixed, paraffin-embedded (FFPE) reference and patient tissue samples, as well as wild-type DNA (no KRAS mutations).

Analytical Accuracy and Comparison to Reference Method

QClamp analytical accuracy is verified and validated through testing of well-characterized samples with known mutations verified by sequencing. Three studies were done to demonstrate concordance in mutation status of FFPE samples tested with the ColoScape™ Mutation Test kit relative to sequencing. A set of samples were chosen for evaluation based on mutation status. Samples were chosen blindly to test with the ColoScape™ CRC mutation Test kit for comparison to the mutation status returned from sequencing. The results demonstrated a 100% match between sequencing and the ColoScape™ Mutation Test kit. A total of 200+ samples were tested in the accuracy and reproducibility studies, 64 of them represent unique FFPE patient samples. 8 gDNA mixed positive controls and 3 FFPE controls were tested multiple times. In samples with sufficient amount of DNA positive samples were correctly detected in >97% of tests, negative samples were detected correctly in >96% cases of samples.

Precision and Reproducibility

Precision of the ColoScape™ kit was determined with defined analyte levels of gDNA of known mutational status and allelic frequencies. To establish lot-to-lot variation, a reproducibility study was performed using two different kit lots. Each lot was tested on one wild-type control and two reference samples containing each mutation at 5% and 1% allelic frequency in nine replicates. Tests were performed on three separate dates. Inter-assay %CV was established with results from two different users at two different instruments performing the assay using the same lot of reagents, with tests run one to two times a day for three days. Intra-assay %CV was established through performance of kit on reference samples run in replicates of nine separately for two different qPCR instrument types. All testing was done using sequence-verified samples from Horizon Diagnostics and ATCC.

Table 13. Summary of Reproducibility Results

Variation	%CV
Intra-assay	≤ 6%
Inter-assay	≤ 8%
Lot-to-Lot Variation	≤ 5%
Operator Variability	≤ 3%

Reproducibility is demonstrated based on %CV of Cq values and rate of % correct mutation calls for all assays on two lots and operators for both Roche and BioRad instruments.

Table 14. Summary of Intra-assay Reproducibility Results on LC480.

Target	WT			5% mutant			1% mutant		
	Average Cq	SD	%CV	Average Cq	SD	%CV	Average Cq	SD	%CV
APC1309	50	0	0.0%	30.2	0.2	0.6%	32.6	0.2	0.6%
APC1367	50	0	0.0%	30.1	0.1	0.4%	32.6	0.3	1.0%
APC1450	42.3	4.3	10.1%	31.4	1.6	5.2%	34.6	0.3	0.8%
BTC 41	39.8	2.2	5.5%	27.9	0.1	0.4%	30.4	0.3	0.8%
BTC 45	39.6	1.3	3.3%	30.4	0.2	0.7%	32.9	0.6	1.7%
KRAS 12	42.2	2	4.7%	32.2	0.4	1.1%	35.6	0.9	2.5%
KRAS 13	42.1	1.2	2.9%	33.3	0.3	0.9%	36.6	0.6	1.7%
BRAF 600	39.1	0.9	2.3%	31.9	0.2	0.7%	35	0.6	1.6%
Internal Control	28.7	0.8	0.3%	28.9	0.8	2.1%	28.6	0.8	2.1%

Table 15. Summary of Inter-assay Reproducibility Results on LC480.

Assay	Run 1	Run 2	Run 3	Average	SD	CV%
	WT					
APC1309/1367	50.0	50.0	50.0	50.0	0.0	0.0%
APC1450	42.3	42.5	43.4	42.7	0.6	1.3%
BTC 41	39.8	40.0	39.7	39.8	0.1	0.3%
BTC 45	39.6	40.4	41.8	40.6	1.1	2.7%
KRAS 12	42.2	43.4	45.3	43.6	1.6	3.6%
KRAS 13	42.1	41.6	42.2	42.0	0.3	0.8%
BRAF 600	39.1	38.3	39.4	38.9	0.6	1.6%
Internal Control	28.7	29.1	29.2	29.0	0.3	0.9%
	5% mutation					
APC1309/1367	30.1	30.3	30.2	30.2	0.1	0.3%
APC1450	31.4	32.0	32.7	32.0	0.5	1.6%
BTC 41	27.9	28.2	27.6	27.9	0.3	0.9%
BTC 45	30.4	30.6	29.3	30.1	0.6	1.9%
KRAS 12	32.2	33.0	31.7	32.3	0.5	1.7%
KRAS 13	33.3	33.9	32.0	33.0	0.8	2.4%
BRAF 600	31.9	32.6	32.2	32.3	0.3	0.9%
Internal Control	29.3	29.7	28.4	29.1	0.5	1.8%
	1% mutation					
APC1309/1367	32.6	32.9	32.7	32.8	0.1	0.4%
APC1450	34.6	34.5	34.8	34.6	0.1	0.4%
BTC 41	30.4	30.6	30.3	30.4	0.1	0.4%
BTC 45	32.9	33.0	31.8	32.6	0.6	1.7%
KRAS 12	35.6	36.3	34.6	35.5	0.7	1.9%
KRAS 13	36.6	37.3	33.8	35.9	1.5	4.2%
BRAF 600	35.0	34.9	34.3	34.7	0.3	0.9%
Internal Control	28.8	29.0	28.8	28.9	0.1	0.3%

Table 16. Summary of Lot-to Lot Reproducibility Results on Roche LC480

WT					
Assay	Lot 1	Lot 2	Average	St Deviation	%CV
APC1309	50.0	50.0	50.0	0.0	0.0%
APC1367	44.3	42.4	43.3	1.3	3.0%
APC1450	40.2	39.9	40.0	0.2	0.5%
BTC CD41	41.7	40.0	40.8	1.2	3.0%
BTC CD45	44.5	42.8	43.7	1.2	2.8%
KRAS CD12	42.8	41.9	42.3	0.7	1.6%
KRAS CD13	39.9	38.7	39.3	0.9	2.3%
BRAF V600	50.0	50.0	50.0	0.0	0.0%

Internal Control	29.0	28.9	28.9	0.1	0.3%
5% Mutation					
APC1309	30.6	30.2	30.4	0.3	1.0%
APC1367	33.1	31.7	32.4	1.0	3.2%
APC1450	27.6	28.1	27.8	0.3	1.2%
BTC CD41	29.4	30.5	29.9	0.8	2.6%
BTC CD45	32.3	32.6	32.5	0.2	0.7%
KRAS CD12	32.0	33.6	32.8	1.1	3.4%
KRAS CD13	32.2	32.3	32.2	0.0	0.1%
BRAF V600	30.6	30.2	30.4	0.3	1.0%
Internal Control	28.8	29.5	29.1	0.5	1.6%
1% Mutation					
APC1309	33.0	32.8	32.9	0.1	0.4%
APC1367	35.0	34.5	34.8	0.4	1.0%
APC1450	30.3	30.5	30.4	0.1	0.4%
BTC CD41	31.7	33.0	32.4	0.9	2.7%
BTC CD45	35.4	35.9	35.6	0.4	1.1%
KRAS CD12	34.5	37.0	35.7	1.8	5.0%
KRAS CD13	34.5	34.9	34.7	0.3	0.8%
BRAF V600	33.0	32.8	32.9	0.1	0.4%
Internal Control	28.8	28.9	28.8	0.1	0.2%

Table 17. Summary of Lot-to Lot Reproducibility Results on BioRad CFX384

Assay	Lot 1	Lot 2	Average	St Deviation	%CV
WT					
APC1309	48.5	48.5	48.5	0.0	0.0%
APC1367	47.2	47.3	47.3	0.1	0.1%
APC1450	43.5	39.5	41.5	2.0	4.7%
BTC CD41	40.7	40.1	40.4	0.3	0.7%
BTC CD45	44.3	40.7	42.5	1.8	4.2%
KRAS CD12	45.4	38.5	41.9	3.4	8.2%
KRAS CD13	40.7	39.9	40.3	0.4	1.0%
BRAF V600	37.6	38.0	37.8	0.2	0.4%
Internal Control	29.0	29.0	29.0	0.1	0.2%
5% Mutation					
APC1309	27.5	28.1	27.8	0.3	1.0%
APC1367	27.6	27.9	27.7	0.2	0.6%
APC1450	29.8	29.6	29.7	0.1	0.2%
BTC CD41	26.2	26.5	26.4	0.2	0.7%
BTC CD45	29.3	30.2	29.8	0.4	1.4%
KRAS CD12	29.8	29.4	29.6	0.3	0.9%
KRAS CD13	30.1	31.6	30.8	0.8	2.5%
BRAF V600	31.0	31.2	31.1	0.2	0.5%
Internal Control	28.8	26.8	26.8	0.3	1.3%

1% Mutation					
APC1309	29.79	30.54	30.17	0.38	1.25%
APC1367	30.16	30.99	30.58	0.58	1.90%
APC1450	32.42	32.05	32.24	0.18	0.57%
BTC CD41	29.06	28.97	29.02	0.06	0.22%
BTC CD45	31.59	32.71	32.15	0.56	1.75%
KRAS CD12	32.33	31.78	32.06	0.39	1.21%
KRAS CD13	32.21	34.22	33.22	1.00	3.02%
BRAF V600	32.10	33.30	32.70	0.84	2.58%
Internal Control	27.27	27.27	27.27	0.49	1.83%

Analytic Sensitivity (LOD)

Table 18. LOD Summary: Reference DNA Samples.

Target mutation	DNA Input, ng/well			
		10	5	2.5
APC 1309		% Correct Call	% Correct Call	% Correct Call
	1% mutation	100%	94%	44%
	0.5% mutation	44%	67%	33%
	0.10% mutation	0%	0%	0%
APC 1367	1% mutation	100%	100%	100%
	0.5% mutation	100%	100%	100%
	0.10% mutation	100%	67%	56%
APC 1450	1% mutation	100%	100%	100%
	0.5% mutation	100%	100%	67%
	0.10% mutation	0%	0%	0%
BCT 41	1% mutation	100%	100%	100%
	0.5% mutation	100%	100%	100%
	0.10% mutation	100%	100%	78%
BCT 45	1% mutation	100%	100%	100%
	0.5% mutation	100%	100%	100%
	0.10% mutation	100%	73%	22%
KRAS 12	1% mutation	100%	100%	100%
	0.5% mutation	100%	100%	67%
	0.10% mutation	87%	33%	24%
KRAS 13	1% mutation	100%	100%	100%
	0.5% mutation	100%	100%	100%
	0.10% mutation	56%	56%	56%
BRAF V600	1% mutation	100%	100%	78%
	0.5% mutation	100%	100%	67%
	0.10% mutation	44%	56%	56%

To determine the limit of detection (LOD) for the kit, the ColoScape™ assay was run using a serial dilution of mutant DNA in wild-type background. Mutant samples were sequence verified by Horizon Diagnostics and ATCC. Mutant concentrations tested were 5, 1, 0.5 and 0.1%. Results demonstrated effective clamping of wild type, providing reproducible detection of mutations at concentrations as low as 0.1% for some targets.

DNA isolated from patient FFPE and plasma samples and reference DNA with known mutational status were used to determine the limit of detections (analytical sensitivity) of this assay. Broad-range serial dilutions of DNA extracted from FFPE and plasma reference materials with known mutations controls were made using wild-type human DNA from FFPE samples to yield various percentages of mutant allele over wild-type background. The assay was performed in multiple replicates (9) near the LOD. For all tested DNA inputs from 10ng/well to 2.5ng/well all target mutations were detected with 100% correct calls at 5% allelic frequency. Table 18 below shows data on 1%, 0.5% and 0.1% mutation detection at various DNA input levels.

Based on the data from the tables below, recommended DNA input is 5ng/well. APC1367, 1450, BTC 41, 45 and KRAS12, 13 can be detected accurately at 0.5% allelic frequency with 2.5 ng/well.

Recommended input of FFPE should not be higher than 20 ng/well due to possible PCR inhibition. Optimal FFPE sample input is between 25 and 31 Cq of the Internal Control reaction.

Analytic Specificity

Analytic specificity of the ColoScape™ test was determined as both the correct calling of the samples with no mutation at different concentrations of WT template and as cross-reactivity of the assays within the kit when one or more mutations are present.

Table 19 shows that all the assays in the kit are accurate in detecting the target mutations. Only in case of KRAS 13 mutations present both KRAS 12 and KRAS 13 signals were detected, but the pattern of lower Cq in KRAS 12 assay and higher Cq in KRAS 13 assay allows for discerning the two targets.

For known WT templates specificity of the assay at 5ng per well is over 97%. There were no false positive calls for up to 320ng of gDNA per well and up to 20ng FFPE DNA. Higher amounts of DNA input were not tested.

For the purposes of detecting colorectal cancer in a test sample, presence of any of the mutations detected by ColoScape™ kit is considered positive signal.

Table 19. Analytic specificity: Cross-reactivity Summary

Assay	Expected mutations in tested 50% templates							
	APC 1309	APC 1367	APC 1450, KRAS12	BCT 41	BCT 45, KRAS 13	KRAS 12	KRAS 13	BRAF 600, BCT 45, KRAS 13
APC 1309	+	-	-	-	-	-	-	-
APC 1367	-	+	-	-	-	-	-	-
APC 1450	-	-	+	-	-	-	-	-
BCT 41	-	-	-	+	-	-	-	-
BCT 45	-	-	-	-	+	-	-	+
KRAS 12	-	-	+	-	*	+	*	*
KRAS 13	-	-	-	-	+	-	+	+
BRAF V600	-	-	-	-	-	-	-	+

* 6Cq higher FAM values than KRAS13

Cut-Offs

Along with studies for analytical accuracy, both reference DNA and FFPE samples were tested to establish cut-offs for the assay. Cutoffs are presented in *Tables 9-11* of the manual.

Limit of Blank

Two lots of reagents were run without template (NTC) in multiple runs (50 data points for each target) to assess the level of background noise when no template is present. No internal control or mutant signals below Cq 48 were detected in any of the runs.

Interfering Substances

A study was performed to evaluate the impact of potentially interfering substances on the performance of the ColoScape™ assay. Potentially interfering substances tested were paraffin (FFPE) and ethanol. The impact of each substance on resultant ΔCq and mutation status of test samples was determined via spiking experiments conducted at three different concentrations, 0.1%, 1% and 5%. DNA extracted from FFPE samples was tested and shown to be compatible with the assay within established cutoffs. None of the potentially interfering substances evaluated at concentrations encountered in normal use impacted the ability of the ColoScape™ to distinguish between mutation-positive and mutation-negative control samples.

Multiple Freeze-Thaw Cycles

The effect of 1-8 freeze-thaw cycles were tested in ColoScape™ Mutation Test kit reagents. There is no effect up to 6 freeze-thaw cycles on the ColoScape™ Mutation detection kit to distinguish between mutation positive and mutation negative samples.

Caution: Repeated freeze-thaw cycles may decrease the reliability of test results.

Shelf-Life

Approximately 12 months after production of kit-see product labels for actual expiration date.

Clinical performance of the assay

Clinical sensitivity and specificity was tested on the samples extracted from FFPE and plasma of patients with different stages of CRC from normal to advanced adenomas (AA), to colorectal cancer stages 1 through 4.










A sample was considered positive if at least one of the target mutations tested positive based on the cutoffs presented in *Tables 9-11*. TP-true positive, FP- false positive, TN- True negative, FN - false negative, NC- no call/inconclusive.

Table 20. Clinical sample testing summary

Types of Clinical Samples		Calls						Clinical parameter	
		Number of samples	TP	FP	TN	FN	NC	specificity	sensitivity
clinical stage/sample type	CRC, including plasma and FFPE	22	22	0	0	0	0	N/A	100%
	CRC, tissue	10	9	0	0	1	0	N/A	90%
	Advanced Adenomas	10	6	0	0	4	0	N/A	60%
	Non-malignant	22	0	1	20	0	1	95%	N/A
	Non-malignant, FFPE only	10	0	1	9	0	0	90%	N/A
Total	FFPE	67	41	1	20	4	1	95%	91%
	Plasma	24	12	0	11	0	1	100%	100%
	Excluding adenomas	57	35	1	21	0	0	95%	100%
	CRC FFPE	35	35	0	0	0	0	N/A	100%

Symbols Used in Packaging

Table 11. Symbols Used in Packaging

Symbol	Definition
	In Vitro Diagnostic Device
	CE Mark
	Authorized Representative in the European Community
	Catalog Number
	Manufactured By
	Temperature Limitation
	Batch Code
	Expiration Date
	Authorized Representative in the European Community
2016-11-16	Date Format (year-month-day)
2016-11	Date Format (year-month)

References

1. Fodde R, Smits R, Clevers H: APC, signal transduction and genetic instability in colorectal cancer. *Nature Rev Cancer* 2001, 1(1):55–67.
2. Dong et al., Detecting Colorectal Cancer in Stool With the Use of Multiple Genetic Targets. *Journal of the National Cancer Institute*, June 6, 2001, Vol. 93, No. 11.
3. Chang et al., Fast simultaneous detection of K-RAS mutations in colorectal cancer. *BMC Cancer*. 2009 Jun 11;9:179.
4. Kobunai et al., The frequency of KRAS mutation detection in human colon carcinoma is influenced by the sensitivity of assay methodology : A comparison between direct sequencing and real-time PCR. *Biochem Biophys Res Commun*. 2010 Apr 23;395(1):158–62.
5. Kwon et al., Frequency of KRAS, KRAS, and KRAS mutations in advanced colorectal cancers: Comparison of peptide nucleic acid-mediated PCR and direct sequencing in formalin-fixed, paraffin-embedded tissue. *Pathol Res Pract*. 2011 Dec 15;207(12):762–8.
6. Ørum, Henrik., PCR Clamping.. *Curr. Issues Mol. gene mutation by real-time fluorescence polymerase chain reaction and peptide nucleic acid clamping. Analytical Biochemistry* 1998; 260: 142–8.
7. **MIQE Reference: "The MIQE Guidelines: Minimum Information Biol. 2000; 2(1), 27–30.
8. Powell et. al., Detection of the hereditary hemochromatosis for Publication of Quantitative Real-Time PCR Experiments". Stephen A. Bustin et. al., *Clin Chem*. 55 (4): 611–22 (2009).
9. Choi et al., Frequency of BRAF mutations in advanced colorectal cancers: Comparison of peptide nucleic acid-mediated PCR and direct sequencing in formalin-fixed, paraffin-embedded tissue. *Pathol Res Pract* 207(12):762–8, 2011.
10. Sparks AB, Morin PJ, Vogelstein B, Kinzler KW: Mutational analysis of the APC/beta-catenin/Tcf pathway in colorectal cancer. *Cancer Res* 1998, 58(6):1130–1134.
11. Smith G, Carey FA, Beattie J, Wilkie MJ, Lightfoot TJ, Coxhead J, Garner RC, Steele RJ, Wolf CR: Mutations in APC, Kirsten-ras, and p53--alternative genetic pathways to colorectal cancer. *Proc Natl Acad Sci U S A*. 2002, 99 (14): 9433–9438. 10.1073/pnas.122612899.
12. Lüchtenborg et al., Mutations in *APC*, *CTNNB1* and *K-ras* genes and expression of hMLH1 in sporadic colorectal carcinomas from the Netherlands Cohort Study. *BMC Cancer* 2005 5:160
13. Vogelstein B, Fearon ER, Hamilton SR, Kern SE, Preisinger AC, Leppert M, et al. Genetic alterations during colorectal-tumor development. *N Engl J Med* 1988;319:525–32.
14. Chan TL, Zhao W, Leung SY, Yuen ST. B-RAF and KRAS mutations in colorectal hyperplastic polyps and serrated adenomas. *Cancer Res* 2003;63:4878–81.
15. Scholtka B, Schneider M, Melcher R, Katzenberger T, Friedrich D, Berghof-Jäger K, et al. A gene marker panel covering the Wnt and the Ras-Raf-MEK-MAPK signalling pathways allows to detect gene mutations in 80% of early (UICC I) colon cancer stages in humans. *Cancer Epidemiol* 2009;33:123–9.
16. Schneider M, Scholtka B, Gottschalk U, Faiss S, Schatz D, Berghof-Jäger K, et al. Detection of up to 65% of precancerous lesions of the human colon and rectum by mutation analysis of APC, KRAS, B-Raf and CTNNB1. *Cancers* 2011;3:91–105.
17. Steinberg, P.; Scholtka, B. Method for conducting non-invasive early detection of colon cancer and/or of colon cancer precursor cells. US Patent: 0,172,823 A1, 2007.

DiaCarta Inc.

2600 Hilltop Drive

Richmond, CA 94806

Email: information@diacarta.com

Tel: +1 (510) 878-6662

www.diacarta.com