

# Istruzioni per l'uso

Ver. NLM 8784 - 03/01/2013

Biologia Molecolare

Estrazione AA1001 non compresa

## HPV Alto Rischio Quantitativo (A9,A7+) Real Time

HPV High Risk Screen Real-TM Quant 2x

REF

AA1115/100

CE



100 TEST

CND

W0105040503

IVD



**NUCLEAR LASER MEDICINE S.r.l.**

SEDE LEGALE: Via Cascina Conighetto - UFFICI OPERATIVI: Viale delle Industrie, 3 - 20049 SETTALA (MI), ITALIA

Tel.: (+39) 02/95.24.51 - Fax (+39) 02/95.24.52.37 - (+39) 02/95.24.52.38 WEB: [www.nlm.it](http://www.nlm.it) - E-MAIL: [segreteria@nlm.it](mailto:segreteria@nlm.it)

Organizzazione con Sistema di gestione qualità certificata ISO 9001 e Sistema di gestione qualità settore medicale certificato ISO 13485  
(organismo di certificazione IMQ - Certificazione CSQ).



For *in Vitro* Diagnostic Use



# HPV High Risk Screen Real-TM Quant 2x HANDBOOK

Real Time PCR kit for quantitative detection of *Human Papillomavirus* (16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59)



100

## NAME

### HPV High Risk Screen Real-TM Quant 2X

## INTRODUCTION

Genital infection with HPV is one of the most common sexually transmitted diseases (STDs) of viral etiology worldwide (20% - 46% in different countries in sexually active young women).

Cervical cancer is the second most common cancer in women worldwide, and a compelling body of clinical, epidemiological, molecular, and experimental evidence has established the etiological relationship between some sexually transmitted HPV genotypes and cervical neoplasia throughout the world. Based on the frequency of detection of HPV genotypes from different grades of Cervical Intraepithelial Neoplasia (CIN Grades I – III), HPV genotypes are subdivided into High-risk HPV types (16, 18, 31 and 45), Intermediate-risk types (33, 35, 39, 51, 52, 56, 58, 59, and 68), and Low-risk types (6, 11, 42-44).

Several methods have been used to diagnose clinical or subclinical infection with HPVs including clinical observation, cytological screening by Pap smear, electron microscopy, immunocytochemistry, but these methods have some disadvantages such as non-standardization and subjectivity, insufficient sensitivity and low predictable value. The most perspective way of HPV diagnosis is a direct detection of DNA of the human papilloma virus of high carcinogenic risk by the polymerase chain reaction. While the value of the Pap smear in routine screening for cervical dysplasia is undisputed, it is now known that 99% of cases of cervical carcinoma are caused by infection with twelve genotypes of the human papilloma virus (HPV). Identification of these high-risk genotypes is very valuable in the management of cervical carcinoma, both as a prognostic indicator and as a secondary screening test where results of a Pap smear are inconclusive. Results from the combination of the Pap smear and the HPV DNA test can aid in determining the intervals for screening.

The PCR-based methods have been used successfully for the detection and typing of genital HPV genotypes in clinical specimens such as cervical swabs or scrapes, cervicovaginal lavages, frozen biopsies and formalin-fixed paraffin-embedded tissues.

## INTENDED USE

Kit **HPV High Risk Screen Real-TM Quant 2x** is an *in vitro* Real Time amplification test for quantitative detection of *Human Papillomavirus* (16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59) in the urogenital swabs. It being known, that the parameter of viral load has a prognostic value and the viral load less than  $10^5$  HPV genomic equivalents in the swab or  $10^3$  genomic equivalents for  $10^5$  cells is considered as insignificant and indicates the presence of transitory infection, however such level of load may have a value only in cases of treatment monitoring. Viral load of more than  $10^5$  genomic equivalents for  $10^5$  cells is considered to be important with high significance and indicates the existence of dysplastic changes or high risk of their occurrence. Quantitative detection of viral load allows to evaluate the character of the infection and to make a forecast concerning the stage of the disease.

HPV High Risk Screen Quant detect the most widespread and oncogenic 12 genotypes of human papilloma virus with determination of clinical significance. Since the human papilloma virus is an intracellular agent, there is need to monitor the presence of cellular material in the sample, in order to avoid false-negative results. HPV High Risk Screen Quant kit contains the internal control (human beta-globine gene), which allows to control the presence of cellular material in the sample.

## PRINCIPLE OF ASSAY

Kit **HPV High Risk Screen Real-TM Quant 2X** is based on two major processes: isolation of DNA from specimens and Real Time amplification of the E1-E2 region of gene HPV (channel *Joe*) and  $\beta$ -globine gene used as Internal Control (channel *Fam*). **PCR-mix-1-FRT A9** tube contains primers directed against regions of A9 groups (HPV types 16, 31, 33, 35, 52, 58) and **PCR-mix-1-FRT A7+** tube contains primers directed against regions of A7 group (18, 39, 45, 59), HPV A5 group (51), A6 group (56). If the swab is not correctly prepared (high quantity of mucous or insufficient quantity of epithelial cells) the Internal Control will not be detected. The kit contains the quantitative standards with known concentration of HPV DNA which allows to determinate the viral load. For the calculation of viral load it is used the relation between the obtained HPV DNA concentration and the quantity of genomic DNA which allows to eliminate the possible errors during the sample preparation.

## MATERIALS PROVIDED

### Module No.1: Real Time PCR kit

Part N° 2 – “HPV High Risk Screen Real-TM Quant 2X ”:

- **PCR-mix-1-FRT A9**, 6 x 0,14 ml;
- **PCR-mix-1-FRT A7+** , 6 x 0,14 ml;
- **PCR-mix-2 buffer**, 6 x 0,30 ml;
- **TaqF DNA Polymerase**, 6 x 0,02 ml;
- **QS HPV K1**, 3 x 0,04 ml (mix HPV DNA C+ 16, 18 and human DNA)\*\*;
- **QS HPV K2**, 3 x 0,04 ml (mix HPV DNA C+ 16, 18 and human DNA);
- **QS HPV K3**, 3 x 0,04 ml (mix HPV DNA C+ 16, 18 and human DNA);
- **Negative Control**, 1,2 ml;\*
- **DNA-buffer (C-)**, 0,5 ml;

Contains reagents for 108 samples.

\* *must be used in the isolation procedure as Negative Control of Extraction.*

\*\* *Standards' concentration is specific for every lot (reported on the HPV High Risk Screen Real-TM Quant Data Card)*

## **MATERIALS REQUIRED BUT NOT PROVIDED**

### **Zone 1: sample preparation:**

- Biological cabinet
- Desktop microcentrifuge for “eppendorf” type tubes
- 65°C ± 2°C dry heat block
- Vortex mixer
- Pipettes with sterile, RNase-free filters tips
- 1,5 ml polypropylene sterile tubes (Sarstedt, QSP, Eppendorf)
- Disposable gloves, powderless
- Tube racks

### **Zone 2: Real Time amplification:**

- Real Time Thermalcycler
- PCR Tubes
- Workstation
- Pipettes with sterile, RNase-free filters tips
- Tube racks

## **STORAGE INSTRUCTIONS**

**HPV High Risk Screen Real-TM Quant 2X** must be stored at -20°C. The kit can be shipped at 2-8°C for 3-4 days but should be stored at 2-8°C and -20°C immediately on receipt.

## **STABILITY**

**HPV High Risk Screen Real-TM Quant 2X** is stable up to the expiration date indicated on the kit label. The product will maintain performance through the control date printed on the label. Exposure to light, heat or humidity may affect the shelf life of some of the kit components and should be avoided. Repeated thawing and freezing of these reagents should be avoided, as this may reduce the sensitivity. Components stored under conditions other than those stated on the labels may not perform properly and may adversely affect the assay results.

## WARNINGS AND PRECAUTIONS



***In Vitro* Diagnostic Medical Device**

For *In Vitro* Diagnostic Use Only

The user should always pay attention to the following:

- Use sterile pipette tips with aerosol barriers and use new tip for every procedure.
- Store extracted positive material (samples, controls and amplicons) away from all other reagents and add it to the reaction mix in a separate area.
- Thaw all components thoroughly at room temperature before starting an assay.
- When thawed, mix the components and centrifuge briefly.
- Use disposable gloves, laboratory coats and eye protection when handling specimens and reagents. Thoroughly wash hands afterwards.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use a kit after its expiration date.
- Dispose of all specimens and unused reagents in accordance with local authorities' regulations.
- Specimens should be considered potentially infectious and handled in a biological cabinet in accordance with appropriate biosafety practices.
- Clean and disinfect all sample or reagent spills using a disinfectant such as 0.5% sodium hypochlorite, or other suitable disinfectant.
- Avoid sample or reagent contact with the skin, eyes, and mucous membranes. If skin, eyes, or mucous membranes come into contact, rinse immediately with water and seek medical advice immediately.
- Material Safety Data Sheets (MSDS) are available on request.
- Use of this product should be limited to personnel trained in the techniques of DNA amplification.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.



Some components of this kit contain sodium azide as a preservative. Do not use metal tubing for reagent transfer.

## SAMPLE COLLECTION, STORAGE AND TRANSPORT

### HPV High Risk Screen Real-TM Quant 2X

- *Cervical swabs:*
  - Remove excess mucus from the cervical os and surrounding ectocervix using a cotton or polyester swab. Discard this swab.
  - Insert the Sampling Cervical Brush 1.0-1.5 centimeters into the cervical os until the largest bristles touch the ectocervix. Do not insert brush completely into the cervical canal. Rotate brush 3 full turns in a counterclockwise direction, remove from the canal.
  - Insert brush into the nuclease-free 2,0 ml tube with 0,3 mL of Transport medium
  - Vigorously agitate brush in medium for 15-20 sec.
  - Snap off shaft at scored line, leaving brush end inside tube.
- *Tissue homogenized with mechanical homogenizer and dissolved in PBS sterile.*
- *Liquid-based cytology samples (Cytoscreen, PreservCyt)*

It is recommended to process samples immediately after collection. Store samples at 2–8 °C for no longer than 24 hours, or freeze at –20/80°C. Transportation of clinical specimens must comply with country, federal, state and local regulations for the transport of etiologic agents.

### PCR REAGENTS PREPARATION:

1. Prepare required quantity of tubes.

Prepare <b>Mix</b> for 14 samples:	Prepare Mix for less than 14 samples (N + 4 controls):
<ul style="list-style-type: none"> <li>• add into the tube with <b>PCR-mix-2 buffer 20 µl</b> of <b>TaqF DNA Polymerase</b> . Vortex the tube and centrifuge for 2-3 sec</li> <li>• add <b>160 µl</b> of the prepared Mix in one <b>PCR-mix-1-FRT A9</b> tube and <b>160 µl</b> in one <b>PCR-mix-1-FRT A7+</b> tube.</li> </ul>	<ul style="list-style-type: none"> <li>• Add in the new sterile tube <b>15*(N+5) µl</b> of <b>PCR-mix-2 buffer</b> and <b>1,0*(N+5)</b> of <b>TaqF DNA Polymerase</b>. Vortex the tube and centrifuge for 2-3 sec</li> <li>• Prepare 2 new tubes and pipette in one tube labeled <b>FRT A9 7*(N+5) µl PCR-mix-1-FRT A9</b> and in the second labeled <b>FRT A7+ 7*(N+5) µl PCR-mix-1-FRT A7+</b>.</li> <li>• Add the half of the prepared mix [(8*(N+5) µl)] <b>PCR-mix-2 buffer/TaqF DNA Polymerase</b> in the tube <b>FRT A9</b> and the second half [(8*(N+3) µl)] the tube <b>FRT A7+</b></li> </ul>
Vortex carefully and centrifuge for 2-3 sec	
<ul style="list-style-type: none"> <li>• Prepare 36 PCR tubes</li> </ul>	<ul style="list-style-type: none"> <li>• Prepare (N*2+8) PCR tubes (2 tubes for each sample and 8 for controls).</li> </ul>
Pipette in the half of tubes <b>15 µl</b> of prepared Mix <b>FRT A9</b> and in the half <b>15 µl</b> of prepared Mix <b>FRT A7+</b> .	

**Example of mixes preparation x sample (with calculation of reagents for controls):**

- Add in new sterile tube:

Samples, N	1	2	3	4	5	6	7
PCR-mix-2 buffer	90	105	120	135	150	165	180
TaqF DNA Polymerase	6	7	8	9	10	11	12
Samples, N	8	9	10	11	12	13	14
PCR-mix-2 buffer	195	210	225	240	255	270	Whole tube
TaqF DNA Polymerase	13	14	15	16	17	18	Whole tube

- Prepare in one sterile tube labeled *FRT A9*:

Samples, N	1	2	3	4	5	6	7
Mix PCR-mix-2 buffer/ TaqF DNA Polymerase	48	56	64	72	80	88	96
PCR-mix-1-FRT A9 (blue cap)	42	49	56	63	70	77	84
Samples, N	8	9	10	11	12	13	14
Mix PCR-mix-2 buffer/ TaqF DNA Polymerase	104	112	120	128	136	142	150
PCR-mix-1-FRT A9 (blue cap)	91	98	105	112	119	126	Whole tube

- Prepare in one sterile tube labeled *FRT A7+*:


Samples, N	1	2	3	4	5	6	7
Mix PCR-mix-2 buffer/ TaqF DNA Polymerase	48	56	64	72	80	88	96
PCR-mix-1-FRT A7 (green cap)	42	49	56	63	70	77	84
Samples, N	8	9	10	11	12	13	14
Mix PCR-mix-2 buffer/ TaqF DNA Polymerase	104	112	120	128	136	142	150
PCR-mix-1-FRT A7 (green cap)	91	98	105	112	119	126	Whole tube

2. Add **10 µl** of **extracted DNA\*** to the tubes with mix **FRT A9** and **10 µl** of **extracted DNA** to the tubes with mix **FRT A7+**.

\*⚠ Re-centrifuge all the tubes with extracted DNA for 2 min at maximum speed (12000-16000 g) and take carefully supernatant. N.B. don't disturb the pellet, sorbent inhibit reaction.

3. Prepare for each mix 2 controls (total 4 PCR tubes):
  - add **10 µl** of **DNA-buffer** as Negative PCR control to the tube labeled  $C_{negA9}$  with prepared mix **FRT A9** and to the tube labeled  $C_{negA7}$  with mix **FRT A7+**;
  - add **10 µl** of each standard (**QS HPV K1**, **QS HPV K2**, **QS HPV K3**) to the 3 appropriate tubes labeled **QS1**, **QS2**, **QS3** with prepared mix **FRT A9** and **10 µl** to the 3 appropriate tubes labeled **QS1**, **QS2**, **QS3** with prepared mix **FRT A7+**;
4. Close tubes and transfer them into the Real Time PCR instrument.
5. Program position of the samples and enter the concentrations of the Quantitative Standards (reported in the Quant Data Card) in the Joe (Yellow)/HEX/Cy3 and Fam (Green) channels in order to generate Standard curves. Use name "Unknown" for the wells that contain samples, K1, K2, K3 for "Standards" and "-" for Negative Controls.

Example of tubes position (for iQ iCycler microplate please check table below):





## Amplification

1. Create a temperature profile on your instrument as follows:

	Rotor-type Instruments <sup>1</sup>			Plate- or modular type Instruments <sup>2</sup>		
Step	Temperature, °C	Time	Repeats	Temperature, °C	Time	Repeats
1	95	15 min	1	95	15 min	1
2	95	5 s	5	95	5 s	5
	60	20 s		60	20 s	
	72	15 s		72	15 s	
3	95	5 s	40	95	5 s	40
	60	20 s fluorescent signal detection		60	30 s fluorescent signal detection	
	72	15 s		72	15 s	

<sup>1</sup> For example Rotor-Gene™ 3000/6000/Q (Corbett Research, Qiagen)

<sup>2</sup> For example, SaCycler-96™ (Sacace), CFX/iQ5™ (BioRad); Mx3005P™ (Agilent), ABI® 7300/7500/StepOne Real Time PCR (Applied Biosystems), SmartCycler® (Cepheid), LineGeneK® (Bioer)

The following programs can also be used:

Rotor-Gene™ 3000/6000/Q (Corbett Research, Qiagen)

	t, °C	Time	Fluorescence detection	Cycles
Hold	95°	15 min	–	1
Hold 2	65°	2 min	–	1
Cycling	95°	20 sec	–	5
	64° Touchdown: 1 deg. per cycle	25 sec	–	
	65°	55 sec	–	
Cycling 2	95°	15 sec	–	40
	60°	25 sec	–	
	65°	40 sec	Fam (Green), Joe (Yellow)	

fluorescence detection on the channels Fam (Green) and Joe (Yellow) on the 2-nd Cycling (65°C)

CFX/iQ5™ (BioRad)

Cycle	Temperature, °C	Time	Fluoresc.detection	Repeats
Cycle 1	95	15 min	–	1
Cycle 3	95	15 s	–	6
	65 Touchdown: 1 deg. per cycle	55 s	–	
	65	25 s	–	
Cycle 4	95	15 s	–	41
	60	55 s	Real-time	
	65	25 s	–	

## INSTRUMENT SETTINGS

### Rotor-type instruments

Channel	Calibrate/Gain Optimisation...	Threshold	More Settings/ Outlier Removal	Slope Correct
FAM/Green	from 4 Fl to 8 Fl	0.03	15 %	On
JOE/Yellow	from 4 Fl to 8 Fl	0.03	15 %	On

### Plate-type instruments

The threshold line should cross only sigmoid curves of signal accumulation of positive samples and should not cross the baseline; otherwise, the threshold level should be raised. Set the threshold at a level where fluorescence curves are linear and do not cross curves of the negative samples.

### Data Analysis:

- The experiment may be considered valid if:
  - the Negative Controls haven't any positive fluorescence signal;
  - the standards have positive signals in all channels (*Fam*, *Joe/HEX/Cy3*)
- The result of the sample is considered:
  - Invalid* in case of absence of any fluorescence signal (positive or internal);
  - Negative* if signal is present only in the *Fam(Green)* channel with the conc. of genomic DNA > 5 x 10<sup>3</sup>;
  - Positive*: if in the *Joe/HEX/Cy3* channel the fluorescence signal is present (Ct ≤ 33) at least in one of the 2 tubes. In these samples fluorescence signal in the channel *Fam* can be absent.

Calculate the concentration of HPV DNA using the following formula:

$$\log\left(\frac{HPV\ DNA\ copies/ reaction}{genomic\ DNA\ copies/ reaction} \times 200000\right) = \log(HPV\ DNA\ in\ 100000\ cells)$$

### RESULTS INTERPRETATION:

Result Ig (HPV in 100.000 cells)	Interpretation
< 3	Clinically insignificant
3-5	Clinically important. Present risk of cervical dysplasia
> 5	Clinically very important. High risk of cervical dysplasia

The results can be calculated automatically using the program in Microsoft® Excel format supplied with the kit.

- Open the program "HPV High Risk Screen Quant 2x" and in the window "Security Warning" click on the button "Enable Macros" (Security level of the Microsoft® Excel must be selected as Medium (Tools→Macro→Security→Medium)).
- Copy with the right button of the mouse the names of the samples from the column "Name" and paste them in the column "Name" of the program "HPV High Risk Screen Quant 2x".
- Copy in the same way the Ct values from the channel FAM (Green) and paste them in the correspond column of the program. Repeat the same procedure for the Joe/HEX/Cy3 channel. Standards must be named as K1, K2, K3 and Negative controls must be marked as "-".
- Select the "Quantitative analysis" and choose "Internal Calibration..."
- At the top right of the window insert in the table "Standards" the concentrations of the Quantitative standards reported in the Quant Data Card.
- Click on the buttons "Sign unnamed" and "Results".
- Save the file with a new name.

## PERFORMANCE CHARACTERISTICS

### Analytical specificity

The analytical specificity of the primers and probes was validated with negative samples. They did not generate any signal with the specific *Human Papillomavirus* primers and probes. The specificity of the kit **HPV High Risk Screen Real-TM Quant 2X** was 100%. The potential cross-reactivity of the kit **HPV High Risk Screen Real-TM Quant 2X** was tested against the group control. It was not observed any cross-reactivity with other pathogens.

### Analytical sensitivity

The kit **HPV High Risk Screen Real-TM Quant 2X** allows to detect *Human Papillomavirus* DNA in 100% of the tests with a sensitivity of not less than 1000 copies/ml. The detection was carried out on the control standard and its dilutions by negative sample.

**Target region:** E2, E1

## KEY TO SYMBOLS USED



List Number



Caution!



Lot Number



Contains sufficient  
for <n> tests



For *in Vitro* Diagnostic  
Use



Version



Store at

**NCA**

Negative Control of  
Amplification



Manufacturer

**C-**

Negative control of  
Extraction



Consult instructions for  
use

**C+**

Positive Control of  
Amplification



Expiration Date

**IC**

Internal Control

- \* SaCycler™ is a registered trademark of Sacace Biotechnologies
- \* CFX™ and iQ5™ are registered trademarks of Bio-Rad Laboratories
- \* Rotor-Gene™ is a registered trademark of Qiagen
- \* MX3005P® is a registered trademark of Agilent Technologies
- \* ABI® is a registered trademark of Applied Biosystems
- \* LineGeneK® is a registered trademark of Bioer
- \* SmartCycler® is a registered trademark of Cepheid