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Ver.3 del 01/07/2020

Biologia Molecolare

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HBV DNA REAL TIME QUANTITATIVO 2.0

REF

AA206/48

REF

GA030/2

48 TEST

CND

W0105020216

CE

0459

IVD

NUCLEAR LASER MEDICINE S.r.l.

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BIBLIOGRAFIA

1. Lanini S, Ustianowski A, Pisapia R, Zumla A, Ippolito G. Viral Hepatitis: Etiology, Epidemiology, Transmission, Diagnostics, Treatment, and Prevention. *Infect Dis Clin North Am.* 2019; 33(4):1045-1062
2. Leoni MC, Ustianowski A, Hamzah Farooq H and Arends JE. HIV, HCV and HBV: a review of parallels and differences. *Infect Dis Ther* 2018; 7:407-419
3. Bitton Alaluf M and Shlomai A. New therapies for chronic hepatitis B. *Liver International* 2016; 36: 775 -782
4. Chen CL, Yang JY, Lin SF, Sun CA, Bai CH, You SL, Chen CJ, Kao JH, Chen PJ. Slow decline of hepatitis B burden in general population: Results from a population-based survey and longitudinal follow-up study in Taiwan. *J Hepatol* 2015; 63: 354-636
5. Coppola N, Alessio L, Gualdieri L, Pisaturo M, Sagnelli C, Caprio N, Maffei R, Strace M, Angelillo IF, Pasquale G, Sagnelli E. Hepatitis B virus, hepatitis C virus and human immunodeficiency virus infection in undocumented migrants and refugees in southern Italy, January 2012 to June 2013. *Euro Surveill* 2015;20:30009
6. Kao J and Chen D. Global control of hepatitis B virus infection. *Lancet Infect Dis* 2002; 2:395-403
7. Nguyen DH, Ludgate L and Hu J. Hepatitis B Virus-Cell Interactions and Pathogenesis. *J Cell Physiol.* 2008; 216(2): 289–294
8. Liaw YF. Clinical utility of HBV surface antigen quantification in HBV e antigen-negative chronic HBV infection. *Nat Rev Gastroenterol Hepatol.* 2019;16(10):631-641
9. Lin CL and Kao JH. Hepatitis B Virus Genotypes and Variants. *Cold Spring Harb Perspect Med* 2015;5
10. Bell TG, Yousif M and Kramvis A. Bioinformatic curation and alignment of genotyped hepatitis B virus (HBV) sequence data from the GenBank public database. *SpringerPlus* 2016; 5:1896
11. Sanbul M. Hepatitis B virus genotypes: Global distribution and clinical importance. *World J Gastroenterol* 2014; 20(18): 5427-5434
12. Lin CL and Kao JH. The clinical implications of hepatitis B virus genotype: Recent advances. *Journal of Gastroenterology and Hepatology* 2011;26 Suppl. 1; 123–130

13. Lin CL and Kao JH. Natural history of acute and chronic hepatitis B: The role of HBV genotypes and mutants. *Best Pract Res Clin Gastroenterol.* 2017;31(3): 249-255
14. Saraswat S, Banerjee K, Chaudhury N, Mahant T, khandekar P, Gupta RK, and Naik S. Post-transfusion hepatitis type B following multiple transfusion of HBsAg-negative blood. *J of Hepatol* 1996; 25:639-643
15. Wang T, Cui D, Chen S, Xu X, Sun C, Dai Y, Cheng J. Analysis of clinical characteristics and S gene sequences in chronic asymptomatic HBV carriers with low-level HBsAg. *Clin Res Hepatol Gastroenterol.* 2019; 43(2): 179-189
16. Hayer J, Jadeau F, Deleage G, Kay A, Zoulim F and Combet C. HBVdb: a knowledge database for Hepatitis B Virus. *Nucleic Acids Research*, 2013, Vol. 41
17. Stramer SL, Wend U, Candotti D, Foster GA, B.A., F B Hollinger, Dodd RY, Allain JP and Gerlich W. Nucleic Acid Testing to Detect HBV Infection in Blood Donors. *M.DN Engl J Med* 2011;364:236-47
18. EASL 2017 Clinical practice guidelines on the management of hepatitis B virus infection. *Journal of Hepatology* 2017 vol. 67: 370–398

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IVD

NUCLEAR LASER MEDICINE S.r.l.

HEADQUARTER: Via Cascina Conighetto – BUSINESS OFFICES: Viale delle Industrie, 3 - 20049 SETTALA (MI), ITALY

Phone: (+39) 02/95.24.51 - Fax (+39) 02/95.24.52.37 - (+39) 02/95.24.52.38 WEB: www.nlm.it - E-MAIL: segreteria@nlm.it

ISO 9001 Quality Management System Organization certified and ISO 13485 Medical Sector Quality Management System certified (IMQ certification body - CSQ certification).

⇒ INTENDED USE

HBV DNA REAL TIME QUANTITATIVE 2.0 provides ready to use master mix for quantitative determination of Hepatitis B Virus (HBV) in human plasma samples by Real Time PCR of the viral HBsAg DNA.

The endogenous Internal Control (Coagulation Factor V) is co-extracted and co-amplified in order to monitor all the procedure.

Four panels of four calibrators are included in the device for quantification (12 test X 4 runs).

If more calibrators are necessary, it is possible to use the accessory device "HBV Standard Curve" (NLM code GA030/2 to be ordered separately), that provides other two set.

The device has to be used in association with "QIASymphony® DSP Virus/Pathogen midi kit" (Qiagen/937055, NLM code AA1440/96) and "MagCore® Viral Nucleic Acid Large Volume Extraction Kit" (RBCBioscience/ MVN1200, NLM code AA1507/96) (sold separately).

Extracted DNA is amplified by Real Time PCR, using CFX (BioRad).

Uracil-DNA Glycosylase (UNG) has been introduced in the ready-to-use mix to prevent possible carry-over contamination.

HBV DNA REAL TIME QUANTITATIVE 2.0 is intended for use in association with clinical presentation and other laboratory markers of HBV infection for the clinical management of patients. The test can be used to assess viral response to antiviral treatment as measured by changes in plasma HBV DNA levels.

The device is not intended for use as a screening test for the presence of HBV DNA in blood products or as a diagnostic test to confirm the presence of HBV infection.

The device is for professional use only.

INTRODUCTION

Hepatitis B virus (HBV) is one of the most common and severe chronic viral infection in the world. Hundred million people worldwide are infected with HBV, that during the chronic phase, can cause a significant hepatic injury due to cirrhosis and hepatocellular carcinoma (HCC) ^{1,2}.

The improvement of socio-economical conditions, the extended vaccination programs and the antiviral treatments used, led to the reduction of disease incidence in areas where the infection was considered as endemic^{3,4}; quite the opposite, the more frequent transfers and the migratory flow caused the disease increment in countries where HBV was less common.⁵

HBV is transmitted by percutaneous or mucosal exposure to infectious body fluids, by sexual contact with an infected person, by blood transfusions, and perinatally from an infected mother to her infant⁶.

HBV is a member of the *Hepadnaviridae* family and it is characterized by circular, partially double-stranded DNA of approximately 3200 nucleotides, that replicates through the reverse transcription of an intermediate RNA. The genome is organized in 4 partially or totally overlapping open reading frames (ORFs, S, C, X and P) which encode for: HBsAg (ORF S), surface proteins involved in the host cell entry; the HBcAg (ORF C), the core protein of the nucleocapsid; the X protein (ORF X) which is essential for viral replication and has also transactivating potential; the polymerase (ORF P) that functions as a reverse transcriptase, DNA polymerase and RNase H activity.⁷

The quantification of blood circulating HBV has an important role in HBV-monitoring and in the therapy response evaluation⁸.

Ten genotypes of HBV have been identified worldwide (A, B, C, D, E, F, H, I and J); they have a characteristic geographical distribution and differ from each other by more than 7,5%^{9,10}.

Moreover some viral recombinant forms between different genotypes have been identified. The genotype and the presence of mutations affect the disease progression and the choice of the anti-viral treatment to be used.^{11,12,13}

The presence of HBV antigens or antibodies in patients infected with HBV has led to the development of specific immunoserological tests that reduced, but not completely eliminated, the incidence of post-transfusion hepatitis B infection.^{14,15}

The detection of viral DNA with PCR allows to monitor the virus presence even in cases that may not be detectable by the immunological tests^{16,17}.

New more sensitive molecular tests may act as “companion diagnostic” not only for the conventional treatment, but also for development and application of new generation therapy.¹⁸

REFERENCES

1. Lanini S, Ustianowski A, Pisapia R, Zumla A, Ippolito G. Viral Hepatitis: Etiology, Epidemiology, Transmission, Diagnostics, Treatment, and Prevention. *Infect Dis Clin North Am.* 2019; 33(4):1045-1062
2. Leoni MC, Ustianowski A, Hamzah Farooq H and Arends JE. HIV, HCV and HBV: a review of parallels and differences. *Infect Dis Ther* 2018; 7:407-419
3. Bitton Alaluf M and Shlomai A. New therapies for chronic hepatitis B. *Liver International* 2016; 36: 775 -782
4. Chen CL, Yang JY, Lin SF, Sun CA, Bai CH, You SL, Chen CJ, Kao JH, Chen PJ. Slow decline of hepatitis B burden in general population: Results from a population-based survey and longitudinal follow-up study in Taiwan. *J Hepatol* 2015; 63: 354-636
5. Coppola N, Alessio L, Gualdieri L, Pisaturo M, Sagnelli C, Caprio N, Maffei R, Strace M, Angelillo IF, Pasquale G, Sagnelli E. Hepatitis B virus, hepatitis C virus and human immunodeficiency virus infection in undocumented migrants and refugees in southern Italy, January 2012 to June 2013. *Euro Surveill* 2015;20:30009
6. Kao J and Chen D. Global control of hepatitis B virus infection. *Lancet Infect Dis* 2002; 2:395-403
7. Nguyen DH, Ludgate L and Hu J. Hepatitis B Virus-Cell Interactions and Pathogenesis. *J Cell Physiol.* 2008; 216(2): 289–294
8. Liaw YF. Clinical utility of HBV surface antigen quantification in HBV e antigen-negative chronic HBV infection. *Nat Rev Gastroenterol Hepatol.* 2019;16(10):631-641
9. Lin CL and Kao JH. Hepatitis B Virus Genotypes and Variants. *Cold Spring Harb Perspect Med* 2015;5
10. Bell TG, Yousif M and Kramvis A. Bioinformatic curation and alignment of genotyped hepatitis B virus (HBV) sequence data from the GenBank public database. *SpringerPlus* 2016; 5:1896
11. Sanbul M. Hepatitis B virus genotypes: Global distribution and clinical importance. *World J Gastroenterol* 2014; 20(18): 5427-5434
12. Lin CL and Kao JH. The clinical implications of hepatitis B virus genotype: Recent advances. *Journal of Gastroenterology and Hepatology* 2011;26 Suppl. 1; 123–130
13. Lin CL and Kao JH. Natural history of acute and chronic hepatitis B: The role of HBV genotypes and mutants. *Best Pract Res Clin Gastroenterol.* 2017;31(3): 249-255

14. Saraswat S, Banerjee K, Chaudhury N, Mahant T, khandekar P, Gupta RK, and Naik S. Post-transfusion hepatitis type B following multiple transfusion of HBsAg-negative blood. J of Hepatol 1996; 25:639-643
15. Wang T, Cui D, Chen S, Xu X, Sun C, Dai Y, Cheng J. Analysis of clinical characteristics and S gene sequences in chronic asymptomatic HBV carriers with low-level HBsAg. Clin Res Hepatol Gastroenterol. 2019; 43(2): 179-189
16. Hayer J, Jadeau F, Deleage G, Kay A, Zoulim F and Combet C. HBVdb: a knowledge database for Hepatitis B Virus. Nucleic Acids Research, 2013, Vol. 41
17. Stramer SL, Wend U, Candotti D, Foster GA, B.A., F B Hollinger, Dodd RY, Allain JP and Gerlich W. Nucleic Acid Testing to Detect HBV Infection in Blood Donors. M.DN Engl J Med 2011;364:236-47
18. EASL 2017 Clinical practice guidelines on the management of hepatitis B virus infection. Journal of Hepatology 2017 vol. 67: 370–398

⇒ PRINCIPLE OF THE PROCEDURE

HBV DNA REAL TIME QUANTITATIVE 2.0 is based on two processes:

1. Viral DNA extraction
2. Amplification and detection of target sequence using Real Time PCR

The Internal Control (Coagulation Factor V) is endogenous and therefore is extracted together with the sample in order to monitor the entire procedure (including the integrity of the plasma sample). Four Panels with four calibrators each are included in order to quantify the amount of HBV DNA in clinical specimens.

1. Viral DNA extraction

HBV and IC DNA extraction from human plasma samples has to be performed using QIASymphony DSP Virus/Pathogen midi kit (NLM code AA1440/96) or MagCore Viral nucleic acid extraction kit Large Volume (NLM code AA1507/96).

2. PCR and Detection

Amplification target

HBsAg region has been chosen as target to detect the HBV presence in human plasma samples because it is highly conserved among HBV genotypes.

The Internal Control has an internal sequence similar for length and base composition to HBV target sequence and specific primers and probe to distinguish between IC and HBV.

Calibrators are obtained by serial dilutions of a synthetic DNA of cloned HBsAg HBV region, therefore they need the same primers and probe as HBV DNA target for Real Time PCR.

Amplification

Following extraction step, DNA (HBV and IC) is amplified by Real Time PCR, using appropriate enzymes and buffer conditions.

Detection

The test is based on the Real-Time PCR, a technique that allows to monitor in real time the amplification of the samples of interest, using probes doubly labeled with a donor fluorophore and a quencher ("dual-labeled probe"). The PCR is carried out in the presence of two specific probes respectively for the internal control and for the HBsAg DNA viral region. The donor fluorophore is different for the two probes: FAM for HBV and JOE for IC. In presence of a light source, when the

probes are intact and donor and quencher are close to each other, the fluorescence emitted by the donor fluorophore is absorbed by the proximal quencher. During the amplification each probe matches to its specific sequence and it is targeted of 5' to 3' nuclease activity of the DNA polymerase. When the donor and quencher fluorophores are separated, the donor fluorescence can be detected at its specific wavelength: in this way, the amplification of the viral DNA and the IC can be monitored as the reaction proceeds.

CFX (BioRad) comprises in a single instrument a thermal cycler for target amplification and a fluorometer for the detection of fluorescence during the cycling process. The computer connected with the Real Time instrument collects fluorescent data that are displayed in a graph through a specific software (*Figure 1*).

Following raw data collection, analysis is carried out. Raw data is normalized to correct the background signal, then a threshold level can be set: this is the level at which fluorescence data are analyzed (*Figure 2*).


Figure 1: raw data is plotted on a graph of fluorescence versus cycle number.


Figure 2: the normalized data is plotted on a log scale graph of fluorescence versus cycle number.

The number of cycles it takes for a sample to reach the threshold level is named Ct-value (threshold cycle) and it is related to the initial amount of target DNA: the higher is the target titer, the earlier the fluorescent signal reaches the threshold level.

If standards with known concentrations are processed, the linear regression analysis produces a standard curve from which the concentration of unknown samples can be determined. The standard curve is automatically generated by plotting the threshold cycle versus standard initial concentration (IU/ml) and calculating the best-fit line; the HBV DNA titre of each sample is determined by locating its Ct on the standard curve (*Figure 3*).


Figure 3: standard curve obtained by linear regression analysis of calibrators Ct plotted versus the respective initial concentration.

PCR product contamination Control

PCR amplification is performed in the presence of deoxyuridine triphosphate (dUTP) to generate amplification products that can be selectively degraded by treatment with uracil-N-glycosylase (UNG) and heat.

UNG cleaves the glycosidic bond between the uracil base and the deoxyribose ring, creating an abasic site that is susceptible to heat. An incubation step at +50°C before the amplification removes uracil bases from any contaminating amplification products present in the reaction mixture. The heat-labile UNG enzyme is completely and irreversibly inactivated when the temperature is increased during the next amplification steps, preventing removal of uracil from the newly synthesized DNA.

➔ REAGENTS COMPOSITION (Store at -25/-15°C)

Components	
HBV MIX mix ready-to-use with primers, probes, amplification enzyme and UNG	2x 1000 µl
No Template Control (NTC) Water to be used as NTC starting from purification step	4 x 1300 µl
RNase, DNase and Protease free water Water for the reconstitution of calibrators	2 x 400 µl
Set of calibrators 1, 2, 3, 4 Freeze-dried HBV DNA template	4 set

The amount of reagents indicated in the table is referred to the standard size of the kit. Reduced packaging is available on request for evaluation and/or demonstration of the product.

⇒ STABILITY AND STORAGE

- All the opened and closed reagents are stable up to the expiry date indicated on the label when stored at -15/-25°C.
- Thaw all the reagents on ice or at +2/+8°C.
- **HBV Mix** is stable up to 5 freeze/thaw cycles. Prepare aliquots if more freeze/thawing cycles are expected.
- **HBV Mix** contains FAM and JOE fluorophores that are photosensitive: avoid prolonged exposure to light.
- **HBV Mix** is ready to use; once dispensed in PCR tubes/plate, avoid prolonged exposure to light.
- **Calibrators** have to be resuspended with DNase, RNase and Protease free water (provided with the device) just before use; it is recommended to dispense them into the PCR tubes/plate after DNA samples to avoid cross contamination.

SPECIMEN COLLECTION AND HANDLING

Viral DNA is isolated and purified from human **plasma** specimens prepared within 6 hours from blood collection.

- Collect blood observing universal precautions for venipuncture.
- Collect blood in sterile tubes with EDTA as anticoagulant and centrifuge at 1000-1500g for 10-15 minutes to separate plasma. Use only EDTA as an anticoagulant, other types of anticoagulants may interfere with the performance of the test.
- ⇒ - Aliquot 1200 µl (NLM code AA1507/96) or 700 µl (NLM code AA1440/96) of plasma using vertical down-flow airbox and operating to avoid DNA degradation.
- Store samples at +2/+8°C until extraction step. If they are not processed immediately, store them at -25/-15°C frozen for up to 72 hours or at ≤ -70°C for longer period
- If samples are freezed, thaw at +2/+8°C and centrifuge at 1000-1500 x g for 10-15 minutes to avoid cryoglobuline interference.
- Avoid repeated freeze-thaw cycles of plasma samples.
- Samples must be handled following good laboratory practices.
- Perform the procedure following good laboratory practices and handle samples as biohazardous and capable of transmitting infection.

⇒ PRECAUTIONS

- Only professional and opportunely trained personnel should use this kit. Handle this product according to established good laboratory practices and universal precautions throughout the assay procedure.
 - The use of the following solutions allows to minimize the risk of cross-contamination:
 - ✓ It is recommended to perform the assay in two different areas:
 - Area 1: pre-PCR (samples handling, extraction and **HBV mix** dispensing. Be careful not to contaminate the mix during its handling.
 - Area 2: post-PCR (Real Time PCR)
- Each area should be provided with dedicated equipment and consumables (lab coats, centrifuge, tubes, pipettes, etc).
- ✓ Workspace should be organized to ensure that the workflow occurs in one direction, from clean areas (pre-PCR) to dirty areas (post-PCR).
 - ✓ At the end of the procedure, it's recommended to clean the working surfaces (bench and vertical downflow airbox) with 5-10% bleach (final concentration of sodium

hypochlorite: 0,5% w/v), freshly prepared. Follow the manufacturer recommendations to clean the instruments.

- ✓ All disposable items (tips and tubes) must be DNase, RNase free. Use aerosol-resistant pipette tips to avoid pipettes contamination. Use a new tip every time a volume is dispensed.
 - ✓ Use vertical downflow airbox with UV lamp
 - ✓ Change gloves frequently
 - ✓ The absence of contamination is guaranteed by the analysis in each session of a NTC, which allows to monitor the whole procedure starting from the extraction step.
- Discard all used material in accordance with the existing local and national regulations in force.
 - Do not eat, drink, smoke or apply cosmetics in areas where reagents or specimens are handled.
 - In case of contact with reagents rinse immediately with water and seek medical advice.
 - Do not use device after its expiration date.
 - Do not mix **HBV Mix** from different lots.
 - Do not use the device if the box is damaged; contact the supplier.
 - It is advisable to have constant and uniform laboratory temperature, avoid to place the instruments near heating/cooling sources that may compromise the right working.

⇒ MATERIAL REQUIRED BUT NOT PROVIDED

AREA 1

Nucleic acid extraction system and related disposable

DNA extraction kit

Vertical downflow airbox with UV lamp

Dedicated adjustable volume pipettes set

Aerosol barrier tips, 2 ml and 1,5ml DNase, RNase free tubes

Vortex mix

Centrifuge

DNase, RNase free, strips or plates for Real TimePCR

AREA 2

CFX (BioRad) and software

ASSAY PROCEDURE

VIRAL DNA EXTRACTION

→ Extraction kit to be used:

- AUTOMATIC EXTRACTION QIASymphony **NLM code AA1440/96**
- SEMI-AUTOMATIC EXTRACTION MagCore **NLM code AA1507/96**

For the preparation, use and disposal of the reagents, refer to the specific instructions for use of the extraction system used.

For the preparation of plasma samples it is necessary to follow the protocol concerning the specific method.

Use the volumes indicated in the table below

Extraction system	Sample volume	Elution Volume
MagCore semi-automatic extraction (NLM code AA1507/96).	1200 µl	60 µl
QIASymphony automatic extraction (NLM code AA1440/96 – protocol Cell free 500)	700 µl	60 µl

Purified DNA can be stored at +2/+8°C if immediately used, otherwise keep it frozen at -25/-15°C or at ≤ -70°C for longer periods. It is recommended to thaw DNA at +2/+8°C.

→ REAL TIME PCR

Set Real Time PCR profile before preparing HBV Mix.

- **Carefully mix the HBV mix before use.**
- Mix gently and dispense in the plate/strip, previously marked, **25 µl** of HBV Mix for each sample + 1 NTC + 4 calibrators.
- Add **25 µl** of purified DNA or NTC to each well/tube and **mix pipetting up and down.**
- Resuspend the calibrators with **30 µl** RNase, DNase and Protease free water, mix properly, then transfer **25 µl** in the specific well/tube.

WARNING: dispense and pipette mix and samples very carefully, avoiding the formation of bubbles. Centrifuge the plate/strip before place it into the instrument.

RUN SETUP

- Turn on the PC and Real Time instrument.
- Open the Bio-Rad CFX Manager Software and click on *Create new run*, select the model of instrument used and click *OK*; the *Run Setup* window will open.
- It is possible to create a new protocol, select/modify an existing one or repeat a run. To repeat a run with the same settings, open the desired run, click on *File - Repeat Run*.

Create a new protocol

- Protocol → create New (Figure 4)
- Protocol Editor – new
- Enter the reaction volume (50 µl) in the Sample Volume box
- Add the steps by clicking on *Insert Step* on the left of the window (it is possible to modify temperatures and times directly on the graph or in the test below with a double click on them):
 1. 37°C for 10 min
 2. 50°C for 15 min
 3. 95° for 20 sec
 4. 95°C for 15 sec
 5. 60°C for 1 min; add “plate read” at this step
 6. Insert GOTO 4, 44 more times
- Click OK. Save the new protocol.


Figure 4

➔ PLATE SETUP

Open *plate* in *Run setup* window. It is possible to create a new plate or select/modify an existing one.

Create a new plate

- Plate → create New”
- In *Plate Editor – new* window (Figure 5) → *Scan Mode* select *All channels*
- Select plate type: *Setting* → *Plate Type* → *BR clear*
- Click *Select Fluorophores* button and check *FAM* and *JOE/HEX*
- Select a well
- Choose *Sample Type* → *Unknown*, *Standard* or *Negative Control* depending on the type of sample
- Type *Target Name* (*HBV* for *FAM* and *IC* for *JOE/HEX*) and Select *Load* check boxes to insert fluorophores (*FAM* and *JOE/HEX* for each well except for the standards, where only *FAM* has to be loaded).
- Type *Sample Name* and click *Load*. Check boxes to insert sample Name

- For the *Standards*, type the respective concentration values as indicated on the label inside the box of the device and click the load button in the concentration window.
- Click *OK*. Save the new plate.


Figure 5

Select/modify existing protocol or plate

- To import an existing protocol or plate click *Select Existing* in protocol or plate window.
- To modify an existing plate click *Select Existing* → *Edit Select*.
- Proceed as explained above.

START RUN

- Prepare the plate/strip with samples and mix.
- In *Start Run* tab click *Open Lid*, place the plate/strip in CFX instrument and click *Close Lid*
- Press *Start Run* to begin.
- It's possible to monitor the Raw data plots during the assay in *Run Details* window in *Real Time status* page.

DATA ANALYSIS

At the end of the run the window *Data Analysis* will open automatically.

To exclude a sample/calibrator from the analysis go to *plate setup-view/edit plate*, select the well of interest and click on *Exclude well in analysis*.

Analyze data for FAM and JOE/HEX separately.

a) **FAM**: deselect JOE/HEX in the amplification window. Click on *Settings* and select (Figure 6):

- *Cq determination mode: single threshold*
- *Analysis mode: Baseline subtracted curve fit*
- *Baseline threshold: Single Threshold-User defined-introduce the value 200*. Click *OK*

In “*Quantification*” for the data analysis: the graph on the left shows the amplification data (linear or log scale), the graph on the right shows the standard curve. The R^2 value should be as close as possible to 1. If one of the four calibrators doesn't match with the linear fitting (i.e. in the presence of a Ct value outside the corresponding acceptability RANGE reported on the label inside the box),

it can be excluded from the analysis, in order to obtain a better fitting of the standard curve and a more precise sample quantification without altering the clinical significance of the results.

If more than one calibrator doesn't match with the linear fitting, Standard Curve cannot be accepted and the run must be repeated.

Check FAM results for valid samples: a value of concentration (IU/ml) is given in the result table.

- If given concentration is lower than Detection Limit, report result as HBV DNA < LOD; in this case positive result reproducibility is not guaranteed.
- If sample has positive signal in the JOE/HEX channel but no signal in the FAM channel report as *HBV DNA Target not detected*.
- Results that are not within the linear range of the assay cannot be considered accurate, therefore they should be carefully interpreted.
- If given concentration is above the upper quantification limit, report result as HBV DNA > 1,00E+08 IU/ml. To get a more accurate result, the test should be repeated using the starting plasma sample diluted with negative plasma. The calculated concentration must be multiplied by the dilution factor to obtain the initial concentration of that sample.


Figure 6

b) **JOE/HEX:** deselect FAM in the amplification window. Click on *Settings* and select

- *Cq determination mode: single threshold*
- *Analysis mode: Baseline subtracted curve fit*
- *Baseline threshold: Single Threshold-User defined-introduce the value 100.* Click **OK**

Check samples JOE/HEX Ct and interpret results as follows:

	JOE/HEX Ct	Interpretation
Sample	Present	Valid
	Absent	Not Valid*

(*) *In absence of IC signal sample result can't be confirmed, therefore sample run should be repeated from the extraction. However, the absence of IC signal in HBV DNA with high-titer (>10⁸ IU/ml) could be explained with the competition between HBV specific target and IC. In this case, the positive result for HBV can be considered valid.*

If NTC is run, interpret results as follows:

	FAM Ct	JOE/HEX Ct	Interpretazione
No Template Control	Absent	Absent	Valid
	Present	Absent	Not Valid*
	Absent	Present	
	Present	Present	

(*) Possible contamination

Create a Report (*Tools* → *Reports*). Select/Deselect useful information.

WARNING: for a correct interpretation of results, always pay attention to the graph of the session.

Example 1 (*figure 7a*): the sample highlighted in red in the graph shows a positive Ct value (35), but the graph itself shows it is an instrumental artifact. The found Ct value should not be taken into account and the sample must be repeated.


Figure 7a

Example 2 (*Figure 7b I and II*): if fluorescence signals at the beginning of the run are inhomogeneous or there are artifacts, it is possible to adjust the signal clicking on *Settings* and selecting:

- Apply *Fluorescence Drift Correction* or, alternatively, *Baseline cycles*; in this case select *User defined* and choose from **10** to **20**.
- After the operation described above, check if the calibrators Ct are in range.


Figure 7b I


Figure 7b II

WARNING

- In case of absence of expected JOE/HEX signal it is advisable to repeat the assay.
- Calibrators haven't Internal Control, therefore JOE/HEX signal is not expected.
- When contamination events occur it's important to find out the source (extraction and/or PCR steps). Provided controls can be very useful to ensure the correct test performance and to find out contamination origin.

➔ PERFORMANCE CHARACTERISTICS

Specificity

The specificity of the HBV DNA REAL TIME QUANTITATIVE 2.0 was determined by analysing 104 HBV negative plasma samples from blood donors. No false positive results were obtained. According to the results the specificity of the assay is 100%.

Analytical sensitivity

Sensitivity was determined by analysing four dilution levels of the 4th WHO International Standard HBV (NIBSC code 10/266). For the dilutions a plasma negative for HBV DNA and antibodies was used. The analysis was performed according to instructions for use, DNA extraction included. As shown in table below, the detection limit of HBV DNA REAL TIME QUANTITATIVE 2.0 is **<8 IU/ml** (Probit analysis at 95%, done with XLSTAT 2014).

QIASymphony cod. NLM AA1440/96:

Dilutions	N tot. runs	N tot. positives	% positives
A (31 IU/ml)	28	28	100%
B (10 IU/ml)	28	28	100%
C (5 IU/ml)	28	23	82,14%
D (2,5 IU/ml)	28	20	71,43%
LOD= 7,04 UI/ml			

*Probit Analysis done with XLSTAT 2014

MagCore cod. NLM AA1507/96:

Dilutions	N tot. runs	N tot. positives	% positives
A (31 IU/ml)	28	28	100%
B (10 IU/ml)	28	28	100%
C (5 IU/ml)	28	22	78,57%
D (2,5 IU/ml)	28	18	64,28%
LOD= 7,84 UI/ml			

*Probit Analysis done with XLSTAT 2014

Linear Range

To evaluate the linear range of the device, 14 dilution levels were prepared using high viral load HBV positive samples diluted in HBV negative human plasma. The test was performed according to the instructions for use, with both the extractions, using two different lots of reagents. Each dilution level has been tested at least in 16 replicates.

Dilutions	Titre (IU/ml)
A	1×10^8
B	$3,3 \times 10^7$
C	1×10^7
D	$3,3 \times 10^6$
E	1×10^6
F	$3,3 \times 10^5$

G	1×10^5
H	$3,3 \times 10^4$
I	1×10^4
L	$3,3 \times 10^3$
M	1×10^3
N	$3,3 \times 10^2$
O	1×10^2
P	$3,3 \times 10^1$

As shown in figures below (Figure 8.1 and Figura 8.2), HBV DNA REAL TIME QUANTITATIVE 2.0 showed a linear response from 3,300E+01 IU/ml to at least 1,00E+08 IU/ml, using the accuracy acceptance criterion of $\pm 0,3 \log 10$.


Figure 8.1: Linear range with QIASymphony NLM code AA1440/96


Figure 8.2: Linear range with MagCore NLM code AA1507/96

Precision

Intra-assay precision

Intra-assay precision was evaluated by using 6 positive samples with different viral load. Two different operators tested three replicates per dilution level using the same lot of reagents. Results of intra-assay tests are showed in the table below.

Viral load (IU/ml)	3,5x10 ⁷	3,2x10 ⁶	2,8x10 ⁵	2,2x10 ⁴	1,8x10 ³	1,03x10 ²
CV%	9,12	4,12	11,31	4,80	7,45	10,66

Inter-assay precision

Inter-assay precision was evaluated by using the same samples used in the intra-assay study. Two different operators tested three replicates per dilution level, in different runs and days, using 3 different lots of reagents.

The test was performed according to the instructions for use, using both the extractions; in this way the precision was evaluated for all the steps of the procedure.

Results of inter-assay tests are showed in the table below.

Viral load (IU/ml)	3,5x10 ⁷	3,2x10 ⁶	2,8x10 ⁵	2,2x10 ⁴	1,8x10 ³	1,03x10 ²
CV%	9,86	6,65	11,30	3,74	10,47	11,24

Genotypes

The performance of HBV DNA REAL TIME QUANTITATIVE 2.0 on different HBV genotypes was determined by analysing dilution levels of different HBV positive specimens for the most relevant genotypes (A-H). Samples from commercial panels were used.

Dilution levels analysed:

- Gen A:** from 1,26 x10⁶ to 1,26x10¹ IU/ml
- Gen B:** from 8,32 x 10⁵ to 8,32 x 10¹ IU/ml
- Gen C:** from 1,23 x 10⁶ to 1,23 x 10¹ IU/ml
- Gen D:** from 1,07 x 10⁶ to 1,07 x 10¹ IU/ml
- Gen E:** from 7,24 x 10⁵ to 7,24 x 10¹ IU/ml
- Gen F:** from 5,75 x 10⁴ to 5,74 x 10¹ IU/ml
- Gen G:** from 6,03 x 10³ to 6,03 x 10¹ IU/ml
- Gen H:** from 2,63 x 10³ to 2,63 x 10¹ IU/ml

As showed in the graphic below (*Figure 9*) quantification is comparable for the different genotypes analysed. HBV DNA REAL TIME QUANTITATIVE 2.0 detection efficiency is therefore independent of the genotype analysed (acceptance criteria: maximum difference between obtained and expected concentration $\pm 0,5$ Log IU/ml).


Figure:9

Potential cross reactive markers

To evaluate the potential cross reactivity of HBV DNA REAL TIME QUANTITATIVE 2.0 with other pathogens the following samples (from infected patients) were analysed:

Other pathogens	Nr of samples tested
HHV6	2
HCMV	2
EBV	2
VZV	2
BKV	2
Chlamydia trachomatis	3
Neisseria gonorrhoeae	3
HPV	3

Each non-HBV sample tested gave a negative result. Therefore HBV DNA REAL TIME QUANTITATIVE 2.0 showed no cross reaction with other pathogens.

Potentially interfering substances

Data not available.

Whole system failure rate leading to false-negative results

To evaluate the whole system failure rate leading to false negative results, 110 HBV positive plasma samples were analyzed (~ 21 IU/ml, 3 times positive cut-off concentration).

The starting plasma was obtained by diluting in human negative plasma a positive HBV sample, calibrated with the 4th WHO International Standard HBV (NIBSC code 10/266).

All HBV positive samples tested gave a positive result, therefore HBV DNA REAL TIME QUANTITATIVE 2.0 reported no false negatives.

Positive results	110
Negative Results for Positive Samples	0

Absence of Cross-contamination

Absence of Cross-contamination for the whole work-flow was demonstrated alternating high positive and negative samples in 5 runs; both the DNA purification methods were used, AA1440/96 and AA1507/96.

All the negative samples tested resulted negative as expected.

Therefore HBV DNA REAL TIME QUANTITATIVE device showed no cross-contamination between positive and negative samples.

Robustness

The robustness of HBV DNA REAL TIME QUANTITATIVE 2.0 kit was evaluated adding different volumes of mix and samples ($\pm 10\%$): 10 low positive HBV plasma samples (21 IU/ml) were tested with +10% of mix and sample and -10% of mix and sample.


All samples tested gave a positive result, therefore HBV DNA REAL TIME QUANTITATIVE 2.0 is reliable during normal usage even when subjected to small volume variations.



Comparison with Roche Cobas TaqMan® HBV test in routine samples

The performance of the HBV DNA REAL TIME QUANTITATIVE 2.0 was compared to Roche Cobas TaqMan® HBV test analysing 103 HBV samples.

The analysis was performed according to the instructions for use, including the extraction step.

Results obtained showed a good correlation between these two assays (*Figure 10*).

					
Lot Number	Catalogue number	Temperature Limitation	Use by	Consult Instructions for use	
					
European Conformity	Manufacturer	<i>in vitro</i> diagnostic service	Sufficient for	No Template Control	Modified text compared with the previous version