BIOLOGIA MOLECOLARE

VER 2 - 22/02/2018

REF AA896/24 **24 TEST**

REF AA896/48 **48 TEST**

CND W0105020307

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HCV RNA REAL TIME QUALITATIVO 2.0

CE
0459

IVD

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6. National Institutes of Health, Consensus Conference Statement. Management of hepatitis C: 2002.
7. P.T. Hrabec et al. Comparative analysis of hepatitis C virus phylogenies from coding and non-coding regions: the 5' untranslated region (UTR) fails to classify subtypes. *Virology Journal* 2006, 3:103
8. M. A. Ansari et al. HCV-Core Region: Its Significance in HCV-Genotyping and Type Dependent Genomic Expression. *Macedonian Journal of Medical Sciences*. 2012 Mar 15; 5(1):30-39.
9. G. Colucci Molecular diagnostic and predictive tests in the evolution of chronic hepatitis C anti-viral therapies *BMC Infect Dis*. 2012; 12(Suppl 2): S8. Published online 2012 November 12

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INTENDED USE

HCV RNA REAL TIME QUALITATIVE 2.0 test provides reagents for the qualitative detection of C-hepatitis virus (HCV) RNA in human plasma samples by Real Time RT-PCR of the 5'-untranslated region (5'-UTR) of the viral RNA; the Core region is co-amplified to make the amplicons compatible with NLM HCV genotyping kit GEN-C 2.0 code AC004/24.

The Internal Control is endogenous (mRNA of the Housekeeping gene GAPDH) and is co-extracted and co-amplified in order to monitor the complete procedure.

This device has to be used in association with viral RNA extraction kit NLM code AA1318 "Viral DNA/RNA Extraction", to be ordered separately.

Viral RNA is retrotranscribed and amplified in one step RT-PCR Real Time, using CFX (BioRad).

⇒ The device can also be used to obtain biotinylated amplification products only for viral genotyping in combination with AC004/24 kit. In this case it is possible to amplify samples using the C1000 thermal cycler (BioRad) or equivalent instrument.

HCV RNA REAL TIME QUALITATIVE 2.0 test is intended for use in association with clinical presentation and other laboratory markers (HCV-Ab, ALT, etc) for the diagnosis of HCV infection in patients.

The device is not intended for use as a screening test for the presence of HCV RNA in blood donors. The device is for professional use only.

INTRODUCTION

Hepatitis C virus is the main etiologic agent of acute viral hepatitis that can lead, in most patients, to liver cirrhosis and then to hepatocellular carcinoma¹. HCV can also cause extrahepatic diseases²⁻⁴.

HCV genome is a linear single strand RNA with one ORF (Open Reading Frame); in 5' position there is an untranslated highly conserved region (5'-UTR), followed by the Core region, they both represent the target for viral molecular biology detection. The 5'-UTR and Core genetic heterogeneity among different HCV strains determines the classification into different genotypes^{5,7,8}. HCV has a parenteral transmission through exposure to infected blood. This exposure can occur in the context of blood transfusion before 1992, injection drug use, solid organ transplantation from infected donors, birth from an infected mother, occupational exposure to infected blood, high-risk sexual practice⁶.

After initial exposure, HCV RNA can be detected in patients blood within one to three weeks; acute infection can be severe but rarely is fulminant⁶.

Although generally asymptomatic, about 85% of the acute infections become chronic: persistence of HCV infection is diagnosed by the detection of HCV RNA in the blood for at least six months⁶.

The recent development of directly acting antiviral (DAA)-based regimens and the approval of the first, boceprevir- and telaprevir-based triple-therapy regimen, have brought new hopes for higher virus eradication rate across different disease settings. This calls for even more accurate and sensitive HCV RNA tests able to monitor faster virus kinetics and promptly detect treatment failures.⁹

With new viral load tests acting as "companion diagnostic", not only the conventional treatment but also the development and applications of the new generation DAA can be fully supported for a more effective management of chronic hepatitis C and its long-term complications⁹.

REFERENCES

1. Sarrazin C. *Diagnosis of hepatitis C: update 2004. Journal of Gastroenterology and Hepatology* (2004) 19, S88-S93
2. Johnson RJ, Gretch DR, Yamabe H et al. *Membranoproliferative glomerulonephritis associated with hepatitis C virus infection. N Engl J med* 1993; 328: 465-470.
3. Agnello V, Chung RT, Kaplan RM. *A role for hepatitis C virus infection in Type II cryoglobulinemia. N Engl J Med* 1992; 327: 1490-1495.
4. Andreone P, Zignego AL, Cursaro C et al. *Prevalence of monoclonal gammopathies with hepatitis C virus infection. Ann Intern Med* 1998; 129: 294-298
5. Zein NN. *Clinical Significance of Hepatitis C Virus Genotypes. Clinical Microbiology Reviews*, Apr. 2000, p. 223-235.
6. *National Institutes of Health, Consensus Conference Statement. Management of hepatitis C: 2002.*
7. P.T. Hrabec et al. *Comparative analysis of hepatitis C virus phylogenies from coding and non-coding regions: the 5' untranslated region (UTR) fails to classify subtypes. Virology Journal* 2006, 3:103

8. M. A. Ansari et al. HCV-Core Region: Its Significance in HCV-Genotyping and Type Dependent Genomic Expression. *Macedonian Journal of Medical Sciences*. 2012 Mar 15; 5(1):30-39.
9. G. Colucci Molecular diagnostic and predictive tests in the evolution of chronic hepatitis C anti-viral therapies *BMC Infect Dis*. 2012; 12(Suppl 2): S8. Published online 2012 November 12

PRINCIPLES OF THE PROCEDURE

HCV RNA REAL TIME QUALITATIVE 2.0 test is based on two processes:

1. Viral RNA extraction
2. Reverse transcription, amplification and detection of target sequence using Real Time RT-PCR.

The Internal Control (GAPDH mRNA) is endogenous and therefore extracted together with the sample in order to monitor all the procedure (including the integrity of the plasma sample).

1. Viral RNA extraction

HCV and IC RNA extraction from human plasma samples has to be performed using “Viral DNA/RNA Extraction” (NLM code AA1318).

2. RT-PCR and Detection

RT-PCR target

The HCV RNA 5'-UTR region was chosen as target to detect the HCV presence in human plasma samples because it is highly conserved among HCV genotypes.

Co-amplification of Core region is included in order to discriminate genotypes 1a, 1b and 6 in case of samples that have to be analyzed with NLM AC004/24 Gen-C 2.0 kit.

The Internal Control is the mRNA of the constitutively expressed GAPDH gene that is co-extracted and co-amplified together with target HCV; the IC amplicon is similar to HCV target sequence in terms of length and base composition.

Reverse transcription and Amplification

Following extraction step, RNA (HCV and IC) is retrotranscribed and amplified in one-step RT-PCR, using appropriate enzymes and buffer conditions.

Detection

Real Time PCR allows target sequences detection during amplification itself, using suitable probes labeled with fluorescent dyes. In the amplification mix there are two specific probes, respectively for the Internal Control and for HCV 5'-UTR region, each labeled with a reporter dye and a quencher dye. The fluorescent reporter dye is different for the two probes: FAM for HCV and JOE for IC. In presence of a light source, when the probes are intact and donor and quencher are close to each other, the fluorescence emitted by the donor fluorophore is absorbed by the proximal quencher. During the amplification each probe matches to its specific sequence and it is target of 5' → 3' nuclease activity of the DNA polymerase. When the donor and acceptor fluorophores are separated, the donor fluorescence can be detected at its specific wavelength: in this way, the amplification of the viral RNA and IC can be monitored as the reaction proceeds.

CFX (BioRad) comprises in a single instrument a thermal cycler for target amplification and a fluorometer for the detection of fluorescence during the cycling process.

The computer connected with the Real Time machine collects fluorescent data that are displayed in a graph through a specific software (*Figure. 1*).

Following raw data collection, analysis is carried out. Raw data is normalized to correct the background signal, then a threshold level can be set: this is the level at which fluorescence data are analyzed (*Figure 2*).


Figure 1: raw data is plotted on a graph of fluorescence versus cycle number.


Figure 2: the normalized data is plotted on a log scale graph of fluorescence versus cycle number

The number of cycles it takes for a sample to reach the threshold level is named Ct-value (threshold cycle) and it is related to the initial amount of target RNA: the higher is target titer, the earlier the fluorescent signal reaches the threshold level.

⇒ **REAGENTS COMPOSITION** (Store at -25°/-15° C)

Components	AA896/24 (24 test)	AA896/48 (48 test)
Mix A Probes Primers Tris-EDTA	1 x 26 µl	1 x 47 µl
Mix B Primers Tris-EDTA	1 x 55 µl	1 x 110 µl
Mix C RT-PCR enzymes MgCl ₂ dNTPs	1 x 170 µl	1 x 335 µl
Negative Control* Negative Control is human plasma non reactive for HIV-Ab and RNA, HCV-Ab and RNA, HBS-Ag and DNA. Negative Control has to be extracted.	2 x 600 µl	4 x 600 µl
Positive Control* Synthetic RNA Positive Control is provided in freeze-dried state. Resuspend it in 14 µl of water (provided) before use in RT-PCR step.	2 tubes	4 tubes
RNase, DNase and Protease free water	1 x 600 µl	1 x 600 µl

⇒ *Together with the samples to be tested, 1 Negative Control and 1 Positive Control shall be analysed in each run.

If the device is used only for viral genotyping in combination with NLM code AC004/24, it is not necessary to include the Positive and Negative controls.

STABILITY AND STORAGE

- All the reagents are stable up to the expiry date indicated on the label when stored at -25/-15°C.
- Thaw all the reagents on ice or at +2/+8°C.
- **Mixes** are stable until expiration date if stored at -25/-15°C. These reagents can be refrozen and thawed max 4 times.
- **Mix A** contains FAM and JOE fluorophores that are photosensitive: avoid prolonged exposure to light.
- **Master Mix** has to be used immediately after preparation; once dispensed in PCR tubes, discard the remaining reagent. Avoid prolonged exposure to light.
- The **Positive Control** is a template RNA, supplied pre-aliquoted and lyophilized: it is recommended to resuspend the reagent with RNase, DNase and Protease free water before use. Then add the master mix into the tube and transfer the entire volume in the amplification strip/plate; do this operation after the addition of the samples into the respective tubes to avoid contamination.

SPECIMEN COLLECTION AND HANDLING

Viral RNA is isolated and purified from human plasma specimens prepared within 6 hours from blood collection.

- Collect blood observing universal precautions for venipuncture
- Collect blood in sterile tubes with EDTA as anticoagulant and centrifuge at 1000-1500g for 10-15 minutes to separate plasma. Use only EDTA as anticoagulant, other types of anticoagulants may interfere with the performance of the test
- Aliquot **450 µl** of plasma for RNA extraction using biohazard vertical down-flow airbox and operate to avoid RNA degradation.
- Store samples at +2/+8°C until extraction step. If they are not processed immediately store them at -25/-15°C for up to 72 hours prior to freezing at ≤ -70°C.
- Avoid repeated freeze-thaw cycles of plasma samples
- Perform the procedure using universal precautions and handle samples as if capable of transmitting infection.

⇒ PRECAUTIONS

- Only professional and opportunely trained personnel should use this kit. Handle this product according to established good laboratory practices and universal precautions throughout the assay procedure.
- The use of the following solutions allows to minimize the risk of cross-contamination:
 - ✓ It is recommended to perform the assay in three different areas:
 - Area 1: pre-PCR (samples handling and extraction)
 - Area 2: Master Mix preparation and sample addition to the RT-PCR mix.
 - Area 3: post-PCR (Real Time PCR)

Each area should be provided with dedicated equipment and consumables (lab coats, centrifuge, tubes, pipettes, etc).

Workflow must proceed in a unidirectional manner, beginning in the RNA extraction area (Area 1), moving to the mastermix preparation and RNA addition area (area 2) and finally to amplification and detection area (Area 3). This means that consumables and PPE (lab coats, gloves, etc.) that have been introduced into the post-PCR room should never be placed back to the pre-PCR room without thorough decontamination.

- ✓ It's recommended to clean the workstation with 5-10% bleach (final concentration of sodium hypochlorite: 0,5% w/v) at the end of the procedure. Prepare bleach solution daily.
- ✓ All disposable items (tips and tubes) must be DNase, RNase free. Use aerosol-resistant pipette tips to avoid pipettes contamination. Use a new tip every time a volume is dispensed.
- ✓ Work using vertical downflow airbox with UV lamp
- ✓ Change gloves frequently
- ✓ The absence of contamination is guaranteed by the analysis in each session of a negative control, which allows to monitor the whole procedure starting from the extraction step
- Discard all used material in accordance with the existing local and national regulations in force

- The negative control is a human plasma sample non reactive for HIV-Ab and RNA, HCV-Ab and RNA, HBV HbsAg and DNA. However it should be handled as potentially infectious and should be treated with the necessary safety precautions.
- Do not eat, drink, smoke or apply cosmetics in areas where reagents or specimens are handled.
- In case of contact with reagents rinse immediately with water and seek medical advice.
- Do not use device after its expiration date.
- Do not mix reagents from different lots.
- Do not use the device if the box is damaged; contact the supplier
- It is advisable to have constant and uniform laboratory temperature, avoid to place the instruments near heating/cooling sources that may compromise the right working.

MATERIAL REQUIRED BUT NOT PROVIDED

⇒ **AREA 1**

Vertical downflow airbox with UV lamp

Dedicated adjustable volume pipettes set and aerosol barrier tips 1,5 and 2 ml DNase, RNase free tubes

Magnetic separator

Vortex mixer and Heat block

Centrifuge

⇒ **AREA 2**

Vertical downflow airbox with UV lamp

Dedicated adjustable volume pipettes set and aerosol barrier tips DNase, RNase free
0,2 ml flat cap Real Time PCR, strips or plates DNase, RNase free

⇒ **AREA 3**

CFX (Biorad) and software

C1000 (BioRad) or equivalent thermal cycler: when the device is used only for viral genotyping in combination with AC004/24 device and not for the qualitative determination of the presence of HCV

ASSAY PROCEDURE

VIRAL RNA EXTRACTION

(Extraction method "Viral DNA/RNA Extraction", NLM code AA1318).

Refer to the instructions for use for preparation, handling and disposal of reagents, taking into account **450 µl** starting sample quantity.

This protocol is for manual use and involves the use of a 2 ml reaction tubes and a suitable magnetic separator. For the automation of the procedure with the Janus workstation (Perkin Elmer) contact Nuclear Laser Medicine.

Prepare the necessary amount of 2 ml centrifuge tubes (samples and negative control).

Before starting with the procedure set heat block to +55°C. Ensure that Proteinase K and Poly(A) RNA have been prepared as described.

1. Dispense in the following order **10 µl Proteinase K** and **450 µl** of **sample** into 2,0 ml centrifuge tubes. Add **450 µl Lysis Buffer Viral RNA** and then **6 µl Poly(A) RNA**; proceed pipetting up and down for at least **8 times** and incubate at **55°C** for **10 minutes**, vortexing ones or twice during incubation.

It is also possible to prepare a premix of Lysis Buffer Viral RNA and Poly(A) RNA, calculating the necessary amount for the number of samples + 1 to process. Use this mix within 15 minutes.

2. Following incubation, briefly centrifuge the samples to avoid contamination when opening the tubes. Add **35 µl** of resuspended **Magnetic Beads Viral DNA/RNA (the Magnetic Bead suspension should be mixed vigorously before dispensing)** and **900 µl** of **Binding Buffer Viral DNA/RNA** to each tube, mix by pipetting up and down for **at least 6 times** and incubate for **5 minutes** at RT.

It is possible to prepare a premix of beads and Binding Buffer Viral DNA/RNA.

3. Place the tubes in a magnetic separator for **at least 2 minutes** to let the beads be attracted by the magnets. Remove and discard supernatant by pipetting, **being careful not to take any beads**.
4. Remove the tubes from the magnetic separator. Add **500 µl Wash Buffer A Viral DNA/RNA** to the tube. Resuspend the beads by pipetting until the beads are completely resuspended and incubate for **1 minute** at RT.
5. Place the tube in the magnetic separator to draw the beads to the side of the tube (about 1 minute). Pipette off the supernatant and then remove the tube from the magnet.
6. Repeat the washing procedure (steps 4 and 5) with **500 µl Wash Buffer B Viral DNA/RNA**. **Discard Wash Buffer B Viral DNA/RNA repeating the aspiration with a fine-tipped pipette; be careful to remove all the residual fluid without disturbing the pellet.**
7. Transfer the tubes to the thermo block and let them air-dry at **55°C** for **20 minutes** to remove any traces of ethanol. The beads-pellet should become light brown when dried, otherwise extend the incubation for further 5 minutes or until completely dry

It is important that the beads are completely dried and that there aren't ethanol drops on the tube walls before continuing with the elution step, as residual ethanol can prevent amplification.

8. Add **70 µl** of **Elution Buffer Viral DNA/RNA** to the tubes and briefly resuspend the beads by pipetting. Close the tubes, gently vortex the tubes and incubate the suspension for **5 minutes** at **55°C**.
9. Separate the magnetic beads by placing the tubes in the magnetic separator for at least. **1 minute**. Transfer the supernatant containing the purified viral RNA to the desired storage tubes.

Purified RNA can be stored at +2/+8°C if immediately used, otherwise keep it frozen at -25/-15°C or at ≤ 70°C for longer periods. It is recommended to thaw RNA at +2/+8°C.

REAL TIME RT-PCR

- Set Real Time PCR profile before preparing Master Mix.

Reagents	Volume per sample
Mix A	0,6 µl
Mix B	0,4 µl
Mix C	5 µl

- **Carefully resuspend the mix before use.**

- ⇒
- Prepare Master Mix for the number of purified samples (clinical specimens and negative control) + positive control + 1 volume (+ 2 volumes for number of samples > 15, to have enough Master Mix for all the samples).
 - Mix gently and dispense **6 µl** of Master Mix in the plate/strip previously marked.
 - Add **14 µl** of purified RNA to each well/tube and mix pipetting up and down.
 - At last (to avoid cross contaminations), resuspend the positive control with **14 µl** RNase, DNase and Protease free water (provided), then add the mix and transfer all the volume in the specific well/tube .

WARNING: dispense mix and samples very carefully, in order to avoid creating air bubbles. If possible, centrifuge the plate/strip before placing it into the Real Time PCR instrument.

RUN SETUP

Turn on the PC and CFX.

Open the “Bio-Rad CFX Manager Software” and click on “*Create new run*”, select the model of instrument used and click “OK”, the “*Run Setup*” window will open.

It is possible to create a new protocol, select/modify an existing one or repeat a run. To repeat a run with the same settings, open the desired run and click on “*File - Repeat Run*”.

Create a new protocol

- “*Protocol → create New*” (figure 3)
- “*Protocol Editor – new*”
- Enter the reaction volume (**20 µl**) in the “*Sample Volume*” box
- Add the steps by clicking on “*Insert Step*” on the left of the window (it is possible to modify temperatures and times directly on the graph or in the test below with a double click):
 1. **50°C for 15 min**
 2. **95° for 20 sec**
 3. **95°C for 15 sec**
 4. **60°C for 1 min**; add “*plate read*” in this step
 5. **Insert GOTO 3, 44 more times**
 6. **Optional: for downstream application with GEN-C 2.0 (NLM code AC004/24), add 10°C forever (type 00:00 and press Enter)**
- Click OK. Save the new protocol.

- ⇒ **NOTE:** The thermal profiles for the RT_PCR step were optimized on the CFX and C1000 (BioRad) thermal cycler. It may be necessary to modify the profile according to the manufacturer’s instruction if a different thermal cycler is used.

PLATE SETUP

In “Run setup” window edit “plate” editor window. It is possible to create a new plate or select/modify an existing one.

Create a new plate

- “Plate → create New”
- In “Plate Editor – new” window (figure 4) in “Scan Mode” select “All channels”
- ⇒ - Select plate type: Setting → Plate Type → BR clear
- Click “Select Fluorophores” button and check FAM and JOE
- Select a well and on the right side of the window Choose “Sample Type → Unknown, Standard, Positive or Negative Control” depending on the type of sample
- In Target Name type HCV for FAM and Internal Control for JOE) and Select Load check boxes to insert fluorophores (FAM and JOE for each well except for the positive control where only FAM has to be loaded).
- Insert “Sample Name” and click “Load”. Check boxes to insert sample Name
- Click OK. Save the new plate.

Select/modify existing protocol or plate

- To import an existing protocol or plate click “Select Existing” in protocol or plate window.
- To modify an existing plate click “Select Existing” → “Edit Select”.
- Proceed as explained above.


Figure 3


Figure 4:

START RUN

- Prepare the plate/strip with samples and mix.
- In “Start Run” tab click –“Open Lid”, place the plate/strip in CFX instrument and click “Close Lid”
- Press “Start Run” to begin.
- It’s possible to monitor the Raw data plots during the assay in “Run Details” window, in “Real Time status” page.

DATA ANALYSIS

At the end of the assay click “Stop” to exit from 10°C forever mode. The window “Data Analysis” will open automatically.

To exclude a sample from the analysis go to “plate setup-view/edit plate” (in the upper right), click on the well to omit and then on “Exclude well in analysis” (in the lower right).

Click on “Settings” and select:

- Cq determination mode: single threshold
- Baseline Setting: Baseline subtracted curve fit
- Apply Fluorescence Drift Correction
- Cycles to analyze: from 1 to 45

Analyze data for FAM and JOE separately.

a) **FAM** (figure 5): deselect Joe in the amplification window.

- Baseline threshold: Single Threshold-User defined-introduce the value **25**. Click OK.

Check samples FAM Ct and interpret results as follows:

	FAM Ct	Interpretation
Sample	Present	positive HCV RNA
	Absent	HCV RNA not detected


Figure 5

b) **JOE**: deselect FAM in the amplification window. Click on “Settings” and select:

- Baseline threshold: Single Threshold-User defined-introduce the value **100**. Click OK

Check samples JOE Ct and interpret results as follows:

	JOE Ct	Interpretation
Sample	Present	Valid
	Absent	Not Valid*

(*) In absence of IC signal sample result can't be confirmed, therefore sample run should be repeated.

If negative and positive controls have been tested, interpret results as follows:


	FAM Ct	JOE Ct	Interpretation
Negative Control	Absent	Present	Valid
	Absent	Absent	Not Valid*
	Present	Present	Not Valid**
Positive Control	Present	Absent	Valid
	Absent	Absent	Not Valid*

(*) Hypothesis of RT-PCR inhibition or RNA degradation in one of the procedure steps

(**) Possible contamination

WARNING: for a correct interpretation of results, always pay attention to the graph of the session.

Example 1 (figure 6a): the sample in red in the graph shows a positive Ct value (35), but the graph itself shows it is an instrumental artifact. The Ct value should not be taken into account and the sample has to be repeated.

↓

Well	Fluor	Target	Content	Sample	Cq	SQ
C03	FAM		Unkn		35.03	5.971E+01
C06	FAM		Unkn		36.00	3.097E+01
C07	FAM		Unkn		37.25	1.332E+01
C08	FAM		Std-1		36.10	8.000E+01
C09	FAM		Std-1		33.59	8.000E+01
C10	FAM		Std		30.51	8.000E+02
E03	FAM		Unkn		N/A	N/A
E04	FAM		Std		20.78	1.000E+06

Figure 6a

Example 2 (figure 6b I e II): If fluorescence signals at the beginning of the run are not homogeneous, it is possible to exclude the first cycles; select from "Settings" - Cycles to analyze: from 10 to 45


Figure 6b I


Figure 6b II

Create a Report ("Tools → Reports"). Select/Deselect useful information.

WARNING

- In case of absence of expected JOE signal it is advisable to repeat the assay.
- CP has not an Internal Control, therefore JOE signal is not expected.
- When contamination events occur it's important to find out the source (extraction and/or RT-PCR steps). Controls provided can be very useful to ensure the correct test performance and to find out contamination origin.

PERFORMANCE CHARACTERISTICS

Specificity

The specificity of HCV RNA REAL TIME QUALITATIVE 2.0 test was determined by analysing 501 HCV negative plasma samples from blood donors. No false positive results were obtained. According to the results the specificity of the assay is 100%.

Analytical sensitivity

Sensitivity was determined by analysing dilution levels of an HCV positive clinical specimen calibrated with the HCV WHO 3rd International Standard (NIBSC code 06/100). For the dilutions a plasma negative for HCV RNA and antibodies was used. The analysis was performed according to instructions for use, RNA extraction included. As shown in table below, the detection limit of HCV RNA REAL TIME QUALITATIVE 2.0 test is **18 IU/ml** (Probit analysis at 95%).

Dilutions	N tot. runs	N tot. positives	% positives
A (100 IU/ml)	24	24	100%
B (50 IU/ml)	24	24	100%
C (25 IU/ml)	24	23	96%
D (12,5 IU/ml)	24	22	92%

*Probit Analysis done with *Statplus 2009* software.

Diagnostic sensitivity

Diagnostic sensitivity was evaluated by analysing 10 HCV seroconversion panels, tested with Roche Cobas TaqMan HCV Test as reference system. All the samples gave the expected result.

HCV Genotypes

The performance of HCV RNA REAL TIME QUANTITATIVE 2.0 test on different HCV genotypes was determined by analysing four dilution levels of different HCV positive clinical specimens genotypes (as listed below). The viral load of the samples used was determined by comparing the signal of the original samples with the signal of the HCV WHO 3rd International Standard (NIBSC code 06/100).

Genotypes tested: 1a, 1b, 2a, 2b, 2a/2c, 3a, 4c/4d, 5 and 6 (determined with NLM GEN C code AC004 and sequencing).

Dilutions analyzed:

- Gen 1a:** from $1,5 \times 10^2$ to $2,8 \times 10^6$ IU/ml
- Gen 1b:** from $1,1 \times 10^2$ to $4,3 \times 10^6$ IU/ml
- Gen 2:** from $9,8 \times 10^1$ to $2,7 \times 10^6$ IU/ml
- Gen 3:** from $1,5 \times 10^2$ to $2,8 \times 10^6$ IU/ml
- Gen 4:** from $4,5 \times 10^2$ to $1,8 \times 10^6$ IU/ml
- Gen 5:** 3×10^2
- Gen 6:** 7×10^2

As showed in the graphic below (*figure 7*) detection of HCV RNA is comparable for the different genotypes analysed. Every tested samples gave a positive result, as expected, therefore detection efficiency of HCV RNA REAL TIME QUALITATIVE 2.0 test is independent from the genotype.


Figure 7 Samples genotype was determined using NLM GEN-C cod. AC004

Potential cross reactive markers

To evaluate the potential cross reactivity of HCV RNA REAL TIME QUALITATIVE 2.0 test with other pathogens the following samples (from infected patients) were analysed:

Other pathogens	Nr of samples tested
HGV	2
HCMV	2
HIV	2
HBV	2
HPV	2
Herpes simplex Virus	2
Chlamydia trachomatis	2

Each non-HCV sample tested gave a negative result. Therefore HCV RNA REAL TIME QUALITATIVE 2.0 test showed no cross reaction with other pathogens.

Potentially interfering substances

Data not available.

Whole system failure rate leading to false-negative results

To evaluate the whole system failure rate leading to false negative results 100 HCV positive plasma samples were analysed (3 x 95% positive cut-off concentration ~ 50 IU/ml); in each run one negative control was added. The starting plasma was obtained diluting in human negative plasma a positive clinical specimen calibrated with the HCV WHO 3rd International Standard (NIBSC code 06/100).

All the HCV positive samples tested gave a positive result, therefore HCV RNA REAL TIME QUALITATIVE 2.0 test doesn't lead to false-negative results.

⇒ Absence of Cross-contamination

Absence of Cross-contamination for the whole work-flow was demonstrated alternating high positive and negative samples in at least 5 runs.

All the negative samples tested resulted negative as expected.

Therefore HCV RNA REAL TIME QUALITATIVE 2.0 test showed no cross-contamination between positive and negative samples.

Comparison with ABBOTT REAL TIME HCV (Abbott)

The performance of the REAL TIME HCV QUALITATIVE test was compared to ABBOTT REAL TIME HCV (Abbott) by analysing 102 HCV positive and negative clinical samples.

All the positive samples gave a positive result and the four negative were found negative. Results obtained show a good correlation between these two assays (figure 8).



Figure 8

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LITERATUR

10. Sarrazin C. *Diagnosis of hepatitis C: update 2004. Journal of Gastroenterology and Hepatology* (2004) 19, S88-S93

11. Johnson RJ, Gretch DR, Yamabe H et al. Membranoproliferative glomerulonephritis associated with hepatitis C virus infection. *N Engl J med* 1993; 328: 465-470.
12. Agnello V, Chung RT, Kaplan RM. A role for hepatitis C virus infection in Type II cryoglobulinemia. *N Engl J Med* 1992; 327: 1490-1495.
13. Andreone P, Zignego AL, Cursaro C et al. Prevalence of monoclonal gammopathies with hepatitis C virus infection. *Ann Intern Med* 1998; 129: 294-298
14. Zein NN. Clinical Significance of Hepatitis C Virus Genotypes. *Clinical Microbiology Reviews*, Apr. 2000, p. 223-235.
15. National Institutes of Health, Consensus Conference Statement. Management of hepatitis C: 2002.
16. P.T. Hraber et al. Comparative analysis of hepatitis C virus phylogenies from coding and non-coding regions: the 5' untranslated region (UTR) fails to classify subtypes. *Virology Journal* 2006, 3:103
17. M. A. Ansari et al. HCV-Core Region: Its Significance in HCV-Genotyping and Type Dependent Genomic Expression. *Macedonian Journal of Medical Sciences*. 2012 Mar 15; 5(1):30-39.
18. G. Colucci Molecular diagnostic and predictive tests in the evolution of chronic hepatitis C anti-viral therapies *BMC Infect Dis*. 2012; 12(Suppl 2): S8. Published online 2012 November 12

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