BIOLOGIA MOLECOLARE

VER 4 – 29/11/2017

REF AC004/24

24 TEST

CND W0105020309

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
24 TEST

GMDN 59866

AMPLIFICATION: AA910/48, AA896/24 or AA896/48
EXTRACTION: AA1318 NOT INCLUDED

GEN-C 2.0

Reverse hybridization strip assay

CE

0459

IVD

NUCLEAR LASER MEDICINE S.r.l.
HEAD OFFICE: Viale delle Industrie, 3 – 20090 SETTALA MI (Italy)
Phone (+39) 02/95. 24. 51 - Fax (+39) 02/95. 24. 52. 37
WEB: www.nlm.it – E-MAIL: segreteria@nlm.it

INTENDED USE

Gen-C 2.0 is an *in vitro* line probe assay for the genotyping of hepatitis C virus. The assay discriminates between HCV genotypes on the basis of variations in the 5'-UTR and Core regions. The assay works with amplicons obtained with NLM code AA910/48, AA896/24 or AA896/48. For optimal results it is strictly recommended to use samples with viral load between 2×10^3 and 1×10^7 (if necessary dilute high viral load plasma before the RNA extraction).

⇒ On the basis of the *in silico* sequence analysis, the device is able to distinguish genotypes 1, 2, 3, 4, 5, 6 and 7 and subtypes 1a, 1b, 2a/2c, 2b, 3a, 3b, 3c, 3k, 4a, 4b, 4c/4d, 4e, 4f, 4h, 5a, 6a/6b, 6g, 6m, 6t and 7a.

The result of this assay should be intended as an aid in selecting type and length of antiviral therapy for patients with chronic HCV infection. Therefore the device is not intended for use as a screening test for the presence of HCV RNA or as a diagnostic test to confirm the presence of HCV infection.

⇒ INTRODUCTION

A very high number of HCV-infected patients develop chronic hepatitis which often results in liver cirrhosis and hepatocellular carcinoma (1). HCV genome is a linear single strand RNA consisting of about 9400 nucleotides. The Core, Envelope and non-structural domains form the coding region of the genome, flanked in 5' and 3' by highly conserved untranslated regions (UTR).

Current international nomenclature proposed by Simmonds in 2005 (2), allows the classification in: genotypes, which have nucleotides difference of about 31-33%, and subtypes, that exhibit nucleotide differences of about 20-25%. Nowadays 7 genotypes and 67 subtypes have been identified: 1a to 1l, 2a to 2r, 3a to 3k (previously defined as 10a), 4a to 4w, 5a, 6a-v and 7a (3). The 5'UTR region is characterized by variable regions among the different HCV genotypes that are useful for distinguishing types 2-5, 6a-b and 7a, but is not able to accurately distinguish between genotype 1 and 6 and to properly subtype 1a and 1b. Analysis of the Core region allows to better identify these types and to improve genotype 1 subtyping accuracy (4-5).

IFN- α treatment

Results from many clinical studies regarding treatment of chronic HCV indicated that different genotypes respond in a dissimilar way if treated with IFN- α . In particular, genotype 1b HCV infection could not be efficiently treated with IFN- α (less than 10% of treated patients with a long-term response), while genotypes 1a, 2a, 2b, and 3a responded favorably (between 50 and 80% long-term response). It has been demonstrated that genotype 1b infections proceed much faster to severe forms of chronic hepatitis, cirrhosis, and hepatocellular carcinoma (6). The use of the new drugs (DAA), in association with IFN- α , increased the therapeutic and pharmacological options available for each genotype (7-8). Discrimination of HCV subtypes is of great diagnostic importance since it helps to find more correct, more effective and personalized therapeutic protocols.

⇒ REFERENCES

1. Sarrazin C. *Diagnosis of hepatitis C: update 2004. Journal of Gastroenterology and Hepatology (2004) 19, S88-S93*
2. Simmonds et al. *Consensus proposals for a unified system of nomenclature of hepatitis C virus genotypes. Hepatology, vol. 42, n° 4, 2005.*
3. Smith et al. *Expanded Classification of Hepatitis C Virus Into 7 Genotypes and 67 Subtypes: Updated Criteria and Genotype Assignment Web Resource HEPATOLOGY, Vol. 59, No. 1, 2014*
4. P.T. Hrabec et al. *Comparative analysis of hepatitis C virus phylogenies from coding and non-coding regions: the 5' untranslated region (UTR) fails to classify subtypes. Virology Journal 2006, 3:103*
5. M. A. Ansari et al. *HCV-Core Region: Its Significance in HCV-Genotyping and Type Dependent Genomic Expression. Macedonian Journal of Medical Sciences. 2012 Mar 15; 5(1):30-39.*
6. Nizar N. Zein. *Clinical Significance of Hepatitis C Virus Genotypes. Clinical Microbiology Reviews, Apr. 2000, p.223-235.*
7. *EASL Clinical Practice Guidelines: Recommendation on treatments of hepatitis C. Journal of Hepatology 2016.*
8. Falade-Nwulia et al. *Oral Direct-Acting Agent Therapy for Hepatitis C Virus Infection A Systematic Review Annals of Internal Medicine Vol. 166 No. 8 18 April 2017*

PRINCIPLE OF THE PROCEDURE

The test is based on reverse-hybridization principle. Biotinylated amplicons, generated by RT-PCR of the 5'-UTR and Core regions of HCV RNA, are hybridized to specific probes that are bound to nitrocellulose strip by a poly-T tail.

Biotinylated hybrids are then detected using streptavidin bound to alkaline phosphatase; amplicons that are not complementary are washed out. Then the substrate reacts with the streptavidin-alkaline phosphatase complex forming a purple precipitate and colouring banding pattern on the strip.

PRODUCT COMPOSITION

Reagents (+2/+8° C)		Code	Quantity	N° Vial	
Strips Nitrocellulose membrane coated with oligonucleotides		-	KC029	24	1
DNAT Denaturation solution containing NaOH	 DANGER	H314, P280, P301+P330+P331, P303+P361+P353, P305+P351+P338	KA107	1 ml	1
Hybridization (*) Saline Solution containing preservatives	 WARNING	H319, EUH208, P280, P305+P351+P338, P337+P313	KA102	60ml	1
Stringent Wash Solution Saline solution containing detergent and preservatives	-	-	KA928	175 ml	1
Conjugate Solution containing streptavidin labeled with alkaline phosphatase, stabilizers and preservatives	-	-	KA574	60 ml	1
Wash Solution B Saline Solution containing preservatives	-	EUH208	KA105-1	175 ml	1
Color developer Solution containing 5-Bromo-4-chloro-3-indolyl-phosphate and 4-Nitroblue tetrazolium	-	-	AA564	60 ml	1
Plastic Trays	-	-	DC000	3	
Transparent interpretation card	-	-	-	1	
Collector sheet For strip storage, not suitable for results interpretation and/or as medical report	-	-	-	1	
Interpretation chart Sheet that provides banding patterns for genotypes identification	-	-	-	1	

(*) The exemption from the requirements for labeling of packages is determined by quantity and severity of hazardous reagent. The criteria are: **amounts <125 ml and not serious hazard categories (ref. Regulation (EC) N° 1272/2008 - list 1.5.2.1.1).**

STORAGE AND STABILITY

- All opened or unopened reagents are stable up to the expiry date indicated on the label when stored at the correct temperature.
- At the end of each assay store reagents at the correct temperature.
- The kit should be kept isolated from any source of contaminating DNA especially amplified products.
- Developed dry strips should be stored in the dark at room temperature (+15/30°C).
- The vial containing the DNAT Solution should be closed immediately after use.

PRECAUTIONS

- Only professional and opportunely trained personnel should use this kit. Handle this product according to established good laboratory practices and universal precautions throughout the assay procedure.
- The use of the following solutions allows to minimize the risk of cross-contamination:

- ✓ It is recommended to perform the assay in three different areas:
 - Area 1: pre-PCR (samples handling and extraction)
 - Area 2: Master Mix preparation and sample addition to the RT-PCR mix.
 - Area 3: post-PCR (Real Time PCR)

Each area should be provided with dedicated equipment and consumables (lab coats, centrifuge, tubes, pipettes, etc).

Workflow must proceed in a unidirectional manner, beginning in the RNA extraction area (Area 1), moving to the mastermix preparation and RNA addition area (area 2) and finally to amplification and detection area (Area 3). This means that consumables and PPE (lab coats, gloves, etc.) that have been introduced into the post-PCR room should never be placed back to the pre-PCR room without thorough decontamination.

- ✓ It's recommended to clean the workstation with 5-10% bleach (final concentration of sodium hypochlorite: 0,5% w/v) at the end of the procedure. Prepare bleach solution daily.
- ✓ All disposable items (tips and tubes) must be DNase, RNase free. Use aerosol-resistant pipette tips to avoid pipettes contamination. Use a new tip every time a volume is dispensed.
- ✓ Change gloves frequently
- Discard all used material in accordance with the existing local and national regulations in force
- Do not eat, drink, smoke or apply cosmetics in areas where reagents or specimens are handled.
- In case of contact with reagents rinse immediately with water and seek medical advice.
- Do not use device after its expiration date.
- Do not mix reagents from different lots.
- Do not use the device if the box is damaged; contact the supplier
- Plastic trays can be reused maximum twice, if washed with bleach and then with water at the end of detection.
- Test strips stored at +2/+8°C in the dark are stable until the expiry date.
- It is advisable to have constant and uniform laboratory temperature, avoid to place the instruments near heating/cooling sources that may compromise the right working.
- **Use the transparent interpretation card of the specific lot.**

- **Hazard statements**

- **H314:** Causes severe skin burns and eye damage.
- **H319:** Causes serious eye irritation.

- **Precautionary statements**

- **P280:** Wear protective gloves/protective clothing/eye protection/face protection.
- **P301+P330+P331:** IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
- **P303+P361+P353:** IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
- **P305+P351+P338:** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
- **P337+P313:** If eye irritation persists: Get medical advice/attention.

- **Supplemental hazard informations**

- **EUH208:** It Contains Kathon. May cause allergic reaction (skin).

AMPLICONS STORAGE

Store amplicons obtained with AA910/48, AA896/24 or AA896/48 at +2/+8° C or at -25/-15° C for later use.

MATERIALS REQUIRED BUT NOT PROVIDED

- Dedicated adjustable volume micropipettes set and tips
- Clean tweezers to handle strips

Manual strip processing:

- Water bath with inclined lid, shaking platform (50-100 rpm) and adjustable temperature (46°C ± 0,5° C)
- Thermometer
- Deionised or distilled water
- Vacuum aspiration apparatus
- Timer

Automatic strip processing:

- Tecan ProfiBlot 48 or equivalent instrument. For specific information on the procedure to be applied contact Nuclear Laser Medicine.

PROCEDURE

MANUAL DETECTION

- Adjust the water level of the shaking bath to approx. 2/3 of the height of the Plastic Tray.
 - Heat the water bath to **46° C** and check the temperature ($46 \pm 0,5^\circ \text{C}$) using a calibrated thermometer.
 - Pre-warm the Hybridization and Stringent Wash Solution to **46° C**; all crystals should be completely dissolved.
 - Allow DNAT, Conjugate, Wash Solution B and Color developer to reach room temperature.
 - Take one strip for each sample using clean tweezers (touch strips with gloves only). Label strips with a pencil (do not use ballpoint pen, markers, etc).
- Warning:** do not allow the strips to dry during the whole procedure.

Hybridization (46° C)

- Pipette **10 µl** of **DNAT** into each lane of the Plastic Tray
Warning: do not use DNAT if is not blue.
- For each sample add **10 µl** of **amplification** product to the corresponding drop of DNAT. Mix thoroughly with a pipette (solution will remain blue).
- Incubate **5 min** at room temperature.
- Add **1 ml** of **Hybridization** (pre-warmed to 46° C) to each sample and mix gently: blue color disappears.
- Place strips into the respective lanes with marker lines facing up. Submerge completely.
- Incubate for **30 min** at **46° C** in the shaking water bath (about 50 rpm). Keep the water bath closed to avoid temperature variations.

Stringent Wash solution (46° C)

- After hybridization step remove the tray from the water bath; angle the tray upwards and aspirate the solution from each lane using a pipette or a vacuum aspiration apparatus, without touching the strip surface to avoid its damage.
- Add **1 ml** of **Stringent Wash solution** (pre-warmed at 46° C) to each lane and rinse briefly.
- Remove the solution and add **1 ml** of **Stringent Wash solution** (pre-warmed at 46°C).
- Incubate **15 min** at **46° C** in the shaking water bath (about 50 rpm).
Warning: put Stringent Wash solution back in the water bath during incubation
- Remove the solution and add **1 ml** of **Stringent Wash solution** (pre-warmed at 46°C).
- Incubate **15 min** at **46° C** in the shaking water bath (about 50 rpm).
- Remove the solution from each lane.

Color development (room temperature)

All the next steps have to be performed at room temperature with shaking. Since water takes too much time to reach room temperature, place a tray (for example of polystyrene) on the lid of the water bath to isolate the Plastic Trays from the hot water and use the agitation function of the water bath only.

- Add **1 ml** of **Conjugate**.
- Incubate **15 minutes** while shaking (about 50 rpm) at room temperature.
- Remove the solution, add **1 ml** of **Wash solution B** and rinse briefly.
- Remove the solution and add **1 ml** of **Wash solution B**.
- Incubate for **5 min** while shaking (about 50 rpm) at room temperature.
- Remove the solution and add **1 ml** of **Wash solution B**.
- Incubate for **5 min** while shaking (about 50 rpm) at room temperature.
- Remove the solution and add **1 ml** of **Color developer**.
- Incubate for **15 min** in the dark, while shaking (about 50 rpm) at room temperature.
- Remove color developer and wash strips with distilled water.
- Let strips dry in the dark on absorbent paper before proceeding with the interpretation of results.

⇒ AUTOMATIC DETECTION

Detection can be performed by using automatic instrument (*Tecan ProfiBlot 48 or equivalent instrument*): select “**Tipizza P**” protocol. In this case **10 µl PCR product + 10 µl DNAT** and **1,5 ml** of each reagent have to be used for each sample.

For more information concerning the protocol load on the instrument to be used contact customercare@nlm.it

RESULTS

Validation


A positive reaction of the uppermost Conjugate Control (CC) line indicates the correct function of the detection step. This line has to be always present and its intensity should be the same among different strips.

The colour intensity among the other positive lines on the strips may vary. This doesn't affect results interpretation.

Interpretation of results

- Align the CC line of the strip with the CC line on the transparent interpretation card. This line should always stain positive.
- The Amplification Control (AC) lines (lines I and II) contain universal probes for conserved sequences of the 5'UTR and Core region respectively. Their positive reaction ensures the addition of the amplified material and indicates the amplification of the corresponding region.
- Read the number of all positive lines for each strip on the Transparent interpretation card: the corresponding HCV typing result can be determined by using the Interpretation Chart (read also the section on diagnostic sensitivity for additional information on genotypes).

Strip scheme:

The part A of the interpretation table shows the possible result combinations on the basis of the 5'UTR region; the part B shows the possible result combinations obtainable on the basis of the Core region.

- If the 5'UTR region of the sample falls in the combinations of genotypes **2, 3, 4, 5, 6a-b or 7a** don't continue with the analysis of the Core region. It is possible to observe positivity on lines **14-18** of the Core, but these **must not** be taken into account in the results interpretation.

- If the 5'UTR region of the sample falls in the combinations of genotypes **1 or 6c-t** the analysis has to be continued with the interpretation of the results on the bands 14-18 (Core) and report only the results associated with the type indicated by the Core region (e.g. if a sample is genotype 1b in the 5'UTR region, but **1a in the Core region, it should be classified as 1a**). If core region information are inconclusive or unavailable, it isn't possible to discriminate between genotype 1a, 1b or 6.

In the presence of genotypes 2, 3, 4, 5 and 7, signals on the probes of the Core region can be absent or present; in this case they should not be interpreted in any way:

Core probes	Genotypes that can show signals
14	5a
15	/
16	2c, 4a, 5a
17	3
18	2, (3), 5a

TROUBLESHOOTING

1. False negative or too weak signals

- Weak signals on the strip could be the consequence of:
 - Too high DNA concentration during the PCR (try to dilute it) especially for samples with a very high titer ($>1 \times 10^7$)
 - Too low DNA concentration during the PCR for samples with a very low titer
- If all the strips show weak signals it could be due to too high hybridization and/or Stringent washing steps water-bath temperature. **Check carefully the water-bath temperature.**
- Hybridization or Stringent Wash Solution were not pre-warmed or mixed carefully. In this case the remaining solution should be discarded and a new one should be used.

2. False positive signals

- Unspecific signals can occur if temperature during hybridization and/or Stringent washing steps is too low. **Check carefully the water-bath temperature.**
- Stringent Wash Solution was not pre-warmed and mixed carefully. In this case the remaining buffer should be wasted and a new one should be used.
- If the same non-specific pattern, corresponding to a given genotype, is visible on all strips, a contamination occurred. Repeat the assay with new solutions. If contamination events occur it is important to discover their origin: extraction, RT and PCR reactions or contamination derived from other amplicons. The use of HCV positive and negative controls can be useful to find out the cause of contamination.

3. Single line non homogeneous staining

Shaking rate during detection steps is too low. Check that strips are submerged totally.

4. Pattern different from regular pattern on the typing list

- Typing probes of 2 different types are present: this may be due to mixed infections or cross contaminations.
- In rare cases, combinations of bands are generated as not interpretable. This may be due to the heterogeneity of the HCV genome, contaminations or recombinant strains.

In both cases we suggest to repeat the test (starting from RNA purification) and eventually contact your distributor.

PERFORMANCE CHARACTERISTICS

⇒ Diagnostic specificity

The specificity of the Gen-C 2.0 assay was determined analysing 101 HCV negative plasma samples. Samples tested were negative for HCV antibody and RNA. The test was performed according to instructions for use. No false positive results were observed, therefore Gen-C 2.0 Test Specificity is 100%.

⇒ Diagnostic sensitivity

The diagnostic sensitivity of the kit Gen-C 2.0 was evaluated using plasma samples from HCV-positive subjects and amplified with kits cod NLM. AA910/48 ver. 2.0.

On the basis of the *in silico* sequence analysis, the device is able to distinguish genotypes 1, 2, 3, 4, 5, 6 and 7 and subtypes 1a, 1b, 2a/2c, 2b, 3a, 3b, 3c, 3k, 4a, 4b, 4c/4d, 4e, 4f, 4h, 5a, 6a/6b, 6g, 6m, 6t and 7a.

The device validation has been carried out with real samples of genotype: 1a, 1b, 2a/2c, 3, 3a, 4, 4c/4d, 5a and 6a. Viral loads were in a range of 2×10^3 and 1×10^7 IU/ml.

The Gen-C 2.0 assay ability to accurately genotype HCV positive samples was evaluated by comparing the results with those based on the previous version of the kit (for genotypes 2-5) and sequencing (for genotype 1 and 6) of individual samples; sequencing results were determined comparing them to HCV database and phylogenetic analysis. For genotypes 5a and 6a, samples from reference panels were used.

Analysis has been carried both with regard to the genotype and subtype (1a and 1b).

Of the 217 samples analyzed, all gave positive and correct results. The overall sensitivity for the assay is therefore 100%.

117/117 genotype 1 samples, were correctly subtyped with AC004 2.0 and were concordant with the reference method. Therefore the accuracy for Gen-C 2.0 assay for subtype 1a and 1b is 100%.

Summary of the analyzed samples and concordance:

Genotype	Nr of samples tested	Concordant to reference method (Gen-C, sequencing or panel)
1a	30	30
1b	87	87
2a/2c	37	37
3a	22	22
3	10	10
4c/4d	21	21
4	2	2
5a	6	6
6a/6b	2	2

Reproducibility


- **INTRA-ASSAY REPRODUCIBILITY:**

Intra-assay reproducibility was evaluated using 3 samples of 3 different genotype in 3 replicates. The same 3 samples were tested with two different batches by two different operators.

- **INTER-ASSAY REPRODUCIBILITY:**

Inter-assay reproducibility was evaluated using 3 samples of different genotype by three different operators in different runs and days, using three different batches.

The test was performed according to the instructions for use, the precision was evaluated for all the steps of the procedure. The reproducibility for Gen-C 2.0 assay is 100%.

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⇒ LITERATUR

1. Sarrazin C. *Diagnosis of hepatitis C: update 2004. Journal of Gastroenterology and Hepatology (2004) 19, S88-S93*
2. Simmonds et al. *Consensus proposals for a unified system of nomenclature of hepatitis C virus genotypes. Hepatology, vol. 42, n° 4, 2005.*
3. Smith et al. *Expanded Classification of Hepatitis C Virus Into 7 Genotypes and 67 Subtypes: Updated Criteria and Genotype Assignment Web Resource HEPATOLOGY, Vol. 59, No. 1, 2014*
4. P.T. Hraber et al. *Comparative analysis of hepatitis C virus phylogenies from coding and non-coding regions: the 5' untranslated region (UTR) fails to classify subtypes. Virology Journal 2006, 3:103*
5. M. A. Ansari et al. *HCV-Core Region: Its Significance in HCV-Genotyping and Type Dependent Genomic Expression. Macedonian Journal of Medical Sciences. 2012 Mar 15; 5(1):30-39.*
6. Nizar N. Zein. *Clinical Significance of Hepatitis C Virus Genotypes. Clinical Microbiology Reviews, Apr. 2000, p.223-235.*
7. *EASL Clinical Practice Guidelines: Recommendation on treatments of hepatitis C. Journal of Hepatology 2016.*
8. Falade-Nwulia et al. *Oral Direct-Acting Agent Therapy for Hepatitis C Virus Infection A Systematic Review Annals of Internal Medicine Vol. 166 No. 8 18 April 2017*

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