

BactoReal® Kit

Riemerella anatipestifer



For veterinary use only

BactoReal® Kit *Riemerella anatipestifer*

Order no.	Reactions	Pathogen	Internal positive control
DVEB05913	100	FAM channel	Cy5 channel
DVEB05953	50	FAM channel	Cy5 channel
DVEB05911	100	FAM channel	VIC/HEX channel
DVEB05951	50	FAM channel	VIC/HEX channel

Kit contents:

- Detection assay for *Riemerella anatipestifer*
- Detection assay for internal positive control (control of amplification)
- DNA reaction mix (contains uracil-N glycosylase, UNG)
- Positive control for *Riemerella anatipestifer*
- Water




Background: *Riemerella anatipestifer* is a Gram-negative bacterium that causes acute or chronic septicaemia. It primarily affects young ducks and less frequently turkeys and geese throughout the world. There are 21 known serotypes and infection is spread horizontally between birds. Infection may be referred to as Duck Septicaemia, Goose 'flu, Riemerellosis, New Duck Disease and Polyserositis.

Description: BactoReal® Kit *Riemerella anatipestifer* is based on the amplification and detection of the 16S-23S ribosomal RNA intergenic spacer of *Riemerella anatipestifer* using real-time PCR. It allows the rapid and sensitive detection of *Riemerella anatipestifer* from DNA samples purified from tissues and swabs (e.g. with the QIAamp DNA Mini Kit).

PCR-platforms: BactoReal® Kit *Riemerella anatipestifer* is developed for the ABI PRISM® 7500 instrument (Life Technologies), LightCycler® 480 (Roche) and Mx3005P® QPCR System (Agilent), but is also suitable for other real-time PCR instruments.

Sensitivity and specificity: BactoReal® Kit *Riemerella anatipestifer* detects at least 20 target copies/PCR reaction, which is equal to seven *R. anatipestifer* colony forming units, CFU (PCR target gene is present three times in the genome of *Riemerella anatipestifer*). The kit is specific for *Riemerella anatipestifer*.



BactoReal®, MycoReal, ParoReal and ViroReal® Kits run with the same thermal cycling conditions. RNA and DNA material can be analysed in one PCR run.

For further information on our products please visit our homepage (www.ingenetix.com)