

# HardyCHROM<sup>™</sup> Salmonella

Cat. no. G309	HardyCHROM <sup>™</sup> Salmonella, 15x100mm Plate, 18ml	10 plates/bag
Cat. no. J37	HardyCHROM <sup>™</sup> Salmonella / XLT-4 Agar, 15x100mm Biplate, 10ml/10ml	10 plates/bag

# **INTENDED USE**

HardyCHROM<sup>TM</sup> Salmonella is a chromogenic medium recommended for the selective isolation and differentiation of *Salmonella* spp. from other members of the family Enterobacteriaceae based on colony color.

## **SUMMARY**

HardyCHROM<sup>TM</sup> Salmonella facilitates the isolation and differentiation of *Salmonella* spp. from other members of the family Enterobacteriaceae. Peptones in the medium supply the necessary nutrients. Selective agents inhibit the growth of gram-positive organisms. Artificial substrates (chromogens) are broken down by specific microbial enzymes which release insoluble colored compounds. *Salmonella* species break down only one of the chromogens and will produce deep pink to magenta colored colonies. Bacteria other than *Salmonella* spp. may break down the other chromogenic substrates and produce blue colonies. If none of the substrates are utilized, natural or white colored colonies will be present.

## **FORMULA**

Ingredients per liter of deionized water:\*

Peptones	10.0gm
Yeast Extract	3.0gm
Sodium Deoxycholate	1.0gm
Chromogenic Mixture	0.5gm
Selective Agents	0.02gm
Agar	15.0gm

Final pH 7.2 +/- 0.2 at 25 degrees C.

\* Adjusted and/or supplemented as required to meet performance criteria.

# STORAGE AND SHELF LIFE

Storage: Upon receipt store at 2-8 degrees C. away from direct light. Media should not be used if there are any signs of deterioration (shrinking, cracking, or discoloration), contamination, or if the expiration date has passed. Product is light and temperature sensitive; protect from light, excessive heat, moisture, and freezing.

The expiration date applies to the product in its intact packaging when stored as directed.

This product has the following shelf life from the date of manufacture:

60 Days:	G309	HardyCHROM <sup>™</sup> Salmonella
J37		HardyCHROM™ Salmonella / XLT-4 Agar

Refer to the document "Storage" on the Hardy Diagnostics Technical Document website for more information.

# PRECAUTIONS

This product may contain components of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not guarantee the absence of transmissible pathogenic agents. Therefore, it is recommended that these products be treated as potentially infectious and handled observing the usual universal blood precautions. Do not ingest, inhale, or allow to come into contact with skin.

This product is for *in vitro* diagnostic use only and is to be used only by adequately trained and qualified laboratory personnel. Observe approved biohazard precautions and aseptic techniques. All laboratory specimens should be considered infectious and handled according to "standard precautions". The "Guideline for Isolation Precautions" is available from the Centers for Disease Control and Prevention at <u>www.cdc.gov/ncidod/dhqp/gl\_isolation.html</u>.

For additional information regarding specific precautions for the prevention of the transmission of all infectious agents from laboratory instruments and materials, and for recommendations for the management of exposure to infectious disease, refer to CLSI document M-29: *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline.* 

Sterilize all biohazard waste before disposal.

Refer to the document "<u>Precautions When Using Media</u>" on the Hardy Diagnostics <u>Technical Document</u> website for more information.

Refer to the SDS Search instructions on the Hardy Diagnostics website for more information.

# PROCEDURE

Specimen Collection: Consult listed references for information on specimen collection. <sup>(2-4,7)</sup> Infectious material should be submitted directly to the laboratory without delay and protected from excessive heat and cold. If there is to be a delay in processing, the specimen should be inoculated onto an appropriate transport media and refrigerated until inoculation.

Consult the listed references for information regarding the processing of specimens.<sup>(2-4,7)</sup>

Method of Use: Allow the plates to warm to room temperature. The agar surface should be dry prior to inoculating. Inoculate and streak the specimen as soon as possible after collection. If the specimen to be cultured is on a swab, roll the swab over a small area of the agar surface. Streak for isolation with a sterile loop. Incubate plates in an inverted position, protected from the light, aerobically at 35-37 degrees C. for 24 hours. Examine plates for colonies showing typical morphology and color.

#### INTERPRETATION OF RESULTS

*Salmonella* spp., including *S. typhi* and *S. paratyphi* A, produce magenta colored colonies. Other members of the Enterobacteriaceae (if present) produce blue, blue-green, white, or colorless colonies. Gram-positive bacteria and non-glucose fermenting bacteria will be inhibited.

Organism	Description	Photo	Color
Salmonella spp.	pink to magenta colonies		
E. coli	blue colonies		

# LIMITATIONS

It is necessary to perform biochemical, immunological, molecular, or mass spectrometry testing on colonies from pure culture for complete identification.

Color-blind individuals may encounter difficulty in distinguishing the color differences on HardyCHROM<sup>TM</sup> Salmonella.

Refer to the document "<u>Limitations of Procedures and Warranty</u>" on the Hardy Diagnostics <u>Technical Document</u> website for more information.

# MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as loops, swabs, applicator sticks, other culture media, incinerators, and incubators, etc., as well as serological and biochemical reagents, are not provided.

# **QUALITY CONTROL**

Test Organisms	Inoculation	Incubation			Dogulta
Test Organishis	Method*	Time	Temperature	Atmosphere	Kesuits
Salmonella enterica ATCC <sup>®</sup> 14028	А	24hr	35°C	Aerobic	Growth; pink to magenta colonies
<i>Escherichia coli</i> ATCC <sup>®</sup> 25922	В	24hr	35°C	Aerobic	Variable growth; blue colonies
Proteus mirabilis ATCC <sup>®</sup> 12453	А	24hr	35°C	Aerobic	Partial to complete inhibition; colorless colonies
Pseudomonas aeruginosa ATCC <sup>®</sup> 27853	В	24hr	35°C	Aerobic	Inhibited

The following organisms are routinely used for testing at Hardy Diagnostics:

\* Refer to the document "<u>Inoculation Procedures for Media QC</u>" on the Hardy Diagnostics <u>Technical Document</u> website for more information.

# **USER QUALITY CONTROL**

Check for signs of contamination and deterioration. Users of commercially prepared culture media may be required to perform quality control testing to demonstrate growth or a positive reaction and to demonstrate inhibition or a negative reaction (where applicable). Refer to the following documents on the Hardy Diagnostics <u>Technical Document</u> website for more information on QC: "<u>Introduction to Quality Control</u>," "<u>Finished Product Quality Control Procedures</u>," or "<u>The CLSI Standard and Recommendations for User Quality Control of Media</u>." Also see the following reference for more specific information. <sup>(1-7)</sup>

# PHYSICAL APPEARANCE

HardyCHROM<sup>TM</sup> Salmonella should appear translucent light beige with a fine white precipitate dispersed evenly throughout.



Salmonella enterica (ATCC<sup>®</sup> 14028) colonies growing on HardyCHROM<sup>™</sup> Salmonella (Cat. no. G309). Incubated aerobically for 24 hours at 35 deg. C.



*Escherichia coli* (ATCC<sup>®</sup> 25922) colonies growing on HardyCHROM<sup>TM</sup> Salmonella (Cat. no. G309). Incubated aerobically for 24 hours at 35 deg. C.



Uninoculated plate of HardyCHROM<sup>TM</sup> Salmonella (Cat. no. G309).

## REFERENCES

1. Wallace, J.S., and K. Jones. 1996. The Use of Selective and Differential Agars in the Isolation of *Escherichia coli* O157 from Dairy Herds. *Jour. Appl. Bacteriol.*; 81:663-668.

2. Versalovic, J., et al. Manual of Clinical Microbiology. American Society for Microbiology, Washington, D.C.

3. Tille, P.M., et al. Bailey and Scott's Diagnostic Microbiology, C.V. Mosby Company, St. Louis, MO.

4. Isenberg, H.D. *Clinical Microbiology Procedures Handbook*, Vol. I, II & III. American Society for Microbiology, Washington, D.C.

5. Koneman, E.W., et al. *Color Atlas and Textbook of Diagnostic Microbiology*. J.B. Lippincott Company, Philadelphia, PA.

6. American Public Health Association. *Standard Methods for the Examination of Dairy Products*, APHA, Washington, D.C.

7. APHA Technical Committee on Microbiological Methods for Foods. *Compendium of Methods for the Microbiological Examination of Foods*, APHA, Washington, D.C.

8. U.S. Food and Drug Administration. *Bacteriological Analytical Manual*. Arlington, VA <u>http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm2006949.htm</u>

ATCC is a registered trademark of the American Type Culture Collection.



1430 West McCoy Lane, Santa Maria, CA 93455, USA Phone: (805) 346-2766 ext. 5658 Fax: (805) 346-2760 Website: www.HardyDiagnostics.com Email: TechService@HardyDiagnostics.com

Distribution Centers: California · Washington · Utah · Arizona · Texas · Ohio · Florida · New York · North Carolina

The Hardy Diagnostics manufacturing facility and quality management system is certified to ISO 13485.

Copyright© 1996 - 2015 by Hardy Diagnostics. All rights reserved.