



## Instructions for Use

### CRITERION™ TRYPTIC SOY BROTH (TSB)

<a href="#">Cat. no. C7140</a>	CRITERION™ Tryptic Soy Broth	60gm
<a href="#">Cat. no. C7141</a>	CRITERION™ Tryptic Soy Broth	500gm
<a href="#">Cat. no. C7142</a>	CRITERION™ Tryptic Soy Broth	2kg
<a href="#">Cat. no. C7143</a>	CRITERION™ Tryptic Soy Broth	10kg
<a href="#">Cat. no. C7144</a>	CRITERION™ Tryptic Soy Broth	50kg

### INTENDED USE

Hardy Diagnostics CRITERION™ Tryptic Soy Broth is recommended for use as a general purpose medium for the isolation and cultivation of a wide variety of microorganisms.

This dehydrated culture medium is a raw material intended to be used in the making of prepared media products, which will require further processing, additional ingredients, or supplements.

### SUMMARY

Tryptic Soy Broth is widely used for the cultivation of microorganisms from environmental sources, supporting the growth of the majority of bacteria and fungi. Tubes of this medium may be used for preparing dilutions of organism for colony counts and preparation of standard inocula. Tryptic Soy Broth is also recommended for use in sterility testing for the detection of contamination with low incidence fungi and aerobic bacteria.<sup>(1)</sup>

Tryptic Soy Broth, also known as Soybean-Casein Digest, conforms to the formula given by the U.S. Pharmacopeia.<sup>(1)</sup> This medium contains digests of soybean meal and casein, which provide amino acids and other nitrogenous substances, making it a highly nutritious medium for a variety of organisms. Sodium chloride is added to maintain the osmotic equilibrium. Dextrose is incorporated as an energy source. Dipotassium phosphate is included in the formulation as a buffer to maintain the pH. The medium can also be supplemented with Tween® 80 or other non-inhibitory surface-active agent to assist in the suspension of water insoluble substances for the microbiological examination of non-sterile products.<sup>(1)</sup>

#### FORMULA\*

Gram weight per liter:	30.0gm/L
Pancreatic Digest of Casein	17.0gm
Sodium Chloride	5.0gm
Papaic Digest of Soybean Meal	3.0gm
Dextrose	2.5gm

Dipotassium Phosphate

2.5gm

Final pH 7.3 +/- 0.2 at 25°C.

\* Adjusted and/or supplemented as required to meet performance criteria.

## STORAGE AND SHELF LIFE

Store the sealed bottle(s) containing dehydrated culture medium at 2-30°C. Dehydrated culture medium is very hygroscopic. Keep lid tightly sealed. Protect dehydrated culture media from moisture and light. The dehydrated culture media should be discarded if it is not free-flowing or if the color has changed from its original light beige.

Store the prepared culture media at 2-30°C.

The expiration dating on the product label applies to the product in its intact packaging when stored as directed. The product may be used and tested up to the expiration date on the product label and incubated for the recommended quality control incubation times.

Refer to the document "[Storage](#)" for more information.

## PRECAUTIONS

This product may contain components of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not guarantee the absence of transmissible pathogenic agents. Therefore, it is recommended that these products be treated as potentially infectious, and handle observing the usual universal blood precautions. Do not ingest, inhale, or allow to come into contact with skin.

This product is for laboratory use only. It is to be used only by adequately trained and qualified laboratory personnel. Observe approved biohazard precautions and aseptic techniques. All laboratory specimens should be considered infectious and handled according to "standard precautions." The "Guidelines for Isolation Precautions" is available from the Centers for Disease Control and Prevention at [www.cdc.gov/ncidod/dhqp/gl\\_isolation.html](http://www.cdc.gov/ncidod/dhqp/gl_isolation.html).

For additional information regarding specific precautions for the prevention of the transmission of all infectious agents from laboratory instruments and materials, and for recommendations for the management of exposure to infectious disease, refer to CLSI document M-29: *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline*.

Sterilize all biohazard waste before disposal.

Refer to the document "[Precautions When Using Media](#)" for more information.

Refer to the document [SDS Search](#) instructions on the Hardy Diagnostics' website for more information.

## METHOD OF PREPARATION FOR DEHYDRATED CULTURE MEDIA

1. Suspend 30.0gm of the dehydrated culture media in 1 liter of distilled or deionized water.
2. Heat to boiling and mix to dissolve completely.
3. If necessary, add 1.0gm/L or appropriate concentration of sterile Tween<sup>®</sup> 80 (Cat. no. U180) or other noninhibitory surface-active agent to assist in the suspension of poorly wettable substances.<sup>(1)</sup>
4. Sterilize in the autoclave at 121°C. for 15 minutes.

## PROCEDURE AND INTERPRETATION OF RESULTS

For information on procedures and interpretation of results, consult listed references or refer to the prepared media

Instructions for Use (IFU) for Cat. no. K89 or K82 (for USP).

## LIMITATIONS

It is recommended that biochemical, immunological, molecular, or mass spectrometry testing be performed on colonies from pure culture for complete identification.

Some formulations may require a settling period before pH testing of the prepared medium. If the pH is tested immediately after preparation and is out of specification, retest the medium after 24 hours to obtain final pH results.

Refer to the document "[Limitations of Procedures and Warranty](#)" for more information.

## MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as autoclaves, incinerators, and incubators, etc., are not provided.

## QUALITY CONTROL

Hardy Diagnostics tests each lot of commercially manufactured media using appropriate quality control microorganisms and quality specifications as outlined on the Certificates of Analysis (CofA). The following organisms are routinely used for testing at Hardy Diagnostics:

Test Organisms	Inoculation Method*	Incubation			Results
		Time	Temperature	Atmosphere	
<i>Candida albicans</i> ** ATCC® 10231	J	1-3 days	30-35°C	Aerobic	Growth
<i>Bacillus subtilis</i> ** ATCC® 6633	J	1-3 days	35°C	Aerobic	Growth
<i>Staphylococcus aureus</i> ATCC® 25923	A	24-48hr	35°C	Aerobic	Growth
<i>Staphylococcus aureus</i> ** ATCC® 6538	J	1-3 days	35°C	Aerobic	Growth
<i>Pseudomonas aeruginosa</i> ** ATCC® 9027	J	1-3 days	35°C	Aerobic	Growth
<i>Escherichia coli</i> ATCC® 25922	A	24-48hr	35°C	Aerobic	Growth
<i>Escherichia coli</i> ** ATCC® 8739	J	1-3 days	35°C	Aerobic	Growth
<i>Salmonella enterica</i> ** ATCC® 14028	J	1-3 days	35°C	Aerobic	Growth

\* Refer to the document "[Inoculation Procedures for Media QC](#)" for more information.

## USER QUALITY CONTROL

Users of dehydrated culture media should perform QC testing in accordance with applicable government regulatory agencies, and in compliance with accreditation requirements. Hardy Diagnostics recommends end users check for

signs of contamination and deterioration and, if dictated by laboratory quality control procedures or regulation, perform quality control testing to demonstrate growth or a positive reaction and to demonstrate inhibition or a negative reaction, if applicable. Hardy Diagnostics quality control testing is documented on the certificates of analysis (CofA) available from Hardy Diagnostics [Certificates of Analysis](#) website. In addition, refer to the following document "[Finished Product Quality Control Procedures](#)," for more information on QC or see reference(s) for more specific information.

## PHYSICAL APPEARANCE

CRITERION™ Tryptic Soy Broth powder should appear homogeneous, free-flowing, and light beige in color. The prepared media should appear translucent, and light amber in color.

## REFERENCES

1. *United States Pharmacopoeia and National Formulary* (USP-NF). Rockville, MD: United States Pharmacopoeial Convention.
2. U.S. Food and Drug Administration. *Bacteriological Analytical Manual*. AOAC, Arlington, VA. <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm2006949.htm>.
3. American Public Health Association. *Standard Methods for the Examination of Dairy Products*, APHA, Washington, D.C.
4. APHA Technical Committee on Microbiological Methods for Foods. *Compendium of Methods for the Microbiological Examination of Foods*, APHA, Washington, D.C.
5. American Public Health Association. *Standard Methods for the Examination of Water and Wastewater*, APHA, Washington, D.C.
6. MacFaddin, J.F. 1985. *Media for Isolation, Cultivation, Identification, Maintenance of Bacteria*, Vol. I. Williams & Wilkins, Baltimore, MD.
7. *Quality Assurance for Commercially Prepared Microbiological Culture Media*, M22. Clinical and Laboratory Standards Institute (CLSI - formerly NCCLS), Wayne, PA.
8. Jorgensen., et al. *Manual of Clinical Microbiology*, American Society for Microbiology, Washington, D.C.

ATCC is a registered trademark of the American Type Culture Collection.

Tween is a registered trademark of ICI Americas, Inc.

IFU-10273[A]



1430 West McCoy Lane, Santa Maria, CA 93455, USA

Phone: (805) 346-2766 ext. 5658

Fax: (805) 346-2760

Website: [www.HardyDiagnostics.com](http://www.HardyDiagnostics.com)

Email: [TechService@HardyDiagnostics.com](mailto:TechService@HardyDiagnostics.com)

[Ordering Information](#)

Distribution Centers:

California · Washington · Utah · Arizona · Texas · Ohio · New York · Florida · North Carolina

The Hardy Diagnostics manufacturing facility and quality management system is certified to ISO 13485.

Copyright© 1996 by Hardy Diagnostics. All rights reserved.

HDQA 2207B [C]